



Cosa Active


Cosa Junior


DE Gebrauchsanweisung	6
EN Instructions for use	10

	Artikelnummer	Konfektionsgröße	Hüftumfang (cm)	Polsternummer	Oberschenkelumfang (cm)
	28L100=80	74–86	57–59	1*	26–32
	28L100=92	86–98	59–62		
	28L100=104	98–110	62–64	2*	32–36
	28L100=116	110–122	64–68		
	28L101=128	122–134	68–72	3*	36–44
	28L101=140	134–146	72–78		
	28L101=152	146–158	78–86	4*	44–54
	28L101=S	158–170	86–96		
		♀ 36–38 ♂ 44–46			
28L101=M	♀ 40–42 ♂ 48–50	96–105	5*	54–62	
28L101=L	♀ 44–46 ♂ 52–54	105–109			

Material	PA, EA, PU Schaum / PU foam
-----------------	-----------------------------

INFORMATION

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INFORMATION

Last update: 2014-12-08

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the Cosa active and Cosa junior abduction pants.

2 Product Description

The orthosis is supplied as a starter set in two versions:

- 28L100 Cosa junior (pants without zip)
- 28L101 Cosa active (pants with zip)

Each starter set contains:

- 2 pants (including strap set)
- 1 abduction pad set (including instructions for use)

Article	Article no. Pants	Quantity	Article no. Abduction pad set	Quantity	Article no. Strap set	Quantity
Cosa junior	29L100=80	2	29L102=1	1	33L110=1	2
	29L100=92		29L102=2			
	29L100=104					
	29L100=116					
Cosa active	29L101=128	2	29L102=3	1	33L110=2	2
	29L101=140		29L102=4			
	29L101=152					
	29L101=S					
	29L101=M					
29L101=L	29L102=5					

3 Application

3.1 Indications for use

The orthosis is intended exclusively for orthotic fittings of the hip and exclusively for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Diplegia, spastic
- Hip dysplasia
- Hip joint subluxation, congenital
- Hypotonicity
- Tetraplegia, spastic

3.3 Contraindications

3.3.1 Absolute Contraindications

- Hip joint luxation

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

Consulting a doctor is also required in case of contractures and shortening in the area of the lumbar spine and the lower limbs.

3.4 Effects

In case of neuromuscular weakness, the orthosis prevents excessive adduction while lying down, sitting, standing and walking through the use of pads between the legs. The orthosis enables a firm basis for sitting and standing, thereby improving the torso position and leading to better balance control.

In patients who are able to walk, the orthosis improves pelvic symmetry and increases knee freedom. The scissor gait is reduced and more physiological walking is promoted.

By using the orthosis, patients who are unable to walk do not have to use separate abductors in the wheelchair. Hip abduction with the orthosis improves balance while sitting by establishing a broader seat basis. Using the hands for support is no longer necessary.

When lying down, e.g. while sleeping, the orthosis reduces muscle shortening due to spasticity.

4 Safety

4.1 Explanation of Warning Symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 General Safety Instructions

CAUTION

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.

CAUTION

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Size Selection

Selection according to clothes size

- ▶ Select the orthosis size based on the clothing size (see sizing table).

Selection according to circumferences

- 1) Measure the hip circumference.
- 2) Determine the size of the pants (see size chart).
- 3) Measure the thigh circumference at the height of the pads.
- 4) Determine the size of the pads (see size chart).

5.2 Adaptation

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

Adapting the pads

- 1) Apply the pads to the thigh so that the guide loops face in the proximal direction (see fig. 1).
- 2) Determine the required abduction angle.
- 3) **If necessary:** Use sandpaper to roughen the surfaces of the pads and glue on additional layers (see fig. 2).
- 4) Grind the surfaces of the pads (see fig. 3).

NOTICE! Do not grind the pad surfaces further than the mark.

Adapting the belt length

- 1) Trial fit the pants on the patient.
- 2) Mark the required length of the belt ends.
- 3) Mark the position of the guide loop (see fig. 4) and the hook-and-loop closure (see fig. 5) for sewing.
- 4) Sew the guide loop and hook-and-loop closure for the belt (see fig. 6).
- 5) **If necessary:** Adapt the fit of the pants.

Adjusting the pads

- 1) Insert the pads into the pad pockets on the pants so that the guide loops face in the proximal direction (see fig. 7).
 - 2) Pull the guide loops of the pads out through the openings in the pants (see fig. 8).
 - 3) Close the hook-and-loop closures on the pockets.
 - 4) Attach the cross straps on the guide loops of the pads (see fig. 9).
 - 5) Trial fit the pants on the patient (see fig. 10).
 - 6) Determine the position of the cross straps and temporarily secure the cross straps with double hook-and-loop (see fig. 11).
 - 7) Check the tension of the cross straps (see fig. 13).
 - 8) Mark the position of the cross straps for later sewing (see fig. 12).
- INFORMATION: Make sure that the routing of the straps is symmetrical in case of symmetrical hips (see fig. 14).**
- 9) Sew the cross straps (see fig. 15).

Functional test

- 1) Verify the fit and functionality of the orthosis.
- 2) **If necessary:** Fine-tune the pads using the straps.

5.3 Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.

⚠ CAUTION

Use on skin sensitive to heat

Skin irritation due to overheating

- ▶ Do not use the product in case of a known heat allergy.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

- > The belt is open.
- > **Cosa active only:** The zip is open.
- 1) Put on the pants.
- 2) Fasten the belt and close the zip if applicable.
- 3) Check the positioning of the pads and readjust them with the cross straps as needed.

5.4 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

- 1) Loosen the cross straps.
- 2) Open the hook-and-loop closures of the pockets.
- 3) Remove the pads.
- 4) Fasten all hook-and-loop closures.
- 5) **Recommendation:** Use a laundry bag or net.
- 6) Wash the pants in warm water at **40 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 7) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 8) Put the pads back in the pants.

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria out-

lined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt · Germany
T +49 5527 848-0 · F +49 5527 72330
healthcare@ottobock.de · www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.