











## 28U24 WalkOn Reaction

## 28U25 WalkOn Reaction Plus


<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	5
<b>EN</b> Instructions for use (qualified personnel) .....	10
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	15
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	20
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	25
<b>PT</b> Manual de utilização (Pessoal técnico) .....	30
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	35
<b>SV</b> Bruksanvisning (Fackpersonal) .....	40
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	45
<b>NO</b> Bruksanvisning (Fagpersonell) .....	49
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	54
<b>PL</b> Instrukcja użytkowania (Personel fachowy) .....	59
<b>HU</b> Használati utasítás (szakszemélyzet) .....	64
<b>CS</b> Návod k použití (Odborný personál) .....	69
<b>RO</b> Instrucțiuni de utilizare (Personal de specialitate) .....	74
<b>HR</b> Upute za uporabu (Stručno osoblje) .....	79
<b>SK</b> Návod na používanie (Odborný personál) .....	83
<b>BG</b> Инструкция за употреба (Специалисти) .....	88
<b>TR</b> Kullanma talimatı (Uzman personel) .....	93
<b>RU</b> Руководство по применению (Квалифицированный персонал) .....	98
<b>JA</b> 取扱説明書 (有資格担当者) .....	104
<b>ZH</b> 使用说明书 (专业人员) .....	108

	Artikelnummer / Article number	Seite / Side	Schuhgröße / Shoe size	Höhe / Height	
				cm	inches
	28U24	L/R	36–39	35 cm	13.8 inches
	28U24	L/R	39–42	36 cm	14.2 inches
	28U24	L/R	42–45	37,5 cm	14.8 inches
	28U24	L/R	45–48	39 cm	15.4 inches
	28U25	L/R	36–39	35 cm	13.8 inches
	28U25	L/R	39–42	36 cm	14.2 inches
	28U25	L/R	42–45	37,5 cm	14.8 inches
	28U25	L/R	45–48	39 cm	15.4 inches

<b>Material</b>	Carbonfaserverbundstoff / Carbon fibre composite Polstermaterial / Padding material: Outlast®, Mikrokleb (weiches Klettband) / Micro hook and loop (soft velcro tape)
-----------------	--

## INFORMATION

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### INFORMATION

Last update: 2014-12-08

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 28U24 WalkOn Reaction and 28U25 WalkOn Reaction Plus ankle foot orthoses.

## 2 Application

### 2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

### 2.2 Indications

- Dorsal flexor weakness, with or without mild spasticity (e. g. post CVA, traumatic brain injury, multiple sclerosis, neuromuscular atrophy, or peroneal paralysis)
- Minor impairment of plantarflexors
- Impairment of plantarflexors (only 28U25)
- Flexible foot deformities that can be corrected with the use of a foot orthotic or use of the supination correction strap
- Minor impairment of knee extensors
- Impairment of knee extensors (only 28U25)

Indications must be determined by the physician.

### 2.3 Contraindications

#### 2.3.1 Absolute Contraindications

- Moderate to severe spasticity in the lower leg
- Moderate to severe oedema
- Leg ulcers
- Moderate to severe foot deformities

#### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which

the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

## 2.4 Further Usage Restrictions

The orthosis must not be used for the prevention of contractures (shortening) etc.

The orthosis must not be used for sports activities that involve jumping, sudden movements or rapid step sequences (e.g. basketball, badminton, riding sports). Sports activities should generally be discussed with the patient.

## 2.5 Effects

The orthosis enables a more natural gait pattern for the user. In the swing phase, the foot remains in the neutral position which permits a natural swing-through for the patient without executing compensating movements for ground clearance.

At heel strike, the carbon spring allows the foot to plantarflex in a controlled and harmonious manner preventing foot slap. The force exhibited is used to support the mid-stance phase.

Prior to toe-off, stored energy is released to preserve momentum into swing phase.

## 3 Safety

### 3.1 Explanation of Warning Symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General Safety Instructions

 **CAUTION**

#### Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.

 **CAUTION**

#### Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.

 **CAUTION**

#### Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

 **NOTICE**

#### Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.

## 4 Handling

### **⚠ WARNING**

#### **Operation of motor vehicles**

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

### **INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

### 4.2 Adaptation

### **⚠ CAUTION**

#### **Improper moulding or application**

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

### **⚠ CAUTION**

#### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.

### **NOTICE**

#### **Use without suitable footwear**

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.

#### **Footwear selection:**

- ▶ Choose a sturdy lace-up shoe with solid heel cup to ensure the optimum effectiveness of the orthosis. The effective heel height should be between 0.5 cm and 1.5 cm.

#### **Adapting the orthosis size:**

- 1) **As needed:** If the patient uses footwear with removable insoles, take the insole out of the shoe on the fitting side.

**INFORMATION: Keep the insole in order to use the shoe without the orthosis.**

- 2) Have the patient sit on a chair with the shoe and without the orthosis, and position the foot at a 90° angle to the lower leg (see fig. 1).
- 3) Mark the position of the malleoli in the shoe with a small piece of adhesive tape (see fig. 1).

- 4) Have the patient sit on a chair without the shoe and with the orthosis, and position the foot at a 90° angle to the lower leg.  
**INFORMATION: Apply the orthosis to the patient so that the ankle joint has sufficient freedom of movement and does not impact the hard medial spring.**
- 5) Determine the position of the ankle according to the "WalkOn® Reaction" lettering (see fig. 2).
- 6) Trim the orthosis sole (see fig. 3) so that the mark on the shoe corresponds to the position on the orthosis (see fig. 4).
- 7) Adapt the sole shape according to the removable insole. Note that the required exterior rotation of the foot must be maintained (see fig. 3).

**Adapting the orthosis:**

- 1) Insert the orthosis in the shoe and cover it with an insole if required (see fig. 5).
- 2) Put on the orthosis together with the shoe (see fig. 6).
- 3) Measure the correct AP position of the knee joint (e.g. on the L.A.S.A.R Posture).  
**INFORMATION: The compromised knee rotation point is 15 mm posterior to the weight line according to Nietert.**
- 4) **As needed:** Wedges can be used to provide a comfortable base of support. For example, an extrinsic forefoot wedge can be used to increase knee extension force (see fig. 10). An extrinsic heel wedge can be used to reduce knee extension force (see fig. 11).  
**INFORMATION: Materials other than those supplied may also be used to build up the outer sole if needed.**
- 5) **As needed:** Sand the top edge and wings of the anterior stop and modify the padding accordingly.  
**NOTICE! The functionality of the orthosis must not be impaired. Do not sand in the area of the medial spring or the surfaces!**
- 6) **As needed:** If the patient has a foot deformity, correct it with an insole or a custom moulded foot orthotic.  
**INFORMATION: Included with the device is the supination correction strap (28Z10), which, if needed, can be attached directly to the device to help control lateral movement of the ankle.**
- 7) Close the textile component at calf height (see fig. 6).

**Final inspection:**

- ▶ Conduct trial walking with the patient. Testing on inclines, ramps and stairs is mandatory. Fine tune as required.

**4.3 Application**

**NOTICE**

**Use of a worn or damaged product**

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- 1) Take the insole out of the shoe on the fitting side.
- 2) Open the calf band.
- 3) Slide the orthosis into the shoe.
- 4) Put on the orthosis together with the shoe.
- 5) Close the calf band.

## 4.4 Removal

- 1) Open the calf band.
- 2) Remove the shoe with the orthosis.
- 3) **As needed:** Remove the orthosis from the shoe and put the original insole back in.

## 4.5 Cleaning

### NOTICE

#### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

#### Textile component:

- Remove the textile component from the orthosis.
- Close the hook-and-loop closures.
- Wash the textile component in warm water at 40 °C/104 °F with a standard mild detergent. Rinse thoroughly.
- Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

#### Plastic component:

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 5 Disposal

Return the product to the specialist dealer for disposal.

All components of the product must be disposed of properly in accordance with the respective national environmental regulations.

## 6 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

### 6.3 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

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Otto Bock HealthCare GmbH