## 12K42 ErgoArm plus, 12K44 ErgoArm Hybrid plus, 12K50 ErgoArm Electronic plus

<b>DE</b> Gebrauchsanweisung .....	2
<b>EN</b> Instructions for Use .....	26
<b>FR</b> Instructions d'utilisation.....	50
<b>IT</b> Istruzioni per l'uso .....	74
<b>ES</b> Instrucciones de uso.....	98
<b>SV</b> Bruksanvisning.....	122
<b>NL</b> Gebruiksaanwijzing .....	146
<b>RU</b> Руководство по применению .....	170

DE

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EN

Date of the last update: 2015-01-16

- Please read this document carefully.
- Follow the safety instructions and the precautions specified in this accompanying document.
- If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

Contents	Page
General .....	26
Scope of delivery .....	29
Laminating the humeral segment .....	29
Sealing with 633F30 special grease and 633F11 Silicone Grease .....	34
Modifying and servicing the forearm .....	36
Removing the elbow joint .....	36
Removing the AFB.....	37
Shortening the forearm, positioning the battery box.....	38
Inserting the AFB.....	39
Changing the right-to-left cable orientation .....	39
Shortening the control cable .....	41
Inserting the elbow joint .....	41
Adjustment of preflexion.....	43
Control harness cable attachment.....	45
Function of the AFB.....	46
AFB compensation adjustment .....	46
Functional range.....	47
Humeral rotation feature.....	47
12K50 - ErgoArm Electronic Plus.....	48
Annex .....	49

 Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

Only put the product into use in accordance with the information contained in the accompanying documents supplied. Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

## 1 Intended use

### 1.1 Medical purpose

The ErgoArm 12K42/12K44/12K50 models of elbow-forearm assembly are **solely** intended for use in upper limb prostheses.

### 1.2 Application

The ErgoArm 12K42/12K44/12K50 models of elbow-forearm assembly can be used for unilateral or bilateral transfemoral amputees. The ErgoArm models are **solely** intended for adults.

### 1.3 Application conditions

The ErgoArm 12K42/12K44/12K50 models of elbow-forearm assembly have been developed for everyday life activities and

must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

### 1.4 Prosthetist qualification

The fitting of patients with the ErgoArm 12K42/12K44/12K50 models of elbow-forearm assembly may only be carried out by an authorized prosthetist.

## 2 Definition of terms

### AFB (Automatic Forearm Balance):

**12K42**  
ErgoArm Plus

**12K44**  
ErgoArm Hybrid Plus

**12K50**  
ErgoArm Electronic Plus


When the prosthetic arm is extended, the AFB stores the energy generated to make flexion easier later on. The complex mechanism of the AFB controls the dynamic movement of the prosthetic forearm, counterbalancing the weight of the forearm, wrist, and terminal device via compensation gears. The force that the patient must generate to flex the elbow is dramatically reduced. The amount of compensation can be easily adjusted by the patient whenever necessary according to the prosthetic forearm weight and clothing type.

### Easy plug:

**12K44**  
ErgoArm Hybrid Plus

**12K50**  
ErgoArm Electronic Plus

1 Battery connection

2 CLOSE electrode

3 OPEN electrode

4 Switch

5 Electronics housing

The Easy Plug feature makes it simple to wire externally powered hybrid prostheses. The electrode, switch, and battery cables are plugged into the internal elbow cover. Since there are no external cables, this enhances both the durability and cosmetic appearance of the ErgoArm.

## Forearm connection cables

**12K44**

ErgoArm Hybrid Plus

**12K50**

ErgoArm Electronic Plus

Forearm cable to clip onto the AFB. Establishes the electrical connection between the elbow joint and system electric hand or system electric gripper.

**Notice!**

Due to different wall thicknesses of the ErgoArm socket, the cut-out must be beveled to allow insertion of the battery box (see Fig.).

**Locking and slip-stop:**

The **internal infinite-position lock** can be unlocked or locked under load in any position. When locked, the ErgoArm can be loaded up to 230 N with a forearm up to 305 mm in length.

**Mechanical:**

The integrated compensation mechanism makes it possible to lock and unlock the elbow while under load by generating a constant low force of 10 N. The **integrated strain-relief device** reduces the risk of possible damage that could result from pulling the cable too aggressively.

A strong pull on the cable (clicking sound) locks, or respectively, unlocks the joint permanently.

A gentle pull on the cable unlocks the joint (without click). After relieving the pull, the joint is locked again (**Slip-Stop Function**).

**Electronic:****12K50**

ErgoArm Electronic Plus

The locking and unlocking of the elbow is controlled by myoelectric signals via electrodes or by a switch depending on the selected switch position (see Attachment).

In some switching modes, the joint can be unlocked by electrode signals so as to lock it again instantly when the signal is turned off (**Slip-Stop Function**).

**3 Description and function of the individual ErgoArm models**

The ErgoArm models have following different features:

	AFB	Easy Plug	Lock	Slip-Stop
<b>12K42</b> ErgoArm plus	✓		mechanical	mechanical
<b>12K44</b> ErgoArm Hybrid plus	✓	✓	mechanical	mechanical
<b>12K50</b> ErgoArm Electronic plus	✓	✓	electronic and mechanical	electronic and mechanical

**12K42**  
ErgoArm plus

Mechanical elbow with **Automatic forearm balance, cable-operated mechanical lock, and Slip-Stop Function**.

**12K44**  
ErgoArm Hybrid plus

Along with the 12K42 ErgoArm's other functions, the Easy Plug internally integrated wiring harness has also been provided.

The joint is therefore particularly well-suited for use in hybrid prostheses.

**12K50**  
ErgoArm Electronic plus

An additional feature of the ErgoArm Electronic Plus 12K50, the **electronic lock**, enables the user to control the elbow's unlocking and locking with myoelectric signals, which are picked up by the electrodes.

To guarantee the optimal adaptation to the patient, eight switching modes are available. The appropriate switching mode should be set by the orthopedic technician by selecting the proper color coding plug.

**Service:**

Because all movable parts are prone to wear, it is necessary to perform yearly service. This service includes a complete inspection and possible adjustment of the ErgoArm by Ottobock Myo-Service. Additionally, the expendable parts are replaced when necessary.

**4 Safety instructions**

Failure to follow the safety instructions mentioned below can lead to faulty control or malfunction of the ErgoArm and result in risk of injuries for the patient.

- Before interrupting or establishing electrical connections (e.g. when removing the hand from the prosthesis), the system must be disconnected from the power supply. For this purpose, either remove the battery from the battery receptacle or turn the prosthesis off by pressing the button in the charging receptacle.
- Educate your patient in the appropriate use of the ErgoArm.
- The ErgoArm has been developed for everyday life activities and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). A careful treatment of the prosthesis and its components is essential not only for prolongation of its life but above all for the patient's safety! Should any extreme loads be applied on the prosthesis (e.g. in case of the fall etc.), then have immediately inspected the prosthesis for possible damages by the prosthetist. If necessary, the prosthet-

ist will pass the prosthesis on to the Ottobock Myo-Service.

- Take special care when lifting heavy loads! Be careful when unlocking the lock under such conditions due to the possibility of injury.
- Ensure that neither solid particles nor liquids can penetrate the ErgoArm.
- The ErgoArm should not be subjected to intense smoke, dust, mechanical vibrations or shocks, or extreme heat.
- Inform the patient that the ErgoArm is to be cleaned with a damp cloth and mild soap only. Do not use solvents such as acetone or benzene as these can make the plastic brittle.
- The electronic components may malfunction in the vicinity of high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores). To reduce this risk, the electrodes should be set to as low a sensitivity as possible.
- In case of malfunction, which could endanger the patient's safety, the system should be turned off by removing the battery. The ErgoArm ElectronicPlus can be locked and unlocked when it is turned off by pulling the control cable mechanically.
- The ErgoArm as well as damaged components of the ErgoArm may be opened or repaired only by the certified Ottobock Myo-Service.
- Should the malfunction persist, please contact the certified Ottobock Myo-Service.
- Risk of pinching where the elbow joint bends. Ensure that fingers and other body parts are not in this area when bending the elbow joint.

 These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

 **Attention!**  
**Car driving capability**  
 The question of the ability of the amputee with an upper limb prosthesis to drive a car does not have a simple answer: It depends on the type of prosthetic fitting (amputation height, uni- or bilateral, residual limb conditions, prosthesis type) as well as the individual capabilities of the upper limb prosthesis user. All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center. In general, Ottobock recommends that the vehicle is adapted by a specialist-manufacturer to fit specific needs of the user (e.g. palm grip) so as to ensure risk-free driving. It is indispensable to ensure that the driver is able to operate the vehicle without any risk.

 **Caution!**  
 Before an ErgoArm with a quick-disconnect wrist can be used, the terminal device must be positioned in such a way that a slight turn of it cannot disconnect the terminal device from the prosthesis!

 **Notice**  
**Distance to communication devices.** If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices) is too small, interference with internal data communications can result in unexpected prosthesis behaviour.

Therefore, keeping the following minimum distances to these HF communication devices is recommended:

- Mobile phone GSM 850/GSM 900: 0.99 m
- Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
- DECT cordless phones incl. base station: 0.35 m
- WiFi (routers, access points,...): 0.22 m
- Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

**Notice**

When operating the product in the immediate vicinity of active, implantable systems (e.g. pacemakers, defibrillators, etc.), ensure that the minimum distances stipulated by the manufacturer of the implant are observed. Interference with active, implantable systems may occur due to the electromagnetic radiation emitted by the product. Furthermore, make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

**12K44**

*ErgoArm Hybrid plus*

**12K50**

*ErgoArm Electronic plus*

The battery can be located either in the forearm socket (connected to the forearm cable) or in the humeral socket (connection via Easy Plug). It is forbidden to use **both** connection options **simultaneously** due to safety reasons.

**5 Technical data**

Maximum load capacity:	Approx. 60 N
Maximum flexion angle:	Approx. 15° – 145°
Storage temperature:	-10° – 60°C
Operating temperature:	5° – 40°C
Storage and transport temperature in packaging:	-25° – +70 °C
Storage and transport temperature without packaging:	-25° – +70 °C
Service life of the elbow joint:	5 years

**6 Symbols on the product**

Declaration of conformity according to the European Directive for Medical Devices 93/42/EEC and 1999/5/EC.

**7 Liability**

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

**8 CE Conformity**

This product meets the requirements of the 93/42/EEC guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

**12K50**  
*ErgoArm Electronic plus*

**Scope of delivery**

- ) Fixation fork
- ) Lamination protection cover (disk with arrow)

Additional parts required:  
... lamination ring


- ) Fixation fork
- ) Lamination protection cover (hollow ball)
- ) Lamination protection cover (disk with arrow)

Additional parts required:  
... Battery Mounting Set  
... lamination ring

- ) Fixation fork
- ) Lamination protection cover (hollow ball)
- ) Lamination protection cover (disk with arrow)
- ) Conical cord stopper

Additional parts required:  
... Battery Mounting Set  
... lamination ring  
... coding plugs / (13E182)

**Laminating the humeral segment**

To dismantle and prepare for lamination, flex the arm fully. Turn the "finger wheel" to minimum compensation and then extend the arm completely.

**Insert the fixation fork 1 into the band cam to secure the band cam to the extension stop.**

Take out the countersunk screw.  
Remove the band clamp.

Remove the lamination ring.

Remove the lamination ring.  
Push the electronics housing out of lamination ring.

**12K42**  
*ErgoArm plus*


**12K44**  
*ErgoArm Hybrid plus*

Lubricate the groove and inner edge of lamination ring with 633F30 special grease.

**12K50**  
*ErgoArm Electronic plus*

Lubricate the elbow cover with 633F30 special grease.

Slide the lubricated elbow cover into the lubricated lamination ring and press **firmly**. **You must hear a loud snapping sound to be certain the cover is fully inserted!**

**Attention:**  
**This illustration shows an elbow cover [1] that is not fully inserted!**  
The visible plastic edge is NOT uniform, so there is risk of lamination resin penetrating the elbow!

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

This illustration shows an elbow cover [1] that is fully inserted. Note that the plastic edge is uniformly visible around the entire circumference.

**12K50**  
*ErgoArm Electronic plus*

Lubricate the lamination protection cover with 633F30 special grease.

Push the lubricated lamination protection cover onto lamination ring.

**You must hear a loud snapping sound to be certain the cover has locked into the groove.**

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

Position the lamination ring.  
**Remember to heed the “posterior” and “anterior” positions!**

**12K50**  
*ErgoArm Electronic plus*

**Note:**  
The lamination ring has three grooves into which the cable can be put when taking off the perlon stockinette. Place the cable for the first perlon stockinette into groove 1. To take off further layers of the perlon stockinette, place the cables into the groove 2 and groove 3 respectively.

Cut off the perlon stockinette at the point where it is twice as long as the socket. Fold the inside of the perlon stockinette outward, pull it over the socket and roll it inward.

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

Place the cable into groove 1 and take off the perlon stockinette. Follow the same procedure for further layers.

**12K50**  
*ErgoArm Electronic plus*

Fold back the perlon stockinette and pull tightly over socket. Repeat procedure for next layers, placing the perlon stockinette into groove 1 and 2 respectively.

**Note:** If the prosthesis is to bear an increased load, strengthen it with carbon.

Laminate the socket. Stretch rubber ring and place around the ungelled resin beneath the lamination protection cap.

**Note:** Do not twist rubber ring!

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

This will ensure a smooth transition to the lamination protection cap.

**12K50**  
*ErgoArm Electronic plus*

Drill four 6mm dia. holes in the humeral section just above the lamination ring to provide ventilation for any moisture that may develop.

**Sealing with 633F30 Special Grease and 633F11 Silicone Grease**


After lamination, clean the lamination ring thoroughly.  
Lubricate the groove, inner edge and thread of the lamination ring with 633F30 special grease.

Remove all dummy plugs from the electronics housing.  
Lubricate the cable connectors with 633F11 silicone grease.

<p><b>12K42</b> <i>ErgoArm plus</i></p>	<p><b>12K44</b> <i>ErgoArm Hybrid plus</i></p>	<p><b>12K50</b> <i>ErgoArm Electronic plus</i></p>
	 <p>Plug in the cables, being certain they are fully inserted into the connectors.</p>	
	 <p>Protect these connections from moisture damage by applying 633F11 silicone grease liberally to the top of the connectors.</p>	
	 <p>Push the electronics housing, including the installed wires, into the lubricated lamination ring and press <b>firmly</b>. <b>You must hear a loud popping sound to be sure the cover is fully inserted!</b></p>	
	 <p><b>Attention: The plate and the brim of the lamination ring must be level with each other!</b></p>	

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

The illustration shows that the plate and brim of the lamination ring are level with each other.


**The electronics housing is fully inserted.**

**12K50**  
*ErgoArm Electronic plus*

**Modifying and servicing the forearm**


**Removing the elbow joint**

Flex the arm fully. Turn the "finger wheel" to minimum compensation and then extend the arm fully.

Insert the fixation fork [1] to hold the cam against the extension stop.

**Lock the elbow joint in this position before performing additional procedures.**


Gently pop off both joint caps.

Carefully pry out the cable plugs.

<p><b>12K42</b> <i>ErgoArm plus</i></p>	<p><b>12K44</b> <i>ErgoArm Hybrid plus</i></p>	<p><b>12K50</b> <i>ErgoArm Electronic plus</i></p>
	 <p>Pull out the cable slightly.</p>	
	 <p>Unscrew the two elbow joint attachment screws.</p>	
	 <p>Pull out the band. Carefully pull elbow joint off the forearm.</p>	
	 <p>Screw the elbow joint completely into the humeral section.</p>	

**Removing the AFB**

	 <p>Unscrew the two attachment screws in the forearm.</p>	
---	---	---

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

Pull the AFB out.

**12K50**  
*ErgoArm Electronic plus*

**Shortening the forearm, positioning the battery box**

Cut or grind the distal forearm to the proper length for the individual.


Insert and bond the 10S1=45/50 lamination ring inside the forearm.

Cut or grind the distal forearm to the proper length for the individual.

Insert and bond the 10S1=45/50 lamination ring inside the forearm.

Wiring will be completed during final assembly of the prosthesis.

**Attention: When mounting the battery box, be sure it does not interfere with the AFB or electric wrist rotator mechanism.**

After mounting the battery box, plug the battery connection cable 13E51=2/4 into the circuit board [1] of the forearm cable. (The circuit board is attached to the AFB.)

**12K42**  
ErgoArm plus


**12K44**  
ErgoArm Hybrid plus

**12K50**  
ErgoArm Electronic plus

Inserting the AFB

Lock the band cam to the extension stop using the fixation fork [1].  
Insert the AFB mechanism into the forearm.

Lock the band cam to the extension stop using the fixation fork [1].  
Insert the AFB mechanism into the forearm.  
**Attention: Be careful not to pinch the cable beneath the AFB !**

Tighten the AFB mounting screws, located in the forearm.

Changing the right-to-left cable orientation

Remove the elbow cap.

Gently disengage telescopic spring.

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

Take off the cord clamp (by cutting at the terminal knot) and remove the telescopic spring.

**12K50**  
*ErgoArm Electronic plus*

Remove the dummy plug and insert it into the hole on the other side.


Thread the telescopic spring into the other side...

...and snap it in onto the projecting metal axle [1].  
Put the ball cap back on.

<p><b>12K42</b> <i>ErgoArm plus</i></p>	<p><b>12K44</b> <i>ErgoArm Hybrid plus</i></p>	<p><b>12K50</b> <i>ErgoArm Electronic plus</i></p>
<p><b>Shortening the control cable</b></p>		
	 <p>Reapply the ladder lock buckle: Thread the control cable as shown in the picture and secure with terminal knot.</p>	
		<p>The ErgoArm 12K50 control cable can be shortened because locking and unlocking is controlled by electrode or switch signals. Thread the supplied ball cord stopper on the control cable end once it is shortened.</p> <p><b>Attention!</b> Complete removal of the control cable is not allowed due to safety reasons. Should the ErgoArm not respond to electrical signals, the control cable may be necessary to set its position.</p> 
<p><b>Inserting the elbow joint</b></p>		
	 <p>Carefully push the elbow joint back into the forearm.</p>	
	 <p><b>Attention:</b> Be careful not to damage the cables during re-assembly!</p>	

**12K42**  
*ErgoArm plus*


**12K44**  
*ErgoArm Hybrid plus*

Apply one drop of Loctite 243 or equal to the two elbow joint attachment screws and tighten them **evenly**.


**12K50**  
*ErgoArm Electronic plus*

Lubricate the **wiring connectors** with 633F11 silicone grease.  
Plug both connectors back into the elbow.

**Note:**  
**The arm must be fully extended!**  
**Press joint covers back into place.**

**Note:**  
**The arm must be fully extended!**  
Push the bevelled joint cap over the cable [1].  
Slide along the cable and then press the joint cap into place.  
This will prevent any cable damage.


Thread band under axis.

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*


**12K50**  
*ErgoArm Electronic plus*

**Pay attention to correct positioning of lamination ring! (See picture!)**

Leave space for the band clamp between the band clamp guide and the lamination ring.

Be certain the band is straight and not twisted.

The illustration shows *incorrect* positioning of the lamination ring:

The band clamp cannot be inserted.

Push the band clamp under the lamination ring and attach with screw.

**Adjustment of preflexion**

Flex the elbow and minimize the compensation force using the "finger wheel".

Extend the forearm and fix the band cam against the extension stop with the fixation fork [1].

Flex the forearm until the band forms a loop.

**12K42**  
*ErgoArm plus*


**12K44**  
*ErgoArm Hybrid plus*


Lock the elbow joint with cable lock.  
Loosen the screw on the band clamp.  
Shorten or lengthen the band.  
Tighten the screw on the band clamp.  
Remove the fixation fork to free the band cam.  
**To decrease preflexion = lengthen the band**  
**To increase preflexion = shorten the band**

**12K50**  
*ErgoArm Electronic plus*

**A 3 mm change in the length of the band causes a 5° change in the flexion angle.**  
Make sure that the forearm swings freely and that the band cam is still fastened to the extension stop with the fixation fork.  
The preflexion angle is initially set at about 10 degrees.

Changing the band length also affects the free swing range.  
When the preflexion angle is increased 10° by tightening the band, the forearm will enter the second compensation range 10° later in flexion as well.  
(See also Compensation Ranges p. 47)

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*


Note: For proper functioning of the AFB flexion compensation mechanism, the band cam must reach its extension stop when the elbow is fully extended.

**12K50**  
*ErgoArm Electronic plus*

**Control harness cable attachment**

The bowing harness should be attached in the marked area.

**Important!**

Take care during the attachment so as to avoid the penetration of the control harness cable into the socket interior and possible damaging of system components.

We recommend following one of the two fixing examples seen in the pictures above.

When using the clamp stopple, use the 21A207 Clamp Stopple Set and proceed according to the 647G231 Instructions for Use.

21A207

**12K42**  
*ErgoArm plus*

**12K44**  
*ErgoArm Hybrid plus*

**12K50**  
*ErgoArm Electronic plus*

## Function of the AFB

### AFB compensation adjustment

The amount of AFB compensation can be adjusted with a "finger wheel" and adapted to individual weight and various clothing.

**Adjustments are easier to make when the forearm is flexed!**

The range of adjustments is limited by definite stops.

The mechanism is also protected by a slip coupling.

The hydrodynamic damper of the AFB mechanism is adjusted to counter-balance the weight of the forearm plus terminal device and electric wrist.

Please be sure to stabilize the humeral section of the prosthesis while testing the joint functions.

Always keep the prosthesis in the fully flexed position when putting it on or taking it off, and when not in use.

**Caution!**

If the AFB compensation is set high and the lock is released while the arm is in a horizontal position, the arm can flex suddenly.

A sudden flexion of either the forearm or humeral section may also occur while working with the prosthesis when it is taken off.

This may result in injury!


**Instruct the patient about this danger!**

**12K42**  
ErgoArm plus

**12K44**  
ErgoArm Hybrid plus

**12K50**  
ErgoArm Electronic plus

### Functional range

**Position 1** Minimal compensation that allows natural free swing of the arm while walking.

**Position 2** Compensation increases progressively with flexion of the arm, and decreases automatically during extension.

**Position 3** Compensation remains constant. If set correctly, the weight of the forearm is balanced by the compensation, and the elbow keeps flexed.

**Important:**

Compensation can be adjusted so that the arm will remain in a flexed position despite the force of gravity. If the patient wants to extend the arm again, extension can be initiated by a quick forward movement of the shoulder or humerus.

**Position 4** Reduced compensation until flexion stop.

### Humeral rotation feature

The humeral rotation turntable has integral stops ( $\pm 80^\circ$ ) to block extreme movements. Turntable friction is easily adjusted by turning the external adjustment screw, using a coin or similar object.

## 12K50 - ErgoArm Electronic Plus

### Setting-up Electrodes

Set the electrodes to the control system of the System Electrohand, or 13E205 MyoRotronic, as applicable. Set-up is easy with 757M11 MyoBoy. For detailed information see Operational Manual of the System Electrohand, or respectively, 13E205 MyoRotronic.

The electrode contacts should be placed flat on intact skin, touching it with as much of their surface space as possible. Should any strong interference be observed due to the electric devices, then the electrode position should be checked and changed, if necessary. If the interference cannot be eliminated, please contact Ottobock Myo-Service.

#### Attention!

Be sure to include breaks during the adjustment of the electrodes, otherwise muscle fatigue may cause irregular results, which could lead the prosthetist to set the electrodes too sensitive a setting.

### Safety management

If the defective function of electronically controlled lock of ErgoArms Electronic plus 12K50 occurs, the user is advised by warning signal lasting about 1 second. The electronically controlled lock will subsequently be disconnected. In such a case remove the battery for a while and then replace it again. Should the electronic lock be disabled again by the Safety Management, than the ErgoArm Electronic plus 12K50 must be sent for repair to Ottobock Myo-Service. By the way the locking and unlocking can be operated manually at any time via control cable even when the lock is disabled.

### Switching mode selection

The ErgoArm Electronic plus 12K50 is provided with several switching modes for electronic lock control. They include various options for electrode signal switching between the hand and the elbow. The required switching mode is selected by inserting a colored coding plug into the electronics housing.

For those patients incapable of generating suitable myoelectric signals, the control system can be set up by Ottobock Myo-Service. (This service is not available in all countries!)

### Switching

#### Switch type:

All Ottobock switches can be used which are designed for myoelectric prostheses. See MYOBOCK Arm Components Catalogue.

#### Electrodes:

Simultaneous contraction of both muscles - short and quick (co-contraction).

#### Battery use

The 12K50 ErgoArm Electronic plus can be powered with the 757B20, 757B21, 757B25=\*, 757B35=\* EnergyPack or with the 757B15 Ottobock Interchangeable Battery. For detailed instructions on the use of Energy Packs see the information attached to EnergyPack.

#### Attention!

For safety reasons as well as system reliability, use only Ottobock Batteries 757B20, 757B21, 757B25=\*, 757B35=\* or 757B15.

### Technical data 12K50

Current demand:		3.5 mA
Power supply:		
EnergyPack	757B20	7.2 V (Li-Ion technology)
EnergyPack	757B21	7.2 V (Li-Ion technology)
Interchangeable Battery	757B15	6 V (NiMH technology)
MyoEnergy Integral	757B25=* /	
	757B35=*	7.4 V (Li-Ion technology)

## Annex

# Switching modes 12K50

When the battery is inserted, the current switching mode is indicated by vibration feedback signal.

Switching mode		Lock	SLIP-STOP	Four Channel Processor II	Hand
1	white	Press the switch → release = lock Press the switch → release = unlock	NO	All switching modes	All versions
2	red	Co-contraction = lock Co-contraction = unlock	NO	Only Program 1 (white) Four-channel Control or Program 9 Adjustment Cap	All versions with 2 electrodes
3	green	Press and hold down switch = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock Release the switch = Hand mode	YES	All switching modes	
4	blue	Press the switch → release = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock Press the switch → release = hand mode	YES		
5	yellow	Press the switch → release = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock 10 s no electrode signal = hand mode or press the switch → release = hand mode	YES		
6	magenta	Co-contraction = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock Co-contraction = hand mode	YES	Only Program 1 (white) Four-channel Control or Program 9 Adjustment Cap	<i>Not recommended for hands with Digital or Digital-Twin Control System</i>
7	orange	Co-contraction = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock 10 s no electrode signal = hand mode or co-contraction = hand mode	YES		

Vibration feedback signal for successful switching over between hand and elbow is activated in Programs 3–7.

**1x vibration = hand mode (the hand is controlled by electrode signals)**  
**2x vibration = elbow mode (the elbow is controlled by electrode signals)**

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