




## 13E205

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	15
<b>FR</b> Instructions d'utilisation .....	27
<b>IT</b> Istruzioni per l'uso .....	39
<b>ES</b> Instrucciones de uso .....	51
<b>PT</b> Manual de utilização .....	64
<b>NL</b> Gebruiksaanwijzing .....	76
<b>SV</b> Bruksanvisning .....	88

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Date of the last update: 2015-02-11

- Please read this document carefully.
- Follow the safety instructions and the precautions specified in this accompanying document.
- If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

## 1 Intended use

### 1.1 Medical purpose

The 13E205 MyoRotronic is to be used **exclusively** for the exoprosthetic fitting of the upper extremities.

Only start up the product in accordance with the information contained in the accompanying documents supplied.

Instruct the patient in the proper use and care of the product. Handover to the patient without prior instruction is not permitted.

### 1.2 Trademarks

All product names mentioned in this accompanying document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners. Should trademarks used in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 1.3 Application

The 13E205 MyoRotronic can be used for unilateral or bilateral amputees.

### 1.4 Application conditions

The 13E205 MyoRotronic has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

### 1.5 Prosthetist qualifications

The fitting of a patient with the 13E205 MyoRotronic may only be carried out by prosthetists who have been trained and certified by Ottobock for the use of MyoBock components.

#### **Attention!**

Failure to follow the safety instructions mentioned below can lead to faulty control or malfunction of the MyoRotronic and result in risk of injuries for the patient.

## 2 Safety instructions

- When operating the product in the immediate vicinity of active, implantable systems (e.g. pacemakers, defibrillators, etc.), ensure that the minimum distances stipulated by the manufacturer of the implant are observed. Interference with active, implantable systems may occur due to the electromagnetic radiation emitted by the product. Furthermore, make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.
- The 13E205 MyoRotronic must not be combined with the 12K100 DynamicArm, since this can destroy the MyoRotronic.
- Do not allow debris or liquids to get into the MyoRotronic.
- The MyoRotronic should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- Wearing the MyoRotronic close to high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores) may cause the MyoRotronic or the connected components to malfunction. Setting the electrode gains as low as possible will minimize this risk.
- Before interrupting or establishing electrical connections (such as when removing the MyoRotronic from the Electric Wrist Rotator), the system must be disconnected from the power supply. For this purpose, either remove the battery from the battery receptacle or turn the prosthesis off by pressing the button in the charging receptacle.
- In case of malfunction, which could endanger the patient's safety, the system should be turned off by removing the battery.
- If malfunctions persist, contact an Ottobock Service facility.
- Modifying component parts and software without the written consent of the manufacturer is strictly prohibited.
- Educate your patient in the appropriate use of the MyoRotronic.
- If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices) is too small, interference with internal data communications can result in malfunctions of the electric wrist rotator. Maintaining the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.99 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points,...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

This product may not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

### 3 Components (Fig. 1)

Delivery includes:

1 piece of 13E205 MyoRotronic (Fig. 1, item 4)

1 piece of 13E183 Buzzer

1 piece of Protective Cover (Fig. 1, item 1)

1 piece of 647G361 Instructions for Use

For pronation and supination, the following is required:

1 piece of 10S17 Electric Wrist Rotator (Fig. 1, item 3)

### 4 Function

In combination with the 10S17 Electric Wrist Rotator, the MyoRotronic enables electromotive pronation and supination as well as the opening and closing of the hand with one or two 13E200/13E202 MyoBock Electrodes or a combination of one 13E200/13E202 Electrode and a 9X50/9X52 Linear Transducer.

To guarantee optimum adaptation to the patient, five different programs are available. Four of the five programs also allow proportional pronation and supination in addition to the digital rotation control. The desired program is set using the 757T13 MyoSelect.

The MyoRotronic can be operated in combination with all Ottobock System Electric Hands and System Electric Greifers that are equipped with Quick Disconnect Wrist. In operation every switchover between System Electric Hand or System Electric Greifer and Electric Wrist Rotator is indicated by brief vibration of the Electric Wrist Rotator (except the programs 1 and 4).

#### 4.1 Training phase

Connecting the 13E183 Buzzer for acoustic feedback on the switchover will increase training success.

In addition, when inserting the battery an acoustic signal tone will indicate the selected program (number of the selected program = number of signal tones).

### 5 Adjusting the electrodes

The MyoRotronic will perform best if the electrodes have been optimally adjusted. This is easy, thanks to the 757M11 MyoBoy. In the first step, using the MyoBoy the electrodes are adjusted to the corresponding hand control. In the next step, the MyoBoy will be used to determine the switching mode for the MyoRotronic that is best suitable for the patient. The electrodes should be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interferences from electric devices, the position of the electrodes should be checked and, if necessary, repositioned. If the interference cannot be eliminated, please contact the Ottobock Myo-Service.

#### **Attention!**

Allow the patient to rest during the adjustment of the electrodes to avoid fatigue and thus incorrect electrode settings (i.e. excess sensitivity).

Signal sources			Strength of myo-signal <sup>1</sup>						Program	
			Pronation			Supination				No.
Switch	Two electrodes	One electrode	ON <sup>2</sup>	LOW <sup>2</sup>	HIGH	ON <sup>2</sup>	LOW <sup>2</sup>	HIGH		
	●			●	●		●	●	Four-channel control	1
	●		●		●	●		●	Cocontraction	2
	●		●		●	●		●	Safety Cocontraction	3
●	●		●			●			Switch	4
		●		●	●		●	●	One-electrode control	5

**Note:** ON, LOW and HIGH correspond to the scale of: 757M11 MyoBoy, MyoSoft/ PAULA (Fig. 4, 5 and 6).

<sup>1</sup> The patient should be able to maintain the signal for two seconds above the respective threshold.

<sup>2</sup> Not suitable for the proportional control of an Electric Wrist Rotator.

**5.1 Setting the programs of the MyoRotronic**

The MyoRotronic provides various modes for switching the electrode signals between the System Electric Hand/ System Electric Greifer and the Electric Wrist Rotator. The programs are set using the 757T13 MyoSelect. The MyoRotronic is delivered with the program 1 preset.

**5.1.1 Selecting a program with 757T13 MyoSelect**

Connect the 757T13 MyoSelect with connected MyoRotronic to the Electric Wrist Rotator and select a program. To do this, read the 647G131 Instructions for Use for the 757T13 MyoSelect.

**Note:**  
Remove the battery from the prosthesis before connecting the MyoSelect!

**5.1.2 Selection of digital and proportional control of the Electric Wrist Rotator**


For the programs 1, 2, 3, and 5, MyoSelect allows to select between digital and proportional control. The setting is made via the MyoSelect. The presently selected control is shown on the display.

**Example, program 1**

**Switch symbol closed:** proportional control of the Electric Wrist Rotator.

Rotation speed of the Electric Wrist Rotator is determined by the speed and strength of muscle contraction.

**Switch symbol opened:** digital control of the Electric Wrist Rotator.

The Electric Wrist Rotator always rotates at maximum speed.

**Program overview**

No.	Method	Proportional control	Digital control
1	Four-channel control with adjustable thresholds	●	●
2	Cocontraction	●	●
3	Safety Cocontraction	●	●
4	Switch		●
5	One-electrode control	●	●

**5.2 Switching**

**Cocontraction (Prog. 2)**

Simultaneous contraction of both muscles – short and quick.

**TIP** e. g. balling a fictitious fist or spreading fictitious fingers. Avoid cramps!

**Safety cocontraction (Prog. 3)**

Normal Cocontraction and subsequent brief contraction of one of the two muscles within 2 seconds (after the end of Cocontraction).

**TIP** Both muscles must be relaxed completely between Cocontraction and the following single contraction to ensure that the switching signal is recognized.

## Impulse switching in case of one-electrode control (Prog. 5)

Contraction of one muscle - short and quick.

**TIP** e.g. balling a fictitious fist or spreading fictitious fingers. Avoid cramps!

MyoRotronic

Switching

Hand control

## 6 Switching modes

Signal sound


Signal sound + vibration

**Program 1:** Four-channel control  
Control with 2 electrodes

*Opening and closing:* By slowly beginning muscle contraction via the electrode: It takes more than 80 ms for the electrode signal to reach the respective upper threshold after exceeding the ON-threshold. The prosthetic component opens or closes until the electrode signal subsides.


*Pronation and supination:* By quickly and strongly beginning muscle contraction via the electrode: The electrode signal reaches the respective upper threshold within 80 ms after exceeding the ON-threshold.

Rotation can be controlled both digitally and proportionally (see Section 5.1.2). The optional 9X24 Switch can be used to turn off pronation and supination.


	<b>A</b>	<b>B</b>
	<b>Upper threshold</b>	<b>Upper threshold</b>
<b>Left arm prosthesis</b>	Rotation: Pronation Hand: CLOSE MyoSoft/PAULA: blue line	Rotation: Supination Hand: OPEN MyoSoft/PAULA: red line
<b>Right arm prosthesis</b>	Rotation: Supination Hand: OPEN MyoSoft/PAULA: red line	Rotation: Pronation Hand: CLOSE MyoSoft/PAULA: blue line

**Program 2:** Cocontraction Control with 2 electrodes


Switching between gripping function and rotation is performed with Cocontraction. See section 5.2 "Switching". In addition, a time control feature can be used to switch back to the gripping function. The switch-back time is set using the 757T13 MyoSelect. The respective setting can be read on scale C on the MyoSelect's display.

Rotation can be controlled both digitally and proportionally (see Section 5.1.2).

MyoSelect scale C	Switch-back time
0	2 sec.
1	4 sec.
2	6 sec.
3	8 sec.
4	10 sec.
5	12 sec.
5,9	14 sec.
6	Automatic switch-back deactivated; switch-back exclusively by cocontraction.

Upon delivery, the automatic switch-back is deactivated. The scale C in the MyoSelect's display is therefore set to 6.

**Program 3:** Safety Cocontraction  
Control with 2 electrodes

Switching between gripping function and rotation is performed with safety Cocontraction. See section 5.2 "Switching". In addition, a time control feature can be used to switch back to the gripping function. The switch-back time is set using the 757T13 MyoSelect. The respective setting can be read on scale C on the MyoSelect's display.



Rotation can be controlled both digitally and proportionally (see Section 5.1.2).

MyoSelect scale C	Switch-back time	
0	2 sec.	
1	4 sec.	
2	6 sec.	
3	8 sec.	
4	10 sec.	
5,9	14 sec.	
6	Automatic switch-back deactivated; switch-back exclusively by safety Cocontraction.	


Upon delivery, the automatic switch-back is deactivated. The scale C in the MyoSelect's display is therefore set to 6.

**Program 4:** Control with 2 electrodes and a MyoBock switch

Switching from gripping function to rotation is performed by actuating and holding the switch. Switch-back to the gripping function takes place as soon as you let go of the switch.

**Program 5:** One-electrode control (1 electrode or 1 linear transducer)

Switching between gripping function and rotation is performed by impulse switching, which means that the control signal must quickly exceed the upper threshold and then quickly subside. See section 5.2 "Switching". In addition, a time control feature can be used to switch back to the gripping function. The switch-back time is set using the 757T13 MyoSelect. The respective setting can be read on scale C on the MyoSelect's display.

Rotation can be controlled both digitally and proportionally (see Section 5.1.2).

The upper threshold is set using MyoSelect (scale D).

Functions		
<b>Opening and closing</b> according to the control mode of the System Electric Greifer or System Electric Hand suitable for one-electrode control (e.g. SensorHand Speed control mode 3).	<b>Supination</b> Electrode: By quickly and strongly beginning muscle contraction  Linear transducer: By quickly beginning to pull on the linear transducer	<b>Pronation</b> Electrode: By slowly beginning muscle contraction  Linear transducer: By slowly beginning to pull on the linear transducer
	<b>Note:</b> To avoid damaging the transducer, the linear transducer must not be pulled beyond the stop.	

MyoSelect scale C	Switch-back time	
0	2 sec.	
1	4 sec.	
2	6 sec.	
3	8 sec.	
4	10 sec.	
5	12 sec.	
5,9	14 sec.	
6	Automatic switch-back deactivated; switch-back exclusively by impulse switching.	

Upon delivery, the automatic switch-back is deactivated. The scale C in the MyoSelect's display is therefore set to 6.

## 7 Assembly into a prosthetic component


Apply 633F11 Silicone Grease to all contact bushings of the MyoRotronic, the red protective plug and the connection cables. The grease will protect the connection from humidity.

- Remove the protective cover from the motor of the Electric Wrist Rotator (Fig. 2, item 1).
- Cut the lock ring (Fig. 1, item 2) where it is notched and remove it.
- Remove the red protective plug (Fig. 2, item 3).
- Remove the protective cover (Fig. 1, item 1) from the MyoRotronic.


Slide the MyoRotronic onto the Electric Wrist Rotator so that the marking 1 or 3 remains free (as illustrated).

### a) Setting for 2 electrodes (Programs 1, 2 and 3)

#### In left arm prostheses


#### In right arm prostheses


### b) Setting for 1 switch (Program 4)

#### In left arm prostheses


#### In right arm prostheses


### c) Setting for 1 electrode (Program 5)

#### In left arm prostheses

#### In right arm prostheses

- Connect the motor, battery and electrode connections as illustrated (depending on the switching mode, also connect the switch/buzzer).
- Lay the motor cable in loops to avoid damaging the cable and carefully put it through the groove (Fig. 3).
- Put on the protective cover (see Fig. 1, item 1). Slide the red protective plug onto the vacant pins (Fig. 1, item 2).
- Push the Electric Wrist Rotator into the socket until it engages in the lamination ring.
- Secure the Electric Wrist Rotator with the 11S4 Lock Ring.

## 8 Handling the batteries

Only fully charged batteries should be used to operate the MyoRotronic. An intelligent "battery manager" alerts the patient to low voltage conditions: the Electric Wrist Rotator becomes slower and slower. At the same time the battery is protected against harmful total discharge.

For more information regarding the batteries please consult the Operating Manual enclosed with the batteries.

### Attention:

As soon as the MyoRotronic recognizes a fully charged 757B20/757B21 EnergyPack or 757B25=\* / 757B35=\* MyoEnergy Integral, it automatically switches to Li-Ion battery technology. Once this has happened, only this technology should be used to operate the system. If a 757B15 Ottobock X-ChangePack is employed after the control has switched to Li-Ion battery technology, the battery can not be used to its full capacity. The 757T13 Ottobock MyoSelect can be used, however, to reset the product for use with X-ChangePack.

**For operational security and reliability, use only 757B20, 757B21, 757B25=\*, 757B35=\* or 757B15 Ottobock Batteries.**

Upon delivery, the default setting is optimized for NiMH Batteries.

## 9 Cleaning and care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Ensure that no liquid penetrates into the system component.

Then, dry the prosthesis component with a soft cloth.

## 10 Technical data

Power supply	757B20	Ottobock EnergyPack (7.2 V)
	757B21	Ottobock EnergyPack (7.2 V)
	757B25=*	MyoEnergy Integral (7.4 V)
	757B35=*	MyoEnergy Integral (7.4 V)
	757B15	Ottobock X-ChangePack (6 V)
Static current		1 mA
Operating temperature range		0 °C – +70 °C
Storage and transport in original packaging		-20 °C – +40 °C
Storage and transport without packaging		-20 °C – +40 °C
Relative humidity		80% non-condensing
Power off		load dependent between 30 ms and 10 s
Weight		10 g
Service life		5 years

## 11 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

## 12 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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