
# ottobock.



## 4R57/4R57=ST

<b>DE</b>	Gebrauchsanweisung .....	4
<b>EN</b>	Instructions for use .....	11
<b>FR</b>	Instructions d'utilisation .....	18
<b>IT</b>	Istruzioni per l'uso .....	25
<b>ES</b>	Instrucciones de uso .....	32
<b>PT</b>	Manual de utilização .....	40
<b>NL</b>	Gebruiksaanwijzing .....	47
<b>SV</b>	Bruksanvisning .....	55
<b>DA</b>	Brugsanvisning .....	62
<b>PL</b>	Instrukcja użytkowania .....	69
<b>HU</b>	Használati utasítás .....	76
<b>CS</b>	Návod k použití .....	83
<b>HR</b>	Upute za uporabu .....	90
<b>BS</b>	Upute za uporabu .....	97
<b>TR</b>	Kullanma talimatı .....	103
<b>EL</b>	Οδηγίες χρήσης .....	110
<b>RU</b>	Руководство по применению .....	118
<b>JA</b>	取扱説明書 .....	126
<b>ZH</b>	使用说明书 .....	132

	■	●	▲
DE	Einzelteile	Einzelteile-Pack	Mindestmenge
EN	Single components	Single-Component Pack	Minimum order quantity
FR	Pièces à l'unité	Kit de pièces de rechange	Quantité minimum
IT	Singoli componenti	Imballo dei singoli componenti	Quantità minima
ES	Componentes	Kit componentes	Cantidad mínima
PT	Componentes	Conjunto de peças de reposo	Quantidade mínima
NL	Onderdelen	Onderdelenpakket	Minimikvantitet
SV	Separata delar	Service-Set	Minimum aantal
DA	Komponenter	Komponentpakke	Mindste mængde
NO	Enkle komponenter	Enkeltkomponentpakke	Minste ordremengde
PL	Podzespoły	Zestaw naprawczy	Ilość minimalna
HU	Alkatrészek	Alkatrész-csomag	Minimum mennyiség
CS	Jednotlivé díly	Sada jednotlivých dílů	Minimální množství pro objednání
RO	Componente individuale	Set componente individuale	Cantitate minimă
HR	Pojedinačni dijelovi	Pakiranje pojedinačnih dijelova	Minimalna količina
TR	Münferit parçalar	Münferit parça ambalajı	Minimum miktar
EL	Μεμονωμένα εξαρτήματα	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα
RU	Отдельные детали	Комплект деталей	Минимальное количество
JA	単一成分パック	単一成分パック	最低発注量
ZH	单个部件	维修组件	最低起订量

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These instructions for use were last updated on 16.01.2015.

### INFORMATION

The items included in the delivery are shown on the cover page.

## Explanation of symbols

English

**⚠ DANGER** Warnings regarding directly impending risks of severe accident or injury.

**⚠ WARNING** Warnings regarding possible risks of severe accident or injury.

**⚠ CAUTION** Warnings regarding possible risks of accident or injury.

**NOTICE** Warnings regarding possible technical damage.

## 1 Safety instructions

### CAUTION

**Non-observance of safety instructions.** Non-observance of the following safety instructions can lead to malfunction of the 4R57 Rotation Adapter, resulting in a risk of injury to the patient and in destruction of the 4R57 Rotation Adapter.

Observe the safety information in these instructions for use and inform the patient.

### INFORMATION

The 4R57 Rotation Adapter cannot be combined with Adapters 2R49, 2R50, 4R69, 4R98 and 4R95, because the connectors have different dimensions.

### CAUTION

**Environment and usage conditions.** The environment and usage conditions may compromise the functionality of the 4R57 Rotation Adapter. This may lead to defects and resulting malfunction, including loss of function due to structural failure. This can cause the patient to fall.

- Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water and acids.
- Do not let foreign particles or liquids get into the rotation adapter.
- The rotation adapter should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.

### CAUTION

**Changes to system components.** Any changes or modifications to system components which are made on your own initiative may result in defects and malfunction of the 4R57 Rotation Adapter, including loss of function due to structural failure. This can cause the patient to fall.

- Any changes or modifications to the 4R57 Rotation Adapter may result in a loss of functionality.

- The rotation adapter and damaged components may only be opened or repaired by certified Otto Bock Service technicians.

### CAUTION

**Mechanical loads.** External mechanical influences and/or loads, such as impacts and vibration, can cause defects and result in malfunction of the 4R57 Rotation Adapter, including loss of function due to structural failure. This can cause the patient to fall.

- Do not subject the system components to mechanical vibrations or impacts.
- Avoid any shock or impact to the release button, as this could impair the rotational function of the rotation adapter.

### INFORMATION

If you notice that the rotation adapter has been damaged in any way, please send it to an Ottobock service centre.

### CAUTION

**Improper use of the 4R57 Rotation Adapter.** Any type of excessive strain or overloading can lead to defects and resulting malfunction of the 4R57 Rotation Adapter, including loss of function due to structural failure. This can cause the patient to fall.

Instruct the patient in the proper use of the 4R57 Rotation Adapter.

## 2 Scope of delivery

- 1 4R57 Rotation Adapter  
or
- 1 4R57=ST Rotation Adapter
- 1 647G258 Instructions for Use

## 3 Single components (Figure 1) ■

### 3.1 Single component

- 4X69=1 Release Button Cover, gray (4)

## 3.2 Single components package ●

–

## 3.3 Minimum order quantity ▲

506G3=M8×12 Set Screws (1)

## 4 Description

### 4.1 Intended use

The 4R57 and 4R57=ST Rotation Adapters are intended **solely** for the exoprosthetic fitting of amputations of the lower limb above the prosthetic knee joint.

Amputation level: transfemoral amputation, hip disarticulation.

### 4.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System**:

#### 4R57, 4R57=ST:

Approved for a **patient weight of up to 150 kg (330 lbs)**.

The Rotation Adapters have been developed for everyday use and **must not** be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

The fitting of patients with the 4R57 and 4R57=ST Rotation Adapters may only be carried out by a prosthetist.

### 4.3 Function

During daily life with a transfemoral prosthesis “combination movements”, such as for sitting with the knees crossed, are possible only to a limited extent. Position changes of the sitting posture, in order to increase the wearing comfort of the patients, are desired as well. The 4R57 and 4R57=ST Rotation Adapters comply with these additional requirements. They must be installed above the prosthetic knee joint so that the lower leg rotates against the socket when flexed, i.e. so that it can be rotated inward and outward. When driving a car or when changing shoes the advantages become obvious (Fig. 2 and 3).

The release button need only be pressed to allow the unit to rotate. The unit will automatically lock in the neutral position.

## INFORMATION

Position the Rotation Adapter's release button medially!

### Structural features:

The rotation joint of the Rotation Adapter consists of two parts that rotate on bearings, as well as a locking device.

The high forces that are generated by the weight of the prosthesis wearer, are absorbed without play by the two bearings. These bearings require a torque moment of only 0.5 Nm/0.37 ft.lbs./4.42 lbs.in. when installed. In the neutral position the unit locks automatically.

### 5 Technical data

Article number	4R57	4R57=ST
Weight	170 g	185 g
System height	22 mm	42 mm
Material	Stainless steel	Stainless steel
Rotation possibility	max. 360° (without foam cover)	max. 360° (without foam cover)
Max. patient weight	150 kg (330 lbs)	150 kg (330 lbs)
Proximal connection	Adjustment pyramid	Thread
Distal connection	Pyramid Receiver	Pyramid Receiver
Mobility Grade	–	–

### 6 Handling

#### 6.1 4R57 Rotation Adapter with adjustment pyramid (Fig. 4)

The Rotation Adapter has a special adjustment pyramid for connection to a socket attachment block with pyramid receiver or double adapter with pyramid receiver. In the case of replacement or disassembly all four set screws of the Adapter must be removed (Fig. 4, item 1). For static corrections, adjustment possibilities are minimal. Use a torque wrench to tighten the set screws (such as Ottobock art. no. 710D4).

## INFORMATION

The torque of the adjustment screws of the proximal adapter connection is **10 Nm** (Fig. 4). This value differs from the torque specified in the relevant user's manual.

**For trial fitting:** Turn the adjustment screws clockwise to a torque of **10 Nm** (Fig. 4).

**For finished prosthesis:** Secure the set screws with 636K13 Loctite and tighten them clockwise to **10 Nm/7.4 ft.lbs./88.5 lbs.in.** and **15 Nm/11 ft.lbs./133 lbs.in.** respectively (Fig. 4).

## 6.2 4R57=ST Rotation Adapter with threaded connector (Fig. 5)

The threaded connector of the Rotation Adapter provides the connection to the lamination anchors 4R43 and 4R111=N.

## INFORMATION

For both lamination anchors the 4x46=ST Lamination Aid is to be used (not included in the delivery).

**Procedure:** Screw the Rotation Adapter in completely. Then unscrew the Rotation Adapter by one turn, **maximally two turns**, until

- the release button is positioned medially and
- after pressing the release button, the housing of the Rotation Adapter can be rotated 360° without contacting the lamination anchor.

Fix the Rotation Adapter in the lamination anchor by tightening the clamping screw (Fig. 5, item 2).

Use a torque wrench to tighten the clamping screw (such as Ottobock art. no. 710D4).

## INFORMATION

The torque for the clamping screw of the 4R111=N and of the 4R43 is **10 Nm/7.4 ft.lbs./88.5 lbs.in.** (Fig. 5).

## INFORMATION

The torque for the set screws is **15 Nm/11 ft.lbs./133 lbs.in.** (Fig. 5).

**For trial fitting:** Tighten the set screws clockwise to **10 Nm/7.4 ft.lbs./88.5 lbs.in.**

**For finished prosthesis:** Secure the set screws with 636K13 Loctite and pre-tighten to **10 Nm/7.4 ft.lbs./88.5 lbs.in.** Then complete assembly by tightening to **15 Nm/11 ft.lbs./133 lbs.in.**

### 6.3 Pyramid receiver

The pyramid receiver of the Rotation Adapters provides the distal connection to the knee joint. Static corrections are possible at any time during alignment and trial walking and after completion of the prosthesis using the four set screws. To exchange a module or during disassembly the previous position may be retained by loosening only the two most deeply screwed in, adjacent set screws.

Replace set screws that appear too long (to prevent foam cover damage) or too short (which can fail) with the appropriate screw length:

The following dimensions are available:

Article no.	Length
<b>506G3=M8x10</b>	10 mm
<b>506G3=M8x12-V</b>	12 mm
<b>506G3=M8x14</b>	14 mm
<b>506G3=M8x16</b>	16 mm

### 6.4 Maintenance instructions

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

## 7 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock

Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

## 8 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock Healthcare Products GmbH