


8E33=7 / 8E34=7

DE Gebrauchsanweisung.....	3
EN Instructions for use.....	8
FR Instructions d'utilisation.....	14
IT Istruzioni per l'uso.....	20
ES Instrucciones de uso.....	25
PT Instrução de utilização.....	31
NL Gebruiksaanwijzing.....	37
SV Bruksanvisning	42

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Last update: 04.03.2015

- **Please read this document carefully.**
- **Follow the safety instructions and the precautions specified in this accompanying document.**
- **If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).**

Only put the product into use in accordance with the information contained in the accompanying documents supplied. Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

1 Application

English

The 8E33=7 and 8E34=7 Ottobock System Electric Greifers are to be used solely for the fitting of upper limb prostheses.

Electric Greifers offer greater function and efficiency to users of the MyoBock System. These Greifers are designed to handle a variety of specialized tasks which are unsuitable with System Electric Hands. They offer the following advantages:

- ✓ **Quick and precise gripping**
i.e. used with precision mechanics, optics, electronics
- ✓ **Handling with heavy objects**
i.e. in industry or farming
- ✓ **Special gripping functions**
by using various fingertips

For all of these situations, the Greifer has the ideal features:

- Universal adjustment of the grip surface by rotating the fingertips
- Passive wrist joint allowing flexion, pronation and supination which permits functional operation in all three planes.

1.1 Qualifications of the orthopaedic technician

Fitting a patient with the product may only be carried out by orthopaedic technicians who have been authorised by Ottobock after completion of a corresponding training course.

2 Safety tips

Before interrupting or establishing any electrical connections (e.g. when removing the Electric Greifer from the Prosthesis) it is absolutely necessary to remove the battery from the battery receptacle and thus to separate the system from the power supply.

If the battery becomes fully discharged for any reason, you can release the Greifer using the safety release lever (fig. 1, item 4). The Greifer can be opened by forcefully prying the fingers apart.

If the System Electric Greifer will not be used for an extended period of time, it should be left in the open position to protect the mechanics.

Do not allow foreign particles or liquids to get into the Greifer mechanism. The System Electric Greifer should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.

Due to strong electromagnetic radiation the System Electric Greifer might open or close unexpectedly when near high-tension power lines, transmitters, or transformers, or when going through safety gates. Setting the electrode gains as low as possible will minimize this risk of inadvertent operation.

If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices) is too small, interference with internal data communications can result in unexpected prosthesis behaviour. Therefore, keeping the following minimum distances to these HF communication devices is recommended:

- Mobile phone GSM 850/GSM 900: 0.99 m
- Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
- DECT cordless phones incl. base station: 0.35 m
- WiFi (routers, access points,...): 0.22 m
- Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

The product and any damaged components may only be opened and/or repaired by certified Ottobock Myo-Service technicians.

When operating the product in the immediate vicinity of active, implantable systems (e.g. pacemakers, defibrillators, etc.), ensure that the minimum distances stipulated by the manufacturer of the implant are observed. Interference with active, implantable systems may occur due to the electromagnetic radiation emitted by the product. Furthermore, make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

Attention:

Driving Vehicles

The ability of arm prosthesis wearers to operate a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the individual abilities of the arm prosthesis wearer. Observe your country's national driving laws when operating a vehicle. For insurance purposes, have your driving ability examined and approved by an authorised test centre. In general, Ottobock recommends that the vehicle be professionally retrofitted (e.g. with a steering fork) to meet each user's individual needs. Ensure that the driver is able to operate the vehicle without any risk with the System Electric Greifer turned off. Driving with the System Electric Greifer ON may present a risk if the System Electric Greifer inadvertently moves due to unintentional muscle contraction or other causes.

Caution!

The 8E33=7 Electric Greifer is equipped with a quick-disconnect adapter. Before using it to drive a vehicle, it must be positioned in such a way that slight rotation of the Greifer- which might occur when steering a car - will not disconnect the Greifer from the prosthesis.

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

3 Scope of delivery

- 1 pc. 8E33=7 System Electric Greifer Digital Twin with system wrist unit and quick-disconnect or
- 1 pc. 8E34=7 System Electric Greifer Digital Twin with friction joint and lamination ring
- 1 pc. 647H382 instructions for use

4 Description and function

- 8E33=7: Electric Greifer with System Hand Quick-Disconnect. Electrical connection to the forearm via coaxial bushing.
- 8E34=7: Electric Greifer with Friction Joint and Lamination Ring. Electrical connection to the forearm via cable (not illustrated).

Both types of Electric Greifers offer Digital Twin Control.

The Digital Twin Control system provides both a digital control mode (dual site) and a double channel control mode (single site) in the same terminal device. The desired control mode can be selected by inserting or removing the 13E189 Function Plug. The digital mode (dual site) is controlled by two electrodes and the double channel mode (single site) is controlled by only one electrode – either 13E68 or 13E125.

The Greifer's coaxial bushing or the cable which exits out of the Greifer are marked with a **white ring** for identification purposes.

The System Electric Greifer is delivered in digital mode (dual site) with the function plug inserted.

4.1 Quick-disconnect adapter (only with 8E33=7) (fig. 1, item 1)

The 8E33=7 Electric Greifer includes a quick-disconnect adapter that fits the MyoBock System. To remove the Electric Greifer from the prosthetic socket and replace it by a System Electric Hand with a quick-disconnect adapter, rotate the hand approximately 360°.

4.2 On/Off switch, flexion joint (fig. 1, item 2)

The electric on/off switch simultaneously activates the wrist flexion joint.

When in position B, the Electric Greifer is ready for operation and the flexion brake is released. This makes it possible for the Greifer to be passively moved approximately 45° in both directions.

When in position A, the Electric Greifer is switched off and the flexion joint is locked.

4.3 Fingertips (fig. 1, item 3)

To improve the angle of approach, the fingertips can be rotated by $\pm 60^\circ$, or they can be replaced with other fingertips having other functions (see "MyoBock Arm Components" Catalog for options).

To adjust the fingertips, open the Greifer as far as possible. Loosen the center screw until the fingertips can be removed from the gear teeth. Replace or rotate the fingertip, then make sure the gear teeth match exactly before tightening the center screw. The gripping surfaces will remain parallel until they have reached a maximum opening of 95mm.

4.4 Safety release (fig. 1, item 4)

In the event of a system malfunction, the safety release lever will disengage the Greifer.

Safety release lever positioned towards control wheel:

The Greifer will operate normally with electronic or manual control (fig. 1, item 5) of the gripping fingers.

Safety release lever positioned away from the control wheel:

The drive train will be disengaged to allow the fingers to move freely.

4.5 Manual grip force control wheel (fig. 1, item 5)

Clockwise rotation will close and counterclockwise rotation will open the fingers. In addition, the Manual Wheel provides visual feedback during the gripping process (i.e. you can see the pause in the time delayed automatic transmission).

5 Control modes

The prosthetist may choose from two different control modes by inserting or removing the 13E189 Function Plug. To reset the system so that it will recognize the new control mode, remove the battery, then re-insert it. To change the control mode, tilt the wrist flexion until the function plug can be attached or removed using tweezers (fig. 4).

5.1 Coding plug inserted: Digital Mode (Dual Site) (fig. 2, as delivered)

After the corresponding electrode signal has exceeded a certain threshold value, the Greifer opens and closes at a constant speed.

5.2 Coding plug removed: Double Channel Mode (Single Site) (fig. 3)

The Greifer is opened and closed by only one electrode.

A slow and gentle muscle signal will close the Greifer.

A quick and strong muscle signal will open the Greifer.

6 Adjusting the 13E68 (round) or 13E125 (rectangle) electrodes

The best possible function of the Digital Twin Greifer is achieved through optimal adjustment of the electrodes. This is easily possible using the 757M11=X-CHANGE MyoBoy together with the 757P23 Test Adapter.

6.1 Procedure with digital mode (Dual Site) (with function plug)

(For detailed information on the procedure, refer to MyoBoy Instructions for Use)

Select the Electric Greifer with digital control (dual site) in the MyoSoft program. Adjust each electrode so that the patient can maintain the corresponding muscle signal for approximately two seconds at a point slightly above the value ON (fig. 5).

6.2 Procedure with double channel mode (Single Site) (without function plug)

(For detailed information on the procedure, refer to MyoBoy Instructions for Use)

Select the Electric Greifer with double channel control (single site) in the MyoSoft program. Adjust the slider on the scale to 3.5. Dial 2 on the electronic relay switch (fig. 3) is already set to the corresponding value when delivered: The marks on dial 2 and on the housing of the electronic relay are identical.

Adjust the electrode so that the muscle signal reaches a value just below the threshold marked in MyoSoft by the broken line at 3.5 (fig. 6): The Greifer closes.

Overcome the threshold with a quick, strong muscle signal to open the Greifer.

If the muscle signal is too slow or too weak, adjust the threshold using Dial 2 (fig. 3):

- + for raising the threshold
- for reducing the threshold

Attention:

Allow the patient to rest during training to avoid fatigue and potentially incorrect electrode settings (i.e. excess sensitivity).

Please take care to place the electrodes so that the entire contact surfaces lie directly on intact skin. If strong interference by electrical appliances still occurs, please check the positioning of the electrode and – if need be – reposition them. If the interference cannot be eliminated, please contact an Ottobock Service Facility.

7 Battery management

Use only fully charged batteries to operate the System Electric Greifer with Digital Twin Control. An intelligent “battery manager” alerts the patient when the voltage is low: The System Electric Greifer will become slower and the gripping force will decrease. At the same time the battery is protected against harmful total discharge.

The Digital Twin Control can be operated with the 757B15 Ottobock X-ChangePack or with the 757B20/757B21 EnergyPack (Li-Ion technology). For more detailed instructions regarding the EnergyPacks, please refer to the information provided with the EnergyPack.

Attention:

As soon as the Digital Twin Control recognizes that a **fully charged** 757B20 or 757B21 Energy-Pack has been connected, it automatically switches to a special mode optimized for this state-of-the-art battery technology. Once this has happened, only EnergyPacks should be used. If a 757B15 Interchangeable Battery is used after the control has automatically switched to the new battery technology. If the user wishes to return to using Interchangeable Batteries, the control may be reset to the factory default for use with this type of batteries at an Ottobock service facility.

For operational security and reliability **only** 757B15 Ottobock Interchangeable Batteries or 757B20/757B21 Energy Packs should be used.

As delivered, the default setting is optimized for NiCad batteries.

8 Cleaning and care

Clean the product with a damp cloth and mild soap when needed. Dry the product with a soft, lint-free cloth afterwards.

Then, dry the prosthetic component with a soft cloth.

9 Technical data

Static current	mA	1
Operating temperature	°C	0-70
Opening width	mm	95
Speed	mm/s	180
Grip force	N	160
Weight (8E33=7)	g	540
Weight (8E34=7)	g	520
Power source	Ottobock X-ChangePack 757B15 - 6 V / 7,2 V 757B20 Ottobock EnergyPack - 7.2 V 757B21 Ottobock EnergyPack - 7.2 V	
Service life of the battery	2 years	
Service life of the gripping component	com-5 years	

10 Service

Since all moving parts are subject to some wear and tear, Ottobock recommends yearly service check-ups by your prosthetist. During the service visit, the entire Greifer will be inspected and any necessary adjustments will be made.

11 Warranty

The Ottobock Warranty applies only when the 8E33=7 and 8E34=7 Ottobock System Electric Greifers with Digital Twin Control are used according to the specific conditions and for the intended purposes, following all manufacturer's recommendations. The System Electric Greifer with Digital Twin Control should not be disassembled or repaired by **non-authorized personnel**. Use with component combinations that have not been tested by Ottobock will void the warranty.

12 CE conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock Healthcare Products GmbH