






## 50K13N Genu Arexa

|  |     |
|--|-----|
| <b>DE</b> Gebrauchsanweisung .....         | 5   |
| <b>EN</b> Instructions for use .....       | 10  |
| <b>FR</b> Instructions d'utilisation ..... | 15  |
| <b>IT</b> Istruzioni per l'uso .....       | 20  |
| <b>ES</b> Instrucciones de uso .....       | 25  |
| <b>PT</b> Manual de utilização .....       | 30  |
| <b>NL</b> Gebruiksaanwijzing .....         | 35  |
| <b>SV</b> Bruksanvisning .....             | 40  |
| <b>DA</b> Brugsanvisning .....             | 45  |
| <b>FI</b> Käyttöohje .....                 | 50  |
| <b>PL</b> Instrukcja użytkowania .....     | 55  |
| <b>HU</b> Használati utasítás .....        | 60  |
| <b>CS</b> Návod k použití .....            | 65  |
| <b>SK</b> Návod na používanie .....        | 70  |
| <b>BG</b> Инструкция за употреба .....     | 75  |
| <b>TR</b> Kullanma talimatı .....          | 81  |
| <b>EL</b> Οδηγίες χρήσης .....             | 86  |
| <b>RU</b> Руководство по применению .....  | 91  |
| <b>JA</b> 取扱説明書 .....                      | 97  |
| <b>ZH</b> 使用说明书 .....                      | 101 |

|  | Seite/<br>Side | Größe/<br>Size | Umfang/Circumference |           |        |           |
|--|----------------|----------------|----------------------|-----------|--------|-----------|
|  |                |                | A (cm)               | A (inch)  | B (cm) | B (inch)  |
|  | L/R            | S              | 32-37                | 12,6-14,6 | 39-46  | 15,4-18,1 |
|  | L/R            | M              | 37-41                | 14,6-16,1 | 46-53  | 18,1-20,9 |
|  | L/R            | L              | 41-46                | 16,1-18,1 | 53-61  | 20,9-24,0 |
|  | L/R            | XL             | 46-50                | 18,1-19,7 | 61-69  | 24,0-27,2 |
|  | L/R            | XXL            | 50-56                | 19,7-22,1 | 69-79  | 27,2-31,1 |

|                 |   |
|-----------------|---|
| <b>Material</b> | AL, SS, klettbarer Velours aus PA/velastic velour of PA, PU-Schaum/PU-foam, PE-Schaum/PE-foam, PA, POM, TPU, PP |
|-----------------|---|

**INFORMATION**

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### INFORMATION

Last update: 2015-03-10

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 50K13N Genu Arexa knee joint brace.

## 2 Product Description

| Item | Component (see fig. 1) | Item | Component (see fig. 1)       |
|------|------------------------|------|------------------------------|
| 1    | Strap guide loop       | 5    | Technogel® condyle pad       |
| 2    | Pad button             | 6    | Pad retainer                 |
| 3    | Clip rivet             | 7    | Extension stop (square grip) |
| 4    | Tibia pad              | 8    | Flexion stop (round grip)    |

## 3 Application

### 3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 3.2 Indications

- After ruptures of the anterior cruciate ligament (ACL) and posterior cruciate ligament (PCL), collateral ligament injuries (MCL, LCL) or after capsular ligament surgery
- Conservative, pre-operative and post-operative joint stabilisation (e.g. after meniscal fixation or meniscus replacement) with limited range of motion
- Conservative therapy of chronic complex knee joint instability

The indication must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

None known.

#### 3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries; inflammation; prominent, swollen scars; reddening and hyperthermia of the treated limb; pronounced vari-

cose veins, especially with impaired return flow; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area where the device will be applied; sensory and circulatory disorders in the legs, e. g. associated with diabetic neuropathy.

### 3.4 Mechanism of Action

The brace relieves and stabilises the knee joint.

## 4 Safety

### 4.1 Explanation of Warning Symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 4.2 General Safety Instructions

**CAUTION**

#### Contact with heat, embers, or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers, and other sources of heat.

**CAUTION**

#### Reuse on other patients and insufficient cleaning

Skin irritation, formation of eczema, or infections due to contamination with germs

- ▶ Use the product only for a single patient.
- ▶ Clean the product regularly.

**INFORMATION**

- ▶ Observe the information in the section "Information on re-use".

**CAUTION**

#### Improper molding or application

Injuries or damage to the brace due to overloading of the material and improper fit of the brace due to breakage of load-bearing components

- ▶ Do not make any improper changes to the brace.
- ▶ Always apply the brace according to the information in the instructions.

**NOTICE**

#### Contact with ointments, lotions, or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to ointments, lotions, or other products that contain oils or acids.

## 5 Handling

**INFORMATION**

- ▶ The duration of daily wear and overall period of use are generally determined by the physician.

- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any unusual changes occur (e.g. worsening of the complaint).

## 5.1 Selecting the Size

1. Measure the circumference of the knee joint (circumference A).
2. Measure the circumference of the thigh (circumference B) 15 cm above the center of the patella.
3. Determine the size of the brace (see size table).

## 5.2 Fitting and Application

### 5.2.1 Adjusting the Range of Motion (Flexion/Extension)

#### INFORMATION

The joint stops (medial and lateral) must have identical degree settings.

|  |                                       |
|--|---------------------------------------|
| Available extension stops (square grip): | 0°, 10°, 20°, 30°, 45°                |
| Available flexion stops (round grip):    | 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90° |

- 1) Remove the pad retainer (and the spacer rings if applicable) (see fig. 4).
- 2) Insert the desired extension and flexion stops into the joint (see fig. 5).
  - The stops will engage with slight resistance.
- 3) **Optional:** Snap the required number of spacer rings onto the joint (see the section "Adaptation to the knee width").
- 4) Attach the pad retainer.

### 5.2.2 Adaptation to the Knee Width

The scope of delivery includes spacer rings (thickness = **5 mm**).

This adjustment can be carried out independently on the medial and lateral joints.

- 1) Loosen the pad retainer from the joint.
- 2) Snap the required number of spacer rings onto the joint (see fig. 6).
- 3) Attach the pad retainer on the spacer rings.

### 5.2.3 Adjusting the strap height

The height of the front/rear top strap for the lower leg can be adjusted for knee joint relief in case of post-operative therapy and/or to adapt the position of the straps to the anatomical shape of the leg.

- 1) Remove the clip rivet from the front/rear, top strap for the lower leg.
- 2) Move the strap.
- 3) Attach the strap with a new clip rivet.

### 5.2.4 Applying the ACL Version

#### ⚠ CAUTION

#### Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the brace is applied properly and fits correctly.

## NOTICE

### Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear, and damage prior to each use.
- ▶ Advise the patient that the product may not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- > The brace is in the delivery condition (ACL/CI version).
- > All straps are opened.
- > The patient is sitting on the front edge of a chair. The knee is angled at 45°.
- 1) Position the brace with the middle of the joints approximately even with the top of the patella (see fig. 7).
  - Although this puts the brace higher than its final position, it settles into the correct position on its own.
- 2) Fasten the rear, top strap for the lower leg (see fig. 8).
- 3) Fasten the rear, bottom strap for the thigh (see fig. 9).
- 4) Engage the rear, top strap for the thigh on the pad button with the strap guide loop and pull it to the back.
  - The strap guide loop engages audibly.
- 5) Tighten the rear, top strap for the thigh (see fig. 10).
- 6) Engage the rear, bottom strap for the lower leg on the pad button with the strap guide loop and pull it to the back.
  - The strap guide loop engages audibly.
- 7) Tighten the rear, bottom strap for the lower leg (see fig. 11).
- 8) Push the brace towards the thigh and fasten the front strap for the thigh (see fig. 12).
- 9) Fasten the front strap for the lower leg (see fig. 13).
- 10) Tighten the brace straps so that the brace fits firmly but is comfortable to wear.
- 11) After **15 min** of active movement (e.g. slight knee bends): slightly retighten the straps in the specified sequence (see fig. 2, straps 1-6).

### 5.2.5 Conversion to CI Version

The brace is in the delivery condition (ACL or CI version). No adaptation is required.

### 5.2.6 Conversion to PCL Version

- > The brace is in the delivery condition (ACL/CI version).
- 1) Switch the Y-hook-and-loop of the front strap for the thigh (strap 2) with the Y-hook-and-loop of the rear, bottom strap for the thigh (strap 5).
- 2) Position the brace correctly (see the section "Applying the ACL version", see fig. 7).
- 3) Fasten the straps in the modified sequence (see fig. 3, straps 1-6).

## 5.3 Delivery

## NOTICE

### Incorrect setting

Damage to the brace due to overloading of the material and improper fit of the brace due to breakage of load-bearing components

- ▶ The brace may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.
- ▶ **Inform the patient.**

- ▶ Check that the brace fits the leg properly at the time of delivery.
- The centre of the patella and the centre of the joints are approximately aligned.

## 5.4 Cleaning

### INFORMATION

- ▶ Ensure proper alignment when inserting the tibia pad. The wide side of the pad must face up.

- 1) Remove the straps and pads from the brace.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at 30 °C (86°F) using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 5) If the brace comes into contact with salt water or dirt, rinse with clear water and allow to air-dry.

## 5.5 Information on re-use (applies for Germany only)

The product is approved for re-use **exclusively** in Germany. Products that are being re-used are subject to increased strain. Re-use is only permitted up to a maximum total period of use of **2 years**.

### Reconditioning

Hygienic reconditioning of the product (according to the Robert-Koch Institute guidelines) is required prior to re-use.

All components subject to pronounced wear or that come into direct contact with skin are not suitable for re-use. These components are contained in the 29K110N replacement kit.

All components suitable for re-use must be checked for proper functionality and damage prior to re-use.

### 5.5.1 Replacing the Straps and Pads

#### Straps:

- 1) The upper and lower straps can be released from the frame using the quick-release locks.
- 2) For the straps closest to the knee, remove the strap attachment caps on the outside of the frame.
- 3) Release the plastic pin from the frame, insert the new strap and secure with the cap.

#### Strap and Frame Pads

The strap and frame pads are attached with hook-and-loop closures and are easy to replace.

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### 7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorized modification of the product.

### 7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria out-

lined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

## 7.3 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

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Otto Bock HealthCare GmbH