
28U90

DE Gebrauchsanweisung	4
EN Instructions for use	7
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	15
ES Instrucciones de uso	18
PT Manual de utilização	22
NL Gebruiksaanwijzing	26
SV Bruksanvisning	29
DA Brugsanvisning	33
PL Instrukcja użytkowania	37
HU Használati utasítás	40
CS Návod k použití	44
RO Instrucțiuni de utilizare	48
BG Инструкции за употреба	52
EL Οδηγίες χρήσης	55
RU Руководство по применению	59

	Artikelnummer / Article Number	Seite / Side	Schuhgröße/ Shoe Size	Farbe/ Colour	Farbe Polster/ Colour Pad
	28U90=L35-37-0	links/left	L35-37	Transparent/Transparent	Hautfarben/Skin colour
	28U90=L35-37-7	links/left	L35-37	Schwarz/black	Schwarz/black
	28U90=R35-37-0	rechts/right	R35-37	Transparent/Transparent	Hautfarben/Skin colour
	28U90=R35-37-7	rechts/right	R35-37	Schwarz/black	Schwarz/black
	28U90=L37-39-0	links/left	L37-39	Transparent/Transparent	Hautfarben/Skin colour
	28U90=L37-39-7	links/left	L37-39	Schwarz/black	Schwarz/black
	28U90=R37-39-0	rechts/right	R37-39	Transparent/Transparent	Hautfarben/Skin colour
	28U90=R37-39-7	rechts/right	R37-39	Schwarz/black	Schwarz/black
	28U90=L39-41-0	links/left	L39-41	Transparent/Transparent	Hautfarben/Skin colour
	28U90=L39-41-7	links/left	L39-41	Schwarz/black	Schwarz/black
	28U90=R39-41-0	rechts/right	R39-41	Transparent/Transparent	Hautfarben/Skin colour
	28U90=R39-41-7	rechts/right	R39-41	Schwarz/black	Schwarz/black
	28U90=L41-44-0	links/left	L41-44	Transparent/Transparent	Hautfarben/Skin colour
	28U90=L41-44-7	links/left	L41-44	Schwarz/black	Schwarz/black
	28U90=R41-44-0	rechts/right	R41-44	Transparent/Transparent	Hautfarben/Skin colour
	28U90=R41-44-7	rechts/right	R41-44	Schwarz/black	Schwarz/black
	28U90=L44-47-0	links/left	L44-47	Transparent/Transparent	Hautfarben/Skin colour
	28U90=L44-47-7	links/left	L44-47	Schwarz/black	Schwarz/black
28U90=R44-47-0	rechts/right	R44-47	Transparent/Transparent	Hautfarben/Skin colour	
28U90=R44-47-7	rechts/right	R44-47	Schwarz/black	Schwarz/black	

Material	PA, PP
-----------------	--------

INFORMATION

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INFORMATION

Last update: 2015-03-17

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use contain important information for fitting and setting up your 28U90 ankle foot orthosis.

2 Application

2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Dorsal flexor weakness, (e.g. after stroke, traumatic brain injury, in multiple sclerosis, neuromuscular atrophy, or peroneal paralysis)

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

- Moderate to severe spasticity in the lower leg
- Moderate to severe oedema
- Leg ulcers
- Moderate to severe foot deformities

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

2.4 Further Usage Restrictions

The orthosis must not be used for the prevention of contractures (shortening) etc.

The orthosis must not be used for sports activities that involve jumping, sudden movements or rapid step sequences (e.g. basketball, badminton, riding sports). Sports activities should generally be discussed with the patient.

2.5 Effects

The ankle foot orthosis provides lifting of the foot and limits plantar flexion.

3 Safety

3.1 Explanation of Warning Symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General Safety Instructions

 CAUTION
Reuse on other patients and improper cleaning Skin irritation, formation of eczema or infections due to contamination with germs <ul style="list-style-type: none">▶ Only use the product for a single patient.▶ Clean the product regularly.
 CAUTION
Contact with heat, embers or fire Injuries (e.g. burns) caused by melting of the material <ul style="list-style-type: none">▶ Keep the product away from open flame, embers and other sources of heat.

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.

4 Handling

⚠ WARNING

Operation of motor vehicles

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

4.2 Adaptation

⚠ CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.

NOTICE

Use without suitable footwear

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.

Adapting the orthosis size:

- 1) **If necessary:** Cut or sand down the orthosis to size according to shoe size.
- 2) Then, deburr the worked edges and polish out any notches.

Adapting the orthosis:

- 1) **If necessary:** Thermoplastically re-form the orthosis by applying gentle heat to a limited extent.

NOTICE! Re-forming of the orthosis in cold conditions or at temperatures above 80°C is not permitted.

- 2) Attach the calf pad and the hook-and-loop strap (see fig. 2).

INFORMATION: The hook-and-loop strap can be fastened medially or laterally.

- 3) Once fitting is complete, reattach the hook-and-loop strap to the calf pad (see fig. 3).

Final inspection:

- ▶ Conduct trial walking with the patient. Testing on inclines, ramps and stairs is mandatory. Fine tune as required.

4.3 Application

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

1. Open the hook-and-loop strap.
2. Insert the orthosis into the shoe.
3. Put the orthosis on together with the shoe.
4. Tighten and close the hook-and-loop strap.

4.4 Cleaning

Clean the orthosis regularly:

Textile component:

- 1) Remove the textile component from the orthosis.
- 2) Close all hook-and-loop closures.
- 3) Hand wash the textile component in warm water at 30 °C with a standard mild detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

Composite component:

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH