








6A40 MagnoFlex Lock

| | |
|--|-----|
| DE Gebrauchsanweisung | 5 |
| EN Instructions for use | 15 |
| FR Instructions d'utilisation | 25 |
| IT Istruzioni per l'uso | 37 |
| ES Instrucciones de uso | 48 |
| PT Manual de utilização | 59 |
| NL Gebruiksaanwijzing | 70 |
| SV Bruksanvisning | 80 |
| DA Brugsanvisning | 90 |
| NO Bruksanvisning | 100 |
| FI Käyttöohje | 110 |
| PL Instrukcja użytkowania | 120 |
| HU Használati utasítás | 131 |
| CS Návod k použití | 142 |
| RO Instrucțiuni de utilizare | 152 |
| HR Upute za uporabu | 162 |
| SL Navodila za uporabo | 173 |
| SK Návod na používanie | 183 |
| BG Инструкция за употреба | 193 |
| TR Kullanma talimatı | 205 |
| EL Οδηγίες χρήσης | 215 |
| RU Руководство по применению | 227 |
| JA 取扱説明書 | 239 |
| ZH 使用说明书 | 249 |
| KO 사용 설명서 | 257 |

INFORMATION

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INFORMATION

Last update: 2015-04-16

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 6A40 MagnoFlex Lock is a shuttle lock with a flexible pin and a magnet integrated in the housing. It serves to secure a suitable liner in the prosthetic socket. The shuttle lock is installed in the check socket by vacuum forming and in the definitive socket by lamination. To connect to the distal prosthetic components, it is equipped with a four-hole connector. The flexible pin is installed on the liner. When the patient slides into the prosthetic socket, the magnet in the MagnoFlex Lock aligns the flexible pin with the opening. This makes it easier to don the prosthesis.

1.2 Combination Possibilities

| 6A40 MagnoFlex Lock | |
|---------------------|------------------|
| Designation | Reference number |
| Socket adapter | 4R54, 4R77 |
| Slider Plate | 6A41 |
| Exhaust Valve | 6A42 |

For additional combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

2 Application

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of Application

Approved for a body weight of up to **125 kg**.

2.3 Environmental Conditions

| Allowable environmental conditions |
|---|
| Temperature range for use: -10 °C to +60°C |
| Allowable relative humidity 0 % to 90 %, non-condensing |

| Unallowable environmental conditions |
|--|
| Mechanical vibrations or impacts |
| Perspiration, urine, fresh water, salt water, acids |
| Dust, sand, highly hygroscopic particles (e.g. talcum) |

2.4 Service Life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of Warning Symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General Safety Instructions

 **CAUTION**

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Use the prosthesis components according to the MOBIS classification (see section "Area of Application").
- ▶ **Inform the patient.**

 **CAUTION**

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

 **CAUTION**

Exceeding the service life and reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service Life").
- ▶ Only use the product for a single patient.
- ▶ **Inform the patient.**

 **CAUTION**

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.

- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).
- ▶ **Inform the patient.**

⚠ CAUTION

Product contains a magnet

Influence on or damage to devices or items due to a strong magnetic field

- ▶ Keep the product away from devices and items that are sensitive to magnetic fields.

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of Delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering according to their respective symbols as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

| 6A40 MagnoFlex Lock | | | | |
|---------------------|------|----------|----------------------|------------------|
| Fig. | Item | Quantity | Designation | Reference number |
| - | - | 1 (■) | Instructions for Use | 647G931 |
| 1 | 1 | 1 | MagnoFlex Lock | - |

| 6A40 MagnoFlex Lock | | | | |
|---------------------|------|----------|------------------------------|-------------------|
| Fig. | Item | Quantity | Designation | Reference number |
| | | 1 (▲) | Set Screw (for ratchet unit) | 506G3=M3X5-"Niro" |
| 1 | 2 | 1 (■) | Lamination Adapter | 6A43 |
| 1 | 3 | 1 (■) | Flexible Pin | 6Y13=F1 |
| 1 | 4 | 1 (▲) | Screw | 503S3 |
| 1 | 5 | 1 (■) | Vacuum Forming Dummy | – |
| 1 | 6 | 4 (▲) | Snap Bushing | – |
| 1 | 7 | 4 (▲) | Cap Screw | 501T28=M6X8 |
| 1 | 8 | 4 (▲) | Flat Head Screw | 501S41=M6X45 |
| 1 | 8 | 4 (▲) | Flat Head Screw | 501S41=M6X50 |
| 1 | 9 | 1 (■) | Ratchet Unit | 5A52=K |

5 Preparation for Use

CAUTION

Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- ▶ Observe the alignment and assembly instructions.

INFORMATION

Not all of the materials listed may be available in your country.

- ▶ If the materials are not available, please contact the local branch of the manufacturer to obtain information on alternative materials.

5.1 Fabricating the Prosthetic Socket

5.1.1 Preparing for Socket Fabrication

> **Required tools and materials:**

6A43 Lamination Adapter, 501T28=M6X8 Cap Screws, Vacuum forming dummy, 503S3 Screw, 633W8 Wax, 636K8=* Plastaband

- 1) Isolate the cap screws (thread and head from the outside) with wax.
- 2) Screw the cap screws into the lamination adapter (see fig. 2).
- 3) Fill the heads of the cap screws with Plastaband (see fig. 3).
- 4) Position the lamination adapter on the plaster positive along the longitudinal axis of the residual limb.
- 5) Place the vacuum forming dummy into the round opening of the lamination adapter and secure the two parts to the plaster positive with the flat head screw (see fig. 4).
- 6) Isolate the head of the flat head screw with Plastaband.

- 7) **If there is an undercut between the lamination adapter and plaster positive:** Fill the undercut with wet plaster.

5.1.2 Fabricating a Check Socket

> **Required tools and materials:**

755T4=360 Outer Ring, 755X104=360 Vacuum Tube (with vacuum sealing ring), 755E9 Vacuum Pump, 756D2, 756B12=* or 756B20=* Oscillating Saw, 616T52=* ThermoLyn rigid or 616T83=* ThermoLyn clear, string

- 1) Prepare for fabrication of the prosthetic socket (see Page 19).
- 2) Position two pieces of string crosswise on the lamination adapter (see fig. 5). The pieces of string form air channels that make it easier to mould the vacuum forming material to the contours.
- 3) Complete the vacuum forming process with ThermoLyn.
- 4) Finish the prosthetic socket (see Page 22).

5.1.3 Laminating a Transtibial Socket

> **Required tools and materials:**

99B81=70X19X5 and 99B81=100X19X5 PVA Bag, 623T3=8 or 623T3=10 Perlon Stockinette, 81A1=8 or 81A1=10 Tube Sock, 616B1=25x* Carbon Fibre Cloth Strap, 616G2 Carbon UD Hose, 616G15 Carbon Fibre Woven Hose, 617H119 Orthocryl Lamination Resin 80:20 PRO, string

- 1) Prepare for fabrication of the prosthetic socket (see Page 19).
- 2) Soak the shorter PVA bag and pull it over the plaster positive.
- 3) Pull a tube sock over the plaster positive.
- 4) Place a layer of carbon fibre cloth strap circularly around the MPT (middle patellar tendon) point (see fig. 6).
- 5) Position a layer of carbon fibre cloth strap medially and laterally from the shuttle lock up to the condyles (see fig. 7).
- 6) Place a layer of carbon fibre cloth strap anteriorly and posteriorly from the shuttle lock to the circular carbon fibre cloth strap (see fig. 8).
- 7) Pull a tube sock over the plaster positive.
- 8) Loop a piece of string around the lamination adapter in the groove and tie tightly (see fig. 9).
- 9) Cut a piece of carbon UD hose (1.5 times the length of the plaster positive).
- 10) Pull the carbon UD hose over the plaster positive up to the socket edge. Tie the upper section of the carbon UD hose and fold the remaining length over the plaster positive.

- 11) Loop a piece of string around the lamination adapter in the groove and tie tightly (see fig. 10).
- 12) Pull a tube sock over the plaster positive.
- 13) Cut a piece of carbon fibre woven hose (1.3 times the length of the plaster positive).
- 14) Pull the carbon fibre woven hose over the plaster positive to the socket edge. Tie the upper section of the carbon fibre woven hose and fold the remaining length over the plaster positive (see fig. 11).
- 15) Cut a piece of Perlon stockinette (twice the length of the plaster positive).
- 16) Pull the Perlon stockinette over the plaster positive up to the socket edge. Tie the upper section of the Perlon stockinette and fold the remaining length over the plaster positive.
- 17) Soak the longer PVA bag and pull it over the plaster positive.
- 18) Perform the casting process with Orthocryl.
- 19) Finish the prosthetic socket (see Page 22).

5.1.4 Laminating a Transfemoral Socket

> Required tools and materials:

99B81=100X26X5 and 99B81=130X26X5 PVA Bag, 623T3=12 or 623T3=15 Perlon Stockinette, 81A1=12 or 81A1=15 Tube Sock, 616B1=50x* Carbon Fibre Cloth Strap, 616G2 Carbon UD Hose, 616G15 Carbon Fibre Woven Hose, 617H119 Orthocryl Lamination Resin 80:20 PRO, string

- 1) Prepare for fabrication of the prosthetic socket (see Page 19).
- 2) Soak the shorter PVA bag and pull it over the plaster positive.
- 3) Pull a layer of Perlon stockinette over the plaster positive.
- 4) Place a layer of carbon fibre cloth strap circularly 3 cm below the perineum.
- 5) Place a layer of carbon fibre cloth strap medially and laterally as well as anteriorly and posteriorly from the shuttle lock up to the circular carbon fibre cloth strap.
- 6) Place a layer of carbon fibre cloth strap circularly 3 cm below the perineum.
- 7) Cut two pieces of carbon fibre cloth strap (length approx. 20 cm).
- 8) Form the carbon fibre cloth straps into Vs by folding them in the middle and position them in the area of the tuberosity support.
- 9) Pull a tube sock over the plaster positive.
- 10) Loop a piece of string around the lamination adapter in the groove and tie tightly.

- 11) Cut a piece of carbon UD hose (1.5 times the length of the plaster positive).
- 12) Pull the carbon UD hose over the plaster positive up to the socket edge. Tie the upper section of the carbon UD hose and fold the remaining length over the plaster positive.
- 13) Loop a piece of string around the lamination adapter in the groove and tie tightly.
- 14) Pull a tube sock over the plaster positive.
- 15) Cut a piece of carbon fibre woven hose (1.3 times the length of the plaster positive).
- 16) Pull the carbon fibre woven hose over the plaster positive to the socket edge. Tie the upper section of the carbon fibre woven tube and fold the remaining length over the plaster positive.
- 17) Cut a piece of Perlon stockinette (twice the length of the plaster positive).
- 18) Pull the Perlon stockinette over the plaster positive up to the socket edge. Tie the upper section of the Perlon stockinette and fold the remaining length over the plaster positive.
- 19) Soak the longer PVA bag and pull it over the plaster positive.
- 20) Perform the casting process with Orthocryl.
- 21) Finish the prosthetic socket (see Page 22).

5.1.5 Finishing the Prosthetic Socket

> Required materials:

Snap Bushings

- 1) Mark and trim the contour of the prosthetic socket.
- 2) Sand the distal end of the socket down to the vacuum forming dummy and the heads of the cap screws (see fig. 12).
- 3) On a level surface, check whether the sanded surface is flat. Rework as needed.
- 4) Remove the vacuum forming dummy and the cap screws from the prosthetic socket.
- 5) Remove the prosthetic socket from the plaster positive.
- 6) Sand the contour of the prosthetic socket.
- 7) Insert the snap bushings into the openings from which the cap screws were removed (see fig. 13).

5.2 Installing the Shuttle Lock and Socket Adapter

CAUTION

Using liquid thread lock

Falling due to screw connections ripping out

- ▶ Do not use liquid thread lock (e.g. Loctite®). The thread lock damages the plastic material.

INFORMATION

The product contains a magnet

- ▶ Maintain a safety distance of at least **12 cm** from a pacemaker.

> Required tools and materials:

Socket Adapter, 501S41=M6X45 Flat Head Screws, 506G3=M3X5-"Niro" Set Screw, 710D4 Torque Wrench, Torque Wrench (adjustable to 0.5 Nm)

- 1) Align the shuttle lock such that the magnet is aligned proximally and the ratchet unit medially or laterally.
- 2) Set the shuttle lock on to the prosthetic socket (see fig. 14).
- 3) **Optional:** Position the 6A41 Slide Plate on the shuttle lock.
- 4) Position the chosen socket adapter on the shuttle lock.
- 5) Select the flat head screws for the socket adapter:
For 4R54: 501S41=M6x45
For 4R77: 501S41=M6x50
- 6) Insert the flat head screws through the socket adapter and the shuttle lock into the lamination adapter and tighten them (see fig. 15) (tightening torque: **12 Nm**).
- 7) Screw the ratchet unit (see fig. 16, item 1) into the shuttle lock and tighten with the torque wrench (tightening torque: **10 Nm**).
- 8) Secure the ratchet unit with the set screw (see fig. 16, item 2) (tightening torque: **0.5 Nm**).

5.3 Mounting the pin

CAUTION

Installation of the pin in a liner that is not approved

Falling due to loosening of the screw connection

- ▶ Only install the pin in liners with a metal thread and observe the specified tightening torque.
- ▶ **Inform the patient.**

The pin to connect the liner to the shuttle lock is screwed into the thread at the distal end of the liner.

> **Required tools and materials:**

Torque Wrench (adjustable to 3 Nm), 636K13 Loctite®

1) Apply Loctite® to the thread.

2) **NOTICE! Do not shorten the pin.**

Screw the pin into the thread and tighten it (see fig. 17) (tightening torque: **3 Nm**).

6 Use

CAUTION

Pin not fully inserted into the shuttle lock

Risk of falling if the pin is not engaged in the shuttle lock

- ▶ Guide the pin into the shuttle lock according to the instructions.
- ▶ Always verify that the pin is engaged in the shuttle lock before using the prosthesis.
- ▶ **Inform the patient.**

The liner and shuttle lock are connected to each other with the pin. The pin engages in the shuttle lock and holds the liner in place.

Putting on the Liner

- 1) Check the pin and liner for damage. The pin must be flexible and always spring back to the original position.
- 2) **NOTICE! Misaligning the pin leads to damage during insertion in the shuttle lock.**
Apply the liner on the end of the residual limb so that the pin is aligned with the longitudinal axis of the residual limb.
- 3) Roll the liner over the residual limb.
- 4) Check the fit of the liner and the alignment of the pin.

Donning the Prosthetic Socket

- 1) With the liner, slide into the prosthetic socket until the pin slides into the opening of the lamination adapter.
- 2) **INFORMATION: While the pin is being inserted, clicking noises caused by the magnet and the ratchet unit can be heard. These noises do not indicate whether the pin is fully engaged in the ratchet unit.**
Insert the pin fully into the shuttle lock.
- 3) Before using the prosthesis, verify that the pin is fully engaged in the shuttle lock.

Doffing the Prosthetic Socket

- ▶ Press and hold the release button on the ratchet unit and pull the residual limb with the liner out of the prosthetic socket.

7 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9 Technical data

| | |
|------------------------------|--------------------------------------|
| Reference number | 6A40 |
| Weight [g] | 225 |
| System height [mm] | 37 |
| Material | High-performance polyamide and metal |
| Max. body weight [kg] | 125 |

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Otto Bock HealthCare GmbH