## 50S1 Dyna Ankle

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for Use .....	6
<b>FR</b> Instructions d'utilisation .....	9
<b>IT</b> Istruzioni per l'uso .....	12
<b>ES</b> Instrucciones de uso .....	16
<b>PT</b> Manual de utilização .....	19
<b>NL</b> Gebruiksaanwijzing .....	22
<b>SV</b> Bruksanvisning .....	26
<b>FI</b> Käyttöohje .....	29
<b>PL</b> Instrukcja użytkowania .....	32
<b>HU</b> Használati utasítás .....	35
<b>CS</b> Návod k použití .....	39
<b>EL</b> Οδηγίες χρήσης .....	42
<b>RU</b> Руководство по применению .....	46

## INFORMATION

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## INFORMATION

Last update: 2015-04-23

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

## 1 Indications for use

The dynamic ankle orthosis 50S1 Dyna Ankle is to be used **exclusively** for the orthotic treatment of acute ankle sprains and ruptures of lateral collateral ligament.

## 2 Indications and function

Acute sprains and ruptures of lateral collateral ligament of the ankle. Indication is given by the physician.

The Dyna Ankle orthosis places the foot in moderate pronation and dorsiflexion. The dynamic restraint limits plantar flexion and supination while preserving the range of motion in pronation and dorsiflexion. Positioning the orthosis with the dynamic strap provides relief of the anterior talofibular ligament.

The trimlines and straps are positioned to leave the outer ankle free, avoiding pressure points, thus allowing immediate application of the orthosis following injury. Heat or ice may be applied with the orthosis in place. It should serve its initial function as a post surgical immobilizer during the first two weeks. Then, subject to physician approval, the orthosis may be worn with a shoe during weight-bearing activities.

The Dyna Ankle is generally worn both day and night for stability for 6 weeks and is removed only for bathing.

The orthosis is waterproof and easy to clean and care for.

## 3 Size selection of the orthosis

The orthosis is selected according to the shoe size:

Article no.	Side	Shoe size	American Shoe Size	
			Men's	Women's
50S1=L35 – 37	left	35 – 37	4 – 5	4 – 6
50S1=R35 – 37	right	35 – 37	4 – 5	4 – 6
50S1=L37 – 39	left	37 – 39	5 – 7	6 – 8
50S1=R37 – 39	right	37 – 39	5 – 7	6 – 8
50S1=L39 – 41	left	39 – 41	7 – 9	8 – 10
50S1=R39 – 41	right	39 – 41	7 – 9	8 – 10
50S1=L41 – 44	left	41 – 44	9 – 12	10 – 13
50S1=R41 – 44	right	41 – 44	9 – 12	10 – 13

## 4 Material and construction

The 50S1 dynamic ankle orthosis consists of polypropylene and is ready for immediate trial fitting. An elastic strap runs medially from the lateral foot edge to the calf and across the top of the foot. The torsion cut at the distal end of the brace leaves the ball of the big toe free, thus lifting the lateral forefoot section.

## 5 Trial fitting, application and delivery of the orthosis

### *Apply the orthosis as follows*

Place the foot in the plastic shell, sliding the heel to the back.

Pull the lower hook & loop strap from the inner side of the foot across the top of the foot through the outer plastic guide, then tighten the strap so that the heel will snugly fit to the shell.

Pull the upper hook & loop strap through the plastic loop and fasten it (fig. 1).

Fasten the hook & loop strap on the forefoot.

Slightly stretch the elastic functional strap and fasten it (fig. 2).

### **Attention!**

The orthosis may be worn in the shoe during weight-bearing activities subject to physician approval. Advise the patient of the risk of slipping when using the orthosis without shoe – especially on smooth surface. An anti-slip pimple rubber is included and can be mounted if need be.

### *Adapting by the orthopedic technician/orthotist*

The injection molded polypropylene sheet may be moderately modified (do not heat to temperatures exceeding 194 °F/90 °C). If it is necessary to modify the trimlines, the edges must be smoothed carefully since the material has minimal tear resistance when notched. The supplied straps are extra long and should be shortened to an appropriate length during fitting.

If pressure develops at the ankle or shin, the pads that are supplied with the orthosis may be glued to the strap (figs. 3 and 4).

If necessary, you can provide the anti-slip properties by attaching the included pimple rubber in the heel area of the sole of the orthosis. For this purpose, clean the adhesion surface area with a fine sanding paper and attach the pimple rubber with double-sided adhesive tape included.

To ensure the effectiveness of this orthosis, it must always be worn snugly.

In addition to proper application and regular check-ups with adjustments, patient cooperation plays an important role. The orthosis must be must worn as prescribed.

Since the orthosis is water-resistant, showers are possible while wearing it.

It is very important that the patient understand the important function of dynamic restraints.

A Patient Information sheet (646D36) is available for this purpose.

## 6 Side effects

### Attention!

All aids – supports and orthoses – applied too tightly to the surface of the body may result in local pressure spots or, in some cases, even restrict the adjacent vessels or nerves.

## 7 Contraindications

In patients with history of deep venous thrombosis in lower limb, the consultation with the treating physician is required. As the calf muscle pump is not functioning due to the use of the orthosis, a medicamentous thrombosis prophylaxis may be required.

### Attention!

In case of the following clinical pictures, the application and use of this type of aid is indicated only after consultation with the physician.

1. Skin diseases/injuries in body regions with fittings, especially in case of inflammatory symptoms. This also applies to prominent scars that are swollen, reddened and feel excessively warm when touched.
2. Lymphatic flow disorders – as well as indistinct soft-tissue swellings not in contact with the applied aid.
3. Sensory disorders - lower leg/foot

## 8 Instructions for use and care

The orthosis must be used in strict conformity with the indication. The patient's acceptance and cooperation play an important part. It is essential to instruct the patient in the proper use of the orthosis! The physician's instructions must be strictly observed.

The daily duration of use and total period of application must be determined by the attending physician.

Patients should be instructed to consult a physician immediately if they experience unusual changes (such as an increase in pain).

### Attention!

The orthosis is not flame-resistant. Patients should be instructed to keep the product away from flames or other heat sources.

Before each use, the orthosis must be checked for functional reliability and for possible wear or damage. If a hook and loop closure or any other part of the orthosis show any signs of wear, the orthosis should no longer be used.

### Attention!

The patient should be instructed not to wear the orthosis during nuclear spin tomography examination.

The product should not come into contact with grease or acidic agents, unguents and lotions.

Improper changes to the product are not permitted.

You can remove the pad and Velcro straps and wash the orthosis by hand at 30°C using a standard mild detergent. We recommend closing the hook & loop closures before washing to prevent any damage to them. To avoid skin irritations, rinse the orthosis thoroughly so that all soap

residue is removed. Do not expose to the sun and do not use oven or radiator heat for drying. Let the support air dry.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH