

1E56, 2F20


DE Gebrauchsanweisung	5
EN Instructions for use	10
FR Instructions d'utilisation	15
IT Istruzioni per l'uso	21
ES Instrucciones de uso	26
PT Manual de utilização	32
NL Gebruiksaanwijzing	37
SV Bruksanvisning	43
DA Brugsanvisning	48
NO Bruksanvisning	53
FI Käyttöohje	58
PL Instrukcja użytkowania	63
HU Használati utasítás	69
CS Návod k použití	74
RO Instrucțiuni de utilizare	79
HR Upute za uporabu	85
TR Kullanma talimatı	90
EL Οδηγίες χρήσης	95
RU Руководство по применению	101
JA 取扱説明書	107
ZH 使用说明书	112

	■	●	▲
DE	Einzelteile	Einzelteile-Pack	Mindestmenge
EN	Single components	Single-Component Pack	Minimum order quantity
FR	Pièces à l'unité	Kit de pièces de rechange	Quantité minimum
IT	Singoli componenti	Imballo dei singoli componenti	Quantità minima
ES	Componentes	Kit componentes	Cantidad mínima
PT	Componentes	Conjunto de peças de reposo	Quantidade mínima
NL	Onderdelen	Onderdelenpakket	Minimikvantitet
SV	Separata delar	Service-Set	Minimum aantal
DA	Komponenter	Komponentpakke	Mindste mængde
NO	Enkle komponenter	Enkeltkomponentpakke	Minste ordremengde
PL	Podzespoły	Zestaw naprawczy	Ilość minimalna
HU	Alkatrészek	Alkatrész-csomag	Minimum mennyiség
CS	Jednotlivé díly	Sada jednotlivých dílů	Minimální množství pro objednání
RO	Componente individuale	Set componente individuale	Cantitate minimă
HR	Pojedinačni dijelovi	Pakiranje pojedinačnih dijelova	Minimalna količina
TR	Münferit parçalar	Münferit parça ambalajı	Minimum miktar
EL	Μεμονωμένα εξαρτήματα	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα
RU	Отдельные детали	Комплект деталей	Минимальное количество
JA	単一成分パック	単一成分パック	最低発注量
ZH	单个部件	维修组件	最低起订量

2

Fußmitte
Middle of Foot
Milleu du pied
Centro del piede
Centro de pie
Meio do Pé
Midden voet
Fotens mitt
Fodmitte
Środek stopy
Jalkaterän keskikohta
Lábfej-közép
Mijlocul labei
Sredina stopala
Střed chodidla
Μέσο πέλματος
Середина стопы
Lütfen çeviriniz
足部の中心
最低発注量

Aufbaulinie
Alignment Reference Line
Ligne de montage
Linea di allineamento
Linea de alineamiento
Linha de montagem
Referenlinie
Opbouwlijn
Opbygningslinie
Linia osiowania
Asennuslinja
Felépítővonal
Stavebni linie
Linia de referință a aliniamentului
Linija postavljanja
Linija uravnoteživanja
Γραμμή ευθυγράμμισης
Ось сборки
対線参考線
アライメント
基準線

15 Nm / 11 ft.lbs.

13 mm

22-25: 25-30 mm; 26-28: 30 mm; 29-31: 30-35 mm

3

INFORMATION

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INFORMATION

Last update: 2015-04-17

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

Explanation of symbols

⚠ CAUTION Warnings regarding possible risks of accident or injury.

📄 NOTICE Warnings regarding possible technical damage.

ℹ INFORMATION Additional information on the fitting / use.

1 Single components (Fig. 1)

1.1 Single component ■ (for reordering)

(1) SL = Spectra sock Protective sock

1.2 Single-component pack ● (for reordering)

(2) 2F20=* Heel wedges

For foot size 22–25 cm: 2F20 = 22–25

For foot size 26–31 cm: 2F20 = 26–31

1.3 Accessories

INFORMATION

The 1E56 Axtion must be used in conjunction with the 2C5 Footshell (inclusive connection cap). The 2C5 Footshell is accessories and not included. It must be ordered separately as an option.

(3) 2C5=Footshell (inclusive connection cap)

Ordering examples: 2C5=L25/4, 2C5=R27/15

Sides	left (L), right (R)									
Sizes (cm)	22	23	24	25	26	27	28	29	30	31
Colors	beige (4), light brown (15)									

(4) 2C100 Tool for Removing the Footshell

2 Description

2.1 Intended use

The 1E56 Axtion Foot is to be **exclusively** used for the prosthetic fitting of amputations of the lower limb.

2.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System**:

Recommended for amputees with **mobility level 3 and 4**

(unrestricted outdoor walkers and unrestricted outdoor walker with particularly high demands)

Approved for a patient weight of max. 125 kg (275 lbs).

2.3 Safety instructions

Please forward the following safety instructions to your patients:

CAUTION

Risk of injury as a result of excessive strain. Excessive impact loading, such as running or jumping activities, may weaken the spring elements.

To avoid the risk of accident and possible injury to the patient, the foot must no longer be used after a noticeable loss of function has occurred.

This noticeable function loss may manifest itself in reduced forefoot resistance.

Measure: Have the prosthetic foot examined and, if necessary, replaced at a specialist workshop.

NOTICE

Damage due to incorrect environmental conditions. Prosthetic components must not be exposed to

- corrosive elements (such as fresh water, salt water, and acids) or
- abrasive substances (e. g. sand).

The footshell material is not UV- or water-resistant.

Do not grind the surface of this prosthetic foot, as this can lead to premature wear.

Using this medical product under such conditions will render all claims against Otto Bock HealthCare null and void.

2.4 Function

This 1E56 Axtion prosthetic foot has been designed for use in endoskeletal lower limb prostheses for moderately to highly active individuals.

It is especially suitable for amputees with very long residual limb (low clearance). Due to the special combination of carbon and polyurethane, the prosthesis foot absorbs the heel impact. Supported by the dynamic heel element, this construction enables a soft transition from heel impact to full ground contact.

The optional use of heel wedges allows the heel characteristics of the foot to be adjusted and customized to the patient's requirements.

The 1E56 Axition mainly consists of the following components (Fig. 2):

- (1) Carbon base spring
- (2) Polyurethane connecting element
- (3) Dynamic heel element of the carbon base spring
- (4) Upper carbon spring
- (5) Pyramid adapter with elastomer wedge
- (6) Footshell (accessory)
- (7) Threaded plate
- (8) Heel wedge (optional)

3 Technical data

Heel height (mm)	13 +/- 5									
Sizes (cm)	22	23	24	25	26	27	28	29	30	31
System height (mm)	34	35	35	35	35	36	36	36	36	37
Weight* (~ g)	310	315	320	325	355	360	370	380	395	410
max. body weight (kg)	125 (275 lbs)									
Mobility grade	3 + 4									

* without footshell

4 Handling

4.1 Initial set-up

4.1.1 Footshell application/removal

NOTICE

Risk of damage as a result of improper use. The foot is designed to be used with the footshell. Please avoid frequent removal of the foot from, and insertion into, the footshell as this can damage the carbon fiber elements.

NOTICE

Damage caused by improper footshell assembly. Please do not use a screwdriver to insert/remove the carbon foot into/from the footshell as this can cause damage to the carbon elements and/or foot component. We recommend using e.g. the 2C100 tool (accessories).

NOTICE

Risk of damage as a result of improper use. The footshell must not be dusted by talcum powder as it would increase the abrasion. Use silicone spray if need be. Replace the footshell if it becomes damaged or worn.

Inserting the foot into the footshell:

- Before inserting the foot into the footshell, attach the tube adapter to the prosthetic foot first. This facilitates the procedure.
- Pull a Spectra sock over the prosthetic foot.
- Insert the prosthetic foot into the footshell using e.g. the 2C100 tool.

Removing the footshell:

- To remove the footshell, insert the 2C100 tool under the heel and carefully pry out the footshell. Now the foot can be easily removed from the footshell.

4.2 Alignment

4.2.1 Bench alignment

Recommendations on bench alignment (foot with footshell):

Foot size (cm)	Placement of mid-foot anterior to alignment line (mm)	Heel height (mm)
22–25	25–30	13 +/- 5
26–28	30	
29–31	30–35	

4.2.2 Static alignment optimization

- If necessary, the alignment recommendations (TF modular leg prostheses: 646F219*, TT modular leg prostheses: 646F336*) may be requested from Ottobock.
- Use the L.A.S.A.R. Posture for static alignment optimisation of the prosthesis, if available.

4.2.3 Dynamic alignment optimization

Adapt the prosthetic alignment in the frontal plane (ML) and in the sagittal plane (AP) by angle adjustments and/or sliding adjustments, in order to ensure correct heel contact, easy roll-over, and optimal shifting of weight onto the contralateral side. For transtibial amputees please make sure that there is physiological knee flexion during stance phase.

4.3 Cosmetic foam cover

The Cosmesis Connection Cap 2C10 serves as a detachable connection between the cosmetic foam cover and the footshell. It snaps into the proximal brim of the footshell.

Please observe the following recommendations:

- Cut the foam cover to length, allowing for distal compression about 0.4 in (10 mm) for transtibial and about 1.6 in (40 mm) for transfemoral covers.
- Pull the foam cover over the prosthesis, mount the foot and mark the outer contour of the connection cap which is attached to the footshell.
- Remove the connection cap and clean with 634A58 Isopropyl Alcohol.
- Apply 636N9 Contact Adhesive or 636W17 Synthetic Adhesive to the distal cut surface of the foam cover and to the connection cap and glue the parts together.
- Finish the outer shape of the foam cover in usual way with the foot installed. Any compression of the foam cover by cosmetic stockings or SuperSkin needs to be taken into account.

4.4 Using the optional heel wedges

CAUTION

Risk of falling due to change of the knee dynamics when using heel wedges. The use of heel wedges can facilitate the stance phase flexion when bearing weight on the prosthesis. Therefore, particular caution is advised during the dynamic alignment optimization – especially when the knee joints without stance phase control is used.

4.4.1 Heel wedge application

Two heel wedges are included to each 1E56 Axtion prosthetic foot:

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Should the patient feel the heel strike too soft or the knee be in hyperextension during the mid-stance phase, the heel can be stiffened by using one of the two heel wedges.

4.4.2 Inserting a heel wedge (Fig. 4)

INFORMATION

Ottobock recommends to test the foot first without using a heel wedge.

- Select a heel wedge of a proper stiffness.
- Position the guide tongue proximally and so it is pointing to the heel.
- Apply a few drops of clean water on the heel wedge to facilitate the insertion.
- Slide the heel wedge from the posterior side into the gap between the carbon base spring (1) and threaded plate (7) until the guide tongue locks into the recess of the threaded plate (7).
- The heel wedge is fully seated when the convex side of the heel wedge hits the concave surface of the connecting element (2) and the guide tongue locks into the recess of the threaded plate (7).
- The heel wedges are designed so that no glue fixation of heel wedge is required.

4.4.3 Removing a heel wedge

- Apply a few drops of clean water on the heel wedge to facilitate the removal.
- Carefully remove the heel wedge from its position from the medial or lateral side. The loosened heel wedge can be removed with fingers.

4.5 Adjustments and final assembly

Because of the high forces in the adapter section, it is indispensable to tighten the adjustment screws of the tube adapter with a torque of **15 Nm**

Use 710D1 Torque Wrench. Once you have attained the final setting for the prosthesis, use 636K13 Loctite® to secure the set screws.

4.6 Maintenance instructions

- Foot assembly should be inspected after first 30 days (4 weeks) of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Footshell may require replacement if wear is excessive.

INFORMATION

This component has been tested according to ISO 10328 standard for two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years.

We recommend carrying out regular yearly safety checks.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Loctite is a registered trademark of Loctite Corp.
MOBIS is a registered trademark of Otto Bock HealthCare GmbH.
Axtion is a registered trademark of Otto Bock HealthCare LP.

The 1E56 Axtion is covered by the following patents:

China: ZL200580031552.8
Canada: Patented CA 2580108 (2011)
USA: Patent US 8,092,550
European Patent: EP 1 510 191 in AT, DE, ES, FR, IT, NL, SE

Patent pending as European Patent Application

Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt/Germany
T +49 5527 848-0 · F +49 5527 72330
healthcare@ottobock.de · www.ottobock.com

Otto Bock HealthCare LP
3820 West Great Lakes Drive
Salt Lake City, UT 84120 · USA
T +1 800 328 4058 · F +1 800 655 4963

Ottobock has a certified Quality Management System in accordance with ISO 13485.