## Probedo 3R31, 3R31=ST

<b>DE</b> Gebrauchsanweisung .....	25
<b>EN</b> Instructions for use .....	37
<b>FR</b> Instructions d'utilisation .....	49
<b>IT</b> Istruzioni per l'uso .....	61
<b>ES</b> Instrucciones de uso .....	73
<b>PT</b> Manual de utilização .....	85
<b>NL</b> Gebruiksaanwijzing .....	97
<b>SV</b> Bruksanvisning .....	109
<b>DA</b> Brugsanvisning .....	121
<b>NO</b> Bruksanvisning .....	133
<b>FI</b> Käyttöohje .....	145
<b>PL</b> Instrukcja użytkowania .....	157
<b>HU</b> Használati utasítás .....	169
<b>CS</b> Návod k použití .....	181
<b>HR</b> Upute za uporabu .....	193
<b>RU</b> Руководство по применению .....	205
<b>ZH</b> 使用说明 .....	219

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## Table of figures

EN

1	1 Notes regarding the document .....	4
2	1.2 Meanings of pictograms .....	4
3	2 Product description .....	5
4	3.2 Area of application (MOBIS) .....	5
5	3.2 Area of application (Amputation height/Side of amputation).....	6
6	3.3 Combination possibilities.....	6
7	3.3 Combination possibilities (MOBIS) .....	7
8	5 Scope of delivery (1/2) .....	8
9	5 Scope of delivery (2/2) .....	9
10	6.2 Bench Alignment of the prosthesis (Sagittal plane).....	10
11	6.2 Bench Alignment of the prosthesis (Frontal plane).....	10
12	6.2 Bench Alignment of the prosthesis (Alignment apparatus) .....	11
13	6.2.1 Adapting the product to the side of amputation .....	11
14	6.2.2 Positioning the prosthetic foot .....	12
15	6.2.3 Installing the distal adapters .....	12
16	6.2.4 Installing the prosthetic knee joint .....	13
17	6.2.5 Installing the proximal adapters.....	13
18	6.2.6 Preparing the prosthetic socket for the lock slide .....	14
19	6.2.7 Marking the prosthetic socket.....	14
20	6.2.8 Positioning the prosthetic socket.....	15
21	6.2.9 Shortening the tube adapter .....	16
22	6.2.10 Installing the lock slide and the lock release cable .....	17
23	6.2.11 Preventing damage in maximum flexion .....	17
24	6.3.1 Checking the static alignment.....	18
25	6.3.2 Optimizing the static alignment (Sagittal plane).....	18
26	6.3.2 Optimizing the static alignment (Frontal plane).....	19
27	6.4.1 Resetting to factory settings .....	19
28	6.4.2 Training to sit down and to stand up .....	20
29	6.4.3 Training to walk .....	21
30	6.5 Attaching a cosmetic cover .....	22
31	6.6 Finishing the prosthesis.....	22
32	12 Technical data .....	23

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Please read this document carefully before using the product.

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<b>1</b>	<b>Notes regarding the document</b> .....	<b>39</b>
1.1	Explanation of Warning Symbols .....	39
1.2	Meanings of pictograms .....	39
<b>2</b>	<b>Product Description</b> .....	<b>39</b>
<b>3</b>	<b>Application</b> .....	<b>39</b>
3.1	Indications for use .....	39
3.2	Area of Application .....	40
3.3	Combination Possibilities .....	40
3.4	Environmental Conditions .....	40
3.5	Reuse and service life .....	41
<b>4</b>	<b>General Safety Instructions</b> .....	<b>41</b>
<b>5</b>	<b>Scope of Delivery</b> .....	<b>41</b>
<b>6</b>	<b>Preparation for Use</b> .....	<b>42</b>
6.1	Information on Fabrication of a Prosthesis .....	42
6.2	Bench alignment of the prosthesis .....	42
6.2.1	Adapting the product to the amputation side .....	42
6.2.2	Positioning the prosthetic foot .....	42
6.2.3	Mounting the distal adapters .....	43
6.2.4	Mounting the prosthetic knee joint .....	43
6.2.5	Mounting the proximal adapters.....	43
6.2.6	Preparing the prosthetic socket for the lock slide.....	43
6.2.7	Marking the prosthetic socket .....	43
6.2.8	Positioning the prosthetic socket .....	43
6.2.9	Shortening the tube adapter.....	44
6.2.10	Mounting the lock slide and lock release cable.....	44
6.2.11	Preventing damage at maximum flexion .....	44
6.3	Optimising the static alignment .....	44
6.3.1	Checking the static alignment .....	44
6.3.2	Optimising the static alignment .....	44
6.4	Optimising during dynamic trial fitting .....	44
6.4.1	Resetting to factory defaults .....	44
6.4.2	Practicing sitting down and standing up .....	45
6.4.3	Practicing walking .....	45
6.5	Attaching the cosmetic cover .....	45
6.6	Finishing the prosthesis .....	45
<b>7</b>	<b>Use</b> .....	<b>46</b>
7.1	Information for Use .....	46
7.2	Using the Lock Slide.....	46
<b>8</b>	<b>Cleaning</b> .....	<b>46</b>
<b>9</b>	<b>Maintenance</b> .....	<b>47</b>

**10 Disposal.....47**

**11 Legal Information .....47**

11.1 Liability .....47

11.2 CE Conformity .....47

**12 Technical data .....47**

# 1 Notes regarding the document

## INFORMATION

Date of last update: 2015-08-06

- ▶ Please read this document carefully before using the product (see Fig. 1 on page 4).
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

This document and the product are intended exclusively for the fabrication of a prosthesis by orthopaedic technicians with technical knowledge of lower limb prosthetic fittings. Ottobock seminars, information materials and service are available for continuing education and questions (for contacts, see the manufacturer's address at the end of this document).

### 1.1 Explanation of Warning Symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 1.2 Meanings of pictograms

Meanings of pictograms in the illustrations (see Fig. 2 on page 4)			
	Section with information about illustration		Read and observe this document
	Read safety notice in section		Cleaning instructions in section
	Right		Wrong
	Numbering for parts of the illustration		Numbering for sequences
	Replace		Note time specification
	Joint locked		Joint unlocked
	Movement against the stop		Number of rotations
	No movement resistance		Movement resistance too low
	Optimum movement resistance		Movement resistance too high
	Actual position		Nominal position

## 2 Product Description

The product 3R31\* distinguishes itself with the following key features (see Fig. 3 on page 5):

- Monocentric prosthetic knee joint ①
- Product versions with the following proximal connectors:
  - Pyramid adapter: 3R31 ②
  - Threaded connector: 3R31=ST ③
- Lock for extension stop of the prosthetic knee joint ④
  - Operating the lock slide ⑤ unlocks the joint ⑦
  - Operating the unlock button ⑥ unlocks the joint ⑦
  - Moving to the extension stop ⑧ unlocks the joint ⑨
- Flexion damping for sitting down ⑩
  - Adjustable using the piston rod of the hydraulic cylinder ⑪

## 3 Application

### 3.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

### 3.2 Area of Application

#### **⚠ CAUTION**

##### **Excessive strain on the product**

Falling due to breakage of load-bearing components

- ▶ Use the product exclusively according to its MOBIS classification (see Fig. [4] on page 5).
- ▶ Use the product exclusively for the approved amputation levels ① and amputation sides ② (see Fig. [5] on page 6).

#### **MG 1: indoor walker**

The patient has the ability or the potential to use the prosthesis for the purpose of moving or proceeding slowly on level floors. The amount of time and the distance he/she can walk are seriously limited due to his/her condition.

### 3.3 Combination Possibilities

#### **⚠ CAUTION**

##### **Unallowable combination of prosthetic components**

Fall due to breakage or deformation of the product

- ▶ Note that the illustrated combinations with the pictogram ✓ are permitted and those with the pictogram ✗ are not permitted (see Fig. [6] on page 6).
- ▶ Only combine the product with prosthetic components that are approved for the patient based on the combination table according to MOBIS (see Fig. [7] on page 7: ① prosthetic hip joints, ② proximal adapters, ③ prosthetic knee joints, ④ distal adapters, ⑤ prosthetic feet, ⑥ weight class: see instructions for use).
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

For additional combination possibilities, please see the catalogue 646K2\* or consult the manufacturer.

### 3.4 Environmental Conditions

#### **⚠ CAUTION**

##### **Use under unallowable environmental conditions**

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

##### **Allowable environmental conditions**

Temperature range for use: -10 °C to +45 °C

Allowable relative humidity 0 % to 90 %, non-condensing

##### **Unallowable environmental conditions**

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

### 3.5 Reuse and service life

#### CAUTION

##### **Reuse on another patient**

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

#### CAUTION

##### **Exceeding the service life**

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life defined in this section is not exceeded.

This product has been tested by the manufacturer as follows:

- **Maximum service life in years: 5**  
(Number of load cycles: 1.000.000, mobility grade: 1)

## 4 General Safety Instructions

#### CAUTION

##### **Use of lubricants**

Fall due to breakage of load-bearing components, material damage to the product caused by lubricants

- ▶ Do not use lubricants on the product.

#### CAUTION

##### **Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

#### CAUTION

##### **Mechanical damage to the product**


Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### **Signs of changes in or loss of functionality during use**

Noticeable changes in functionality may include poor response, abnormal noises, etc.

## 5 Scope of Delivery

The scope of delivery (product components and quantities ①) is shown in the illustrations  8 and  9 on pages 8 and 9. Only product components with article numbers ② can be ordered separately.

## 6 Preparation for Use

### 6.1 Information on Fabrication of a Prosthesis

#### ⚠ CAUTION

##### **Incorrect alignment, assembly or adjustment**

Injuries due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

#### ⚠ CAUTION

##### **Initial use of the prosthesis by the patient**

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail, wheeled walker) during initial standing and walking.

The alignment and adjustments are mainly completed in the following steps:

1. Bench alignment of the prosthesis (see Page 42)
2. Optimising the static alignment (see Page 44)
3. Optimising during dynamic trial fitting (see Page 44)

These steps are first carried out with a test prosthesis in order to determine the best possible combination and positioning of the prosthetic components relative to each other. Once all settings have been tailored to the patient, the definitive prosthesis is produced. The same steps are repeated for this purpose.

In dynamic trial fitting, the patient must be given sufficient time to become familiar with the functions of the prosthesis, thereby learning safe use for everyday life.

### 6.2 Bench alignment of the prosthesis

#### INFORMATION

The prosthetic components and prosthesis shown in the illustrations are used as examples to illustrate the general process. The instructions for use of the prosthetic components selected for the patient contain detailed information and have to be used for fabrication of the prosthesis.

- ▶ **CAUTION!** To enable the patient to stand safely, align the prosthesis using the alignment images for the sagittal plane (see Fig. [10](#) on page 10 – ① alignment line, ② prosthetic foot, ③ prosthetic knee, ④ prosthetic socket) and for the frontal plane (see Fig. [11](#) on page 10 – ① alignment line, ② prosthetic foot, ③ prosthetic knee, ④ prosthetic socket) and the following subsections.

**INFORMATION:** Using an alignment apparatus ① with lasers ② is recommended for bench alignment. The retaining bits ③ are used to clamp the prosthetic knee joint in the knee bracket ④ (see Fig. [12](#) on page 11).

#### 6.2.1 Adapting the product to the amputation side

- ▶ Adapt the product to the amputation side of the patient as shown in Fig. [13](#) on page 11.

#### 6.2.2 Positioning the prosthetic foot

- ▶ Position the prosthetic foot as shown in Fig. [14](#) on page 12.
- ▶ **NOTICE!** To prevent rapid wear, adhere to the tolerance range for exterior rotation and do not adapt to the exterior rotation of the sound foot.

### 6.2.3 Mounting the distal adapters

#### **⚠ CAUTION**

##### **Improper assembly of the screw connections**

Fall due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified installation torque values.
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

- ▶ Mount the distal adapters of the prosthetic knee joint on the prosthetic foot (see Fig. [15](#) on page 12).

### 6.2.4 Mounting the prosthetic knee joint

#### **⚠ CAUTION**

##### **Incorrect mounting on the tube**

Falling due to breakage of load-bearing components

- ▶ Do not use the insertion range of the prosthetic knee joint for length compensation.

- ▶ **CAUTION!** Clean the contact surfaces of the prosthetic knee joint and tube using a degreasing cleaning agent compatible with the material.
- ▶ Mount the prosthetic knee joint as shown in Fig. [16](#) on page 13.

### 6.2.5 Mounting the proximal adapters

- ▶ Mount the proximal adapters as shown in Fig. [17](#) on page 13.

### 6.2.6 Preparing the prosthetic socket for the lock slide

- ▶ Prepare the prosthetic socket for mounting the lock slide as shown in Fig. [18](#) on page 14 and note the following points:
  - Positioning at the optimum height for the patient's hand ①
  - Laying the perlon cable in the prosthetic socket, in a straight line from the prosthetic knee joint to the lock slide ②
  - Lamination plate secured against rotation in the laminate ③
  - In and out holes for the perlon cable with sufficient spacing from the prosthetic knee joint and lock slide for reliable lock function ④
  - Plastaband as protection for the thread on the lamination plate while laminating ⑤
  - Greased perlon cable dummy for making a channel centred in the laminate ⑥

### 6.2.7 Marking the prosthetic socket

- ▶ Mark the prosthetic socket as shown in Fig. [19](#) on page 14 for positioning.

### 6.2.8 Positioning the prosthetic socket

- ▶ For positioning the prosthetic socket, note Fig. [20](#) on page 15 and the instructions that follow.

#### **Thomas test**

- ▶ Place the flat hand ② under the lumbar spine.
- ▶ Instruct the recumbent patient ① to flex the contralateral knee joint and hip joint.
- ▶ Bring the residual limb to hip flexion until the hollow back disappears.
- ▶ Extend the residual limb slowly by pushing it until the pressure of the spine on the back of your hand diminishes.
- ▶ In this residual limb position, measure angle  $\alpha$  of the flexion malposition.

#### **Socket flexion**

- ▶ Based on angle  $\alpha$ , determine the optimum socket flexion  $\beta$  ③ according to the illustration.

- ▶ Set the socket flexion  $\beta$  so that the lateral socket reference point (4) lies on the alignment reference line.

**INFORMATION:** If the socket reference point (5) cannot be positioned on the alignment reference line with the adapters that are installed, special adapters have to be used (6 examples of adapters).

### Socket adduction

- ▶ Based on the socket length, determine the required angle  $\gamma$  (7) for socket adduction, and reconcile this with a possible residual limb malposition.
- ▶ Set the socket adduction  $\gamma$  so that the frontal socket reference point (8) lies on the alignment reference line.

**INFORMATION:** If the socket reference point (9) cannot be positioned on the alignment reference line with the adapters that are installed, special adapters have to be used (10 examples of adapters).

## 6.2.9 Shortening the tube adapter

### **CAUTION**

#### **Incorrect processing of tube**

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ For shortening the tube, use only a tube cutter.

**INFORMATION:** The overall system height of all prosthetic components corresponds to the socket-floor measurement (see Fig. [21] on page 16).

- ▶ Shorten and deburr the tube adapter.

## 6.2.10 Mounting the lock slide and lock release cable

- ▶ **CAUTION!** Use the lock slide only as shown in Fig. [22] on page 17 and mount it together with the lock release cable.

## 6.2.11 Preventing damage at maximum flexion

- ▶ **NOTICE!** Perform the steps shown in Fig. [23] on page 17 to prevent damage at maximum flexion.

## 6.3 Optimising the static alignment

### 6.3.1 Checking the static alignment

- ▶ To check the static alignment, position the patient on the measuring device as shown in Fig. [24] on page 18.

### 6.3.2 Optimising the static alignment

- ▶ In the sagittal plane, check the course of the load line (1) relative to the alignment reference point (2) and correct this as shown if needed (see Fig. [25] on page 18).
- ▶ In the frontal plane, check the course of the load line (1) relative to the marked points (2 prosthetic foot, 3 prosthetic knee joint, 4 Spina iliaca anterior superior) and correct if needed (see Fig. [26] on page 19).

## 6.4 Optimising during dynamic trial fitting

### 6.4.1 Resetting to factory defaults

- ▶ Before performing the exercises in the sections that follow, reset settings that have been changed to factory defaults according to Fig. [27] on page 19.

## 6.4.2 Practicing sitting down and standing up

### **⚠ CAUTION**

#### **Adjusting the Settings**

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
  - ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.
- ▶ The patient has to learn safe operation of the lock for all exercises (lock by loading at full extension/unlock using the lock slide).
- ▶ **CAUTION!** With the patient, practice sitting down ① and standing up ② (see Fig. 28 on page 20).
- ▶ ③ If the movement resistance is too low (flexion damping), the prosthetic knee joint has to be adjusted as illustrated.
- ▶ ④ If the movement resistance is too high (flexion damping), the prosthetic knee joint has to be adjusted as illustrated.
- ▶ **CAUTION!** With the patient, practice standing up ⑤ with a highly flexed joint. While standing up, it is necessary to wait until the joint moves to extension again.

## 6.4.3 Practicing walking

- ▶ **CAUTION!** Using parallel bars or a handrail ①, practice walking ② with the patient while the prosthetic knee joint is locked (see Fig. 29 on page 21).
- ▶ **NOTICE!** Be sure to establish functional exterior foot rotation ③.
- ▶ Check the exterior rotation ④ of the prosthetic knee joint.
- ▶ Check socket flexion ⑤.
- ▶ Check socket adduction ⑥.

## 6.5 Attaching the cosmetic cover

### **⚠ CAUTION**

#### **Use of talcum**

Falling, damage to the prosthetic knee joint due to lack of lubrication

- ▶ Do not use talcum on the prosthetic knee joint or other prosthetic components.
- ▶ Attach the cosmetic cover as shown in Fig. 30 on page 22, noting the following points:
- ① To maintain the positioning, only unscrew the two set screws that are screwed in the furthest.
  - ② The cosmetic cover must not interfere with the movement of the lock release cable.
  - ③ After completion, verify functioning of the lock.
  - ④ Check the functioning of flexion damping.
- ▶ To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.

## 6.6 Finishing the prosthesis

- ▶ **CAUTION!** To avoid product damage and the risk of falling, finish the prosthesis by replacing set screws that are too short or too long, and by tightening the screw connections according to Fig. 31 on page 22.

## 7 Use

### 7.1 Information for Use

#### CAUTION

##### **Mechanical overload**

Falling due to breakage of load-bearing components, impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### 7.2 Using the Lock Slide

#### CAUTION

##### **Malfunctions and improper use of the lock**

Falling due to unexpected deactivation of the prosthetic knee joint lock

- ▶ Fully extend the prosthetic knee joint to activate the lock. Then carefully check whether the lock has actually engaged.
- ▶ Only deactivate the lock when you are able to control flexion of the prosthetic knee joint.

## 8 Cleaning

#### CAUTION

##### **Use of unsuitable cleaning agents or disinfectants**

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Clean the product only according to the instructions given in this section.
- ▶ Remove dirt as soon as possible.
- ▶ Test the cleaning agents and disinfectants for material compatibility on non-critical parts of the product.
- ▶ For the prosthesis, observe the cleaning instructions for all prosthetic components.

##### **Cleaning light soiling**

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

##### **Cleaning heavier soiling**

- > **Required materials:** cleaning cloth, 634A58 isopropyl alcohol, soft cloth
- ▶ **NOTICE!** Verify material compatibility. Do not remove lubricant from bearing points.  
Clean the product with a cleaning cloth and isopropyl alcohol.
- ▶ Dry the product with a cloth.

##### **Cleaning with Disinfectants**

- > **Required materials:** colourless, alcohol-free disinfectant (verify material compatibility!), soft cloth
- ▶ Disinfect the product with the disinfectant.
- ▶ Dry the product with a cloth.
- ▶ Allow to air dry in order to remove residual moisture.

## 9 Maintenance

### **⚠ CAUTION**

#### **Failure to follow the maintenance instructions**

Risk of injuries due to changes in or loss of functionality and damage to the product

► Observe the following maintenance instructions.

- **CAUTION!** For trial fittings, inspect the prosthesis several times according to the trial period and use.
- **CAUTION!** For all maintenance and before each use, inspect the plastic components of the product for deformation, damage and cracks.
- **NOTICE!** Do not lubricate and grease the prosthetic joint.
- **NOTICE!** Repair work must be performed exclusively by manufacturer service.
- The prosthetic components should be inspected after the first 30 days of use.
- Arrange regular maintenance intervals with the patient depending on the level of use.
- Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.
- As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

## 10 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

## 11 Legal Information

### 11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 11.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

## 12 Technical data

For the technical data, see Figure [32](#) on page 23.

Abbreviations in the technical data					
①	$h_{s-prox}$	Proximal system height	②	$h_{s-dist}$	Distal system height

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Otto Bock HealthCare GmbH