1E80*, 1E81*, 1E82*, 1E87*

DE Gebrauchsanweisung	3
EN Instructions for use	13
FR Instructions d'utilisation	22
IT Istruzioni per l'uso	32
ES Instrucciones de uso	42
PT Manual de utilização	52
NL Gebruiksaanwijzing	62
SV Bruksanvisning	72
DA Brugsanvisning	81
NO Bruksanvisning	90
FI Käyttöohje	99
PL Instrukcja użytkowania	108
HU Használati utasítás	118
CS Návod k použití	127
RO Instrucțiuni de utilizare	136
HR Upute za uporabu	146
SL Navodila za uporabo	156
SK Návod na používanie	165
BG Инструкция за употреба	175
TR Kullanma talimatı	185
EL Οδηγίες χρήσης	194
RU Руководство по применению	205
JA 取扱説明書	216
ZH 使用说明书	225
KO 사용 설명서	233

1

1

2

3

4

5

INFORMATION

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INFORMATION

Last update: 2015-08-11

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 1E80*, 1E81*, 1E82* und 1E87* Chopart Footplates are used together with the 2C5* and 2E3* Footshells as prosthetic feet in prostheses for Chopart, Pirogoff or Syme amputations.

The workflow for fabricating a Chopart prosthesis is described in detail and illustrated with images in the 1.1.3 Technical Information (646T113*).

The fabrication of prostheses for Chopart, Pirogoff and Syme amputations is essentially the same. The distance between the prosthetic socket and the Chopart footplate must be the same for all prostheses. For this reason, the distal area of the prosthetic socket must be designed to compensate for the amputation level.

1.2 Combination Possibilities

Designation	Reference number
Chopart Glue Kit	SL=P078
Footshell	2C5 (for 1E80*, 1E81*, 1E82*) 2E3 (for 1E87*)
Fill Foam	SL=P071

For combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

2 Application


2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of Application

Area of application according to the MOBIS mobility system:

1E80*, 1E81*, 1E82*

Recommended for mobility grades **3 and 4** (unrestricted outdoor walker and unrestricted outdoor walker with particularly high demands). Approved for a body weight of up to **136 kg**.

1E87*

Approved for a body weight of **max. 50 kg**.

The product has been developed for use in a children's prosthesis.

2.3 Environmental Conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand

3 Safety

3.1 Explanation of Warning Symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions (see the section "Environmental conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

CAUTION

Reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification (see section "Area of application").

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved according to the section "Combination possibilities".

- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Reduced spring effect or decreased forefoot resistance are indications of loss of functionality.

4 Scope of Delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

1E80*, 1E81*, 1E82*, 1E87*					
Fig.	Item no.		Quantity	Designation	Reference number
1	①	■	1	Chopart Footplate	–
–	–	■	1	Instructions for Use	647G1076
–	–	■	1	CD with 1.1.3 Technical Information	646C62

Accessories					
Fig.	Item no.		Quantity	Designation	Reference number
1	②	■	0	Chopart Glue Kit	SL=P078
1	③	■	0	Fill Foam	SL=P071
only for 1E80*, 1E81*, 1E82*					
1	④	■	0	Footshell	2C5
only for 1E87*					
1	⑤	■	0	Footshell for Children	2E3

5 Preparation for Use

5.1 Information on Fabrication of a Prosthesis

⚠ CAUTION

Initial use of the prosthesis by the patient

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail) during initial standing and walking.
- ▶ Do not use devices with insufficient stability (e.g. crutches or walking sticks).

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

NOTICE

Grinding the footshell

Premature wear resulting from damage to the outer footshell

- ▶ Do not damage the outer surface of the footshell when grinding the interior.
- ▶ **Inform the patient.**

INFORMATION

Not all of the materials listed may be available in your country.

- ▶ If the materials are not available, please contact the local branch of the manufacturer to obtain information on alternative materials.

INFORMATION

The legal warranty covers only product defects of the Chopart footplate. It does not include the adhesive bonding of the prosthetic socket and the Chopart footplate, which is included in the prosthetist's area of responsibility. A clean workspace and compliance with all instructions found in the instructions for use are essential for a stable adhesive bond.

Fabricating a prosthesis that is optimally tailored to the patient demands technical knowledge of orthopaedic technology. The alignment and adjustments may only be carried out by a prosthetist. The workflow for fabricating a Chopart prosthesis is described in detail and illustrated with images in the 1.1.3 Technical Information (646T113*).

The key points of the workflow are summarised in these instructions for use in the following sections.

The technical information is provided as a PDF file on the included CD. The file can be viewed on a computer using a PDF display program and printed if necessary.

5.2 Fabricating the Prosthetic Socket

Fabricating the prosthesis is included in the prosthetist's area of responsibility. The design is based on the needs of the patient. The prosthetic socket must be very stiff in the distal area for subsequent gluing with polyurethane. In case of laminated sockets, at least 2 layers of carbon must be applied over the distal area in addition to the reinforcement.

In the course of subsequent processing, the prosthetic socket lies above the footplate at a distance of 3 mm (0.12"). This gap is needed so that polyurethane can penetrate fully underneath the prosthetic socket.

In case of special physical features of the patient (e.g. growth deficits), the distal area of the prosthetic socket must be fabricated so that the height difference is compensated.

5.3 Measuring plumb lines on prosthetic socket

Preparing the prosthetic socket (photos see 646T113*, Section 3.1.1)

- 1) Apply white adhesive tape to the sagittal side of the prosthetic socket.
- 2) Apply white adhesive tape to the frontal side of the prosthetic socket.

Preparing the patient (photos see 646T113*, Section 3.1.2)

- 1) Remove the shoe.
- 2) Put on the protective sock.
- 3) Apply the prosthetic socket.

- 4) Use adhesive tape or strap material to secure the prosthetic socket on the leg.
- 5) **CAUTION! In order to prevent the prosthetic socket from sliding on smooth surfaces, apply adhesive tape to the distal side of the prosthetic socket.**
- 6) **CAUTION! To prevent falling and for load relief, ensure the patient stands between parallel bars.**
- 7) Posture of the patient:
 - standing upright
 - looking straight ahead
- 8) Use foot blocks to compensate for any height difference.
- 9) Check the upright pelvis position.
- 10) Based on the needs and anatomy of the patient, use foot blocks to compensate for the length difference.

Marking the frontal plumb line (photos see 646T113*, Section 3.1.3)

- ▶ Project the frontal plumb line onto the prosthetic socket with the laser and mark it.

Marking the sagittal plumb line (photos see 646T113*, Section 3.1.4)

- ▶ Project the sagittal plumb line onto the prosthetic socket with the laser and mark it.

5.4 Preparing the Prosthesis for the Trial Fitting

Preparing the prosthetic socket and footplate (photos see 646T113*, Section 3.2.1)

- 1) Roughen the distal side of the prosthetic socket.
- 2) Roughen the footplate to the marked line.
- 3) Glue a rubber ring from the Chopart glue kit onto the centre of the prosthetic socket using superglue.
- 4) Draw crosshairs onto an underlay.
- 5) Position the footshell with the inserted footplate so that the lines on the footplate line up with the crosshairs.
- 6) Mark the outside contour of the footshell on the underlay.
- 7) Apply double-sided adhesive tape to the lower side of the footplate.
- 8) Lay the crepe sole onto the marked contour.
- 9) Glue the footplate onto the crepe sole so that the lines of the footplate line up with the crosshairs.
- 10) Put the footplate with crepe sole into the shoe.
- 11) Mark the spacer washer with chalk.
- 12) Ensuring it is plumb, press the prosthetic socket onto the footplate in order to transfer the marking. Make sure the plumb lines are positioned vertically.
- 13) The marking helps adjust an initial a – p position of the footplate to the prosthetic socket in the alignment apparatus.
- 14) Remove the footplate from the shoe.

Clamping in the alignment apparatus (photos see 646T113*, Section 3.2.2)

- 1) Position the footplate in the foot adapter unit on the alignment table.
- 2) Set the effective heel height of the shoe.
- 3) Set the external rotation.
- 4) Transfer the lateral and medial markings onto the edges of the footplate at 90 degrees to the walking direction.
- 5) Clean the footplate and the roughened area of the prosthetic socket with isopropyl alcohol.
- 6) Clamp the prosthetic socket in the alignment apparatus and align it according to the frontal and sagittal plumb lines.
- 7) Lower the prosthetic socket.
- 8) Adjust the footplate adapter with the footplate in the a – p direction so that the lateral markings ② line up with the spacer washer ①.

Gluing the prosthetic socket to the footplate (photos see 646T113*, Section 3.2.3)

- 1) Mix Orthocryl sealing resin and talcum powder to make a firm putty.
- 2) Add the hardener.
- 3) Apply a dot of putty onto the footplate.
- 4) Lower the prosthetic socket onto the footplate.
- 5) Allow the putty to harden.
- 6) Remove the joined components from the alignment table.
- 7) Use a marker to transfer the outline of the footplate onto the crepe sole.
- 8) Before applying more putty, pull the crepe sole off the footplate.
- 9) Add putty to the joint between the prosthetic socket and the footplate so that the unit can withstand a dynamic trial fitting.
- 10) Allow the putty to cool.
- 11) If desired tear-resistant adhesive tape may be used to secure the glued joint between the footplate and prosthetic socket.
- 12) Glue the crepe sole back onto the footplate in the marked position.

5.5 Trial Fitting of the Prosthesis

Donning the prosthesis (photos see 646T113*, Section 3.3.1)

- 1) Apply the prosthetic socket and secure it with adhesive tape or strap material.
- 2) Use filler (e.g. foam) to stabilise the footplate in the shoe.
- 3) Put on the shoes.

Checking the static alignment (photos see 646T113*, Section 3.3.2)

- ▶ Use the L.A.S.A.R. Posture to optimise the static alignment based on the available technical information.

Checking the alignment in the dynamic trial fitting (photos see 646T113*, Section 3.3.3)

- 1) Only test the prosthesis during normal walking since the glued joint is not designed for greater stress.
- 2) Conduct the dynamic trial fitting.
- 3) If required, modify the alignment in regards to plantar flexion and dorsal flexion as well as pronation and supination by inserting small spacers into the shoe.
- 4) If the positioning of the prosthetic socket to the footplate has to be corrected, follow the instructions for separating the components in the section "Definitive Gluing of the Prosthetic Components" (see Page 18) – heading "Preparing the prosthetic socket and the footplate for gluing" and then continue with the instructions in the section "Preparing the Prosthesis for the Trial Fitting" (see Page 17), – heading "Clamping in the alignment apparatus", taking into account the corrections.
- 5) Repeat the dynamic trial fitting and correct the footplate position until the prosthesis wearer exhibits a harmonious gait pattern and can walk physiologically and without stress.

5.6 Definitive Gluing of the Prosthetic Components

Creating the mould (photos see 646T113*, Section 3.4.1)

- 1) Open but do not remove the cover on component A in order to prevent the formation of over-pressure during heating.
- 2) Warm component A in the heating oven in order to decrystallise it (min. 3 h at 65 °C / 149 °F).
- 3) The processing time for the polyurethane can be extended by cooling it (min. > 17 °C / > 62.6 °F).
- 4) Remove any pads that are positioned and/or secured on the prosthetic socket.
- 5) Pull the crepe sole off the footplate.
- 6) Remove glue residue from the footplate.
- 7) Remove the white adhesive tape from the prosthetic socket.
- 8) Grind the putty until it is straight to slightly concave.

- 9) The shape that is achieved by grinding will be the shape of the PU glued joint.
- 10) Grinding must result in a shape that permits the fitting of a footshell.
- 11) Clamp the prosthesis so it faces up at an angle of approximately 40°, with the toes pointing down.
- 12) Use a cloth to practise the application of the polyethylene plate to the heel area of the prosthesis.
- 13) Use isopropyl alcohol to clean both sides of the polyethylene plate.
- 14) Heat the polyethylene plate at 130 °C (266 °F) until it becomes transparent and can be moulded.
- 15) Remove it from the heating oven wearing temperature-resistant gloves.
- 16) Apply the polyethylene plate over the heel of the prosthesis (1/2 over the prosthetic socket, 1/2 over the footplate).
- 17) Mould the heated material. In doing so, avoid excessive undercutting and close the dorsal region by butting the edges up against each other rather than overlapping them.
- 18) Do not squeeze the polyethylene tight together in the dorsal area.
- 19) Apply the polyethylene plate carefully in the region with the putty.
- 20) Mould the polyethylene plate on the sides, in the anterior area of the footplate.
- 21) The mould has to firmly encompass the edges of the footplate and the edge in the heel area of the prosthesis, so that the footplate and the prosthesis are kept in their correct position to each other later during casting.
- 22) Allow the mould to cool.
- 23) Mark the trimming edge so that the edges of the mould hold the footplate and the prosthetic socket in place.
- 24) Expose the dorsal area of the prosthesis so the polyurethane can be poured in later.
- 25) Trim the mould so that the prosthesis can be removed.
- 26) Remove the prosthesis from the mould.
- 27) Trim the mould so that the prosthetic socket and the footplate can be reinserted and are held firmly in place.
- 28) Deburr and clean the mould.

Preparing the prosthetic socket and the footplate for gluing (photos see 646T113*, Section 3.4.2)

- 1) Use a band saw to cut the glued joint.
- 2) Cut within the 3 mm gap between the prosthetic socket and footplate in order to avoid damage.
- 3) Completely remove the putty from the prosthetic socket.
- 4) Avoid grinding the resin of the laminate as far as possible during this process, so that fibre exposure is minimised.
- 5) Remove the putty from the footplate by grinding.
Avoid grinding or damaging the epoxy facing and the carbon.
- 6) If fibres are exposed during grinding or too much acrylic resin has been removed from the acrylic resin prosthetic socket, clean the affected area thoroughly and then apply a thin coating of sealing resin.
Allow the material to dry fully for at least 24 h prior to further processing.
- 7) If too much epoxy resin has been removed from the epoxy resin prosthetic socket or the footplate, or if fibres are exposed, clean the affected area thoroughly and then apply a thin coating of epoxy resin.
Allow the material to dry fully for at least 24 h prior to further processing.
- 8) **CAUTION! Only for acrylic resin prosthetic sockets, use the primer as a bonding agent.**
NOTICE! Shake the container with the primer thoroughly so that all components that may have settled are dissolved and mixed.

- 9) Clean the prosthetic socket with isopropyl alcohol and allow it to dry.
- 10) **CAUTION! Risk of falling – Do not apply the primer to the Chopart footplate and an epoxy resin prosthetic socket, as the polyurethane bond will not be durable enough and will break.**

Apply a thin coat of primer to the acrylic resin prosthetic socket on the surfaces being glued.

- 11) The primer has to dry for 15 - 20 min before gluing.
- 12) Cut a piece (approx. 1/4) from the rubber ring.
- 13) Clean the mould with isopropyl alcohol.
- 14) Roughen the footplate with the sandpaper.
- 15) Clean the footplate and sides with isopropyl alcohol in the area being glued, and do not touch it again.
- 16) Insert the footplate in the mould.
- 17) Insert the prosthetic socket into the mould without touching the areas being glued.
- 18) Push the piece of the rubber ring under the centre of the prosthetic socket using a clean stir stick.

INFORMATION: This ensures that the original 3 mm gap is maintained between the footplate and the prosthetic socket.

- 19) Secure the position of the prosthetic socket relative to the footplate in the mould using adhesive tape.
- 20) Clamp the assembled components facing down at an angle of approximately 45°, with the dorsal area (toes) pointing up.

Casting the polyurethane (photos see 646T113*, Section 3.4.3)

- 1) Add the entire contents of component B to the container with component A. Slowly and thoroughly stir the polyurethane to mix it (2 – 3 min).
- 2) Trapped air bubbles must be avoided.
- 3) In order to ensure that the polyurethane is thoroughly mixed, scrape the sides and bottom of the container while stirring.
- 4) Slowly pour the polyurethane into the mould to the side of the spacer ring. Trapped air bubbles must be avoided.
- 5) **NOTICE! Do not exceed the marked line on the footplate when pouring.**
- 6) The height of the joint to the prosthetic socket can be influenced by changing the angle of the prosthesis.
- 7) If the polyurethane settles as it hardens, add more until the level of the marked line on the footplate is reached.
- 8) After casting, allow the glued joint to harden undisturbed (min. 3 h / min. 20 °C / 68 °F).

Tempering the polyurethane (photos see 646T113*, Section 3.4.4)

- 1) Place the prosthesis with the mould into the heating oven on its side and temper it to its final strength (16 h at 80 °C / 176 °F).
- 2) Then remove the prosthesis from the heating oven and allow it to cool.
- 3) Remove the mould.
- 4) Remove excess polyurethane from the bottom of the footplate.
- 5) Trim excess polyurethane on the socket edge.
- 6) Clean the prosthesis with isopropyl alcohol.

5.7 Finishing the Prosthesis

Preparing the footshell (photos see 646T113*, Section 3.5.1)

- 1) The footshell is adapted step-by-step according to the instructions below, so that the prosthesis fits exactly in the footshell and the edge of the footshell fully contacts the prosthetic socket.
- 2) Grind the proximal edge of the footshell to 0 on the inside.
- 3) Remove the wedges in the heel area.

- 4) Insert the prosthesis into the footshell as far as possible.
- 5) Use a washable marker to mark an initial trim line.
- 6) Trim the edge step-by-step.
- 7) Grind the edge from the inside.
- 8) Check the fit by inserting the prosthesis.
- 9) Repeat the process until the prosthesis fits exactly into the footshell and the edge of the footshell fully contacts the prosthetic socket.

Foaming the footshell (photos see 646T113*, Section 3.5.2)

- 1) Practise inserting the prosthesis into the footshell quickly and confidently.
- 2) Place the shoehorn within reach.
- 3) Apply plaster parting agent cream to the outside of the footshell in order to prevent the adhesion of fill foam.
- 4) **INFORMATION: Work quickly when processing the foam.**
- 5) Add component A to the container with component B.
- 6) Mix the components by stirring quickly with a stir stick (stirring time approx. 10 s).
- 7) Pour most of the fill foam into the footshell.
- 8) Distribute the fill foam in the footshell by tilting it back and forth.
- 9) Insert the prosthesis into the footshell.
- 10) If required, raise the upper edge of the footshell and add more fill foam.
- 11) In order to avoid changing the heel height due to excessive fill foam between the footplate and footshell, push down on the prosthesis with moderate pressure until the foam has set.
- 12) Remove foam residue from the prosthesis so that the edge of the footshell fully contacts the prosthetic socket.
- 13) Add buckles or straps to the prosthesis. Position and secure pads if applicable.
- 14) Clean the prosthesis.

6 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

7 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

8 Legal Information

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH