


17LK3

DE	Gebrauchsanweisung	4
EN	Instructions for use	11
FR	Instructions d'utilisation	18
IT	Istruzioni per l'uso	26
ES	Instrucciones de uso	33
PT	Manual de utilização	41
NL	Gebruiksaanwijzing	49
SV	Bruksanvisning	56
DA	Brugsanvisning	63
NO	Bruksanvisning	70
FI	Käyttöohje	77
CS	Návod k použití	84
SK	Návod na používanie	91
RU	Руководство по применению	98
ZH	使用说明书	107


1

22

2

INFORMATION

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INFORMATION

Last update: 2015-10-20

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the processing of the 17LK3 System Knee Joints.

2 Product description

2.1 Available sizes

The unilateral system knee joint is available in four system sizes. It can be used either on the medial or lateral side, or on both the medial and the lateral sides of the knee. The choice of system size is based on the patient's body weight and conditions of use.

	Maximum patient weight [kg]			
	17LK3=12	17LK3=14	17LK3=16	17LK3=20
Unilateral use (lateral or medial)	20	50	85	110
Bilateral use (lateral and medial)	40	80	120	160

INFORMATION

- When using the product unilaterally, and in case of flexion contracture in the knee or hip > 10° or distinct torsion or valgus/varus instabilities or valgus/varus malpositions > 10° or increased bodily activity, the next higher size must be used.
- For tubercle seat the product must be fitted bilaterally.

2.2 Components

Item	Component (see Fig. 1)	Joint	Article no.
1	Knee lever		30Y265=*
2	Knee joint upper part		
3	Lock washer		507S96=2.3
4	Washer		30Y267=3.3
5	Washer		30Y267=4.2
6	Bearing nut	(17LK3=12/=14) (17LK3=16) (17LK3=20)	17Y93=9x9.2xM6 17Y93=9x10.75xM6 17Y93=9x11.75xM6
7	Lock washer	(17LK3=12/=14) (17LK3=16/=20)	507S96=3.2 507S96=4
8	Compression spring	(17LK3=12) (17LK3=14/=16/=20)	513D83=0.75x4.3x19 513D83=0.75x4.6x20
9	Lock wedge	(17LK3=12/=14) (17LK3=16) (17LK3=20)	17Y156=* 17Y157=* 17Y158=*
10	Bushing	(17LK3=12/=14)	30Y87=3
11	Bearing washer	(17LK3=16/=20)	170Z93=1
12	Bushing	(17LK3=16/=20)	30Y87=1
13	Knee joint lower part		
14	Axle	(17LK3=12/=14) (17LK3=16) (17LK3=20)	4A101=5x10.6 4A101=6x12.7 4A101=6x13.7
15	Screw with flattened half-round head	(17LK3=12/=14) (17LK3=16/=20)	501F7=M6x8 501F7=M6x10
16	Axle	(17LK3=12/=14) (17LK3=16) (17LK3=20)	4A101=3.2x15.1 4A101=3.2x17.1 4A101=3.2x18.1
17	Axle	(17LK3=12/=14) (17LK3=16) (17LK3=20)	4A101=4x15.1 4A101=4x17.1 4A101=4x18.1
18	Protective plug	(17LK3=12/=14) (17LK3=16/=20)	30Y91=3 30Y91=1

Item	Component (see Fig. 1)	Joint	Article no.
19	Temporary switch		170Y162=*
20	Joint dummy		30Y268=*
22	Countersunk head screw (Allen screw)	(17LK3=12/=14) (17LK3=16/=20)	501S41=M4x10 501S41=M5x12

The scope of delivery does NOT include:

Item	Component (see Fig. 1)	Article no.
21	Shoulder screw (part of 743R6 adjustment Set)	30Y89
Item	Component (see Fig. 2)	Article no.
23	Joint bar	17LS3=*

3 Application

3.1 Indications for use

The system knee joints are to be used **exclusively** for the orthotic fitting of the lower limbs on a patient.

In case of unilateral use, the manufacturer recommends that the system knee joint be processed by means of the lamination resin technique or carbon prepreg technology. When used bilaterally, there is no limitation on the processing technology used.

3.2 Indications

Partial or total paralysis of the leg muscles.

4 Safety

4.1 Explanation of warning symbols

WARNING	Warning regarding possible serious risks of accident or injury.
CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

4.2 General safety instructions

WARNING

Operation of motor vehicles

Risk of accidents due to restricted body function

- Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

► **Inform the patient.**

△ CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- Only use the product for the defined area of application.
- If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).

△ CAUTION

Incorrect alignment, assembly or adjustment

Injuries may result from incorrectly assembled, incorrectly adjusted or damaged components

- Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- Observe the alignment, assembly and adjustment instructions.
- **Inform the patient.**

△ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- Use caution when working with the product.
- Check the product for proper function and readiness for use.
- In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

△ CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- Do not reach into the joint mechanism during daily use.
- Close attention is required during assembly and adjustment tasks.
- **Inform the patient.**

△ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at a temperature above 300 °C (570 °F).

NOTICE

Thermal overloading of product components

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Processing

5.1 Laminating the orthosis

INFORMATION

The instructions for manufacturing an orthosis with unilateral joint systems are included in the 646T5=9.1 Technical Information. The laminate structure described here is a recommendation of Ottobock and can be adjusted to the needs of the patient.

- > **Required tools and materials:** 646G15 woven carbon stockinette, 17Y106 PVC profile material, 616G12 carbon fibre cloth, 99B81=* PVA bag, 623T3=* Perlon stockinette, 636K8 plastaband, 616G6 Dacron felt, 633W8 wax, 17LS3=* joint bar, 617H55 C-Orthocryl lamination resin, 743R6 adjustment set for orthotic joints.
 - > **Optional:** 30Y268=* lamination dummy
- 1) Place cut pads or shim layers on plaster model.
 - 2) Pull a layer of perlon stockinette over the plaster model.
 - 3) Align the bar parts on the plaster model.
 - 4) **When using a lamination dummy:** Fix the joint bars to the alignment adapter (743Y56=4) using the lamination dummy and the shoulder screw.
 - 5) **When not using a lamination dummy:** Fix the joint bars to the alignment adapter (743Y56=4) using the base body (no illustration).
 - 6) Pull a PVA bag over the plaster model to isolate it.
 - 7) Pull a layer of perlon stockinette over the plaster model.
 - 8) Fix the finished reinforced bar parts to the plaster model. Carry out the reinforcement in accordance with patient's needs.
 - 9) Seal the lamination dummy or base body of the joint with Plastaband and fill hollows with Dacron felt.
 - 10) Pull another layer of perlon stockinette over the plaster model.
 - 11) Affix further reinforcement layers for the support surfaces as described in Section 8.
 - 12) Complete reinforcement with a layer of Perlon stockinette and designer fabric.
 - 13) Pull on a PVA bag and vacuum.
 - 14) Then laminate with C-Orthocryl lamination resin.

6 Assembly

	Torque [Nm]			
	17LK3=12	17LK3=14	17LK3=16	17LK3=20
Screw with flattened half-round head (item 15)	8	8	8	8
Countersunk head screw (item 22)	3.6	3.6	7.1	7.1

- 1) Demould the model and cut free the lamination dummy (20) or the system knee joint.
- 2) Insert the knee joint lower part (13) with the bushings (10) (**only 17LK3=12/=14**) or with the bushings (12) and washers (11) (**only 17LK3=16/=20**) into the knee joint upper part (2).

- 3) Fasten the screw (15) and bearing nut (6) and secure the screw connection with Loctite® 241 or tighten to the specified torque value (see table).
- 4) Fasten the joint on the joint bars (23) using the countersunk head screws (22) and secure the screw connections with Loctite® 241 and tighten to the specified torque value (see table).

6.1 Mounting the accessories

CAUTION

Improper handling

Risk of falling due to accidentally opening or engaging the lock

- Before loading the orthosis, verify that the temporary switch is in the desired position.
- Avoid accidentally opening or engaging the lock.
- Instruct the patient in the correct usage of the system joints.

A temporary switch (17Y162=*) which can be subsequently installed in the system knee joint is included as an accessory for the system knee joint. The temporary switch prevents the lock from engaging, e.g. while cycling. To install the temporary switch, please follow the assembly instruction (647G103).

7 Cleaning

After contact with water containing salt, chlorine or soap, or when they get dirty, the system joints must be cleaned promptly.

- 1) Rinse the system joints with clean fresh water.
- 2) Dry the system joints with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

8 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH