ottobock.

Liner

DE	Gebrauchsanweisung	2
EN	Instructions for use	8
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	20
ES	Instrucciones de uso	27
PT	Manual de utilização	33
NL	Gebruiksaanwijzing	39
SV	Bruksanvisning	46
DA	Brugsanvisning	51
NO	Bruksanvisning	57
FI	Käyttöohje	63
PL	Instrukcja użytkowania	69
HU	Használati utasítás	75
CS	Návod k použití	81
RO	Instrucțiuni de utilizare	87
HR	Upute za uporabu	93
SL	Navodila za uporabo	99
SK	Návod na používanie	105
BG	Инструкция за употреба	111
TR	Kullanma talimatı	117
EL	Οδηγίες χρήσης	123
RU	Руководство по применению	130
JA	取扱説明書	137
ZH	使用说明书	142
KO	사용 설명서	147

INFORMATION

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INFORMATION

Last update: 2015-10-30

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

INFORMATION

- ▶ The patient must be fitted with the liner no later than the date printed on the inside of the packaging or on the label (MHD/BBD).

1.1 Function

In combination with an appropriate suspension system, liners increase wear comfort and improve the suspension of prostheses. Together, they reduce the forces that affect the residual limb and, depending on the type, may also serve further functions.

These instructions for use apply to the following liners:

TF fitting	
Material	Reference number
Polyurethane	6Y430, 6Y540

TT fitting	
Material	Reference number
Polyurethane	6Y400, 6Y416, 6Y510, 6Y512, 6Y520, 6Y521, 6Y522, 6Y523

Polyurethane

Polyurethane liners are made of a polyurethane gel. The material is viscoelastic. It adapts very well to the form and structure of the residual limb and protects very well against impact forces.

These liners are appropriate for sensitive residual limbs, for instance with scar tissue and bony structures.

Ottobock recommends a total surface weight-bearing socket (the "Harmony" principle) with a harmony system or an exhaust valve.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system. As a result of its design, it may be combined with other prosthetic components of the modular system. You can find possible restrictions in this section. Contact the manufacturer with questions.

2 Application

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Service Life

The product is a wear part, which means it is subject to normal wear and tear.

3 Safety

3.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

3.2 General Safety Instructions

 WARNING
Product is pulled over mouth and nose
Danger of suffocation due to lack of oxygen supply
► Store the product away from the reach of children.

⚠ CAUTION

Contact with cold, heat, embers or fire

Injuries (e.g. burns) caused by melting of the material and damage to the product

- ▶ Keep the product away from open flames, embers and other sources of heat.
- ▶ Do not expose the product to temperatures above +80 °C or below -20 °C.

⚠ CAUTION

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.

⚠ CAUTION

Incorrect or excessively tight application of the product to the body

Risk of pressure points and constriction of blood vessels and nerves due to improper application

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ Do not continue to use the product if there are problems with the correct fit.

⚠ CAUTION

Use of the product on damaged skin

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product on healthy skin.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Liner	–

5 Preparation for use

The following technical information is available if needed:

5.1 Adapting the liner

- > **Tools needed for liners with textile cover:** 719S7 trimming scissors or 719S20 special scissors, cyanoacrylate superglue (e.g. 636K36)
- > **Tools needed for liners without textile cover:** 756L9 trimming device or 756L10 liner trimmer
- 1) **NOTICE! Do not use wave edge scissors.**
Cut the liner horizontally at the required length.
- 2) **For liners with textile cover:** secure open stitch end with superglue.

5.2 Optional: attaching the pin

INFORMATION

The shape of the prosthetic socket must allow load on the end of the residual limb.

The pin to connect the liner to the shuttle lock is screwed into the thread at the distal end of the liner.

> **Required tools and materials:**

Torque wrench (adjustable to 3 Nm), 636K13 Loctite®

- 1) Apply Loctite® to the thread.
- 2) **NOTICE! Do not shorten flexible pins and pins with exhaust valve.**
Screw the pin into the thread and tighten it (installation torque: **3 Nm**).

6 Use

⚠ CAUTION

Use of powder (e.g. baby powder or talcum)

Skin irritation on the residual limb and loss of function of prosthetic components due to clogging with particles or deprivation of lubricant

- ▶ Keep the product away from powders.
- ▶ **Inform the patient.**

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

INFORMATION

Some liners are treated with the Skinguard Technology additive. Skinguard Technology has an antibacterial effect and prevents the formation of unpleasant odours.

Substances in sweat containing sulphur may cause discolouration. This discolouration has no effect on the function or durability of the liner but cannot be removed. Sulphur may enter sweat through food, for example, protein products, cabbage, nuts, legumes and sulphured products (e.g. dried fruit).

Discolouration of the liner does not justify customer complaint claims.

INFORMATION

- ▶ Liners that have already been worn generally feel softer than new liners.
- ▶ Because of the relatively long drying times of polyurethane liners, a spare liner is recommended.
- ▶ Clean new liners before using for the first time.

INFORMATION

Derma Prevent may damage the outer and inner coating of polyurethane liners.

6.1 Putting on and taking off the liner

Putting on the liner

- 1) Check the liner for damage. If a pin is attached, also check the pin for damage.
- 2) Roll up the liner and place it on the end of the residual limb.
NOTICE! If a pin is attached, then align the pin along the axis of the residual limb.
- 3) Unroll the liner onto the residual limb so that it has no wrinkles, without shifting the soft tissue.
- 4) Check the fit of the liner and the alignment of the liner.
INFORMATION: Do not wrap the liner over the prosthetic socket (exception: 6Y100).

Taking off the liner

- 1) Roll up the liner from the open end to the closed end and remove it from the residual limb.
- 2) **Recommendation:** clean the liner according to the "Cleaning and care" section.

6.2 Cleaning and Care

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- Only clean the product with the approved cleaning agents.

INFORMATION

Clean the liner and skin more frequently in case of increased perspiration. Wearing a liner requires daily intensive skin care. Ottobock recommends the use of the Derma skin care system.

- 1) Turn the liner inside out.
- 2) Hand wash the liner in warm water (**30 °C**) using neutral detergent (e. g. Ottobock 453H10=1 Derma Clean).
- 3) Thoroughly rinse the liner with clear water.
- 4) Put a hand towel into the liner and leave the liner to dry lying flat on a drying rack. Alternatively, use a drying stand (TT: 5Y100=1 TF: 5Y100=2).

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive.

The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

7.3 Warranty

The manufacturer warrants this device from the date of fitting. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH