











## 743A220 PROS.A. Assembly

|  |     |
|--|-----|
| <b>DE</b> Gebrauchsanleitung (Benutzer / Techniker) .....                              | 23  |
| <b>EN</b> Instructions for use (user / technician) .....                               | 49  |
| <b>FR</b> Mode d'emploi (Utilisateur / Technicien) .....                               | 75  |
| <b>IT</b> Istruzioni d'uso (utilizzatore / tecnici) .....                              | 101 |
| <b>ES</b> Instrucciones de uso (usuario / técnico) .....                               | 127 |
| <b>TR</b> Kullanım kılavuzu (Kullanıcı / Teknisyen) .....                              | 153 |
| <b>RU</b> Руководство по применению (для пользователей / технического персонала) ..... | 179 |

**13**

**7.1**

**646F219=D** ottobock.

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# Table of contents

EN

|           |   |           |
|-----------|---|-----------|
| <b>1</b>  | <b>Important information on this document .....</b>                       | <b>51</b> |
| 1.1       | Scope of this document .....  | 51        |
| 1.2       | Document release date .....   | 51        |
| 1.3       | Structure of this document .....  | 52        |
| 1.4       | Meanings of safety instructions and pictograms .....                      | 53        |
| 1.4.1     | Meanings of signal words in safety instructions.....                      | 53        |
| 1.4.1.1   | Signal word WARNING .....   | 53        |
| 1.4.1.2   | Signal word CAUTION .....   | 53        |
| 1.4.1.3   | Signal word NOTICE .....  | 53        |
| 1.4.1.4   | Structure of the safety instructions .....                                | 53        |
| 1.4.2     | Meanings of pictograms.....   | 53        |
| <b>2</b>  | <b>Product Description.....</b>   | <b>54</b> |
| 2.1       | Intended Use.....   | 54        |
| 2.2       | Product overview .....  | 54        |
| 2.3       | Accessories .....   | 54        |
| 2.3.1     | Mounting frame with laser holder .....                                    | 54        |
| 2.3.2     | Spacer .....  | 54        |
| 2.3.3     | Casters.....  | 54        |
| <b>3</b>  | <b>Safety.....</b>  | <b>55</b> |
| <b>4</b>  | <b>Technical data .....</b>   | <b>55</b> |
| 4.1       | Dimensions and weight.....  | 55        |
| 4.2       | Environmental conditions .....  | 55        |
| <b>5</b>  | <b>Transport.....</b>   | <b>56</b> |
| <b>6</b>  | <b>Storage.....</b>   | <b>56</b> |
| <b>7</b>  | <b>Unpacking.....</b>   | <b>56</b> |
| 7.1       | Checking the delivered items.....   | 57        |
| <b>8</b>  | <b>Positioning.....</b>   | <b>57</b> |
| <b>9</b>  | <b>Installation .....</b>   | <b>58</b> |
| 9.1       | Move the crank handle to the working position.....                        | 58        |
| 9.2       | Adjusting the feet of the product .....                                   | 58        |
| 9.3       | Installing the mounting frame with laser holder.....                      | 58        |
| 9.4       | Installing the casters.....   | 58        |
| 9.5       | Using the mounting plate .....  | 58        |
| <b>10</b> | <b>Initial operation.....</b>   | <b>59</b> |
| 10.1      | Installing the laser gauge .....  | 59        |
| 10.2      | Adjusting the lasers.....   | 59        |
| 10.3      | Removing the laser gauge.....   | 59        |
| <b>11</b> | <b>Use.....</b>   | <b>60</b> |
| 11.1      | Information for use .....   | 60        |
| 11.2      | Aligning a prosthesis .....   | 60        |
| 11.2.1    | Important information about aligning a prosthesis.....                    | 60        |
| 11.2.2    | Positioning the prosthetic foot.....                                      | 61        |
| 11.2.2.1  | Determining the required measurements of the shoe .....                   | 61        |
| 11.2.2.2  | Determining the required measurements of the prosthetic foot.....         | 61        |
| 11.2.2.3  | Adjusting the height of the holder for the footplate.....                 | 61        |
| 11.2.2.4  | Setting the heel height .....   | 61        |
| 11.2.2.5  | Positioning the footplate in a-p direction .....                          | 62        |
| 11.2.2.6  | Positioning the prosthetic foot.....                                      | 62        |
| 11.2.2.7  | Setting the external rotation of the foot .....                           | 62        |
| 11.2.3    | Positioning the prosthetic knee joint .....                               | 63        |
| 11.2.3.1  | Determining the required measurements for the prosthetic knee joint ..... | 63        |

---

|           |   |           |
|-----------|---|-----------|
| 11.2.3.2  | Adjusting the height of the holder for the knee bracket.....              | 63        |
| 11.2.3.3  | Mounting the knee bracket.....  | 63        |
| 11.2.3.4  | Mounting the retaining bits.....  | 64        |
| 11.2.3.5  | Clamping the prosthetic knee joint.....                                   | 64        |
| 11.2.4    | Installing the prosthetic components proximal to the prosthetic foot..... | 64        |
| 11.2.5    | Installing the prosthetic components above the prosthetic knee joint..... | 64        |
| 11.2.6    | Positioning of the prosthetic socket.....                                 | 65        |
| 11.2.6.1  | Selecting and installing socket clamping fixture.....                     | 65        |
| 11.2.6.2  | Configuring and installing the centring unit.....                         | 65        |
| 11.2.6.3  | Marking the prosthetic socket and determining important measurements..... | 66        |
| 11.2.6.4  | Securing the prosthesis socket on the socket clamping fixture.....        | 66        |
| 11.2.6.5  | Positioning the prosthetic socket to the alignment reference line.....    | 67        |
| 11.2.6.6  | Adjusting the height of the holder for the socket clamping fixture.....   | 67        |
| 11.2.7    | Documenting the settings.....   | 67        |
| 11.3      | Measuring a prosthesis.....   | 70        |
| 11.4      | Cleaning.....   | 70        |
| <b>12</b> | <b>Maintenance.....</b>   | <b>71</b> |
| 12.1      | Inspection.....   | 71        |
| 12.1.1    | Checking the lasers.....  | 71        |
| 12.2      | Maintenance.....  | 71        |
| 12.3      | Repair.....   | 71        |
| <b>13</b> | <b>Disposal.....</b>  | <b>71</b> |
| <b>14</b> | <b>Legal Information.....</b>   | <b>71</b> |
| 14.1      | Liability.....  | 71        |
| <b>15</b> | <b>Appendices.....</b>  | <b>72</b> |
| 15.1      | List of pictograms used.....  | 72        |
| 15.1.1    | Warning symbols.....  | 72        |
| 15.1.2    | Mandatory symbols.....  | 72        |
| 15.1.3    | Prohibition symbols.....  | 72        |
| 15.1.4    | Transport and storage pictograms.....                                     | 72        |

# 1 Important information on this document

- see fig. 1 (see page 2)

|   |   |
|---|---|
|  |  <b>WARNING</b>  |
|   | <p><b>Working with the product without observing the instructions for use</b><br/>         Serious injury due to reckless manner of working</p> <ul style="list-style-type: none"> <li>▶ Observe the instructions for use when working with the product.</li> <li>▶ Familiarise yourself with the structure of the instructions for use by reading the chapter <b>Structure of this document</b> (see Page 52) and the table of contents.</li> <li>▶ In particular, observe the chapter <b>Safety</b> (see Page 55) and the safety instructions in the chapters relating to the work to be carried out.</li> <li>▶ Please always keep these instructions for use by the product.</li> </ul> |

## Relationship between the figure pages and the chapters

The figure pages at the beginning of this document contain figures about the product and its use. The figures are related to the chapters.

For easier navigation, a chapter reference (chapter number below the symbol of the title page of the document) is located below the figure number of each figure, and a figure reference (• see fig. X (see page Y)) is included in the chapters.

The figure reference is always located above the related text.

## 1.1 Scope of this document

- see fig. 2 (see page 2)

This document is valid for the following product:

| Reference number | Designation      | From year of manufacture (year-calendar week)<br>From factory no./from lot no. |
|------------------|------------------|--|
| 743A220          | PROS.A. Assembly | 2014-27<br>0001201427  |

| Nameplate legend |                    |   |                     |
|------------------|--------------------|---|---------------------|
| ①                | Product name       | ② | Reference number    |
| ③                | Fabric no./lot no. | ④ | Year of manufacture |

## 1.2 Document release date

This document is valid as of the following release date:

| Printing reference number | Version | Release date [year-month-day] |
|---------------------------|---------|-------------------------------|
| 647G1000                  | 01      | 2016-04-12                    |

### 1.3 Structure of this document

The chapters of this document contain all important product information.

The document structure is essentially oriented towards the various phases in the product life cycle and the associated tasks.

The product can only be used safely if all chapters relevant to the given work have been read.

The following list shows the chapters important for all work, those important for specific tasks and those containing general information:

#### Chapters important for all work

- **Important information on this document** (see Page 51)
- **Product description** (see Page 54)
- **Safety** (see Page 55)
- **Legal information** (see Page 71)

#### Chapters important for specific tasks

- **Technical data** (see Page 55)
- **Transport** (see Page 56)
- **Storage** (see Page 56)
- **Unpacking** (see Page 56)
- **Positioning** (see Page 57)
- **Installation** (see Page 58)
- **Initial Operation** (see Page 59)
- **Use** (see Page 60)
- **Servicing** (see Page 71)
- **Disposal** (see Page 71)

#### Chapter containing general information

- **Appendices** (see Page 72)

## 1.4 Meanings of safety instructions and pictograms

### 1.4.1 Meanings of signal words in safety instructions

This document contains safety instructions that provide information about hazards connected to the product.

The hazard severity level and consequences are highlighted and classified in the safety instructions by means of signal words.

In the following subchapters the signal words and the structure of the safety instructions are explained.

#### 1.4.1.1 Signal word WARNING

 **WARNING** The signal word "WARNING" draws attention to dangerous actions or situations which could possibly have the severe consequences specified in the safety instruction (e.g. death or serious injury) in case of noncompliance.

#### 1.4.1.2 Signal word CAUTION

 **CAUTION** The signal word "CAUTION" draws attention to dangerous actions or situations which will have the minor consequences specified in the safety instruction (minor injury) in case of noncompliance.

#### 1.4.1.3 Signal word NOTICE

 **NOTICE** The signal word "NOTICE" draws attention to dangerous actions or situations which could possibly lead to the material damage specified in the safety instruction (e.g. damage to the product) in case of noncompliance.

#### 1.4.1.4 Structure of the safety instructions

In this document the following types of safety instructions are used:

- Safety instructions
- Embedded safety instructions

The following examples illustrate the manner in which these types of safety instructions are structured. The signal word "CAUTION" is used as an example.

#### Structure of a safety instruction

|   |  |
|---|--|
|  | <p> <b>CAUTION</b></p> <p><b>In the first paragraph under the signal word the type and source of the hazard are indicated. These are dangerous activities or situations connected with the product. The optional pictogram represents either the hazard type and source or the precautionary measure to be taken.</b></p> <p>In the second paragraph the consequences of disregarding the hazard are described.</p> <p>► From the third paragraph on handling instructions for avoiding the hazard are given.</p> |
|---|--|

#### Structure of an embedded safety instruction

The embedded safety instruction is only used in handling instructions.

- **CAUTION! The first part of the handling instruction is the embedded safety instruction. After the signal word a handling instruction is given as a measure for avoiding the hazard.**

The second part contains the handling instruction which is directly connected with the embedded safety instruction.

### 1.4.2 Meanings of pictograms

In this document symbols and pictograms which are either self-explanatory or whose meaning can be deduced from the context are used.

The chapter "Appendix" contains an explanation of the symbols and pictograms in the section "List of pictograms used" (see Page 72).

## 2 Product Description

### 2.1 Intended Use

- see fig. 3 (see page 2)

The 743A220 PROS.A. Assembly is exclusively intended for the alignment of modular leg prostheses by trained prosthetists (required technical knowledge: prosthetics lower extremities).

Using the optional mounting frame with 743A210 laser holder, the various adjustment options and scales of the product make it possible to adhere to the specifications which are provided in the associated documentation of the prosthetic components that are required for the alignment of the prosthesis (e. g. instructions for use, alignment recommendations, quick reference guide and technical information).

The exact documentation of the various settings simplifies the reproduction of a prosthesis.

The alignment measurements that have been used can be determined when a finished prosthesis is clamped into the alignment apparatus.

Prostheses for the following types of fittings can be aligned according to the amputation level of the prosthesis wearer:

- ① Transtibial fitting – TT (transtibial)
- ② Knee disarticulation fitting (figure shows KD prosthetic knee joint)
- ③ Long residual limb fitting (figure shows ST prosthetic knee joint)
- ④ Transfemoral prosthetic fitting – TF (transfemoral)
- ⑤ Hip disarticulation fitting (figure shows HD prosthetic knee joint)

### 2.2 Product overview

- see fig. 4 (see page 2)

|   |   |   |   |
|---|---|---|---|
| ① | Base frame  | ② | Retaining bits for prosthetic knee joints                                       |
| ③ | Height-adjustable holder for the socket clamping fixtures                       | ④ | 743Y594 socket clamping fixture for prosthetic socket diameter of 145 to 200 mm |
| ⑤ | 743Y594=1 socket clamping fixture for prosthetic socket diameter of 90 to 145mm | ⑥ | 743Y627 centring unit   |
| ⑦ | Height-adjustable holder for the knee bracket                                   | ⑧ | Knee bracket for the prosthetic knee joint                                      |
| ⑨ | Height-adjustable holder for the footplate                                      | ⑩ | Footplate   |

### 2.3 Accessories

|   |
|---|
| <b>INFORMATION</b>  |
| The accessories described in the following subsections are not included in the scope of delivery. |

#### 2.3.1 Mounting frame with laser holder

- see fig. 5 (see page 3)

The 743A210 mounting frame with laser holder has 3 lasers. It is installed on the PROS.A. Assembly. The lasers project the alignment reference line on to the prosthesis and thereby enable an exact positioning of the prosthetic components to the alignment reference line.

#### 2.3.2 Spacer

- see fig. 6 (see page 3)

The 743Y621=1 spacer is used if the 743Y621 standard spacer is not sufficient to centre the prosthetic socket. It is mounted on the 743Y627 centring unit.

#### 2.3.3 Casters

- see fig. 7 (see page 3)

The casters are intended as a replacement for the feet of the product. They make it possible to easily transport the PROS.A. Assembly to a different installation location in the workshop. The 743Y741 set includes the casters together with the required screws and slot nuts.

### 3 Safety

#### **⚠ WARNING**

##### **Damage to the housing or the optic of the lasers**

Severe injury to eyes through laser beam

- ▶ Do not look directly into the laser beam.
- ▶ Do not switch on the laser if the housing or the optic of the laser are damaged.
- ▶ Switch off the laser immediately if you notice any changes to the laser beam.
- ▶ Replace a damaged laser immediately.

#### **NOTICE**

##### **Severe vibration**

Inaccurate alignment of the prosthesis through distorted or damaged product components

- ▶ Do not expose the product to severe vibration.

#### **NOTICE**

##### **Unallowable environmental conditions**

Inaccurate alignment of the prosthesis through distorted product components, damage to product through unallowable environmental conditions (e. g. moisture)

- ▶ Do not expose the product to unallowable environmental conditions (see chapter "Environmental conditions" - see Page 55; chapter "Storage" - see Page 56).

### 4 Technical data

#### 4.1 Dimensions and weight

- see fig. 8 (see page 3)

|  |      |
|--|------|
| ① Height (with socket clamping fixture) [mm] | 1959 |
| ② Height (to mounting plate) [mm]            | 1659 |
| ③ Width (with crank handle) [mm]             | 724  |
| ④ Depth [mm]                                 | 684  |
| ⑤ Width (base frame) [mm]                    | 585  |
| Weight (without packaging) [kg]              | 45   |

#### 4.2 Environmental conditions

|   |  |
|---|--|
| Admissible temperature range (storage)              | -10 – +60 °C   |
| Admissible temperature range (operation)            | +5 – +40 °C  |
| Recommended maximum temperature (operation)         | +22 °C   |
| Allowable relative humidity (storage and operation) | 70%, non-condensing                                  |
| Unallowable environmental conditions                | Potentially explosive atmospheres, humid environment |

## 5 Transport

- see fig. 9 (see page 4)

### **⚠ WARNING**

#### **Standing under suspended loads**

Death or severe injury due to falling or tipping product

- ▶ Make sure there is no one in the region of suspended loads during transport.

- ▶ **WARNING! The product must be transported carefully to avoid it tipping over. Several persons are required for transport. Only use transport aids with sufficient load carrying capacity.**

Transport the transport crate ① with a lift truck ② or a forklift ③.

## 6 Storage

- see fig. 10 (see page 4)
- ▶ For storage observe the technical data in the chapter "Environmental conditions" (see Page 55).
- ▶ Protect the product from moisture and extreme temperature fluctuations.
- ▶ To prevent the product being contaminated, store the product in its packaging or under a protective film.

## 7 Unpacking

- see fig. 11 (see page 4)
- ▶ **NOTICE! Do not damage the product while unpacking it.**
- ▶ Unscrew the Phillips head screws using a cordless screwdriver or a screwdriver.
- ▶ Open the transport crate.
- ▶ Remove the cardboard box containing the accessories ②.
- ▶ Remove the posters ③.
- ▶ Remove the base frame ① with the aid of another person.

## 7.1 Checking the delivered items

► Check the delivery for completeness.

The following products are included in the delivery:

- see fig. 12 (see page 4)
- Base frame
- ① 743Y594 socket clamping fixture
- ② Knee bracket for the prosthetic knee joint
- see fig. 13 (see page 5)
- 646F219=D TF alignment poster
- 646F219=GB TF alignment poster
- 646F336=D TT alignment poster
- 646F336=GB TT alignment poster
- see fig. 14 (see page 5)
- ① 647G1000 instructions for use
- ② 743Y627 centring unit
- ③ 743Y594=1 socket clamping fixture
- see fig. 15 (see page 5)
- ① 743Y579=R retaining bit
- ② 743Y579=L retaining bit
- ③ 743Y580=R retaining bit
- ④ 743Y580=L retaining bit
- see fig. 16 (see page 5)
- ① 743Y581=R retaining bit
- ② 743Y581=L retaining bit
- ③ 743Y582=R retaining bit
- ④ 743Y582=L retaining bit
- see fig. 17 (see page 6)
- ① 743Y583=R retaining bit
- ② 743Y583=L retaining bit
- see fig. 18 (see page 6)
- ① 743Y584=R retaining bit
- ② 743Y584=L retaining bit
- see fig. 19 (see page 6)
- ① 743Y608=R retaining bit
- ② 743Y608=L retaining bit
- ③ 743Y629=R retaining bit
- ④ 743Y629=L retaining bit
- see fig. 20 (see page 6)
- ① 743Y681 retaining bit
- ② 743Y681 retaining bit

## 8 Positioning

- see fig. 21 (see page 7)
- **INFORMATION: The recommended distances to the walls ensure safe and comfortable working conditions.**  
Set up the product at a suitable location in the workshop.

## 9 Installation

### 9.1 Move the crank handle to the working position

- see fig. 22 (see page 7)
- ▶ Fold up the crank handle from the transport position.
- ▶ Allow the crank handle to engage in the working position.

### 9.2 Adjusting the feet of the product

- see fig. 23 (see page 7)
- ▶ Align the PROS.A. Assembly with a spirit level and by adjusting the feet of the product.

### 9.3 Installing the mounting frame with laser holder

- see fig. 24 (see page 7)
- ▶ Have the required materials to hand (6x AG13).
- ▶ ❶ Unscrew the battery compartment.
- ▶ ❷ Insert the batteries (2x AG13) in the battery compartment.
- ▶ ❸ Screw the battery compartment onto the laser housing.
- see fig. 25 (see page 8)
- ▶ Have the required materials to hand to install the mounting frame on the base frame.
- see fig. 26 (see page 8)
- ▶ Thread the slot nuts into the profiles from the side.
- ▶ Ensure the correct alignment of the slot nuts in the profile.
- ▶ Position the slot nuts in the profiles by moving them so that they are aligned with the positions of the screws in the mounting frame.
- see fig. 27 (see page 8)
- ▶ Insert the screws in the mounting frame.
- ▶ Position the mounting frame on the base frame.
- ▶ Screw the screws into the slot nuts.
- ▶ Correctly position the arrows on the mounting frame to the lower edge of the mounting plate.
- ▶ **NOTICE! The base frame must stand on a level surface and be aligned with a spirit level.**  
Align the mounting frame with a spirit level.
- ▶ Tighten the screws in the slot nuts.

### 9.4 Installing the casters

- see fig. 28 (see page 8)
- ▶ Have the materials to hand that are required to install the casters on the base frame (743Y741).
- see fig. 29 (see page 9)
- ▶ ❶ Remove the cover caps from the base frame.
- ▶ ❷ Unscrew the screws.
- ▶ ❸ Remove the feet of the product from the base frame.
- ▶ ❹ Remove the slot nuts from the base frame.
- see fig. 30 (see page 9)
- ▶ ❶ Insert the slot nuts into the base frame and position them.
- ▶ **INFORMATION: The best position for the swivel casters is at the front and for the fixed casters at the back.**
- ▶ ❷ Position the casters at the profiles.
- ▶ ❸ Tighten the screws in the slot nuts.
- ▶ ❹ Push the cover caps onto the base frame.

### 9.5 Using the mounting plate

- see fig. 31 (see page 9)
- ▶ **INFORMATION: The figure shows the recommended positioning.**  
Insert the retaining bits, the centring unit and the socket clamping fixture into the mounting plate.

## 10 Initial operation

### 10.1 Installing the laser gauge

- see fig. 32 (see page 9)
- ▶ ❶ Loosen the locking mechanism of the socket clamping fixture with the clamping lever.
- ▶ ❷ Remove the socket clamping fixture from the holder.
- see fig. 33 (see page 10)
- ▶ Remove the knee bracket from the holder.
- see fig. 34 (see page 10)
- ▶ Unscrew the star handle.
- ▶ Remove the star handle from the footplate holder.
- see fig. 35 (see page 10)
- ▶ Remove the footplate.
- ▶ Insert the 743Y689 laser gauge in the footplate holder.
- ▶ Align the laser gauge according to the markings.

### 10.2 Adjusting the lasers

- see fig. 36 (see page 10)
- ▶ Switch on the anterior laser.
- see fig. 37 (see page 11)
- ▶ ❶ Loosen the knurled screw.
- ▶ ❷ Adjust the anterior laser so that the laser line runs through the 0-line of the scale and the marking of the laser gauge.
- ▶ ❸ Tighten the knurled screw.
- ▶ ❹ Switch off the laser.
- see fig. 38 (see page 11)
- ▶ Switch on the right laser.
- see fig. 39 (see page 11)
- ▶ ❶ Loosen the knurled screw.
- ▶ ❷ Adjust the right laser so that the laser line runs through the marking of the laser gauge.
- ▶ ❸ Tighten the knurled screw.
- ▶ ❹ Switch off the laser.
- see fig. 40 (see page 11)
- ▶ Switch on the left laser.
- see fig. 41 (see page 12)
- ▶ ❶ Loosen the knurled screw.
- ▶ ❷ Adjust the left laser so that the laser line runs through the marking of the laser gauge.
- ▶ ❸ Tighten the knurled screw.
- ▶ ❹ Switch off the laser.

### 10.3 Removing the laser gauge

- see fig. 42 (see page 12)
- ▶ Remove the laser gauge from the footplate holder.
- ▶ Place the footplate on the footplate holder.
- see fig. 43 (see page 12)
- ▶ Place the star handle in the footplate holder.
- ▶ Tighten the star handle.
- see fig. 44 (see page 12)
- ▶ Insert the knee bracket in the holder.
- see fig. 45 (see page 13)
- ▶ ❶ Place the socket clamping fixture into the holder.
- ▶ ❷ Tighten the locking mechanism of the socket clamping fixture with the clamping lever.

## 11 Use

### 11.1 Information for use

- ▶ Check the lasers daily prior to use (chapter "Checking the lasers" - see Page 71)
- ▶ For use of the product, also observe the instructions for use of the 743A210 mounting frame (647G999).
- ▶ For use of the product, also observe the quick reference guide of the line laser (647G999 instructions for use – chapter "Quick reference guide line laser").
- ▶ If required, replace weak or empty batteries and dispose of them correctly.
- ▶ Remove the batteries during longer periods of non-use.
- ▶ Switch off the laser after use to avoid unnecessary discharge of the batteries.

### 11.2 Aligning a prosthesis

#### 11.2.1 Important information about aligning a prosthesis

#### ⚠ CAUTION

##### Deviation from alignment recommendations

Patient has a fall due to wrongly aligned prosthesis

- ▶ Align the prosthesis according to the alignment recommendations (alignment poster, instructions for use of the prosthetic components).
- ▶ Apply the most recent alignment recommendations, by using the instructions for use provided with the prosthetic component.
- ▶ Take all safety measures required for the patient (e. g. use of parallel bars).
- ▶ Take note that any deviation from the alignment instructions takes place at your own risk. Deviations must be based on technical expertise and medical necessity.

Fabricating a prosthesis that is optimally tailored to the patient demands technical knowledge of orthopaedic technology. The alignment and adjustments may only be carried out by a prosthetist.

The alignment and adjustments are mainly completed in the following steps:

1. Bench alignment
2. Static alignment
3. Dynamic trial fitting

These steps are first carried out with a test prosthesis in order to determine the best possible combination and positioning of the prosthetic components relative to each other. Once all settings have been tailored to the patient, the definitive prosthesis is produced. The same steps are repeated for that purpose.

In a dynamic trial fitting, the patient must be given sufficient time to become familiar with the functions of the prosthesis and thereby learn how to use it safely in everyday life.

## 11.2.2 Positioning the prosthetic foot

### 11.2.2.1 Determining the required measurements of the shoe

- see fig. 46 (see page 13)
- > **Recommended tools:**  
743S12 heel height measuring device
- ▶ ❶ Measure the posterior heel height  $h_{h1}$  of the shoe.
- ▶ ❷ Measure the anterior heel height  $h_{h2}$  of the shoe.
- ▶ ❸ Calculate the effective heel height  $h_{eh}$  ( $h_{eh} = h_{h1} - h_{h2}$ ).
- ▶ **INFORMATION: The 5 mm are intended as a safety margin for the patient.**
- ▶ ❹ Calculate the value for the total heel height  $h_{th}$ , which must be set on the PROS.A. Assembly ( $h_{th} = h_{eh} + 5$  mm).

### 11.2.2.2 Determining the required measurements of the prosthetic foot

- see fig. 47 (see page 13)

#### Prosthetic feet with markings (a–p, m–l) on the footshell ❶ and the connection cap ❷

- ▶ Continue with chapter "Adjusting the height of the footplate" (see Page 61), because the calculation is not required due to the m-l marking.

#### Prosthetic feet without markings (a–p, m–l) on the footshell and the connection cap

- see fig. 48 (see page 13)
- ▶ ❶ Calculate the middle of the prosthetic foot  $f_{s/2}$  ( $f_{s/2} = f_s / 2$ ).
- ▶ ❷ Use the instructions for use of the prosthetic foot to check if the default value (30 mm) has to be used for the anterior placement ( $f_{m-al}$ ) of the middle of the prosthetic foot ❷ to the alignment reference line ❶, or if a different value has to be used.
- ▶ ❸ Calculate the setting value  $f_{al-h}$  (distance of alignment reference line to heel) for the footplate at the PROS.A. Assembly ( $f_{al-h} = f_{s/2} - f_{m-al}$ ).

### 11.2.2.3 Adjusting the height of the holder for the footplate

- see fig. 49 (see page 14)

|   |
|---|
| <b>INFORMATION</b>  |
| The holder for the footplate is height-adjustable. In many cases, sliding it upwards can result in a more comfortable working position. |

- ▶ ❶ Loosen the left locking mechanism with the left clamping lever.
- ▶ ❷ Loosen the right locking mechanism with the right clamping lever.
- ▶ ❸ Position the holder for the footplate at the required working height.
- ▶ ❹ Tighten the left locking mechanism with the left clamping lever.
- ▶ ❺ Tighten the right locking mechanism with the right clamping lever.

### 11.2.2.4 Setting the heel height

- see fig. 50 (see page 14)
- ▶ ❶ Loosen the thumb screw.
- ▶ **INFORMATION: The heel height adjustment column is grooved at 5 mm intervals for easier positioning.**
- ▶ ❷ Position the adjustment plate so that the distance between the adjustment plate and the footplate corresponds to the total heel height  $h_{th}$ .
- ▶ ❸ Tighten the thumb screw.

### 11.2.2.5 Positioning the footplate in a-p direction

- see fig. 51 (see page 14)
- ▶ Leave the external rotation of the foot at 0° during alignment in the a-p direction.
- ▶ ❶ Loosen the star handle.
- ▶ **Positioning of prosthetic feet with markings (a-p, m-l):**  
Continue with operating instructions in chapter "Positioning the prosthetic foot" (see Page 62).
- ▶ **Positioning of prosthetic feet without markings (a-p, m-l):**  
❷ Position the footplate in a-p direction according to the setting value  $f_{al-h}$  (Distance of alignment reference line to heel).

### 11.2.2.6 Positioning the prosthetic foot

- see fig. 52 (see page 14)
- ▶ Switch on the anterior laser to project the frontal alignment reference line.
- see fig. 53 (see page 15)
- ▶ For a right prosthetic foot, switch on the left laser to project the sagittal alignment reference line.
- see fig. 54 (see page 15)
- ▶ For a left prosthetic foot, switch on the right laser to project the sagittal alignment reference line.
- see fig. 55 (see page 15)
- ▶ Position the prosthetic foot on the footplate so that the heel is flush with the adjustment plate.
- see fig. 56 (see page 15)
- ▶ **Positioning of prosthetic feet with markings (a-p, m-l):**  
Leave the external rotation of the foot at 0° when positioning the prosthetic foot.
- ▶ ❶ Position the prosthetic foot on the footplate, so that the frontal alignment reference line ❷ runs through the a-p markings of the prosthetic foot.
- ▶ ❷ Move the footplate in the a-p direction, so that the sagittal alignment reference line ❶ runs through the m-l markings on the prosthetic foot.
- see fig. 57 (see page 16)
- ▶ **Positioning of prosthetic feet without markings (a-p, m-l):**  
Leave the external rotation of the foot at 0° when positioning the prosthetic foot.
- ▶ Position the prosthetic foot on the footplate, so that the frontal alignment reference line ❷ runs through the middle of the longitudinal axis of the prosthetic foot (middle of the heel and the pyramid adapter).

### 11.2.2.7 Setting the external rotation of the foot

- see fig. 58 (see page 16)
- ▶ **INFORMATION: The default value for the external rotation of the foot is 7°.**  
Adjust the footplate to the required external rotation of the foot, depending on the side of the prosthetic foot.
- see fig. 59 (see page 16)
- ▶ Tighten the star handle.

## 11.2.3 Positioning the prosthetic knee joint

### 11.2.3.1 Determining the required measurements for the prosthetic knee joint

- see fig. 60 (see page 16)
- ▶ ❶ Measure the MTP-floor measurement  $h_{k-f}$  at the patient.
- ▶ ❷ Measure the socket end-to-floor measurement  $h_{s-f}$  at the patient.
- ▶ ❸ The distance  $h_{r-k}$  between the medial tibial plateau and the alignment reference point (compromise knee rotation point) is 20 mm.
- ▶ ❹ Determine the proximal system height  $h_{ps}$  according to the prosthetic knee instructions for use.
- ▶ ❺ Determine the distal system height  $h_{ds}$  according to the prosthetic knee instructions for use.
- ▶ ❻ Calculate the distance  $h_{s-r}$  between the socket end and the alignment reference point ( $h_{s-r} = h_{k-f} - h_{r-k}$ ).
- ▶ ❼ Calculate the distance  $h_{r-f}$  between the alignment reference point and the floor ( $h_{r-f} = h_{s-f} - h_{s-r}$ ).

#### INFORMATION

When dealing with knee disarticulation fittings and long residual limb fittings, it is not possible to position the alignment reference point at the calculated height. In such cases, it is necessary to position it low enough to leave enough room for the adapter above the prosthetic knee.

Using the sum  $h_{s-r}$  of the system heights of the adapters that are used proximally, and the proximal system height of the prosthetic joint, the alignment reference point is calculated using the following formula:

$$h_{r-f} = h_{s-f} - h_{s-r}$$

### 11.2.3.2 Adjusting the height of the holder for the knee bracket

- see fig. 61 (see page 17)
- ▶ ❶ Calculate the required distance  $h_{kf}$  between the height  $H_f$  of the holder for the footplate and the height  $H_k$  of the holder for the knee bracket of the prosthetic knee joint ( $h_{kf} = h_{r-f} + h_{th}$ ).
- ▶ ❷ Read off the height  $H_f$  of the holder for the footplate on the scale of the PROS.A. Assembly.
- ▶ ❸ Calculate the required height  $H_k$  of the holder for the knee bracket of the prosthetic knee ( $H_k = H_f + h_{kf}$ ).
- see fig. 62 (see page 17)
- ▶ ❶ Loosen the left locking mechanism with the left clamping lever.
- ▶ ❷ Loosen the right locking mechanism with the right clamping lever.
- ▶ ❸ Position the holder for the knee bracket at the required height  $H_k$ .
- ▶ ❹ Tighten the left locking mechanism with the left clamping lever.
- ▶ ❺ Tighten the right locking mechanism with the right clamping lever.

### 11.2.3.3 Mounting the knee bracket

- see fig. 63 (see page 17)
- ▶ **INFORMATION: The holder is configured for a fixed value of 5° external rotation of the prosthetic knee joint.**  
Insert the knee bracket into the holder depending on the prosthesis side of the prosthetic foot.

#### 11.2.3.4 Mounting the retaining bits

- see fig. 64 (see page 17)
- ▶ **INFORMATION: The retaining bits are marked with the reference numbers of the prosthetic knee joints for which they are suitable. The instructions for use of newer prosthetic knee joints contain appropriate information on the retaining bits that must be used.**  
Select the retaining bits intended for the prosthetic knee joint.
- see fig. 65 (see page 18)
- ▶ **CAUTION! Risk of fingers and hands being crushed when inserting retaining bits**  
**INFORMATION: The arms of the knee bracket and the retaining bits are marked with one or two points to provide a positioning aid. Only the 743Y681 retaining bits are not marked and can be used on the left and the right side.**  
Insert the retaining bits in the knee bracket.

#### 11.2.3.5 Clamping the prosthetic knee joint

- see fig. 66 (see page 18)
- ▶ **CAUTION! Risk of fingers and hands being crushed when clamping the prosthetic knee joint**  
Clamp the prosthetic knee joint into the retaining bits of the knee bracket.

#### 11.2.4 Installing the prosthetic components proximal to the prosthetic foot

- see fig. 67 (see page 18)
- ▶ Select the prosthetic components to be installed proximal to the prosthetic foot according to the measured and calculated values (see also 646K2\* catalogue) and assemble them according to their instructions for use.

##### Determine the required length of the tube adapter

- ▶ Install all adapters between the prosthetic foot and the prosthetic knee joint.
- ▶ **CAUTION! The tube adapter must be inserted fully into the tube clamp adapter or the prosthetic knee joint.**  
Determine the distance between the height of the holder for the footplate  $H_f$  and the height of the holder for the knee bracket  $H_k$  and in relation to this, calculate the difference to the distance value  $j_{kf}$ .
- ▶ Shorten the tube adapter by the difference calculated.

##### Transtibial fittings (TT)

- ▶ Select the prosthetic components to be installed proximal to the prosthetic foot according to the measured and calculated values (MTP-floor measurement  $h_{k-f}$ , socket end-to-floor measurement  $h_{s-f}$  and MTP-socket end measurement  $h_{k-s}$ ) (see also 646K2\* catalogue) and assemble them according to their instructions for use.

#### 11.2.5 Installing the prosthetic components above the prosthetic knee joint

- see fig. 68 (see page 18)
- ▶ Select the prosthetic components to be installed proximal to the prosthetic knee joint according to the measured and calculated values (see also 646K2\* catalogue) and assemble them according to their instructions for use.

## 11.2.6 Positioning of the prosthetic socket

### 11.2.6.1 Selecting and installing socket clamping fixture

- see fig. 69 (see page 19)
- ▶ Adjust the height of the holder for the socket clamping fixture with the help of the crank handle in such a way that the prosthetic socket can be positioned without any difficulty later.
- see fig. 70 (see page 19)

#### INFORMATION

The socket clamping fixtures supplied are particularly suited for transfemoral fittings. The 743Y594 socket clamping fixture is used for a socket diameter of 145 to 200 mm. The 743Y594=1 socket clamping fixture is used for a socket diameter of 90 to 145 mm. It is important that the prosthetic socket is always securely fixed in position. If required, additional measures must be taken to securely position and hold the prosthetic socket in the alignment apparatus.

- ▶ Select the socket clamping fixture that is suited to the specific application.
- ▶ Place the socket clamping fixture into the holder.
- ▶ Tighten the locking mechanism of the socket clamping fixture with the clamping lever.

### 11.2.6.2 Configuring and installing the centring unit

- see fig. 71 (see page 19)
- ▶ Configure the centring unit according to the length and diameter of the prosthetic socket.
- ▶ Install the centring unit on the socket clamping fixture.

### 11.2.6.3 Marking the prosthetic socket and determining important measurements

- see fig. 72 (see page 19)

#### All prosthetic sockets

- ▶ Mark two points both frontally as well as laterally (① - distal, ② - proximal).
- ▶ Draw a line ③ through both points.

#### TF prosthetic socket only

- ▶ Draw a third point ④ on the lateral line at the height of the tuberosity.
- ▶ Draw the lateral socket reference point ⑤ 30 mm above this point.
- ▶ Draw the frontal socket reference point ⑥ on the intersection of the frontal line with the prosthetic socket edge.
- see fig. 73 (see page 20)

#### Hip flexion angle $\alpha$ in transfemoral fittings (TF)

- ▶ Flex the contralateral knee joint and hip joint.
- ▶ With supine patients, place a flat hand under the lumbar spine.
- ▶ Flex the residual limb until the hollow back disappears.
- ▶ Carefully bring the residual limb into extension until the lumbar spine starts to flex back in the direction of the hollow back.
- ▶ Determine the flexion angle  $\alpha$  in this position.

#### Individual socket flexion $\alpha$ with transtibial fittings (TF)

- ▶ Measure the socket flexion  $\alpha$  when patient is standing.

### 11.2.6.4 Securing the prosthesis socket on the socket clamping fixture

- see fig. 74 (see page 20)
- ▶ ① Loosen the retaining screw ①.
- ▶ ② Loosen the clamping lever for the a-p positioning of the socket clamping fixture.
- ▶ ③ Loosen the clamping lever for the m-l positioning of the socket clamping fixture.
- ▶ ④ Position the socket clamping fixture so that all positioning is in the starting position (0 ° or 0 mm).
- ▶ ⑤ Tighten the clamping lever for the m-l positioning of the socket clamping fixture.
- ▶ ⑥ Tighten the clamping lever for the a-p positioning of the socket clamping fixture.
- ▶ ⑦ Tighten the retaining screw ①.
- see fig. 75 (see page 20)
- ▶ ① Push the prosthetic socket onto the socket clamping fixture.
- ▶ ② Use the bellows to pump up the socket clamping fixture so that the prosthetic socket is securely fixed on the socket clamping fixture.
- ▶ ③ Position the prosthetic socket on the socket clamping fixture so that the frontal alignment reference line ① runs along or parallel to the frontal middle line ② of the prosthetic socket.
- ▶ ④ Position the prosthetic socket on the socket clamping fixture so that the sagittal alignment reference line ① runs along or parallel to the lateral middle line ② of the prosthetic socket.

### 11.2.6.5 Positioning the prosthetic socket to the alignment reference line

- see fig. 76 (see page 20)

#### INFORMATION

The rotation axes of the socket clamping devices in a-p direction ① and m-l direction ② lie at the height of the top ③ of the upper display plate.

- ▶ Loosen the retaining screw ④.
- ▶ Follow the alignment instructions to adjust the socket clamping fixture to the required socket flexion in relation to the alignment reference line ⑤.
- ▶ Tighten the retaining screw ④.
- see fig. 77 (see page 21)

#### INFORMATION

The scale ① provides a rough overview of the angle setting. The scale ② is used to set the angle in degrees.

- ▶ ① Loosen the clamping lever for the a-p positioning of the socket clamping fixture.
- ▶ ② Adjust the socket clamping fixture in relation to abduction and adduction according to the patient measurements and the alignment instructions.
- ▶ ③ Tighten the clamping lever for the a-p positioning of the socket clamping fixture.
- see fig. 78 (see page 21)
- ▶ ① Loosen the clamping lever for the a-p positioning of the socket clamping fixture.
- ▶ ② Move the socket clamping fixture in the a-p direction, until the prosthetic socket is positioned according to the alignment instructions in relation to the sagittal alignment reference line.
- ▶ ③ Tighten the clamping lever for the a-p positioning of the socket clamping fixture.
- see fig. 79 (see page 21)
- ▶ ① Loosen the clamping lever for the m-l positioning of the socket clamping fixture.
- ▶ ② Move the socket clamping fixture in the m-l direction, until the prosthetic socket is positioned according to the alignment instructions in relation to the frontal alignment reference line.
- ▶ ③ Tighten the clamping lever for the m-l positioning of the socket clamping fixture.


### 11.2.6.6 Adjusting the height of the holder for the socket clamping fixture

- see fig. 80 (see page 21)
- ▶ Adjust the height of the holder for the socket clamping fixture with the help of the crank handle in such a way that the prosthetic socket is positioned on the prosthesis according to the alignment instructions and patient measurements.

### 11.2.7 Documenting the settings





- ▶ Document the following measurements and scale values in relation to the prosthesis produced:



#### Settings at the PROS.A. Assembly, prosthesis measurements, patient measurements

|   |   |
|---|---|
|  | <p>Prosthesis side (left/right)</p> <hr/> |
|---|---|

| Settings at the PROS.A. Assembly, prosthesis measurements, patient measurements                 |  |
|---|--|
| <p> <math>h_{eh} = h_{h1} - h_{h2}</math><br/> <math>h_{th} = h_{eh} + 5 \text{ mm}</math> </p> | <p>Posterior heel height (<math>h_{h1}</math>) [mm]</p> <hr/> <p>Anterior heel height (<math>h_{h2}</math>) [mm]</p> <hr/> <p>Effective heel height (<math>h_{eh}</math>) [mm]</p> <hr/> <p>Footplate - total heel height (<math>h_{th}</math>) [mm]</p> <hr/> |
|   | <p>Foot size (<math>f_s</math>) [mm]</p> <hr/> <p>Distance of the alignment reference line to the heel (<math>f_{al-h}</math>) [cm]</p> <hr/>  |
|   | <p>Footplate - external rotation [°]</p> <hr/>   |
|   | <p>Footplate - height (<math>H_f</math>) [cm]</p> <hr/>  |

**Settings at the PROS.A. Assembly, prosthesis measurements, patient measurements**

|   |  |
|---|--|
|    | <p>Knee bracket - height (<math>H_k</math>) [cm]</p> <hr/>                             |
|   | <p>Socket clamping fixture - socket flexion [°]</p>                                    |
|  | <p>Socket clamping fixture - abduction and adduction [°]</p>                           |
|  | <p>Socket clamping fixture - position in a-p direction (<math>s_{a-p}</math>) [mm]</p> |

| Settings at the PROS.A. Assembly, prosthesis measurements, patient measurements    |  |
|--|--|
|   | Socket clamping fixture - position in m-l direction [mm] |
|  | Socket clamping fixture - height ( $H_s$ ) [cm]          |

### 11.3 Measuring a prosthesis

- ▶ Clamp the prosthesis into the PROS.A. Assembly. While doing so, observe the operating instructions in chapter "Aligning a prosthesis" (see Page 60) and its subsections.
- ▶ Document the settings according to the chapter "Documenting the settings" (see Page 67).

### 11.4 Cleaning

#### Cleaning the laser optic

- 1) In case of dust, use an optical brush with bellows.
- 2) In case of light soiling, use a dry cleaning cloth for lenses or glasses.
- 3) In case of heavier soiling, use a moist cleaning cloth for lenses or glasses.
- 4) Clean the laser optic.

#### Cleaning light soiling

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

#### Cleaning heavier soiling

- > **Required materials:** Cleaning cloth, 634A58 Isopropyl Alcohol, soft cloth
- 1) **NOTICE! Verify material compatibility! Do not remove lubricant from bearing points!**  
Clean the product with a cleaning cloth and isopropyl alcohol.
- 2) Dry the product with a cloth.

## 12 Maintenance

Product maintenance encompasses the following areas:

- **Inspection**  
During the inspection the product is checked to determine whether it is in good condition.
- **Maintenance**  
During maintenance actions are taken to maintain the good condition of the product.
- **Repair**  
Repair refers to repair measures taken to restore the product to good condition.

### 12.1 Inspection

#### 12.1.1 Checking the lasers

- > **Inspection interval:** daily prior to use
- ▶ Check the lasers individually.
- ▶ Observe the quick reference guide of the line laser (chapter "Quick reference guide line laser" - Quick reference guide line laser).
- ▶ Check the laser optic for soiling and clean if required (chapter "Cleaning" - see Page 70).
- ▶ Check the course of the laser beam (marking on the footplate and for the anterior laser using the scale of the height-adjustable holder for the socket clamping fixtures - detailed information in chapter "Adjusting the lasers - see Page 59).

### 12.2 Maintenance

The product requires no maintenance.

### 12.3 Repair

|  |  |
|--|--|
|  | <p><b>⚠ WARNING</b></p> <p><b>Repair by inappropriate personnel</b><br/>Serious injury as well as damage to the product due to insufficient expertise</p> <p>▶ Repair work may only be carried out by a service technician authorised by the manufacturer.</p> |
|--|--|

|                    |   |
|--------------------|---|
| <b>INFORMATION</b> | <p>The 743Y643 repair set is available for the repair of the silicone cuffs of the socket clamping fixture.</p> |
|--------------------|---|

## 13 Disposal

|   |   |
|---|---|
|  | <p><b>INFORMATION</b></p> <p>This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.</p> |
|---|---|

## 14 Legal Information

### 14.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document (Instructions for Use). The manufacturer will not assume liability for damage caused by disregard of this document.

## 15 Appendices

### 15.1 List of pictograms used

In the following subchapters the symbols and pictograms used in this document along with their meanings are given.

#### 15.1.1 Warning symbols

|  |   |
|--|---|
|   | <ul style="list-style-type: none"> <li>• Tip hazard warning</li> </ul>                                |
|   | <ul style="list-style-type: none"> <li>• Lower extremity crush hazard warning</li> </ul>              |
|   | <ul style="list-style-type: none"> <li>• Upper extremity (e.g. hands) crush hazard warning</li> </ul> |
|   | <ul style="list-style-type: none"> <li>• Upper extremity cut hazard warning</li> </ul>                |
|  | <ul style="list-style-type: none"> <li>• Laser radiation warning</li> </ul>                           |

#### 15.1.2 Mandatory symbols

|   |   |
|---|---|
|  | <ul style="list-style-type: none"> <li>• Read and observe this document</li> </ul>                |
|  | <ul style="list-style-type: none"> <li>• Work in pairs or with a second person present</li> </ul> |

#### 15.1.3 Prohibition symbols

|   |  |
|---|--|
|  | <ul style="list-style-type: none"> <li>• Standing in the vicinity or under suspended loads prohibited</li> </ul> |
|---|--|

#### 15.1.4 Transport and storage pictograms

|   |   |
|---|---|
|  | <ul style="list-style-type: none"> <li>• Forklift truck fork insertion point</li> </ul> |
|---|---|

|   |  |
|---|--|
|  | <ul style="list-style-type: none"> <li>• This side up for transport and storage</li> </ul> |
|  | <ul style="list-style-type: none"> <li>• Protect from moisture</li> </ul>                  |

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Otto Bock HealthCare GmbH