


## 28F10N

|           |                                  |    |
|-----------|----------------------------------|----|
| <b>DE</b> | Gebrauchsanweisung .....         | 6  |
| <b>EN</b> | Instructions for use .....       | 11 |
| <b>FR</b> | Instructions d'utilisation ..... | 17 |
| <b>IT</b> | Istruzioni per l'uso .....       | 22 |
| <b>ES</b> | Instrucciones de uso .....       | 28 |
| <b>PT</b> | Manual de utilização .....       | 34 |
| <b>NL</b> | Gebruiksaanwijzing .....         | 39 |
| <b>SV</b> | Bruksanvisning .....             | 45 |
| <b>DA</b> | Brugsanvisning .....             | 50 |
| <b>PL</b> | Instrukcja użytkowania .....     | 55 |
| <b>HU</b> | Használati utasítás .....        | 61 |
| <b>CS</b> | Návod k použití .....            | 66 |
| <b>SK</b> | Návod na používanie .....        | 71 |
| <b>TR</b> | Kullanma talimatı .....          | 77 |
| <b>EL</b> | Οδηγίες χρήσης .....             | 82 |
| <b>RU</b> | Руководство по применению .....  | 88 |

|                 |                                    |
|-----------------|------------------------------------|
| <b>Material</b> | PP, PA, PUR, PE, NBR/SBR, EVA, PES |
|-----------------|------------------------------------|

**INFORMATION**

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### **INFORMATION**

Last update: 2016-06-14

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

These instructions for use contain important information on the fitting and application of the 28F10N heel relief brace designed by Dr. Settner and CPO Münch.

## **2 Product description**

### **2.1 Components**

| Item | Component                            | Item | Component                    |
|------|--------------------------------------|------|------------------------------|
| 1    | Longitudinal arch piece (left/right) | 9    | Compensation sole (Lunasoft) |
| 2    | Forefoot cap                         | 10   | Heel wedge (Lunasoft)        |
| 3    | Adhesive strips for calf pad         | 11   | Compression pads             |
| 4    | Mounting rivets                      | 12   | Plantar pad (foam)           |

| Item | Component                                                   | Item | Component                          |
|------|-------------------------------------------------------------|------|------------------------------------|
| 5    | Adjusting nuts and fastening screws (hook-and-loop closure) | 13   | Heel container                     |
| 6    | Fastening screws (forefoot cap)                             | 14   | Tailorable calf pad                |
| 7    | Hook-and-loop closure (instep)                              | 15   | Hook-and-loop closure (lower shin) |
| 8    | Self-adhesive hook-and-loop strap                           | 16   | Hook-and-loop closure (upper shin) |

## 2.2 Design

The heel relief brace is made of glass fiber-reinforced polypropylene.

The heel container (13) is designed to enclose the dorsal calf and support the longitudinal arch (1) of the foot, thus enabling suspended positioning of the calcaneus in the brace. Plantar foot support is provided by the molded foot piece for the longitudinal arch and the compensation sole (9), which can be affixed to the heel container. The heel wedge (10) can be used to compensate for the height of the shoe on the unaffected side.

To facilitate check-ups and for improved air circulation, the brace has an opening at the back of the heel container. The weight-bearing pad set for the heel, which includes a plantar pad (12) and compression pads (11), is applied according to the therapy plan and allows the weight borne by the calcaneus to be built up gradually.

The hook-and-loop closures (7, 15, and 16) and the forefoot cap (2) hold the foot in place. The zip fastener on the forefoot cap makes it easy to apply and remove the brace and also protects the toes. The enclosed fastening screws and adjusting nuts (5) are used to affix the hook-and-loop closure (7) to the brace. Optionally, this hook-and-loop closure can also be riveted. The hook-and-loop closures at shin level (15 and 16) can also be riveted or sewn together to increase stability. The tailorable calf pad (14) establishes contact with the heel container.

The underside of the brace features a non-skid profile for increased safety when walking.

## 3 Intended use

### 3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 3.2 Indications

- Functional follow-up treatment for calcaneal fractures (unilateral or bilateral), regardless of the type of fracture and course of primary treatment.
- Weight-bearing relief following arthrodesis of the lower ankle

The indication must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

- External fixator (wire osteosynthesis)

#### 3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

### 3.4 Mechanism of Action

The brace relieves the calcaneus by providing targeted support to the longitudinal arch beneath the os naviculare as well as a form-fitting enclosure for the metatarsals and the calf. Accelerated bone healing is achieved by applying the basic principle of the Ilizarov technique (compression and distraction of the fractured areas). An almost natural rollover of the injured foot with the brace exerts a positive influence on the neuromuscular structures and calcareous salt content of the bone. The functioning muscle pump may also contribute to thrombosis prophylaxis.

### 3.5 Treatment and therapy plan

The therapy plan is independent of the type of fracture and treatment approach.

After subsidence of swelling/removal of stitches: fitting of the brace (4 hrs)

|           |                                                                                                                                                   |
|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Days 8-12 | Full weight-bearing without crutches                                                                                                              |
| Week 4    | 1st x-ray evaluation (lateral view of upper ankle, special/lateral view of calcaneus, upper ankle under full weight bearing/Hallgrímsson imaging) |
| Week 6    | 1st compression pad                                                                                                                               |
| Week 8    | 2nd X-ray evaluation<br>2nd compression pad                                                                                                       |
| Week 10   | 3rd compression pad                                                                                                                               |
| Week 11   | Scan of the foot and construction of orthopedic shoes within 5-6 days                                                                             |
| Week 12   | Conclusion of medical treatment/test for fitness for work for worker's compensation cases                                                         |

(The time indications are from the BG Unfallklinik (hospital for accident cases Duisburg), Germany and are calculated from the date of the injury)

## 4 Safety

### 4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 General safety instructions

 **CAUTION**

#### Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

 **CAUTION**

#### Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

 **CAUTION**

#### Use of a product with neoprene

Allergic reactions due to neoprene or heat allergies

- ▶ Do not wear products with neoprene for more than 3 – 4 hours without interruption (unless explicitly instructed otherwise by the doctor).

- Do **not** use the product in case of a known neoprene or heat allergy.

#### NOTICE

#### Contact with ointments, lotions, or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- Do not expose the product to ointments, lotions, or other products that contain oils or acids.

#### NOTICE

#### Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- Avoid storage in condensing ambient humidity.
- Avoid contact with abrasive substances (e.g. sand, dust).
- Do not expose the product to temperatures below -10 °C (14 °F) or above +50 °C (122 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

## 5 Handling

#### INFORMATION

- The daily duration of use and period of application are generally determined by the physician.
- The initial fitting and application of the product must be carried out by qualified personnel.
- Instruct the patient in the handling and care of the product.
- Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 5.1 Selecting the Size

- Select the brace size based on the shoe size (see sizing table).

| Size | Shoe size | Size<br>Heel container      | Size<br>Longitudinal arch<br>piece        | Size<br>Forefoot cap           |
|------|-----------|-----------------------------|-------------------------------------------|--------------------------------|
| S    | 36 - 41   | Container size 1<br>(36-41) | Longitudinal arch<br>piece size S (36-41) | Forefoot cap size 1<br>(36-41) |
| M    | 42 - 46   | Container size 2<br>(42-50) | Longitudinal arch<br>piece size M (42-46) | Forefoot cap size 2<br>(42-46) |
| L    | 47 - 50   | Container size 2<br>(42-50) | Longitudinal arch<br>piece size L (47-50) | Forefoot cap size 3<br>(47-50) |

### 5.2 Adaptation

- 1) Place the longitudinal arch piece (1) in the heel container (13).

**INFORMATION: To determine the placement, affix the longitudinal arch piece to a piece of paper. Position the patient on the longitudinal arch piece and mark the contour of the foot using a pencil (see fig. 2). Then cut out the contour (see fig. 3) and place it in the heel container (see fig. 4) so the end of the heel meets the dot in the heel container.**

- 2) Mark the position of the longitudinal arch piece in the heel container (see fig. 5).
- 3) Fix the longitudinal arch piece along the mark inside the heel container.
- 4) Insert the plantar pad (12) into the heel container.

- 5) Carry out an initial weight bearing test with the longitudinal arch piece and the forefoot cap (2) in place.

**INFORMATION: When the heel is under the maximum weight, use the opening cut in the heel container to check the degree of relief provided by the brace.**

- 6) **Optional:** Adapt the heel container to the individual patient using milling (see fig. 9) or heat (see fig. 10), then widening it as necessary.
- 7) **Optional:** Secure the compensation sole (9) with the hook-and-loop closure (8) to further increase the height of the longitudinal arch.
- 8) Position and then secure the plantar pad (12) in the heel area so that it maintains continual contact with the heel.

**INFORMATION: During the course of therapy, compression pads (11) are inserted under the plantar pad according to the therapy plan and serve the gradual build-up of weight borne by the calcaneus.**

- 9) Mount the hook-and-loop closure and the guide loop (7) at the instep using the supplied adjusting nuts and screws (5) (see fig. 11).

**INFORMATION: The longitudinal arch piece must be removed temporarily before mounting the hook-and-loop closure.**

- 10) **Optional:** Cut the calf pad (14) to the required size (see fig. 12).
- 11) Attach the adhesive strips (3) and secure the calf pad in the heel container.
- 12) **Optional:** Cut the sides of the forefoot cap. The two pieces can then be riveted or sewn together.
- 13) Push the forefoot cap over the heel container until the desired shoe size is set (scaling displayed on the side of the heel container).
- 14) Fix the forefoot cap to the heel container using screws (6) (see fig. 16).
- 15) Carry out a weight bearing test with the fully assembled brace.

**INFORMATION: Use a strip of paper in the opening cut in the heel container to check the degree of relief provided by the brace (see fig. 17). Place the strip of paper between the heel and the compression pad. When the calcaneus is bearing the user's full weight, pull the paper strip out of the brace in one even stroke. There should be slight resistance as you pull the paper strip out of the brace, but it should not tear.**

### 5.3 Application

#### CAUTION

##### **Direct skin contact with the product**

Skin irritation due to chafing or sweat formation

- Do not wear this product directly on the skin.

#### CAUTION

##### **Improper molding or application**

Injuries or damage to the brace due to overloading of the material and improper fit of the brace due to breakage of load-bearing components

- Do not make any improper changes to the brace.
- Always apply the brace according to the information in the instructions.

#### CAUTION

##### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- Ensure that the brace is applied properly and fits correctly.

#### NOTICE

##### **Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Open the hook-and-loop closures and the zip fastener.
- 2) First, place the toes in the brace (see fig. 18).
- 3) Close the hook-and-loop closures, ensuring they exert even pressure (see fig. 19).
- 4) Close the zip fastener (see fig. 20).

## **5.4 Cleaning**

#### NOTICE

##### **Use of improper cleaning agents**

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

##### **Brace frame:**

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

##### **Straps and pads:**

- 1) Remove the straps and pads from the brace.
- 2) Close all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at 86 °F using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## **6 Disposal**

Dispose of the product in accordance with national regulations.

## **7 Legal Information**

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### **7.1 Liability**

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorized modification of the product.

### **7.2 CE Conformity**

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

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Otto Bock HealthCare GmbH