

**8E38, 8E39, 8E41, 8E44**

<b>DE</b> Patienteninformation .....	3
<b>EN</b> Patient information .....	8
<b>FR</b> Informations aux patients.....	12
<b>IT</b> Informazioni per il paziente .....	17
<b>ES</b> Información para el paciente.....	22
<b>PT</b> Informação para o paciente .....	27
<b>NL</b> Patiënteninformatie.....	32
<b>SV</b> Brukarinformation .....	37
<b>DA</b> Patientinformation .....	42
<b>PL</b> Informacja dla pacjentów .....	46
<b>HU</b> Betegtájékoztató .....	51
<b>CS</b> Informace pro pacienty .....	56
<b>RO</b> Informații pentru pacienți.....	60
<b>TR</b> Hasta bilgisi .....	66
<b>EL</b> Πληροφορίες για τον ασθενή.....	70

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**Date of the last update: 2016-08-22**

- **Please read this document carefully.**
- **Follow the safety instructions and the stated precautions in this accompanying document.**
- **Please contact the responsible prosthetist if assistance is required during the start-up, use or maintenance of the product, or in the event of unexpected operating behaviour or circumstances.**

With their high-power motors and microprocessor control electronics, the System Electric Hands offer versatile functional possibilities and exceptional wearer comfort. The following explanations highlight some of the technical features in order to optimise their use.

Please follow your prosthetist's directions and the safety instructions.

## Safety Instructions

**Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the stated precautions in this accompanying document.**

- Before interrupting or establishing electrical connections (such as when removing the Hand from the prosthesis), the system must be disconnected from the power supply. For this purpose, either remove the battery from the battery receptacle or turn the prosthesis off by pressing the button in the charging receptacle.
- Do not let foreign particles or liquids enter the Electric Hand. The Electric Hand should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- The System Electric Hand and the connected components might malfunction **when near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores)**. Setting the electrode gains as low as possible will minimise this risk. If such malfunctions occur repeatedly, please have the electrode settings checked by your prosthetist.
- The System Electric Hand has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).
- Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures your safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Myo-Service.
- If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WLAN devices) is too small, interference with internal data communications can result in malfunctions of the product. Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.99 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WLAN (routers, access points,...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m
- The product and any damaged components may only be opened and/or repaired by certified Ottobock Myo-Service technicians.
- When operating the product in the immediate vicinity of active, implantable systems (e.g. pacemakers, defibrillators, etc.), ensure that the minimum distances stipulated by the manufacturer of the implant are observed. Interference with active, implantable systems may occur due to the electromagnetic radiation emitted by the product. Furthermore, make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

- **Water and Humidity**

The electrical and mechanical systems of your System Electric Hand are not water-resistant. You must prevent water from entering the Electric Hand. You may wash the Electric Hand, do the dishes, and similar activities only when wearing a prosthetic glove free of tears and punctures. Be careful not to let water run over the top of the glove and enter the hand mechanism. If water enters the inside of the Hand for any reason, please contact your prosthetist immediately to avoid further damage.

**Attention!**

**Driving Vehicles**

An upper extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test center. For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork). It is indispensable that the car can be driven safely with the Hand turned OFF. Driving with the Electric Hand turned ON may present a risk if the Hand inadvertently opens due to unintentional muscle contraction or other causes.

**CAUTION!**

Before a terminal device with a quick-disconnect wrist can be used, the terminal device must be positioned in such a way that a slight turn of it cannot disconnect the terminal device from the prosthesis.

## 1 Instructions for use

### 1.1 Hand switch

Your System Electric Hand has an inconspicuous on/off switch inside the inner hand.

**System Electric Hand, sizes 7¼, 7¾, 8¼**

To turn on the Hand, press the switch on the back of the Hand (Fig. 1).

To turn off the Hand, press the switch in the thumb region (Fig. 2).

**System Electric Hand, size 7** (Fig. 1)

**Transcarpal hand** (Fig. 3)

To turn the Hand on and off, press the button on the back of the Hand. To prevent inadvertent activation or deactivation, the ON/OFF button must be held down for approx. 1 second to turn the prosthesis on or off.

Turn the Hand off whenever you would like to maintain a certain hand position for an extended period of time, such as while carrying a bag, holding silverware, while writing or during longer periods of passive use, e.g. air or train journey, visit to a theater or cinema.

Turning the Hand off prevents unintentional opening of the Hand due to accidental muscle tension or electrical interference. It also extends battery life.

The hand switch can be pressed easily with the other hand, as well as on a table edge or against your leg.

### 1.2 Quick-disconnect wrist

The quick-disconnect wrist allows a quick change from the System Electric Hand to a System Electric Greifer and vice versa.

If your System Electric Hand is supplied with a quick-disconnect wrist (allowing passive pronation and supination) you can rotate it to a better position by turning the Hand. This makes it easier to grasp many objects. You can also easily separate the Hand from the socket (to replace it with a System Electric Greifer) by rotating the Hand in either direction (approx. 360°) until you feel slight resistance (Fig. 4). Once you overcome the resistance, the Hand will disconnect from the wrist.

### 1.3 Locking the quick-disconnect wrist

Insert the Hand's quick-disconnect wrist into the lamination ring and firmly press it in. Next, turn the System Electric Hand slightly to the left or right. The pressure activates the locking mechanism and anchors the Hand to the lamination ring.

**Make sure the Hand or Greifer is properly attached before picking anything up.**

### 1.4 Rechargeable batteries

Only a fully charged Ottobock battery should be used to operate the System Electric Hand.

An intelligent "battery manager" alerts you to low voltage conditions: The Hand will grow slower and the grip force will decrease.

We recommend that you operate the System Electric Hand with the 757B20 or 757B21 EnergyPack or 757B35=\* MyoEnergy Integral, in order to be able to take advantage of the Hand's full capacity over an extended period of time. The Hand can also be operated with the 757B15 X-ChangePack, but this will result in capacity restrictions.

For detailed instructions on care of the battery, please refer to the instructions for use included with the battery. To recharge the batteries, use only the Ottobock Battery Chargers mentioned in the instructions for use.

#### **Attention!**

As soon as the DMC plus or the Digital Twin electronic system recognises a fully charged 757B20/757B21 EnergyPack or 757B35=\* MyoEnergy Integral, it automatically switches to Li-Ion battery technology. Once this has happened, only this technology should be used to operate the system. If a 757B15 X-ChangePack is employed after the control has switched to Li-Ion battery technology, the battery can not be used to its full capacity. The 757T13 Ottobock MyoSelect can be used, however, to reset the System Hand control for use with X-ChangePack or Interchangeable Battery.

For operational safety and reliability, only 757B20, 757B21, 757B15 or 757B35=\* batteries are to be used.

As delivered, the System Electric Hand's default setting is optimised for NiMH Batteries.

## 2 Repair/maintenance

If the System Electric Hand is not in use, it should be left in the **open** position to protect the sensor and mechanism. Instructions for care of the prosthetic glove are enclosed with the glove. For further information, please contact your prosthetist.

## 3 Additional information for patients with a System Electric Hand with Dynamic Mode Control Plus\*

Your prosthetist may choose from two different control modes by inserting or removing a function plug from the System Electric Hand with DMC plus.

### 3.1 DMC Control Mode: function plug attached.

Grip speed and grip force are determined by the strength of the muscle signal. As soon as the strength of the muscle signal changes, grip speed and grip force immediately and proportionally adapt to the new muscle signal. The gripping process is much more natural, and makes it easier to gently grip small or fragile objects.

### 3.2 DMC plus Control Mode: function plug not attached.

Grip speed and grip force are determined by the strength of the muscle signal. As soon as the strength of the muscle signal changes, grip speed and grip force immediately and proportionally adapt to the new muscle signal. After gripping once with maximum force, the EMG signal required to open the Hand will be set to a higher value. This reduces the risk of the Hand opening due to unintentional muscle contraction.

Gripping security for activities such as eating will be increased.

\* *Dynamic Mode Control plus is a proportional control system using microcontroller technology.*

## 4 Disposal

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

## 5 Cleaning and care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Make sure the product is not damaged and that no liquid can penetrate into system components.

Then, dry the prosthesis component with a soft cloth.

## 6 Technical data

Technical data, gripping components				
	System Electric Greifer 8E33=9, 8E34=9	8E38=8*, 8E39=9*, 8E41=8* SensorHand Speed	8E38=9, 8E39=9, 8E41=9 MyoHand Vari-Plus Speed	8E33=7, 8E34=7 System Electric Greifer with Digital Twin control
Static current	1 mA	2 mA	1 mA	1 mA
Operating temperature	0-70 °C	0-70 °C	0-70 °C	0-70 °C
Storage temperature	-20 °C to +60 °C (-4 °F to +140 °F)	-20 °C to +60 °C (-4 °F to +140 °F)	-20 °C to +60 °C (-4 °F to +140 °F)	-20 °C to +60 °C (-4 °F to +140 °F)
Relative humidity	to 80%, non-condensing	to 80%, non-condensing	to 80%, non-condensing	to 80%, non-condensing
Storage and transport temperature in packaging	-20 ° – +60 °C	-20 ° – +60 °C	-20 ° – +60 °C	-20 ° – +60 °C
Storage and transport temperature without packaging	-20 ° – +60 °C	-20 ° – +60 °C	-20 ° – +60 °C	-20 ° – +60 °C
Opening width	95 mm	100 mm	100 mm	95 mm
Proportional speed	8-200 mm/s	15-300 mm/s	15-300 mm/s	180 mm/s
Proportional grip force	0-160 N	0-100 N	0-100 N	0-160 N
Service life	5 years	5 years	5 years	5 years
Technical data, batteries				
	757B20 EnergyPack	757B21 EnergyPack	757B35* MyoEnergy Integral	757B15 X-ChangePack
Power supply	7.2 V	7.2 V	7.4 V	6 V
Service life	2 years	2 years	2 years	2 years

## 7 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise.

The device may only be opened and repaired by authorised Ottobock technicians.

## 8 CE conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

## 9 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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Otto Bock Korea HealthCare Inc.
South East Asia Co., Ltd.
Otto Bock HealthCare GmbH