

50S1 Dyna Ankle

DE Gebrauchsanweisung	4
EN Instructions for use	7
FR Instructions d'utilisation	10
IT Istruzioni per l'uso	13
ES Instrucciones de uso	16
PT Manual de utilização	20
NL Gebruiksaanwijzing	23
SV Bruksanvisning	26
FI Käyttöohje	29
PL Instrukcja użytkowania	32
HU Használati utasítás	35
CS Návod k použití	38
EL Οδηγίες χρήσης	41
RU Руководство по применению	45

	Schuhgröße/Shoe size
	35 – 37 37 – 39 39 – 41 41 – 44

INFORMATION

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INFORMATION

Last update: 2016-10-24

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 50S1 Dyna Ankle orthosis.

2 Intended use

2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

2.2 Indications

- Ruptured lateral ligament
- Acute distortion

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

2.4 Effects

The orthosis sets the foot to slight pronation and dorsal extension. The dynamic restraint limits plantar flexion and supination while pronation and dorsiflexion movements remain unrestricted. The orthosis relaxes the anterior talofibular ligament.

3 Safety

3.1 Explanation of warning symbols

⚠ WARNING Warning regarding possible serious risks of accident or injury.

⚠ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

⚠ WARNING

Operation of motor vehicles

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.

⚠ CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.

NOTICE

Use without suitable footwear

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

4.2 Fitting and Application

- 1) **Optional: thermoforming the orthosis:** Thermoform the orthosis to adjust the fit. In order to do so, remove all hook and loop straps, heat the orthosis to 90°C max. and reshape it.
- 2) **Optional: attach additional pads with hook-and-loop:** In case of pressure sensitivity in the area of the ankle or shin, attach the additional pads with hook and loop straps (see fig. 1 and see fig. 2).
- 3) **Optional: apply anti-slip protection:** Clean the bonding surface in the area of the heel, roughen with fine sandpaper and attach the pimple rubber in the heel area using the supplied double-sided adhesive tape.
- 4) Open all hook and loop straps.
- 5) Set the foot into the orthosis. In doing so, ensure that the heel and forefoot are positioned comfortably.
- 6) Pull the middle hook and loop strap from the inside of the foot across the top of the foot and through the outer loop, and tighten it until the heel fits in the shell. Fasten the hook and loop strap in doing so.
- 7) Put the upper hook and loop strap through the loop and fasten it (see fig. 3).
- 8) Put the lower hook and loop strap through the loop and fasten it.
- 9) Stretch the functional strap and fasten it (see fig. 4).

4.3 Cleaning

- 1) Remove the straps from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH