
## 17LA3N

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	9
<b>FR</b> Instructions d'utilisation .....	16
<b>IT</b> Istruzioni per l'uso .....	23
<b>ES</b> Instrucciones de uso .....	30
<b>PT</b> Manual de utilização .....	38
<b>NL</b> Gebruiksaanwijzing .....	44
<b>SV</b> Bruksanvisning .....	51
<b>DA</b> Brugsanvisning .....	58
<b>NO</b> Bruksanvisning .....	64
<b>FI</b> Käyttöohje .....	71
<b>CS</b> Návod k použití .....	77
<b>SK</b> Návod na používanie .....	84
<b>RU</b> Руководство по применению .....	91
<b>ZH</b> 使用说明书 .....	98

1

## INFORMATION

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### INFORMATION

Last update: 2016-12-14

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide you with important information on the processing of the 17LA3N ankle joints.

## 2 Product description

### 2.1 Available sizes

	Maximum body weight in kg (lbs)	
	Unilateral use (lateral or medial)	Bilateral use (later or medial)
17LA3N=10*	15 (33)	25 (55)
17LA3N=12*	20 (44)	40 (88)
17LA3N=14*	50 (110)	80 (176)
17LA3N=16*	85 (187)	120 (264.5)
17LA3N=20*	110 (242.5)	160 (352.5)

#### INFORMATION

- ▶ When using the product unilaterally, and in case of flexion contracture in the knee or hip  $> 10^\circ$  or distinct torsion or valgus/varus instabilities or increased bodily activity, the next higher size must be used.
- ▶ For tubercle seat or valgus/varus malpositions  $> 10^\circ$  the product must be fitted bilaterally.

### 2.2 Components/design (see Fig. 1)

Item	Component	Reference no. for spare part
1	Set screw	30Y349=12 (size 10, 12) 30Y309=12 (size 14) 30Y309=20 (size 16, 20)
2	Stop bumper	30Y350=12 (size 10, 12) 30Y350=14 (size 14) 30Y350=20 (size 16, 20)
3	Cylinder pin ISO 8734	506A8=4x20 (size 10, 12) 506A8=6x24 (size 14) 506A8=8x26 (size 16, 20)
4	Stop sleeve	30Y346=12 (size 10, 12) 30Y346=14 (size 14) 30Y346=20 (size 16, 20)
5	Screw connection (joint nut, joint screw)	17LA32N=10 (size 10) 17LA32N=12 (size 12) 17LA32N=14 (size 14)

Item	Component	Reference no. for spare part
		17LA32N=16 (size 16) 17LA32N=20 (size 20)
6	Axial washers (set)	17LA33N=10 (size 10) 17LA33N=12 (size 12) 17LA33N=14 (size 14) 17LA33N=16 (size 16) 17LA33N=20 (size 20)
7	Countersunk head screw	501S137=M4x8 (size 10, 12) 501S137=M4x10-2 (size 14) 501S137=M5x10-1 (size 16, 20)
8	Bearing ball	509Y1=5.0 (size 10, 12) 509Y1=7.0 (size 14) 509Y1=9.0 (size 16, 20)
9	Spring guide housing	30Y345=12 (size 10, 12) 30Y345=14 (size 14) 30Y345=20 (size 16, 20)
10	Compression spring set / compression spring	513D87=0.75x2.95x2 (size 10, 12) 513D87=1.1x4.8x2.8 (size 14) 30Y254=16-1 (size 16) 30Y254=20-1 (size 20)

### 3 Intended use

#### 3.1 Indications for use

The system ankle joints are **exclusively** intended to be used for the patient's orthotic fitting of the lower limbs.

In case of unilateral use, the manufacturer recommends that the system ankle joint be processed by means of lamination resin technique or carbon prepreg technology. When used bilaterally, there is no limitation on the processing technology used.

#### 3.2 Indications

Partial or total paralysis of the leg muscles.

Indications must be determined by the physician.

## 4 Safety

### 4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 Safety instructions

 **CAUTION**

**Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).

 **CAUTION**

**Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

 **CAUTION**

**Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

 **CAUTION**

**Excessive strain due to use on more than one patient**

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

#### NOTICE

##### **Thermal overloading of the product**

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at a temperature above 300 °C (570 °F).
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

#### NOTICE

##### **Thermal overloading of product components**

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

#### NOTICE

##### **Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

## 5 Preparation for use

### 5.1 Processing

#### INFORMATION

The instructions for manufacturing an orthosis with unilateral joint systems are included in the 646T5=9.1 Technical Information. The laminate structure described here is a recommendation of Ottobock and can be adjusted to the needs of the patient.

- > **Required materials:** 17LD1N=\* lamination dummy, 17LD1N=\*-P dummy cover, 17LF3N=\* or 17LF31N foot stirrup, 17LS3=\* lamination bar, 636K8 Plastaband, 616G6 Dacron felt, 743Y56=4 alignment adapter
- > The plaster model has been produced in the neutral zero position.
- 1) Mount the lamination bar on the lamination dummy with the screws.

- 2) Insert the foot stirrup into the dummy cover.
- 3) Place the dummy cover and the inserted foot stirrup on the lamination dummy.
- 4) Insert the screw through the dummy cover, the foot stirrup and the lamination dummy and affix it to the alignment adapter.
- 5) Seal the dummy with Plastaband and fill hollows with Dacron felt.
- 6) **NOTICE! Align the marking of the foot stirrup in the neutral zero position with the marking of the counterpart. In doing so, the foot stirrup must be positioned such that the ankle movement in dorsiflexion and plantar flexion permits the full range of motion.**
- 7) Laminate the orthosis.
- 8) Cut free and remove the lamination dummy.

## 5.2 Assembly

### Mounting the foot stirrup

	Torque in Nm (lbf in.)				
	17LA3N=1-0*	17LA3N=1-2*	17LA3N=1-4*	17LA3N=1-6*	17LA3N=2-0*
Joint screw	3 (27)	3 (27)	5 (44.5)	8 (71)	8 (71)

- 1) Insert the foot stirrup with the axial washers into the basic piece. Use one axial washer per side.
- 2) Check whether the foot component can be moved smoothly and without play. Use axial washers of varying thicknesses if necessary.
- 3) Make sure that the axial washers prevent direct contact between the metal basic piece and the foot stirrup.
- 4) Insert the joint screw and the joint nut and secure the screw connection with Loctite® 241 and the corresponding torque.

### Information on axial washers

- 1) Use the axial washers in accordance with the recommendation (see following table).
- 2) **In case of play in the m-l direction or resistance in flexion/extension:** Replace the axial washers with thicker or thinner washers until the orthosis joint moves as desired.

Orthosis joint size	Thickness of axial washers [mm]
17LA3N=10*, 17LA3N=12*	0.70, 0.75, <b>0.80</b> , 0.85, 0.90
17LA3N=14*, 17LA3N=16*, 17LA3N=20*	0.85, 0.90, <b>0.95</b> , 1.0, 1.05, 1.10
<b>Start by inserting 2 axial washers of the thickness indicated in bold.</b>	

## Mounting the lamination bar

	Torque in Nm (lbf in.)				
	17LA3N=1-0*	17LA3N=1-2*	17LA3N=1-4*	17LA3N=1-6*	17LA3N=2-0*
Counter-sunk head screw	3,5 (31)	3,5 (31)	3,5 (31)	7 (62)	7 (62)

- ▶ Fasten the lamination bar to the basic piece using the countersunk head screws and secure the screw connections with Loctite® 241 and the corresponding torque.

### 5.3 Adaptation

The system ankle joint is adapted to the needs of the patient through the individual combination of the cylinder pin and spring with the set screws.

- 1) **Insert the spring:** Insert the ball into the basic piece with the spring guide housing and the spring.
- 2) **Insert the stop:** Insert the stop sleeve into the basic piece with the cylinder pin and the stop bumper. In doing so, position the cylinder pin with the tapered side on the foot stirrup.
- 3) Secure the thread of the set screw with Loctite® 241 and screw the set screw into the basic piece. Observe the angle adjustment of the ankle joint.

## 6 Cleaning

After contact with water containing salt, chlorine or soap, or when they get dirty, the system joints must be cleaned promptly.

- 1) Rinse the system joints with clean fresh water.
- 2) Dry the system joints with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

## 7 Maintenance

### INFORMATION

**The product may be exposed to increased loads by the patient.**

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Only use 633F7 special lubricant.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH