



17CF1

DE	Gebrauchsanweisung	2
EN	Instructions for use	13
FR	Instructions d'utilisation	24
IT	Istruzioni per l'uso	35
ES	Instrucciones de uso	46
PT	Manual de utilização	57
NL	Gebruiksaanwijzing	68
SV	Bruksanvisning	79
DA	Brugsanvisning	90
NO	Bruksanvisning	101
PL	Instrukcja użytkowania	112
HU	Használati utasítás	123
CS	Návod k použití	134
SK	Návod na používanie	145
TR	Kullanma talimatı	156
EL	Οδηγίες χρήσης	167
RU	Руководство по применению	179
JA	取扱説明書	191

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




Explanation of Symbols

NOTICE	Warnings regarding possible technical damage.
INFORMATION	Information regarding operation. Information for service personnel.

INFORMATION

Last update: 2017-02-14

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

INFORMATION

For preparing the plaster cast, please use the enclosed paperboard template.

INFORMATION

Tolerances with respect to spring thickness. Springs of a particular size generally have the same stiffness. Different thicknesses of the spring are due to production tolerances and have no effect on stiffness.

Table of contents

Explanation of Symbols	14
1 Description	15
1.1 Intended purpose.....	15
1.2 Special features of the carbon fiber spring.....	15
1.3 Indications and contraindications.....	15
1.4 Safety instructions	16
1.5 Scope of delivery and accessories.....	16
1.6 Support	17
2 Handling	17
2.1 Fabrication of the plaster model.....	17
2.2 Installation of carbon fiber spring using thermoplastic technique:	17
2.3 Final assembly.....	21
2.4 Trial fitting and delivery of orthosis	22
2.5 Proper functioning and premature wear	22
3 Legal information	22
3.1 Liability	22
3.2 CE conformity.....	22

1 Description

1.1 Intended purpose

The Ottobock Ankle Seven Carbon Fiber Spring facilitates the fabrication of light lower leg orthoses that are subject to high dynamic loads. The special design allows energy to be stored in the carbon fiber spring during heel strike and releases it again at toe-off, so that people who wear orthoses can walk naturally and expend less energy doing so. The orthosis does not have any mechanical orthosis joint and is therefore described as jointless.

1.2 Special features of the carbon fiber spring

The special design of the carbon fiber spring takes into consideration a natural outer foot position of 7° (see figure 1) from the dorsal perspective. This means that once a suitable plaster model has been prepared and constructed, the patient's foot during walking will have something close to a physiological alignment (see sketch 2). The outer position of the carbon fiber spring in the bending area makes it possible to provide support for a more beneficial roll-over.

1.3 Indications and contraindications

The carbon fiber spring is designed exclusively for use with lower extremity orthoses and with the functional defects listed here:

Paralysis or weakness of foot levering and/or foot dropping. Weakness of the knee extensor (muscle tone for the knee extensor >3 according to Janda). For example, spina bifida or other neuromuscular diseases, post-traumatic conditions).

Not illustrated here is a fitting in the case of equinus position with more than a 25 mm reduction, severe spasticity or other causes that seriously compromise the patient's ability to walk.

The right fiber strength of the carbon fiber spring is selected on the basis of body weight and degree of activity by using the classification system below. It is required that the patient can walk. The requirements made of the feather spring's strength are greater for active users of orthoses, i.e. those that run.

INFORMATION

The classification cited applies to AFO orthoses.

The use of carbon springs for KAFO orthoses may require a deviation from the classification.

Choice of carbon springs as per classification for AFOs

	Activity level		
Body weight	Normal activity	High activity	Fiber spring width
up to 100 kg (220 lbs)	17CF1=L/R1	17CF1=L/R1	30 mm
up to 90 kg (198 lbs)	17CF1=L/R2	17CF1=L/R1	30 mm
up to 80 kg (176)	17CF1=L/R3	17CF1=L/R2	30 mm
up to 70 kg (154 lbs)	17CF1=L/R4	17CF1=L/R3	30 mm
up to 60 kg (132 lbs)	17CF1=L/R5	17CF1=L/R4	30 mm
up to 50 kg (110 lbs)	17CF1=L/R6	17CF1=L/R5	25 mm/ 30 mm
up to 40 kg (88 lbs)	17CF1=L/R8	17CF1=L/R7	25 mm
up to 30 kg (66 lbs)	17CF1=L/R9	17CF1=L/R8	25 mm

Choice of carbon springs as per classification for AFOs

	Activity level		
Body weight	Normal activity	High activity	Fiber spring width
up to 20 kg (44 lbs)	17CF1=L/R11	17CF1=L/R10	22 mm
up to 10 kg (22 lbs)	17CF1=L/R12	17CF1=L/R11	22 mm
Activity level	Normal activity	High activity	

1.4 Safety instructions

NOTICE

Damage due to wrong positioning of boreholes. Please make bore holes in the middle of the carbon fiber spring, since otherwise the carbon fiber spring could suffer premature wear and possibly break in the area of the bore holes.

NOTICE

Damage due to wrong tools. Only use sharp tools to avoid ripping out material when drilling. File the bore holes with countersinks, since otherwise the carbon fiber spring could suffer premature wear and possibly break in the area of the bore holes.

NOTICE

Damage due to wrong positioning of boreholes. Please pay attention to the position of the bore holes in the fixed and bending area, since otherwise the carbon fiber spring could suffer premature wear and possibly break (suffer a loss of function).

NOTICE

Damage due to thermal deformation. The carbon spring must not be postformed or reshaped by heat application, because this will result in a destruction of the component.

NOTICE

Damage due to material weakening. The carbon fiber spring may not be ground or cut in the movement area, since otherwise the carbon fiber spring could suffer premature wear and possibly break (suffer a loss of function).

1.5 Scope of delivery and accessories

	for 17CF1 < 60 kg	for 17CF1 > 60 kg
1 Carbon fiber spring		
4 Welding nut	502E3=M4×6	502E3=M5×7.5
4 Rosette washer	507U9=M4	507U9=M5
4 Countersunk head screw	501S86=M4×12	501S86=M5×14

1.6 Support


The Ottobock Team Germany will be pleased to answer your technical questions. Contact addresses and telephone numbers are listed on the last page of this manual.

2 Handling

2.1 Fabrication of the plaster model

The special construction makes it necessary to ensure that the plaster cast is prepared with a slight outer foot position. For better positioning of the carbon fiber spring under the foot, a slight heel elevation of approx. 15-20 mm can be provided in the plaster cast (see sketch 4). The enclosed template may be used for the sagittal positioning of the lower leg (small template for sizes 10-12). The flexion cutout will be added to this afterwards. The flexion cutout runs flat and parallel to the walking direction, while the outer foot position is recognizable (fig. 6).

2.2 Installation of carbon fiber spring using thermoplastic technique:

Pull layer of perlon stockinette over the plaster model and fabricate Pedilin insert (617S3=W5).

INFORMATION

A material thickness of 5 mm can be used for all size models.

Please subdivide the model into two parts to determine the fixed area and the movement area for the carbon fiber spring. Measure the distance from the medial tibial plateau to the floor and mark the middle (see sketch 3).

Trace the carbon fiber spring on the bottom of the foot. The length of the carbon fiber spring extends up to approx. 20 mm in front of the roll-over area when viewed from the calcaneus. Then determine the calf length of the carbon fiber spring. The end of the carbon fiber spring reaches up to 20 mm below the flexion cutout.


The carbon fiber spring may be shortened no further than the marking, since otherwise it is not possible to fix it in the foot part safely. Round off cut edges.

NOTICE

Damage due to wrong tools. Please use the 649Z12 Sanding Brush to smooth the sanded areas, since otherwise the carbon fiber spring could suffer premature wear and possibly break (suffer a loss of function).

To adapt the carbon fiber spring to the lower leg, the insert must be mounted on the appropriate places and sanded. If need be, cavities must be filled in.

The integrated outer position allows the carbon fiber spring to run in the direction of the second metatarsophalangeal joint when it is positioned correctly in the foot part (cf. sketch 1).

Dorsal position of the carbon fiber spring.


Medial view of completely adapted carbon fiber spring. If the contours of the Pedilin construction are not correct, the plastic material could be pulled below the carbon fiber spring during vacuum forming.


Prepare model for vacuum forming. Use the Polyethylene Adhesive Tape (627B4) to attach the carbon fiber spring to the model above the ring/span and below the calf. Pull 1 Perlon stockinette over the model.

INFORMATION

Isolate with 519L5 Silicone Spray. This makes it easier to remove the Perlon stockinette from the cooled polypropylene.

Please use 4 mm for all sizes in order to vacuum form 616T20=2000×4 Polypropylene. The plastic must be shaped on the model under vacuum.

Mark the edges and remove the model. Next separate the calf part from the foot part.

NOTICE

Damage due to wrong cutting. Ensure that when cutting open the model the lining is separated in the middle only.

Mark the bore holes on the calf piece. A ruler can be used to mark the middle.

NOTICE

Damage due to wrong positioning of boreholes. The boreholes in the carbon spring must be positioned centrally, otherwise there may be premature wear and possible breakage at the boreholes.

Mark the bore holes on the foot part. Please do not make proximal bore holes at the highest point of the heel, since this is the flexion area of the carbon fiber spring.

Use a 4 mm drill for M4 to make the hole in the polypropylene. Use a 5 mm drill for M5.

Hold carbon fiber spring below polypropylene and mark the places for the bore holes on it.


Drill a 5 mm hole (for M4) for the 502E3 welding nut. Drill a 6.5 mm hole for M5.

NOTICE

Damage due to wrong positioning of bore-holes. Make sure that no material is ripped out during drilling. Only use sharp tools.

INFORMATION

Always drill one hole after the other and press in the welding nut; screw on the shell component and check the position of the remaining bore holes.

File bore holes by hand or on a machine with countersink (726S9=90x11.5).

NOTICE

Damage due to wrong tools. Do not use a twist drill **for** deflashing. Use a counterbore to deflash boreholes.

Press in welding nuts.

2.3 Final assembly

Screw together carbon fiber spring with the provided 501S86 Niro Countersunk Head Screw and the provided 507U9 Rosette Washers. The rosette washers prevent you from pressing the screws into the polypropylene and ensures that there is a large contact surface.

INFORMATION

Making slotted holes on the lower leg part permits its mounting in a moveable way.

2.4 Trial fitting and delivery of orthosis

Please cut down the subcondylar orthosis edge in accordance with the required level of stability. We recommend shortening the orthosis edge to the height of the medial tibial plateau in cases where the patient has a relatively stable knee joint.

INFORMATION

A condylar support may be indicated in cases of knee joint instability.

The insert does not have to be glued for the trial fitting. To make an assessment of the sagittal fitting, we recommend using the L.A.S.A.R Posture to check the patient with the orthosis in place. The weight bearing line of the laser projected onto the leg should fall 15 mm in front of the compromise pivot point in an ideal scenario, according to Nietert (fig. 20). Appropriate modifications to the shoe can also influence the fit. Check the fit and function of the orthosis. (Fig. 21/22)

Self-adhesive Velcro hook and loop band can be glued to the shell components for the trial fitting. Usually a closure below the tibia head is sufficient. An additional closure can be attached above the back of the foot if necessary. For closing and padding, it is possible to use the 623Z1 Hook and Loop Closure Band, 623P5 Padding Tape, 623P5 Terry Cloth Padding and self-adhesive 616T25 Padding Material. We recommend Ottobock 636W71 CP Quick-Drying Rubber Cement for gluing the insert into the orthotic shell, since any discharged colors must be avoided below the transparent plastic materials. Ottobock 636N9 Quick-Drying Rubber Cement can be used for gluing in lamination shells. We recommend 620P4 Thermoplastic Microcork for building up height and creating a flat surface on the foot part of the orthosis. The sagittal position of the finished orthosis model should be 90° after alignment. Screw connections with Loctite 241 (636K13) secure the orthosis.

Using the orthosis in active life. (Fig. 23)

2.5 Proper functioning and premature wear

We recommend having the orthosis checked for proper functioning and premature wear every six months. This should include an inspection of the carbon fiber spring for delamination or fractures and a check of the screw connections. If need be, the fit should be adapted for growing patients.

3 Legal information


All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.





3.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

3.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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
Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
F +43 1 5267985
service-admin.vienna@ottobock.com · www.ottobock.at

Otto Bock Adria d.o.o. Sarajevo
Ramiza Salčina 85
71000 Sarajevo · Bosnia-Herzegovina
T +387 33 255-405 · F +387 33 255-401
obadria@bih.net.ba · www.ottobockadria.com.ba

Otto Bock Bulgaria Ltd.
41 Tzar Boris IIIrd Blvd. · 1612 Sofia · Bulgaria
T +359 2 80 57 980 · F +359 2 80 57 982
info@ottobock.bg · www.ottobock.bg

Otto Bock Suisse AG
Luzerner Kantonsspital 10 · 6000 Luzern 16 · Suisse
T +41 41 455 61 71 · F +41 41 455 61 70
suisse@ottobock.com · www.ottobock.ch

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Otto Bock South Africa (Pty) Ltd
Building 3 Thornhill Office Park · 94 Bekker Road
Midrand · Johannesburg · South Africa
T +27 11 564 9360
info-southafrica@ottobock.co.za
www.ottobock.co.za

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Otto Bock HealthCare Canada
5470 Harvester Road
Burlington, Ontario, L7L 5N5, Canada
T +1 800 665 3327 · F +1 800 463 3659
CACustomerService@ottobock.com · www.ottobock.ca

Oficina Ottobock Habana
Calle 3ra entre 78 y 80.
Edificio Jerusalem · Oficina 112 · Calle 3ra.
Playa, La Habana. Cuba
T +53 720 430 69 · +53 720 430 81
hector.corcho@ottobock.com.br
www.ottobock.com.br

Otto Bock HealthCare Andina Ltda.
Calle 138 No 53-38 · Bogotá · Colombia
T +57 1 8619988 · F +57 1 8619977
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
Prolongación Calle 18 No. 178-A
Col. San Pedro de los Pinos
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290 · F +52 55 5575 0234
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare LP
11501 Alterra Parkway Suite 600
Austin, TX 78758 · USA
T +1 800 328 4058 · F +1 800 962 2549
USCustomerService@ottobock.com
www.ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Suite 1.01, Century Corporate Centre
62 Norwest Boulevard
Baulkham Hills NSW 2153 · Australia
T +61 2 8818 2800 · F +61 2 8814 4500
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.
B12E, Universal Business Park
10 Juxianqiao Road, Chao Yang District
Beijing, 100015, P.R. China
T +8610 8598 6880 · F +8610 8598 0040
news-service@ottobock.com.cn
www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza
1 Science Museum Road, Tsim Sha Tsui
Kowloon, Hong Kong · China
T +852 2598 9772 · F +852 2598 7886
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India Pvt. Ltd.
20th Floor, Express Towers
Nariman Point, Mumbai 400 021 · India
T +91 22 2274 5500 / 5501 / 5502
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.
Yokogawa Building 8F, 4-4-44 Shibaura
Minato-ku, Tokyo, 108-0023 · Japan
T +81 3 3798-2111 · F +81 3 3798-2112
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.
4F Agaworld Building · 1357-74, Seocho-dong
Seocho-ku, 137-070 Seoul · Korea
T +82 2 577-3831 · F +82 2 577-3828
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.
1741 Phaholyothin Road
Kwaeng Chatuchark · Khet Chatuchark
Bangkok 10900 · Thailand
T +66 2 930 3030 · F +66 2 930 3311
obsea@ottobock.co.th · www.ottobock.co.th

The 17CF1 is covered by the following patents:

China: ZL 200610153119.8

USA: Patent US 7,815,587

European Patent: EP 1795156 in DE, FR, IT, SE, TR, NL,GB.

Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt / Germany
T +49 5527 848-0 · F +49 5527 72330
healthcare@ottobock.de · www.ottobock.com