
# ADI®

A Stealth Product

**Installationshandbuch**

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

# 1.0 2 Point PRO Mounting Hardware

---

## For Aluminum & Carbon Series Back Supports

### INSTALLATION INSTRUCTIONS FOR

**Stealth Products:** BR1103

### Compatible with

**Stealth's:** ADI Aluminum/Carbon Fiber Backs

## 1.1 Customer Satisfaction

Stealth Products strives for 100% customer satisfaction. Your complete satisfaction is important. Please contact us with feedback or suggested changes that will help improve the quality and usability of these products. You may reach us at:

Stealth Products, LLC  
104 John Kelly Drive, Burnet, TX 78611  
Phone: (512) 715-9995      Toll Free: 1(800) 965-9229  
Fax: (512) 715-9954      Toll Free: 1(800) 806-1225  
info@stealthproducts.com      www.stealthproducts.com

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## 1.2 Warning Labels

### 1.2.1 Warning Labels

Warnings are included for the safety of the user, client, operator and property. Please read and understand what the signal words **SAFETY, NOTICE, CAUTION, WARNING** and **DANGER** mean, how they could affect the user, those around the user, and property.

 <b>DANGER</b>	Identifies an <b>imminent</b> situation which (if not avoided) will result in <b>severe injury, death, and property damage.</b>
 <b>WARNING</b>	Identifies a <b>potential</b> situation which (if not avoided) will result in <b>severe injury, death, and property damage.</b>
 <b>CAUTION</b>	Identifies a <b>potential</b> situation which (if not avoided) will result in <b>minor to moderate injury, and property damage.</b>
<b>NOTICE</b>	Identifies important information not related to injury, but possible <b>property damage.</b>
<b>SAFETY</b>	Indicates steps or instructions for safe practices, reminders of safe procedures, or important safety equipment that may be necessary.

### 1.2.2 Limited Liability

Stealth Products, LLC accepts no liability for personal injury or damage to property that may arise from the failure of the user or other persons to follow the recommendations, warnings, and instructions in this manual.


### 1.2.3 Testing

Initial setup and driving should be done in an open area free of obstacles until the user is fully capable of driving safely.

The product should always be tested without any person sitting in the wheelchair until every alteration of the physical installation or adjustment is complete.

<b>1.0 2 Point PRO Mounting Hardware</b>	
<b>1.1 Customer Satisfaction</b>	<b>22</b>
<b>1.2 Warning Labels</b>	<b>22</b>
1.2.1 Warning Labels	22
1.2.2 Limited Liability	22
1.2.3 Testing	22
<b>1.3 Table of Contents</b>	<b>23</b>
<b>1.4 Parts and Accessories</b>	<b>24</b>
1.4.1 2 Point PRO Mounting Hardware	24
1.4.2 Tools Required for Installation & Adjustment	25
1.4.3 Torque Specification	25
<b>1.5 Installation Instructions</b>	<b>26</b>
1.5.1 Installing Stud Plate and U-Block Assembly	26
1.5.2 Fastening 2 Point Assembly to U-Block	26
1.5.3 2 Point PRO Clamp Assembly to Back Post	27
1.5.4 Height Adjustment AL/AFT Series	28
1.5.5 Width Adjustment	28
1.5.6 Angle Adjustment	30
1.5.7 Depth Adjustment	30
1.5.8 Back Installation Onto Clamps	31
<b>2.0 WC20 Approved 2 Point PRO Mounting Hardware</b>	
<b>2.1 WC20 Approved</b>	<b>33</b>
2.1.1 WC20 2 Point PRO Mounting Hardware	33
<b>2.2 WC20 Transit Information</b>	<b>33</b>
2.2.1 Before Transit	33
<b>2.3 WC20 Safety</b>	<b>36</b>
2.3.1 User Safety	36
<b>2.4 WC20 Hardware Installation</b>	<b>37</b>
2.4.1 WC20 2 Point PRO Hardware Installation	37
2.4.2 Installing Label onto Cover	37
<b>2.5 Maintenance</b>	<b>38</b>
2.5.1 Cleaning	38
2.5.2 Disinfection	38
<b>2.6 First Time Use</b>	<b>39</b>
2.6.1 Dealer Assistance	39
2.6.2 User Testing	39
2.6.3 Conditions of Use	39
<b>2.7 Warranty</b>	<b>40</b>
2.7.1 In Case of Product Failure	40

### 1.4.1 2 Point PRO Mounting Hardware


		Description
<b>A</b>	2 Point PRO U-Block w/Stud Plate (2)	<b>A1</b> M6 x 1 Nylon Lock Nut
		<b>A2</b> M6 x 18mm Fender Washer
		<b>A3</b> Stud Plate (2)
<b>B</b>	Left/Right Side Clamp Assembly	<b>B1</b> M6 x 18mm Fender Washer (4)
		<b>B2</b> M6 Split Washer (4)
		<b>B3</b> M6 x 1 x 20mm BHS (4)
<b>C</b>	Back Post Clamp (2)	<b>C1</b> M5 x .8 x 25mm SHS <sup>1</sup> (4)
		<b>C2</b> M5 x .8 x 25mm SHS <sup>2</sup> (4)

1. Used for mounting onto 3/4" and 7/8" canes.
2. Used for mounting onto a 1" cane.

### 1.4.2 Tools Required for Installation & Adjustment

Tools*		Description
1	Wrench	10mm Wrench
2	Allen Wrench/T-Handle	4mm Allen Wrench
3	Torque Wrench	Lb-in

### 1.4.3 Torque Specification

Part Number	Torque Specs
BR5001 (A2)	117in-lbs/13.2N-m
BR1100 (C2)	45.1in-lbs./5.1N-m
BR1101 (B3)	77in-lbs./8.7N-m


\*Tools not included in hardware package.

### 1.5.1 Installing Stud Plate and U-Block Assembly

**Step 1:** Determine the desired location for the stud plate to be placed along the pre-cut slots in the back shell.

**Step 2:** Insert the stud plate onto the back. Align the stud plate to the u-block.

**Step 3:** Tighten the lock nuts (A1, A2) with the torque wrench to firmly secure. See *Torque Specs* in *Section 4.3* for the correct torque specifications for tightening the hardware.

#### NOTICE

Installing the U-Block onto the back will be required if the hardware is purchased separately from the back. Stealth will install hardware if simultaneously ordered with the back.


**CAUTION** Install U-Block before placing the cover over the back.

### 1.5.2 Fastening 2 Point Assembly to U-Block

**Step 1:** Align 2 Point mounting bracket and cane clamp with the assembled and mounted U-Block.

**Step 2:** Insert M6 screw (B3) (x2) into washers (B1, B2) (x2) and attach assembly to left and right side U-Block, using 4mm Allen Wrench.

**Step 3:** After final adjustments are complete, tighten all M6 nuts with torque wrench. See *Section 4.3* for torque specifications.


### 1.5.3 2 Point PRO Clamp Assembly to Back Post

**Step 1:** Visually determine the desired height of your ADI back support. Mark on each back post the needed location to attach clamp assemblies.

**Step 2:** Loosen the M5 screws (C1), and with two fingers, pinch and hold the back of the post clamp open in order to slide clamp assembly over top of back post (receiving post facing inward). Screws do not have to be completely removed from hardware in order to slide clamp over post.

**Step 3:** Slightly tighten each M5 screw with a 4mm Allen Wrench. Align clamp assembly parallel to the wheelchair seat base.

**CAUTION**

For 1" cane, use the provided 25mm socket head screws (C2).

**NOTICE**


The mounting hardware should be installed at equal heights on the left and right sides, and parallel to the seat frame.

**NOTICE**

Tighten clamp screws, but allow for movement so you can adjust the receiving posts.

**Step 4:** To achieve proper alignment, slightly move the mounting hardware in or out and/or up or down. Make adjustments until the backrest smoothly attaches and removes from the top and bottom receiving posts.

**Step 5:** After final adjustment is complete, tighten all screws (C1) with torque wrench to the proper specifications. This will ensure secure placement and safe operation.


### 1.5.4 Height Adjustment AL/AFT Series

**Step 1:** With a 10mm wrench, loosen M6 nuts (A2) on the stud plate.

**Step 2:** Slide stud plate up or down along the pre-cut slot in the shell in order to achieve desired height.

**Step 3:** Tighten and secure M6 nuts with torque wrench. See *Section 4.3* for torque specifications.

#### ⚠ CAUTION


Stealth does not recommend drilling holes into the CF backs to reposition the hardware to a different height.

### 1.5.5 Width Adjustment

**Step 1:** On each side of the hardware, loosen the top and bottom screws (A2) using a 4mm Allen Wrench and a 10mm wrench.

**Step 2:** Adjust hardware inward or outward to achieve required mounting hardware width.

**Step 3:** After final adjustments are complete, tighten all screws to the proper torque specification.


#### NOTICE

Minimum and maximum growth for the width of the hardware is 1".

#### NOTICE

AFT Series backs offer an extra inch of width adjustment (2 inches total) on the shell due to a slot pre-cut in the shell.

Minimum and maximum measurements\* for width adjustment:


\*Measurements shown are for an 16" AL series back.

### 1.5.6 Angle Adjustment

**Step 1:** On each side of the hardware, loosen the top and bottom screws (B3) with a 4mm Allen Wrench.

**Step 2:** Adjust back support to achieve desired user angle.

**Step 3:** After final adjustments are complete, tighten all screws to the proper torque specification.


### 1.5.7 Depth Adjustment

**Step 1:** On each side of the hardware, loosen the top and bottom screws (B3) with a 4mm Allen Wrench.

**Step 2:** Slide back support forward or aft to achieve desired users depth.

**Step 3:** After final adjustments are complete, tighten all screws to the proper torque specification.



### 1.5.8 Back Installation Onto Clamps

**Step 1:** Align mounting block holes with the slots on the clamp bracket. Insert washers (B1, B2) and M6 BHS (B3).

**Step 2:** Tighten screws until the back is secure and not loose. Make any angle adjustments if needed.


**Step 3:** Tighten screws to proper torque specification.

---

## 2.0 WC20 Approved 2 Point PRO Mounting Hardware For All Aluminum Back Supports

### **INSTALLATION INSTRUCTIONS FOR**

**Stealth Products:** BR1103 with BR6006

### **Compatible with**

**Stealth's:** ADI Aluminum/AFT Series Backs

### 2.1.1 WC20 2 Point PRO Mounting Hardware

The 2 Point PRO hardware can be ordered as a WC20 package. This includes the hardware and transit pins.

The 2 Point PRO mounting hardware has met WC20 and ISO16840-4 requirements and provisions.

WC20 establishes design and performance requirements and test methods for complete seating systems that consist of a seat, back support, and attachment hardware.

#### CAUTION

**The WC20 approved 2 Point PRO hardware and AL/AFT backs should only be used as described in the manufacturer's instructions.**

#### NOTICE

If the BR6006 order is placed after a shell has been ordered, the corresponding label **MUST** be permanently affixed to the cover in order for the hardware to be WC20 compliant.

#### SAFETY

ADI's CF series backs are NOT WC20 approved.

### 2.2.1 Before Transit

- Clients' who use wheelchairs should transfer to the vehicle seat and use the vehicle-installed occupant-restraint systems.
- Prior to transporting, always double check the tie-downs and hardware.
- On tie-downs and restraints, look for tears or worn webbing straps. Replace as necessary and never use torn straps.
- Look for damaged hooks or loose tie-down hardware.
- On the seating system, check for loose screws attaching the hardware to the seat. These should be snug and not over tightened.
- The 2 Point PRO hardware should be securely attached to the back, and the transit pins should be properly inserted into the hardware.

#### DANGER

**The WC20 approved 2 Point PRO hardware and the AL/AFT backs are designed specifically for occupancy when used on wheelchairs that are facing forward in motor vehicles.**

#### DANGER

The WC20 approved 2 Point PRO hardware and AL/AFT backs have been designed to accommodate the proper use of vehicle-anchored belt restraints and the ratings of the seating systems with regard to ease of properly positioning vehicle-anchored belt restraints on the wheelchair passenger, and the degree to which proper positioning of belt restraints can be achieved.


#### DANGER

Alterations or substitutions should not be made to the seating system or to its parts and components without consulting the manufacturer.

## Belt Placement

- The pelvic belt restraint should be worn low, across the front of the pelvis near the junctions of the pelvis and the thighs.
- The belt restraint buckle of three-point belt restraints should be placed in contact with the occupant's body and away from wheelchair components.
- The junction of the shoulder belt and pelvic belt of three-point-belt restraints should be located near the hip opposite the shoulder over which the diagonal belt crosses and not near the midline of the occupant.

### Proper Positioning of Belts on Occupant

- Upper torso belt restraints should fit directly over, and in contact with, the middle of the shoulder.

### ⚠ CAUTION

Belt restraints should be adjusted to fit as snugly as possible, consistent with user comfort.

### SAFETY

Belt restraints should be positioned on the wheelchair occupant in accordance with WTORS manufacturer's instructions.

### Improper Positioning of Belts on Occupant


- Belt restraints should not be routed outside of the wheelchair wheels or over the wheelchair arm supports and should not be held away from the body by wheelchair components or parts.

### SAFETY

**Belt restraints should not be worn twisted in a manner that reduces the area of contact of the belt webbing with the occupant.**

Both pelvic-belt and shoulder-belt restraints that comply with RESNA WC-4:2012, Section 18 and/or 4.6,5.2,5.3 or RESNA WC-4:2012, Section 19 should be used to limit occupant movement in a crash and reduce the likelihood of occupant ejection from the vehicle and/or injurious contact of the occupant with the vehicle interior.

The pelvic belt should be angled between 30° to 75° horizontally when viewed from the side, and ideally between 45° to 75° horizontally, as shown in the figure below.

#### **⚠ WARNING**

**Back Supports with adjustable recline angles should not be tilted back more than 30° to the vertical during travel in vehicles unless necessary for the medical and postural needs of the occupant.**

#### **⚠ WARNING**

**If the back support must be adjusted to an angle that is greater than 30° to the vertical during travel, the upper shoulder-belt anchor point should be adjusted rearward to maintain belt contact with the wheelchair occupant's shoulder and chest.**

### **Trays**


Wheelchair-mounted rigid trays that are not specifically designed for use during travel in motor vehicles should be:

- Removed and secured separately in the vehicle.
- Secured to the wheelchair so they will not break free in a crash.
- Be positioned with a gap of at least 75mm (3in.) between the back edge of the tray and the wheelchair occupant's abdomen and/or chest so as not to interfere with proper belt-restraint use.
- Have energy absorbing padding placed between the back edge of the tray and the wheelchair occupant.

### 2.3.1 User Safety

The use of a postural pelvic belt attached to the wheelchair base or seat frame is encouraged during travel, but these belts should be positioned so that they do not interfere with the proper positioning of crashworthy belt restraints, and they should not be relied on for occupant protection in crash situations unless the postural belt has been designed to comply with, and perform to, the requirements 4.6,5.2, and 5.3 of RESNA WC-4:2012, Section 19.

Sufficient forward and rearward clear space should be provided around the wheelchair occupant.

Torque specifications for the hardware are listed in the *Parts and Accessories* section on page 2.

The forward clear-space zone (FXZ) needs to be larger when a shoulder-belt restraint is not used.

When mounting the WC20 2 Point PRO hardware to your AL or AFT back, refer to the installation instructions starting on page 3.

#### **⚠ DANGER**

The WC20 approved 2 Point PRO hardware and AL/AFT backs are designed **ONLY** for use with wheelchair frames that have been successfully crash tested to the frontal-impact requirements of RESNA WC-4:2012, Section 19.

#### **⚠ WARNING**

The WC20 approved 2 Point PRO hardware and AL/AFT backs **MUST** be used on a wheelchair frame that provides four securement points and pelvic-belt anchor points that conform to the requirements of RESNA WC-4:2012, Section 19.


### 2.4.1 WC20 2 Point PRO Hardware Installation

Follow the installation instructions starting on page 3 to install the WC20 approved hardware.

To install the WC20 safety tag onto the hardware, follow the instructions below:


**Step 1:** Attach the WC20 safety tag onto the tab.

**Step 2:** Align the tab, washers and screw with the hardware and secure.

### 2.4.2 Installing Label onto Cover

The WC20 Label will need to be secured to the cover of the ADI back if the WC20 pins are ordered separately from the back and hardware. The label must be sewn to cover in order to be WC20 compliant.

**Warning!**

You must comply with the following for proper transit use!  
 Failure to comply may result in serious injury or death to the user!

Minimum: \_\_\_\_\_

Maximum: \_\_\_\_\_

Stealth Products, LLC

Date: \_\_\_\_\_

Model: \_\_\_\_\_

Serial: \_\_\_\_\_

This seating system has passed crash testing with a midsize male ATD in accordance to **RESNA WC-4:2012/ISO16840** using methods described in Annex A of WC-4:2012, Sec.20, and has met requirements when attached to a wheelchair equipped with tie downs specified in WC-4:2012, Sec. 20.

**DO NOT REMOVE TAG!**

### 2.5.1 Cleaning

- Use a soft cloth to clean the back shell.
- Hand wash the foam pad cover with mild, soapy water or detergent and let dry.
- Foam pad should NOT get wet.

**⚠ WARNING** Do NOT hand wash in water temperatures above 160°F (70°C).

**⚠ WARNING** Do NOT wash or machine dry.

### 2.5.2 Disinfection

- Wipe gently with a cloth dampened with a household disinfectant.
- Allow foam to air dry thoroughly.

**⚠ WARNING** Do NOT immerse the foam in water.

#### SAFETY

Always check all mounting hardware, making sure each fastener is properly tightened before using the back support.

**⚠ WARNING** 250lbs. (113.40 kg) user weight limit.

#### NOTICE

Altering an ADI back in any way and/or improper installation will void the warranty.

## 2.6.1 Dealer Assistance

During first time use by the client, it is advised that the dealer or service technician assists and explains the set-up to the customer (the user and/or the attendant). If needed, the dealer can make final adjustments.

**⚠ CAUTION** Do not attach the back with the user in the chair.

### SAFETY

Do not attach the ADI backs with the user in the chair or use them during transportation unless you are using the 2 Point PRO hardware with the WC20 approved pins.

## 2.6.2 User Testing

It is important that the customer is fully aware of the installation, how to use it, and what can be adjusted in order to gain as much mobility as possible. As a dealer, proceed as follows:

- Explain and show the customer how you have executed the installation, and explain the function of the hardware.
- If needed, adjust the hardware to the proper position.
- Explain to the customer possible problems and how to address them.

## 2.6.3 Conditions of Use

The ADI hardware is intended for use as installed by the dealer, in accordance to the installation instructions in this manual.

- The foreseen conditions of use are communicated by the dealer or service technician to the user and/or attendant during the first time use.
- If the usage conditions change significantly, please contact your dealer or a qualified service technician to avoid excessive wear and tear or unintended damage.

Our products are designed, manufactured, and produced to the highest of standards. If any defect in material or workmanship is found, Stealth Products will repair or replace the product at our discretion. Any implied warranty, including the implied warranties of merchantability and fitness for a particular purpose, shall not extend beyond the duration of this warranty. Stealth Products, LLC does not warrant damage due to, but not limited to:

- Misuse, abuse, or misapplication of products.
- Modification of product without written approval from Stealth Products, LLC.
- Any alteration or lack of serial number, where applicable, will automatically void this warranty.
- Stealth Products, LLC is liable for replacement parts only.
- Stealth Products, LLC is not liable for any incurred labor costs.

No person is authorized to alter, extend, or waive the warranties of Stealth Products, LLC.

Stealth Products warrants against failure due to defective materials or workmanship:

Covers: 180 days

Hardware: 5 years

Electronics: 3 years

### 2.7.1 In Case of Product Failure

In the even of product failure covered by our warranty, please follow the procedures outlined below:


1. Call Stealth at +1 (512) 715-9995 or toll free +1-800-965-9229.
2. Request the Returns Department or obtain an RA from the Returns Department and follow department or documentation instructions.

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- Modificación sin permiso escrito de Stealth Products, LLC.
- Stealth Products, LLC es responsable solo por la sustitución de las piezas.
- Stealth Products, LLC no es responsable de ningún coste de mano de obra en que se incurra.
Ninguna persona está autorizada a alterar, extender o renunciar a las garantías de Stealth Products, LLC.