## 50S3 Malleo Sprint

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	8
<b>FR</b> Instructions d'utilisation .....	12
<b>IT</b> Istruzioni per l'uso .....	16
<b>ES</b> Instrucciones de uso .....	20
<b>PT</b> Manual de utilização .....	25
<b>NL</b> Gebruiksaanwijzing .....	29
<b>SV</b> Bruksanvisning .....	33
<b>DA</b> Brugsanvisning .....	37
<b>NO</b> Bruksanvisning .....	41
<b>PL</b> Instrukcja użytkowania .....	45
<b>HU</b> Használati utasítás .....	49
<b>CS</b> Návod k použití .....	54
<b>HR</b> Upute za uporabu .....	58
<b>SK</b> Návod na používanie .....	62
<b>BG</b> Инструкция за употреба .....	66
<b>TR</b> Kullanma talimatı .....	70
<b>EL</b> Οδηγίες χρήσης .....	74
<b>RU</b> Руководство по применению .....	78
<b>ZH</b> 使用说明书 .....	83

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### INFORMATION

Last update: 2017-11-22

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 50S3 Malleo Sensa ankle orthosis.

## 2 Intended use

### 2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

### 2.2 Indications

- Functional treatment of ligament instability
- Post-operative treatment of ligament ruptures
- Ankle sprains
- Functional, conservative therapy of lateral ligament ruptures

Indications must be determined by the physician.

## 2.3 Contraindications

### 2.3.1 Absolute Contraindications

Not known.

### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the foot area, e.g. in case of diabetic neuropathy.

## 2.4 Effects

This orthosis stabilises the ankle with lateral elements and straps. Eversion and inversion of the ankle joint are restricted without negatively affecting plantar and dorsal flexion.

## 3 Safety

### 3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

 **CAUTION**

**Using the product without following the instructions for use**

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

 **CAUTION**

**Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

**⚠ CAUTION****Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

**NOTICE****Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 4 Handling

**INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Size Selection

- ▶ Select the orthosis size based on the circumference of the ankle (see size chart).

### 4.2 Fitting and Application

**⚠ CAUTION****Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

**NOTICE****Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

#### NOTICE

#### **Improper use and changes**

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

The orthosis is delivered ready for trial fitting. Adjustment is not required.

> The strap on the orthosis is open and the fasteners are loosened.

1) Pull the orthosis over your foot (see fig. 1) and tie it like a normal shoe (see fig. 2).

**Optional:** Use the quick fasteners (see fig. 3).

2) Pass the outer end of the strap over the instep and pull it under the foot from the inner (medial) side.

3) Fasten the outer end of the strap over the ankle on the outer (lateral) side (see fig. 4).

4) Pass the inner end of the strap over the instep and pull it under the foot from the outer (lateral) side.

5) Fasten the inner end of the strap over the ankle on the inner (medial) side (see fig. 5).

6) If the physician has allowed weight-bearing on the foot:  
Loosen the ends of the strap and re-tighten (see fig. 6).

#### **4.3 Cleaning**

- 1) Close all hook-and-loop closures.
- 2) Hand wash the orthosis in warm water at 40°C with mild laundry detergent.
- 3) Rinse well.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

#### **5 Disposal**

Dispose of the product in accordance with national regulations.

## 6 Legal information

### 6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

### 6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH