


17LK3

DE	Gebrauchsanweisung	5
EN	Instructions for use	12
FR	Instructions d'utilisation	20
IT	Istruzioni per l'uso	28
ES	Instrucciones de uso	36
PT	Manual de utilização	44
NL	Gebruiksaanwijzing	52
SV	Bruksanvisning	60
DA	Brugsanvisning	67
NO	Bruksanvisning	75
FI	Käyttöohje	82
CS	Návod k použití	89
SK	Návod na používanie	97
RU	Руководство по применению	105
ZH	使用说明书	114


1


22

2


3


17Y162=*


4


5

6

7

INFORMATION

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INFORMATION

Last update: 2018-08-31

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

These instructions for use provide important information on the processing of the 17LK3 unilateral knee joint.

2 Product description

2.1 Available sizes

The orthosis joint is available in five sizes and can be used unilaterally or bilaterally. The choice of size is based on the patient's body weight and the conditions of use.

	Maximum patient weight [kg]				
	17LK3=10	17LK3=12	17LK3=14	17LK3=16	17LK3=20
Unilateral use (lateral or medial)	15	20	50	85	110
Bilateral use (lateral and medial)	25	40	80	120	160

INFORMATION

- When using the product unilaterally, and in case of flexion contracture in the knee or hip > 10° or distinct torsion or valgus/varus instabilities or valgus/varus malpositions > 10° or increased bodily activity, the next higher size must be used.
- For tubercle seat the product must be fitted bilaterally.

2.2 Components

Item	Component (see Fig. 1)	Joint size 17LK3=L*, 17LK3=R* *=10, 12, 14, 16, 20	Reference number
1	Knee lever		30Y265=*-A
2	Upper joint section		30X94=*
3	Lock washer		507S96=2.3
4	Washer		30Y267=3.3
5	Washer		30Y267=4.2

Item	Component (see Fig. 1)	Joint size 17LK3=L*, 17LK3=R* *=10, 12, 14, 16, 20	Reference number
6	Bearing nut	17LK3=10/=12/=14 17LK3=16 17LK3=20	17Y93=9x9.2xM6 17Y93=9x10.75xM6 17Y93=9x11.75xM6
7	Lock washer	17LK3=10/=12/=14 17LK3=16/=20	507S96=3.2 507S96=4
8	Compression spring	17LK3=10/=12 17LK3=14/=16/=20	513D83=0.75x4.3x19 513D83=0.75x4.6x20
9	Lock wedge	17LK3=10/=12/=14 17LK3=16 17LK3=20	17Y156=* 17Y157=* 17Y158=*
10	Bushing	17LK3=12/=14	30Y87=3
11	Bearing washer	17LK3=16/=20	170Z93=1
12	Bushing	17LK3=16/=20	30Y87=1
13	Lower joint section		30U132=*
14	Axle	17LK3=10/=12/=14 17LK3=16 17LK3=20	4A101=5x10.6 4A101=6x12.7 4A101=6x13.7
15	Joint screw		501F9=M6x8
16	Axle	17LK3=10/=12/=14 17LK3=16 17LK3=20	4A101=3.2x15.1 4A101=3.2x17.1 4A101=3.2x18.1
17	Axle	17LK3=10/=12/=14 17LK3=16 17LK3=20	4A101=4x15.1 4A101=4x17.1 4A101=4x18.1
18	Protective plug	17LK3=10/=12/=14 17LK3=16/=20	30Y91=3 30Y91=1
19	Temporary switch		170Y162=*
20	Lamination dummy		30Y268=*
22	Screw	17LK3=10/=12 17LK3=14 17LK3=16/=20	501S137=M4x8 501S137=M4x10-2 501S137=M5x10-1

The scope of delivery does NOT include:

Item	Component (see Fig. 1)	Reference number
21	Shoulder screw (part of 743R6 adjustment set)	30Y89
Item	Component (see Fig. 2)	Reference number
23	Lamination bar	17LS3=*

3 Intended use

3.1 Indications for use

The manufacturer recommends that the orthosis joint be processed using lamination technology or carbon prepreg technology if used on one side of the body.

The orthosis joint is intended **exclusively** for orthotic fittings of the lower limbs on one patient.

3.2 Indications

Partial or total paralysis of the leg muscles.

Indications must be determined by the physician.

4 Safety

4.1 Explanation of warning symbols

WARNING	Warning regarding possible serious risks of accident or injury.
CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

4.2 General safety instructions

WARNING

Operation of motor vehicles

Risk of accidents due to restricted body function

- Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- Only use the product for the defined area of application.
- If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

⚠ CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

⚠ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Thermal overloading of product components

Loss of function due to improper thermal treatment

- Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- Avoid storage in condensing ambient humidity.
- Avoid contact with abrasive substances (e.g. sand, dust).
- Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Processing

5.1 Laminating the orthosis

INFORMATION

The instructions for fabricating an orthosis with unilateral joints are included in the 646T5=5441 technical information. The laminate structure described here is a recommendation from Ottobock and must be adjusted to the needs of the patient.

- > **Required tools and materials:** 616G15 woven carbon stockinette, 17Y106 PVC profile material, 616G12 carbon fibre cloth, 99B81=* PVA bag, 623T3=* Perlon stockinette, 636K8 Plastaband, 633W8 wax, 17LS3=* lamination bars, 617H55 C-Orthocryl lamination resin, 743R6 adjustment set for orthosis joints.
- > **Optional:** 30Y268=* lamination dummy
- 1) Place cut pads or separating layers on plaster model.
- 2) Pull a layer of Perlon stockinette over the plaster model.
- 3) Align the bars on the plaster model.
- 4) **If using a lamination dummy:** Fix the lamination bars to the 743Y56=4 alignment adapter using the lamination dummy and the shoulder screw.
- 5) **If not using a lamination dummy:** Fix the lamination bars to the 743Y56=4 alignment adapter using the base body (no illustration).
- 6) Pull a PVA bag over the plaster model to isolate it.
- 7) Pull a layer of Perlon stockinette over the plaster model.

- 8) Secure the reinforced bars on the plaster model. Carry out the reinforcement in accordance with the patient's needs.
- 9) Seal the lamination dummy or base body of the joint with Plastaband. Fill in hollow areas.
- 10) Pull another layer of Perlon stockinette over the plaster model.
- 11) Affix further reinforcement layers for the support surfaces as described in Section 8.
- 12) Complete reinforcement with a layer of Perlon stockinette and designer fabric.
- 13) Pull on a PVA bag and apply vacuum.
- 14) Then laminate with C-Orthocryl lamination resin.

6 Assembly

	Torque [Nm]				
	17LK3-10	17LK3-12	17LK3-14	17LK3-16	17LK3-20
Joint screw (item 15)	8	8	8	8	8
Screw (item 22)	3.5	3.5	3.5	7	7

- 1) Demould the model and cut free the lamination dummy (20) or the orthosis joint.
- 2) Insert the lower joint section (13) with the bushings (10) (**17LK3-10=/12=/14**) or with the bushings (12) and washers (11) (**17LK3-16=/20**) into the upper joint section (2).
- 3) Fasten the joint screw (15) and bearing nut (6); secure the screw connection with Loctite® 241 and tighten to the specified tightening torque (see table).
- 4) Fasten the joint on the lamination bars (23) using the screws (22); secure the screw connections with Loctite® 241 and tighten to the specified tightening torque (see table).

6.1 Mounting the accessories

CAUTION

Improper handling

Risk of falling due to accidentally opening or engaging the lock

- Before placing weight on the orthosis, verify that the lock is in the desired position.
- Avoid accidentally opening or engaging the lock.
- Instruct the patient in the correct handling of the system joints.

17Y162=* temporary switch

A temporary switch that can be subsequently installed in the orthosis joint is included with the orthosis joint. The temporary switch prevents the lock from engaging, e.g. while riding a therapy bike (see fig. 7).

- 1) Join the button to the tilt lever and insert it into the orthosis joint (see fig. 3, see fig. 4).
- 2) Insert the dowel pin into the orthosis joint and slide it through the hole in the tilt lever in order to secure the temporary switch (see fig. 5, see fig. 6).

7 Disassembling the lock

- > Disassemble the lock to use the product as a free-moving knee joint.
- 1) Unload the orthosis joint using the knee lever (1).
 - 2) Remove the lock rings (3) and lock washers (4, 5).
 - 3) Remove the axles (16, 17).
 - 4) Pull off the knee lever.
 - 5) Push the lock wedge (9) with the compression spring (8) sideways out of the upper joint section (2).

8 Cleaning

Clean the orthosis joint promptly after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the orthosis joint with clean, fresh water.
- 2) Dry the orthosis joint with a cloth or allow it to air-dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

9 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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