

## 170K1=80 Free Walk, 170K1=120 Free Walk









<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	4
<b>EN</b> Instructions for use (qualified personnel) .....	12
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	20
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	29
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	38
<b>PT</b> Manual de utilização (Pessoal técnico) .....	47
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	55
<b>SV</b> Bruksanvisning (Fackpersonal) .....	64
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	72
<b>NO</b> Bruksanvisning (Fagpersonell) .....	80
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	88
<b>TR</b> Kullanma talimatı (Uzman personel) .....	96

1

2

## INFORMATION

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### INFORMATION

Last update: 2018-10-17

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These Instructions for Use provide you with important information on the processing of the 170K1=80 and 170K1=120 Free Walk orthosis system.

## 2 Product description

### 2.1 Available sizes

The orthosis is available in two sizes. The choice of size is based on the patient's body weight.

Article number	Max. body weight [kg]
170K1=80	80
170K1=120	120

## 2.2 Components

Item	Component (see fig. 1)	Article number
1	Y-hook-and-loop	170D29=*
2	Hook and loop closure	170Z4=600
3	Hook and loop closure	170Z4=400
4	Strap guide loop with tube connection	170D18=80
5	Base tube	170B2=80 170B2=120
6	3-phase switch	170D71=*
7	Base tube	170B2=80 170B2=120
8	Ankle joint, upper part	170C3=80 170C3=120
9	Ankle joint, lower part	170X17
10	Foot stirrup	
11	Pad button	170D69
12	Strap guide loop	170D20
13	Foam pads	170D100=1
14	Foam pads	170D100=2

## 2.3 Design

The Ottobock Free Walk orthoses include an orthosis knee joint with lock. 3 modes can be selected on the orthosis knee joint:

- **Free Walk mode:** The knee joint locks in the stance phase and unlocks in the swing phase.
- **Permanent lock:** The knee joint is permanently locked.
- **Free motion mode:** The knee joint is permanently unlocked.

## 3 Intended use

### 3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

### 3.2 Indications

- Flaccid paresis or paralysis of the lower limb
- Adequate muscle strength of the hip or knee extensors (scale according to Janda: 3 to 5) or hyperextension in the knee joint

Indications must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

- Ankle movement < 10°
- Knee flexion contractures > 10°
- Unstable varus position of the knee when fully extended
- Valgus deformity > 15°

#### 3.3.2 Relative Contraindications

- Valgus deformity 0° to 15° (fitting only with 170D30 medial support)

- Spasticity resulting in involuntary movements in the ankle joint

## 4 Safety

### 4.1 Explanation of warning symbols

**WARNING** Warning regarding possible serious risks of accident or injury.

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 4.2 General safety instructions

#### **CAUTION**

##### **Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

#### **CAUTION**

##### **Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

#### **CAUTION**

##### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

#### **CAUTION**

##### **Improper moulding or application**

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

#### **NOTICE**

##### **Thermal overloading of the product**

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at a temperature above 300 °C (570 °F).
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

#### NOTICE

##### **Thermal overloading of product components**

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

#### NOTICE

##### **Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)**

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

#### NOTICE

##### **Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

#### NOTICE

##### **Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 5 Preparation for use

### 5.1 Assembly

#### **Assembly of the foot stirrup**

> **Required tools:** Torx screwdriver T20.

- 1) Press in the ankle setting nuts on the inside of the upper foot stirrup part (see fig. 3).
- 2) Set the lower foot stirrup part onto the outside of the upper foot stirrup part (see fig. 4).
- 3) Insert the screws and tighten them slightly (see fig. 5).

#### **Assembly of the orthosis**

> **Required tools:** Torx screwdriver T20, T30.

- 1) Shape the stirrup with insert sole or inner sandal.
- 2) Loosen the central screw of the ankle joint lower part and remove the threaded plate.
- 3) Place the foot stirrup onto the ankle joint lower part (see fig. 6).
- 4) Reinsert the threaded plate (see fig. 7).
- 5) Install the foot stirrup on the ankle joint lower part with the central screw (see fig. 8).
- 6) Put the orthosis on together with the shoe and insert.
- 7) Position the foot stirrup so that the Torx screw of the orthosis ankle joint is at the level of the lateral malleolus.
- 8) Tighten the screws to secure the foot stirrup.
- 9) Mark the position of the foot stirrup in the a-p direction and the exterior rotation on the insert.
- 10) Drill the holes for riveting the insert.
- 11) Rivet the insert to the foot stirrup.

## 5.2 Trial fitting

After assembly, the orthosis can be fitted and the fit and function can be checked on the patient.

### 5.2.1 Functional test

To verify technically correct functionality, take the foot section in one hand and the femoral section in the other:

#### Free Walk mode

- > The 3-phase switch is in the upper position.
- ▶ Push the orthosis ankle joint in the dorsiflexion direction and extend the orthosis knee joint.
  - The lock disengages automatically and the orthosis knee joint can be flexed.

#### Permanent locking

- > The knee joint is set to Free Walk mode.
- 1) Slide the 3-phase switch to the lower position with the knee extended.
- 2) Push the orthosis ankle joint in the dorsiflexion direction and extend the orthosis knee joint.
  - The knee joint is permanently locked.

#### Free motion mode

- > The knee joint is set to permanently locked mode.
- 1) Slide the 3-phase switch to the upper position.
- 2) Push the orthosis ankle joint in the dorsiflexion direction and extend the orthosis knee joint.
- 3) Flex the orthosis knee joint.
- 4) Slide the 3-phase switch to the lower position with the knee flexed.
  - The orthosis knee joint can be flexed freely at all times.

In case of problems, please contact your Ottobock Free Walk specialist.

### 5.2.2 Verification of fit

#### Checking the fit in the frontal plane

- > **Required tools:** 170W2 tool set, 170W13 aluminium clamping jaws, 170W11 frontal contour bending tool
- 1) Check the fit of the base tube on the leg contour.
- 2) Mount the 170W13 aluminium clamping jaws on a vice.
- 3) Adapt the base tube in the vice. Use the 170W11 frontal contour bending tool for this purpose.

**NOTICE! Avoid crushing the base tube by excessive shaping.**

#### Verifying the fit of the support tubes

Support tubes that have too small a diameter may result in medial pressure areas. Support tubes with too large a diameter could obstruct sitting, especially in the thigh area.

In either case, the patient's individual wishes should be taken into account while checking the fit of the support tubes.

- > **Required tools:** 170W13 clamping jaws, 170W12\* pad holder bending tool, 170W18 bending iron
- 1) Loosen the screw on the strap guide loop with tube connection by 2 turns.
- 2) Push the screw in until a distinct clicking is heard. Use an Allen key to push it in.
- 3) Unscrew the strap guide loop with tube connection.
- 4) Remove the support tube. Loosen the two screws in order to do so.
- 5) Mount the 170W13 aluminium clamping jaws on a vice.
- 6) Secure the support tube. Use the 170W12\* pad holder bending tool in order to do so (see fig. 9).
- 7) **NOTICE! Avoid crushing the support tube by excessive clamping.**
- 8) Correct the shape of the support tubes. Use the 170W18 bending iron in order to do so.

- 9) Install the support tubes. In order to do so, apply Loctite® 241 to the screws and tighten them alternately until the tensioning wedge is flush with no gaps (**torque value 4 Nm**).

### Checking the fit in the sagittal plane

The fit of the orthosis in the sagittal plane is of crucial importance for the locking and unlocking function of the orthosis knee joint.

- > To verify technical functionality, have the patient seated on a chair. If necessary, help the patient unlock the knee joint to allow sitting with the knee flexed.
- ▶ Take the calf in your hand and bring it slowly into the extended position.
  - If the anatomical limit in the knee joint is reached before the orthosis extension stop, then the orthosis must be aligned with slightly more flexion.

The corrections for the standing position can be performed by re-shaping the femoral support tubes.

- 1) Remove the axis of the orthosis knee joint.
  - The upper and lower parts are separated from each other.
- 2) Correct the fit in the sagittal plane in 2° increments.

**INFORMATION: Corrections of 2° or 4° may result in significant differences in the fit.**

### 5.2.3 Checking freedom of movement in the ankle

#### NOTICE

#### Improper settings

Loss of functionality due to damaged steel cable

- ▶ Do not use the set screw of the steel cable on the ankle joint bottom part to adjust the lock.

- 1) Check the freedom of movement of the orthosis ankle joint.

**INFORMATION: The orthosis ankle joint must have at least 5° mobility up to the dor-siflexion stop in order to ensure automatic unlocking of the orthosis knee joint.**
- 2) **If needed:** Use the central screw of the orthosis ankle joint (see fig. 8, item 1) to adjust the neutral position according to the scale on the foot stirrup (see fig. 8, item 2).

**INFORMATION: In a standing position, ensure that the orthosis ankle joint is not too far in dorsal or plantar flexion since this can result in unintentional unlocking of the orthosis knee joint.**
- 3) **If needed:** Bend the foot stirrup with insert with the help of the 711S7=3.5x3 bending iron.

### 5.3 Finishing the orthosis

#### CAUTION

#### Unallowable use of lubricants

Risk of injury and loss of functionality as well as damage to the product

- ▶ The product is designed to be free of lubricants. Avoid contact with lubricants (e.g. grease, oil).
- ▶ After contact with lubricants, promptly clean the product or parts of the product with a degreasing cleaning agent (e.g. acetone or isopropyl alcohol).

> Required tools: Torque wrench

- 1) Apply Loctite® 241 to the central screw of the orthosis ankle joint and tighten it (**torque value 14 Nm**).
- 2) Apply Loctite® 241 to the screws of the foot stirrup and tighten them (**torque value 2.5 Nm**).
- 3) **Optional:** Trim all hook-and-loop straps of the closure straps so that sufficient material remains for attaching the Y hook-and-loop closure.
- 4) Attach the Y hook-and-loop closure.
- 5) **Optional:** Sew the Y hook-and-loop closure of the orthosis to the loop strap.

## 6 Handling

### **WARNING**

#### **Operation of motor vehicles**

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

### **CAUTION**

#### **Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).

### **CAUTION**

#### **Use on more than one patient**

Skin irritation, formation of eczema or infections

- ▶ Only use the orthosis on one patient.

### **CAUTION**

#### **Improper handling**

Risk of falling due to accidentally opening or engaging the lock

- ▶ Avoid accidentally opening or engaging the lock.
- ▶ Instruct the patient in the correct handling of the orthosis joints.

### **CAUTION**

#### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.

### **CAUTION**

#### **Using orthosis without padding**

Risk of local pressure points, superficial injuries or skin irritation due to skin contact with rigid or sharp orthosis components

- ▶ Only use the orthosis with undamaged padding.

### **6.1 Trial walking**

The Ottobock Free Walk orthosis deviates from conventional orthosis designs. This means that patients who have been using a conventional orthosis for a long time, in particular, will have to learn the principle.

For further information, please consult the document on therapeutic application and gait training (646A214).

## 7 Cleaning

### CAUTION

#### Improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

Clean the orthosis regularly:

- 1) Remove the straps from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 8 Maintenance

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

The wear parts of the knee joint and ankle joint are to be replaced by the maintenance sets 170Z99=1/2 (see fig. 2, item 1), 170X18 (see fig. 2, item 2) and 170X17=L/R (see fig. 2, item 3). Pads and closures for replacement are shown in the component overview (see Page 13) (colour code: -7 = black, -0 = skin colour).

### Replacing the pads

- 1) Loosen the screw on the strap guide loop with tube connection by 2 turns.
- 2) Push the screw in until a distinct clicking is heard. Use an Allen key to push it in.
- 3) Unscrew the strap guide loop with tube connection.
- 4) Remove and replace the pads.
- 5) Remove the strap guide loops from the end of the support tube.

### Replacing the lock

> **Required tools and materials:** Loctite® (636K13), Allen key (**2 mm**), torque wrench (adjustable to **0.5 Nm**).

- 1) Remove the retaining ring from the knee axis in the orthosis knee joint.
- 2) Push the knee axis out to the inside.  
→ The orthosis is disassembled in the orthosis knee joint.
- 3) Loosen the set screw of the steel cable with the Allen key (2 mm).
- 4) Remove the set screw and pressure ball from the orthosis ankle joint.
- 5) Remove the retaining ring of the lock axis.
- 6) Push the lock axis out to the inside.
- 7) Remove the lock with steel cable upwards from the orthosis knee joint and replace the lock.
- 8) Insert the flange bushings into the lock by hand and press them in lightly using the vice.
- 9) After pressing them in, finish the flange bushings with a 6H7 reamer.
- 10) Knot the pull cable at the lock.

**INFORMATION: Note that the knot has to be at least 5 cm away from the lock.**

- 11) Put the orthosis lower part with the strap guide loops down onto the workbench.
- 12) Insert the steel cable into the hole in the orthosis knee joint and push it through the base tube.

**INFORMATION: Note that the steel cable comes out the back of the base tube but before the orthosis ankle joint.**

- 13) With one hand, pull the steel cable further out of the orthosis ankle joint and insert the lock in the orthosis knee joint with the other hand.

**INFORMATION: Verify proper fit in the bearing bushings.**

- 14) Position the lower spring of the lock so that it does not obstruct the steel cable.

- 15) Insert the lock axis from medial into the orthosis knee joint and slide on the retaining ring.
- 16) Bring the orthosis ankle joint forward to the stop and thread the steel cable from above through the orthosis ankle joint.
- 17) Insert the upper part of the orthosis with new bearing washers into the orthosis knee joint.
- 18) Slide the knee axis from the medial side into the orthosis knee joint and slide on the retaining ring.
  - The orthosis knee joint now has to be locked.
- 19) Insert the pressure ball into the orthosis ankle joint.
- 20) Apply Loctite® 241 to the locking screw and install it. Do not tighten the locking screw.
- 21) Keep slight tension on the steel cable by hand and lift the lock by about 1 mm. Tighten the locking screw at the same time (**torque value 0.5 Nm** ).
  - The lock has to unlock due to dorsiflexion in the orthosis ankle joint and extension in the orthosis knee joint.
- 22) Shorten the pull cable. The pull cable has to project 2 cm out of the ankle joint.

## 9 Disposal

Dispose of the product in accordance with national regulations.

## 10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 10.2 CE conformity

This custom product meets the requirements of the European Directive 93/42/EEC for medical devices according to Annex VIII. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The product does not bear the CE marking since it is a custom fabrication.

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