X2 - Bionic Prosthetic System 3B5-X2/3B5-X2=ST

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1	Important information.....	60
2	Foreword	62
3	Product description	62
3.1	Design	62
3.2	Function.....	62
4	Application.....	63
4.1	Indications for use	63
4.2	Conditions of use	63
4.3	Indications.....	63
4.4	Qualification	64
5	Safety.....	64
5.1	Explanation of warning symbols.....	64
5.2	Structure of the safety instructions	64
5.3	General safety instructions.....	64
5.4	Information on the Power Supply/Battery Charging.....	67
5.5	Battery charger information.....	67
5.6	Information on Proximity to Certain Areas	68
5.7	Information on Use	69
5.8	Notes on the safety modes.....	70
5.9	Information on the use of a mobile device with the cockpit app.....	71
6	Scope of Delivery and Accessories.....	72
6.1	Scope of delivery.....	72
6.2	Accessories	72
7	Charging the battery	73
7.1	Connecting the power supply and battery charger	73
7.2	Charging the prosthesis battery.....	73
7.3	Display of the current charge level	74
7.3.1	Display of battery charge level without additional devices	74
7.3.2	Display of the current charge level using the Cockpit app	74
7.3.3	Display of the current charge level during the charging process	75
7.3.4	Display of current charge level using remote control (optional accessory)	75
8	Cockpit app.....	75
8.1	System Requirements	76
8.2	Initial connection between cockpit app and prosthesis	76
8.2.1	Starting the cockpit app for the first time	76
8.3	Control elements for cockpit app.....	77
8.3.1	Cockpit app navigation menu.....	78
8.4	Managing components.....	78
8.4.1	Adding component	78
8.4.2	Deleting a component	79

8.4.3	Connecting component with multiple devices.....	79
9	Use	80
9.1	Movement patterns in basic mode (mode 1)	80
9.1.1	Standing	80
9.1.1.1	Stance function	80
9.1.2	Walking.....	80
9.1.3	Running short distances ("walk-to-run" function)	81
9.1.4	Sitting down	81
9.1.5	Sitting	81
9.1.5.1	Sitting function.....	81
9.1.6	Standing up	82
9.1.7	Walking up stairs step-over-step.....	82
9.1.8	Overcoming obstacles	83
9.1.9	Walking down stairs	83
9.1.10	Walking down a ramp	83
9.2	Changing prosthesis settings.....	83
9.2.1	Changing the prosthesis setting using the cockpit app.....	84
9.2.2	Overview of adjustment parameters in basic mode.....	84
9.2.3	Overview of adjustment parameters in MyModes	85
9.3	Turning Bluetooth on the prosthesis on/off.....	86
9.3.1	Switching Bluetooth off/on using the cockpit app	87
9.4	Querying the prosthesis status	87
9.4.1	Query status through cockpit app.....	87
9.4.2	Status display in the cockpit app	87
9.5	Mute mode (silent mode)	87
9.5.1	Turning mute mode on/off using the Cockpit app	87
9.6	OPG function (Optimised Physiological Gait).....	88
10	MyModes	88
10.1	Running function as configured MyMode.....	89
10.2	Switching MyModes with the cockpit app	89
10.3	Switching MyModes using motion patterns	90
10.4	Switching from a MyMode back to basic mode.....	91
11	Additional operating states (modes)	91
11.1	Empty battery mode	91
11.2	Mode for charging the prosthesis	92
11.3	Safety mode.....	92
11.4	Overheating mode	92
12	Storage and bleeding	92
13	Cleaning	93
14	Maintenance	93
15	Legal information	93
15.1	Liability	93

15.2 Trademarks 93

15.3 CE conformity 93

15.4 Local Legal Information 94

16 Technical data 95

17 Appendices 97

17.1 Symbols Used 97

17.2 Operating states/error signals 99

17.2.1 Signals for operating states 99

17.2.2 Warnings/error signals 99

17.2.3 Error messages while establishing a connection with the cockpit app 102

17.2.4 Status signals 102

17.3 Directives and manufacturer’s declaration 104

17.3.1 Electromagnetic environment 104

1 Important information

The functional scope of products that have been serviced or repaired by an authorised Ottobock Service Centre since October of 2018 has been extended.

Extended functional scope:

- Switching to pre-configured MyModes and changing the prosthesis settings by using the remote control or an app for smartphones (4X441 Cockpit) (see Page 75, see Page 84). If the joint bears the marking "DUAL", the remote control cannot be used anymore (see next section).
- Possibility of using an intentional stance function (see Page 80)
- Displays of the current battery level without additional devices (see Page 74)
- Deactivation of feedback and vibration signals (mute mode – see Page 87)
- Stance, sitting, and stairs function can be switched off (see Page 84)
- Bluetooth function can be switched off (see Page 86)

INFORMATION

Remote control function only possible with version 1.0.0.6 or above

Please note that you can only use this and/or the updated knee joint with the **4X350-X3 remote control in version 1.0.0.6 or higher**.

Therefore, check the version of the remote control as follows:

- ▶ With the remote control switched on, use the ▲ and ▼ keys to select the menu option **RC settings** and confirm with the ■ key.
- ▶ Use the ▲ and ▼ keys to select the menu option **Info** and confirm this with the ■ key.
- ▶ Use the ▲ and ▼ keys to select the **RC version: 1.0.0.4** menu option. If an older version appears than **1.0.0.6**, the remote control must be updated by an authorised Ottobock Service Centre.


The function of the remote control with this or the updated knee joint is only possible with version 1.0.0.6 or higher.

INFORMATION

Adjustment of the knee joint possible only with X-Soft 4X1 version 1.8 and higher

Please note that you need X-Soft 4X1 version 1.8 or higher in order to adjust this knee joint or the updated one. If X-Soft version 1.0 or higher is available, it can be updated (Alignment)

Using the 4X350-X3 remote control is no longer possible on products with the marking "DUAL"

2 Foreword

INFORMATION

Date of last update: 2018-12-21

- ▶ Please read this document carefully before using the product.
- ▶ Have yourself instructed by qualified personnel in the proper and safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). The contact information for the manufacturer is found on the back page.
- ▶ Please keep this document in a safe place.

The product "X2 bionic prosthetic system" is called the product/prosthesis/knee joint/component in the following.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

3 Product description

3.1 Design

The product consists of the following components:

1. Proximal pyramid adapter
2. Optional flexion stops
3. Rechargeable battery
4. Hydraulic unit
5. LED (blue) as indicator for the Bluetooth connection
6. Receiver of the inductive charging unit

3.2 Function

This product features microprocessor control of the stance and swing phase.

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the damping behaviour of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

Thanks to the microprocessor-controlled stance and swing phase, the product can be individually adapted to your needs.

For this purpose, the product is adjusted by qualified personnel using adjustment software.

The product features MyModes for special motion types (e.g. cycling ...). These are pre-configured by the orthotics and prosthetics professional (O&P professional) using the adjustment software and can be activated with special motion patterns and the Cockpit app (see Page 88).

In case of a product malfunction, safety mode makes restricted operation possible. Resistance parameters that are predefined by the product are configured for this purpose (see Page 92).

Empty battery mode permits safe walking when the battery is drained. Resistance parameters that are predefined by the product are configured for this purpose (see Page 91).

The microprocessor-controlled hydraulic unit offers the following advantages

- Approximation of the physiological gait pattern
- Stability while standing and walking
- Adaptation of product characteristics to various surfaces, inclines, gait situations and walking speeds

4 Application

4.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.

4.2 Conditions of use


The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see Page 95).

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.


Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

Knee joint with attached 2R20 AXON tube adapter

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of up to **150 kg (330 lbs)**.

Knee joint with attached 2R21 AXON tube adapter with torsion

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of up to **125 kg (275 lbs)**.

4.3 Indications

- For users with unilateral knee disarticulation and unilateral transfemoral amputation
- Dysmelia patients with residual limb characteristics corresponding to knee disarticulation, transfemoral amputation or hip disarticulation
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

4.4 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

5 Safety

5.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

5.2 Structure of the safety instructions

 WARNING
<p>The heading describes the source and/or the type of hazard</p> <p>The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:</p> <ul style="list-style-type: none">> E.g.: Consequence 1 in the event of failure to observe the hazard> E.g.: Consequence 2 in the event of failure to observe the hazard▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

5.3 General safety instructions

 WARNING
<p>Non-observance of safety notices</p> <p>Personal injury/damage to the product due to using the product in certain situations.</p> <ul style="list-style-type: none">▶ Observe the safety notices and the stated precautions in this accompanying document.

 WARNING
<p>Using the prosthesis while operating a vehicle</p> <p>Accidents due to unexpected behaviour of the prosthesis because of changed damping behaviour.</p> <ul style="list-style-type: none">▶ All users are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.▶ Observe national legal regulations for retrofitting your vehicle in accordance with the type of fitting.

 WARNING
<p>Use of damaged power supply unit, adapter plug or battery charger</p> <p>Risk of electric shock due to contact with exposed, live components.</p> <ul style="list-style-type: none">▶ Do not open the power supply unit, adapter plug or battery charger.▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ CAUTION

Failure to observe warning/error signals

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see Page 99) and corresponding change in damping settings must be observed.

⚠ CAUTION

Failure to observe activated mute mode (silent mode)

Falling due to unexpected product behaviour because of changed damping behaviour.

The following feedback signals are deactivated when mute mode is activated:

- > Long vibration signal if the hydraulic unit overheats.
- > Beep and vibration signal to confirm that the movement pattern has been recognised (switching to a MyMode/basic mode with movement pattern).
- > Beep and vibration signal to indicate successful switching to a MyMode/basic mode.
- > Beep and vibration signal upon successfully switching to deep sleep mode.
- ▶ Before activating mute mode, note that these feedback signals will be deactivated. For more information about mute mode, see the section "Mute mode" (see Page 87).
- ▶ After switching to a MyMode/basic mode, verify the changed damping settings.
- ▶ Ensure that you stand securely during all switching processes.
- ▶ To deactivate mute mode, connect and then disconnect the battery charger.

⚠ CAUTION

Independent manipulation of the product and the components

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by authorised, qualified Ottobock personnel (no replacement by the user).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Mechanical stress on the product

> Falling due to unexpected product behaviour as the result of a malfunction.

- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Use of the product when battery charge level is too low

Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the prosthesis if required.
- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

⚠ CAUTION

Risk of pinching in the joint flexion area

Injuries due to pinching of body parts.

- ▶ Ensure that fingers/body parts or soft tissue of the residual limb are not in this area when bending the joint.

⚠ CAUTION

Penetration of dirt and humidity into the product

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that neither solid particles, foreign objects nor liquids penetrate into the product.
- ▶ Do not expose the product to splashed water.
- ▶ In the rain, thick clothing should be worn over the product as a minimum.
- ▶ If water has penetrated system components, remove the protector and allow the components to dry. The prosthesis must be inspected by an authorised Ottobock Service Centre. The prosthetist is your contact person.
- ▶ If salt water has penetrated the prosthesis, the protector must be removed immediately. The prosthesis must be inspected by an authorised Ottobock Service Centre. The prosthetist is your contact person.

⚠ CAUTION

Use of the product without protector or with damaged protector

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ If the protector has been removed, ensure prior to the next use of the product that the protector has been properly installed.
- ▶ Using the product with a damaged protector or without the protector is not permitted.
- ▶ Using the product with a foam cover is not possible, since the protector would have to be removed for this purpose.

⚠ CAUTION

Signs of wear to the product components

Falling due to damage or malfunction of the product.

- ▶ In the interest of your own safety and in order to maintain operating reliability and protect the warranty, the service intervals must be observed.

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product with a damp cloth only (fresh water).

INFORMATION

Knee joint movement noise

When using exoprosthetic knee joints, servomotor, hydraulic, pneumatic or brake load dependent control functions can cause movement noise. This kind of noise is normal and unavoidable. It generally does not indicate any problems. If movement noise increases noticeably during the lifecycle of the knee joint, the knee joint should be inspected by an O&P professional immediately.

5.4 Information on the Power Supply/Battery Charging

⚠ CAUTION

Charging the product without taking it off

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ For safety reasons, take the product off before charging the battery.

⚠ CAUTION

Charging the product with damaged power supply unit/charger/charger cable

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

5.5 Battery charger information

NOTICE

Penetration of dirt and humidity into the product

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

NOTICE

Mechanical stress on the power supply/battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

NOTICE

Operating the power supply unit/charger outside of the permissible temperature range

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see Page 95).

NOTICE

Independent changes or modifications carried out to the battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

NOTICE

Contact of the battery charger with magnetic data storage devices

Wiping of the data storage device.

- ▶ Do not place the battery charger on credit cards, diskettes, audio or video cassettes.

5.6 Information on Proximity to Certain Areas

CAUTION

Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.50 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.35 m
 - DECT cordless phones incl. base station: 0.18 m
 - WiFi (routers, access points,...): 0.11 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.11 m

CAUTION

Operating the product in very close proximity to other electronic devices

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for people (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations...).
- If this cannot be avoided, make sure to at least have a safeguard when walking or standing (e.g. a handrail or the support of another person).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected changes in the damping behaviour of the product.
- ▶ In general, monitor the product for unexpected changes in the damping behaviour when electronic or magnetic devices are in the immediate vicinity.

CAUTION

Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of strong magnetic fields.
- ▶ Take off the product before entering a room or area with strong magnetic fields and store the product outside this room or area.

- ▶ Damage to the product caused by exposure to strong magnetic fields cannot be repaired.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Falling due to malfunction or the breakage of load-bearing product components.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see Page 95).

5.7 Information on Use

⚠ CAUTION

Walking up stairs

Falling due to foot placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking up stairs and place most of the area of the sole of your foot on the stair surface.
- ▶ Particular caution is required when carrying children up stairs.

⚠ CAUTION

Walking down stairs

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and roll over the edge of the step with the middle of the shoe.
- ▶ Observe the warning/error signals (see Page 99).
- ▶ Be aware that resistance in the flexion and extension direction can change in case of warning and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

⚠ CAUTION

Overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)

- > Falling due to unexpected behaviour of the product because of switching into overheating mode.
- > Burns due to touching overheated components.
- ▶ Be sure to pay attention when pulsating vibration signals start. They indicate the risk of overheating.
- ▶ As soon as these pulsating vibration signals begin, you have to reduce the activity level so the hydraulic unit can cool down.
- ▶ Full activity may be resumed after the pulsating vibration signals stop.
- ▶ If the activity level is not reduced in spite of the pulsating vibration signals, this could lead to the hydraulic element overheating and, in extreme cases, cause damage to the product. In this case, the product should be checked for damage by a prosthetist. If necessary, the prosthetist will forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Overloading due to unusual activities

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by a prosthetist. If necessary, the prosthetist will forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Improper mode switching

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.
- ▶ Take the weight off the product and correct the switching, if required.

⚠ CAUTION

Improper use of the stance function

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that you are standing safely when using the stance function and check the lock of the knee joint before placing your full weight on the prosthesis.
- ▶ Make sure you have been instructed in the proper use of the stance function by a prosthetist or therapist. Information on the stance function see Page 80.

⚠ CAUTION

Quickly pushing the hip forward with the prosthesis extended (e.g. serve while playing tennis)

- > Falling due to unexpected activation of the swing phase.
- ▶ Note that the knee joint may flex unexpectedly when the hip is pushed forward quickly while the prosthesis is extended.
- ▶ Therefore, familiarise yourself with swing phase activation in such situations under secure conditions (e.g. while holding on to parallel bars) and with the instruction of trained, qualified personnel.
- ▶ For sports where this movement pattern may occur, use a corresponding pre-configured MyMode. For further information about the MyModes, see the section 'MyModes' (see Page 88).

5.8 Notes on the safety modes

⚠ CAUTION

Using the product in safety mode

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see Page 99) have to be observed.
- ▶ Particular caution is necessary when using a bicycle without a freewheel (with a fixed gear).

⚠ CAUTION

Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Do not continue using the defective product.
- ▶ Contact your prosthetist immediately.

⚠ CAUTION

Safety mode cannot be deactivated

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ If you cannot deactivate safety mode by recharging the battery, a permanent malfunction has occurred.
- ▶ Do not continue using the defective product.
- ▶ The product must be inspected by an authorised Ottobock Service Centre. The prosthetist is your contact person.

⚠ CAUTION

Safety signal occurs (ongoing vibration)

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Observe the warnings/error signals (see Page 99).
- ▶ Do not continue using the product after the safety signal has been emitted.
- ▶ The product must be inspected by an authorised Ottobock Service Centre. In this case, the contact person is your prosthetist.

5.9 Information on the use of a mobile device with the cockpit app

⚠ CAUTION

Improper use of the device

Falling due to altered damping behaviour as a result of unexpected switching into MyMode.

- ▶ Make sure you have been instructed in the proper use of the device with the cockpit app.

⚠ CAUTION

Improper use of the setting parameters in the MyModes

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Get the O&P professional and/or therapist to instruct you regarding the functionality and adjustment options for **all parameters** of the MyModes.

⚠ CAUTION

Independently applied changes or modifications made to the mobile device

Falling due to altered damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Do not make any independent changes to the hardware of the mobile device on which the app is installed.
- ▶ Do not make any independent changes to the software/firmware of the mobile device that are not included in the update function of the software/firmware.

⚠ CAUTION

Improper mode switching with the device

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter (beeper) and the display on the device.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.

NOTICE

Destruction of the device due to falling or penetration of water

Malfunction of the device.

- ▶ Note the instructions for use for the device.
- ▶ If it should no longer be possible to switch back from a MyMode to basic mode, the component can only be switched back to basic mode by using a movement pattern (see Page 91) or by connecting/disconnecting the inductive battery charger.

NOTICE

Failure to observe the system requirements for the installation of the Cockpit app

Malfunction of the mobile device.

- ▶ Install the Cockpit app only on the operating systems listed in the section "System requirements" (see Page 76).

INFORMATION

The illustrations in these instructions for use are only examples and may deviate from the respective mobile device being used and the version.

6 Scope of Delivery and Accessories

6.1 Scope of delivery

- 1 pc. 3B5-X2=ST; X2 bionic prosthetic system (with threaded connector) or 1 pc. 3B5-X2; X2 bionic prosthetic system (with pyramid connector) with preinstalled 4X900 Protective Cover
- 1 pc. 2R20 AXON tube adapter or 1 pc. 2R21 AXON tube adapter with torsion
- 1 pc. 646D1356 instructions for use (user)
- Android app "Cockpit 4X441-Andr=V*" for download from the website: <http://www.ottobock.com/cockpitapp>
- 1 pc. 757L16* power supply
- 1 pc. 4E60* inductive charger
- 1 pc. cosmetic case for battery charger and power supply
- 1 pc. 646C107 Bluetooth PIN card
- 1 pc. 647F542 prosthesis passport
- 1 pc. card holder for prosthesis passport

Only when "DUAL" marking is present (see Page 60)

- iOS app "4X441-IOS=V* Cockpit" for download from the website: <http://www.ottobock.com/cockpitapp>

6.2 Accessories

The following components are not included in the scope of delivery and may be ordered separately:

- Genium X3 Protective Cover: 4X900

- Genium X3 Protective Cover: 4X193-1

Only when "DUAL" marking is not present (see Page 60)

- 4X350-X3 remote control
The user interface is available in the following languages: German, English, French, Italian, Spanish, Dutch, Swedish
- 4X194 wrist band for 4X350-X3 remote control

7 Charging the battery

The following points must be observed when charging the battery:

- With average use, the capacity of the fully charged battery is sufficient for about 5 days.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.
- Use the 757L16* power supply and 4E60* battery charger to charge the battery.
- Note the permissible temperature range for charging the battery (see Page 95).
- The distance between the battery charger and the receiver on the product must not exceed 2 mm.

7.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
 - 2) Connect the round, **three-pin** plug of the power supply to the receptacle on the battery charger (see fig. 2) so that the plug locks into place.
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
 - 3) Plug the power supply unit into the outlet (see fig. 3).
 - The green LED on the back of the power supply lights up.
 - The LED ring (status indicator) on the rear of the charger lights up green to indicate the correct connection to the power supply.
- If the green LED on the power supply and the LED ring on the battery charger do not light up, there is an error (see Page 99).

7.2 Charging the prosthesis battery

INFORMATION

When the Protective Cover is installed, the battery charger cable has to point to the upper closure. A correct knee joint charging process is only ensured with this alignment.

- 1) Connect the inductive charger to the receiver of the charging unit on the rear of the product. The charger is held in place by a magnet.
 - The LED ring on the rear of the charger pulsates purple (4-second cycle).
 - If the LED ring lights up in a different colour, this indicates an error (see Page 99).
- 2) The charging process starts.
 - Once the product battery is fully charged, all LEDs on the side of the battery charger light up.
- 3) After the charging process is complete, hold the prosthesis still and remove the inductive charger from the receiver.
 - A self-test is performed. The joint is operational only after corresponding feedback (see Page 102).

7.3 Display of the current charge level

7.3.1 Display of battery charge level without additional devices

INFORMATION

The charge level cannot be displayed during the charging process, e.g. by turning the prosthesis over. The product is in charging mode.

- 1) Turn the prosthesis by 180° (the sole of the foot must face up).
- 2) Hold still for 2 seconds and wait for beep signals.

Beep signal	Battery charge level
5x short	more than 80%
4x short	60% to 80%
3x short	40% to 60%
2x short	20% to 40%
1x short	less than 20%

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app (see Page 83) or if mute mode (silent mode) is activated, there are no beep signals.

7.3.2 Display of the current charge level using the Cockpit app

Once the Cockpit app has been started, the current charge level is displayed in the bottom line of the screen:

1. 38% – Charge level of battery for currently connected component

7.3.3 Display of the current charge level during the charging process


During the charging process, the current battery charge level is indicated by the LEDs lit on the side of the charger.

	Quantity	Battery charge level
	0	0%-10%
	1	10%-30%
	2	30%-50%
	3	50%-70%
	4	70%-90%
5	> 90%	

7.3.4 Display of current charge level using remote control (optional accessory)

Only when "DUAL" marking is not present (see Page 60)

The current charge level is displayed in the status bar of the remote control when it is connected and switched on:

1. – Charge level of battery for currently connected prosthesis
2. – Prosthesis is being charged
3. 68% – Charge level of battery for currently connected prosthesis, as a percentage

8 Cockpit app

The cockpit app enables switching from basic mode into the pre-configured MyModes. In addition, information about the product (step counter, battery charge level, etc.) can be called up.

The behaviour of the product can be changed to a certain extent on a day-to-day basis using the app (e.g. while becoming accustomed to the product). At the next appointment, the prosthetist can trace the changes using the adjustment software.

INFORMATION

The Cockpit app can be downloaded free of charge from the respective online store. For more information, please visit the website: <http://www.ottobock.com/cockpitapp>. To download the Cockpit app, the QR code on the supplied Bluetooth PIN card can be read with the mobile device (requirement: QR code reader and camera).

INFORMATION

The language of the Cockpit app user interface can only be changed by the prosthetist using the adjustment software.

INFORMATION

The serial number of the component to be connected has to be registered with Ottobock the first time it is connected. If the registration is not accepted, use of the Cockpit app for this component will be limited.

INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (sole of the foot must point up) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 86).

8.1 System Requirements

The functioning of the Cockpit app is assured on devices that support the following operating systems:

Only when "DUAL" marking is not present (see Page 60)

- **Android:** version 4.0.3 or higher

Only when "DUAL" marking is present (see Page 60)

- The 4X441=* Cockpit app must be installed in version 2.0 or higher
- **Android:** version 4.1 or higher
- **iOS (for iPhone, iPad, iPod):** version 9.3 or higher. The mobile device has to support Bluetooth 4.0 (BT LE).

8.2 Initial connection between cockpit app and prosthesis

The following points should be observed before the initial connection:

- Bluetooth of the component must be switched on (see Page 86).
- Bluetooth on the mobile device must be switched on.
- The device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- **The device must be connected to the internet.**
- The serial number and Bluetooth PIN of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The serial number begins with the letters "SN".

INFORMATION

If the Bluetooth PIN card with the Bluetooth PIN and serial number of the component is lost, contact your O&P professional.

8.2.1 Starting the cockpit app for the first time

- 1) Tap the symbol of the Cockpit app ().
→ The end user license agreement (EULA) is displayed.
- 2) Accept the end user license agreement (EULA) by tapping the **Accept** button. If the end user license agreement (EULA) is not accepted, the Cockpit app cannot be used.
→ The welcome screen appears.
- 3) Hold the prosthesis with the sole of the foot facing up, or connect and then disconnect the battery charger, in order to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
- 4) Tap the **Add component** button.

- The Connection Wizard opens and guides you through the process of establishing a connection.
- 5) Follow the subsequent instructions on the screen.
- 6) After the Bluetooth PIN is entered, a connection to the component is established.
 - While the connection is being established, 3 beep signals sound and the  symbol appears.
 - The  symbol is displayed when the connection has been established.
- Once the connection has been established, the data are read from the component. This process may take up to a minute.
- Then the main menu appears with the name of the connected component.

INFORMATION


After the initial connection to the component has been established successfully, the app will connect automatically each time it is started. No further steps are required.


Only when "DUAL" marking is present (see Page 60)

INFORMATION


After activating the "visibility" of the component (holding the component with the sole of the foot facing up, or connecting and then disconnecting the battery charger), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case, hold the component with the sole of the foot facing up again, or connect and then disconnect the battery charger.

8.3 Control elements for cockpit app

1. ☰ Access the navigation menu (see Page 78)
2. Product
 - The component name can only be changed with the adjustment software.
3. If connections to more than one component have been saved, you can switch between the saved components by tapping the entry **change**.
4. MyModes configured with the adjustment software.
 - Switching the mode by tapping the corresponding symbol and confirming by tapping "OK".
5. Currently selected mode
6. Charge level of the component.
 -  Component battery fully charged
 -  Component battery empty
 -  Component battery is being charged
 - The current charge level is also displayed in %.
7. Display and name of the currently selected mode (e.g. **1. Basic Mode**)
8.  Mute mode is activated
9.  Connection to component has been established
 -  Connection to component has been interrupted. The app is attempting to re-establish the connection automatically.
 -  No existing connection to the component.

8.3.1 Cockpit app navigation menu

Tap the ☰ symbol in the menus to display the navigation menu. Additional settings for the connected component can be configured in this menu.

Product

Name of the connected component

MyModes

Return to the main menu to switch MyModes

Functions

Call up additional functions of the component (e.g. turn off Bluetooth) (see Page 86)

Settings

Change settings of the currently selected mode (see Page 83)

Status

Query status of the connected component (see Page 87)

Manage components

Add or delete components (see Page 78)

Imprint/Info

Display information/legal notices for the cockpit app

8.4 Managing components

Connections with up to four different components can be stored in the app. However, a component can only be connected to one device at a time.

The following points need to be observed before establishing the connection:

- Bluetooth of the component must be switched on (see Page 86).
- Bluetooth on the mobile device must be switched on.
- The device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- **The device must be connected to the internet.**
- The serial number and Bluetooth PIN of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The serial number begins with the letters "SN".

INFORMATION

If the Bluetooth PIN card with the Bluetooth PIN and serial number of the component is lost, contact your O&P professional.

8.4.1 Adding component

- 1) Tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Hold the prosthesis with the sole of the foot facing up, or connect and then disconnect the battery charger, in order to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
- 4) Tap the "+" button.
→ The Connection Wizard opens and guides you through the process of establishing a connection.
- 5) Follow the subsequent instructions on the screen.

- 6) After the Bluetooth PIN is entered, a connection to the component is established.
 - While the connection is being established, 3 beep signals sound and the  symbol appears.
 - The  symbol is displayed when the connection has been established.
 - Once the connection has been established, the data are read from the component. This process may take up to a minute.

Then the main menu appears with the name of the connected component.

INFORMATION

If establishing a connection to a component is not possible, perform the following steps:

- ▶ Delete the component from the Cockpit app if applicable (see the section "Deleting a component")
- ▶ Add the component again in the Cockpit app (see the section "Adding a component")

Only when "DUAL" marking is present (see Page 60)

INFORMATION

After activating the "visibility" of the component (holding the component with the sole of the foot facing up, or connecting and then disconnecting the battery charger), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case, hold the component with the sole of the foot facing up again, or connect and then disconnect the battery charger.


8.4.2 Deleting a component

- 1) Tap the  symbol in the main menu.
 - The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Tap the "**Edit**" button.
- 4) Tap the  symbol under the component you want to delete.
 - The component is deleted.

8.4.3 Connecting component with multiple devices

The connection for a component can be stored on more than one device. However, only one device can be connected to the component at one time.

If there is an existing connection between the component and a different device, the following information appears while the connection is being established with the current device:


- ▶ Tap the **OK** button.
- The connection to the last connected device is broken off and established with the current device.

9 Use

9.1 Movement patterns in basic mode (mode 1)

9.1.1 Standing

Knee control through high hydraulic resistance and static alignment. A stance function can be enabled by the prosthetist. Please see the following section for further information on the stance function.

9.1.1.1 Stance function

INFORMATION

To use this function, it needs to be enabled in the adjustment software. It also has to be activated using the Cockpit app (see Page 84).

The stance function (standing mode) is a functional supplement to the basic mode (mode 1). This function makes it easier, for example, to stand on an inclined surface for a longer time. The joint is fixed in the flexion direction.

The standing function must be enabled by the O&P professional. In addition, the type of joint locking (conscious/intuitive) has to be established by the O&P professional. Changing the locking type using the Cockpit app is not possible.

Intuitive locking of the joint

The intuitive stance function recognises any situation that puts strain on the prosthesis in the flexion direction but where flexion is not permitted. Examples of this include standing on uneven or sloping surfaces. The knee joint is always locked in the flexion direction when the prosthetic leg is not fully extended, under some amount of load and at rest. Upon forward or backward rollover or extension, the level of resistance is immediately reduced to stance phase resistance again.

The knee joint is not locked when the above conditions are met and a sitting position is assumed (for example while driving).

Deliberate locking of the joint


- 1) Assume the desired knee angle.
- 2) Do not remove the entire load from the prosthesis.
- 3) Do not change the knee angle for a brief period (1/8 second). This time period prevents unintentional activation of the standing function.

→ The blocked joint can now be fully loaded in the flexion direction.

Deliberate unlocking of the joint

- By deliberately extending or unloading the knee joint, it is unlocked again.

9.1.2 Walking


Initial attempts at walking with the prosthesis always require the instruction of trained, qualified personnel.

The hydraulics stabilise the knee joint in the stance phase and release the knee joint in the swing phase so that the leg can swing forward freely.

Switching to the swing phase requires that the prosthesis roll over to the front out of the stride position.

9.1.3 Running short distances ("walk-to-run" function)

For covering short distances quickly, the knee joint detects a transition from walking to running in basic mode and automatically changes the following settings:

- The swing phase angle is increased
- Preflexion of 4° at heel strike (PreFlex) is reduced to 0°

The requirements to automatically switch to the running motion are fast forward movement of the prosthetic leg and high dynamic load on the knee joint. When stopping from the running motion, the changed settings are set back to the standard values.

INFORMATION

For covering longer distances, the prosthetist can configure a "Running" MyMode (see Page 89).

9.1.4 Sitting down

The resistance in the prosthetic knee joint while sitting down ensures evenly bending into the sitting position.

The O&P professional can set either a constant resistance or a continuously increasing resistance for sitting down.


- 1) Position both feet side by side at the same level.
- 2) While sitting down, distribute your weight evenly between both legs and use the arm supports where applicable.
- 3) Move the buttocks in the direction of the back support and lean the upper body forward.

INFORMATION: Resistance while sitting down can be changed with the Cockpit app via the parameter "Resistance" (see Page 84).

9.1.5 Sitting

INFORMATION

While sitting, the knee joint also switches to energy saving mode. This energy saving mode is activated regardless of whether the sitting function is activated or not.

If the patient is in a sitting position for more than two seconds (i.e. the thigh is close to horizontal and there is no load on the leg), the knee joint switches the resistance to a minimum in the extension direction.

A sitting function can be enabled by the prosthetist. For more information about the sitting function, see the following section.

9.1.5.1 Sitting function

INFORMATION

To use this function, it needs to be enabled in the adjustment software. It also has to be activated using the Cockpit app (see Page 84).

In the sitting position, the resistance in the flexion direction is reduced in addition to the reduction of resistance in the extension direction. This makes it possible to swing the prosthetic leg freely.

9.1.6 Standing up

Flexion resistance is increased steadily while standing up.


- 1) Place the feet at the same level.
- 2) Lean the upper body forward.
- 3) Put the hands on armrests, if available.
- 4) Stand up with support from the hands while keeping weight evenly distributed on feet.

9.1.7 Walking up stairs step-over-step

INFORMATION

To use this function, it needs to be enabled in the adjustment software. It also has to be activated using the Cockpit app (see Page 84).

Although the knee joint is passive, which means it cannot actively initiate movements, negotiating stairs step-over-step is possible. This function must be practised and executed consciously.

- 1) Lift the extended prosthesis off the floor.
- 2) Immediately after lifting the extended leg off the floor, extend the hip briefly and then abruptly flex it. This requires a sufficiently secure hold in the socket and a certain level of residual limb strength.
→ This whip motion flexes the knee, because the knee joint automatically recognises the movement and sets the flexion resistance to minimum.


INFORMATION: Take note of people behind you before executing the whip motion.

- 3) When sufficient knee flexion has been achieved, the knee joint increases extension resistance so that there is enough time to position the foot on the next step before the knee joint is extended again.
- 4) Set the foot onto the next step.
The support area for the foot on the step must be large enough that the heel does not extend back too far over the edge. With too little support area, the lower leg would extend too early and position the leg too far backwards. In this phase, the knee joint has already set the flexion resistance to maximum (blocked). The knee joint can no longer be flexed but only extended. This ensures that the leg does not buckle if the hip strength is not sufficient for the extending motion.
- 5) Support yourself with your hand on the contralateral side. A smooth wall will also work. This lateral support is intended to prevent the residual limb from twisting in the socket. Twisting can lead to unpleasant surface tension between the skin and the socket. Lateral support also improves balance.
- 6) Bring the knee into extended position. When the knee joint is fully extended, the initial position has been reached.
- 7) You can climb the next step or continue walking normally.

9.1.8 Overcoming obstacles

INFORMATION

To use this function, it needs to be enabled in the adjustment software. It also has to be activated using the Cockpit app (see Page 84).


The stair function can also be used to cross obstacles:

- 1) Lift the extended prosthesis off the floor.
- 2) Briefly extend the hip.
- 3) Quickly flex the hip. This causes the knee to flex.
- 4) With the knee flexed, step over the obstacle.

With sufficient knee flexion, the extension resistance is increased to allow enough time for crossing the obstacle.

9.1.9 Walking down stairs


This function must be practised and executed consciously. Only when the sole is properly positioned can the knee joint react correctly and permit controlled flexion.

- 1) Hold the handrail with one hand.
- 2) Position the leg with the prosthesis on the step so that the foot projects halfway over the edge of the step.
→ This is the only way to ensure a secure rollover.
- 3) Roll the foot over the edge of the step.
→ This flexes the prosthesis slowly and evenly in the knee joint.
- 4) Place the foot of the other leg onto the next step.
- 5) Place the foot of the prosthetic leg on the next step after that.

INFORMATION: The flexion speed of the knee joint can be changed using the Cockpit app via the parameter "Resistance" (see Page 84).

9.1.10 Walking down a ramp

Under increased flexion resistance, permit controlled flexion of the knee joint which lowers the body's centre of gravity.

INFORMATION: The flexion resistance at which the knee bends can be changed using the Cockpit app via the parameter "Resistance" (see Page 84).

9.2 Changing prosthesis settings

Once an active connection to a component has been established, the settings **of the respective active mode** can be changed using the Cockpit app.

INFORMATION


Bluetooth on the prosthesis must be switched on to change the prosthesis settings.

If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. The connection must be established during this period.

Information for changing the prosthesis settings

- Before changing settings, always check the main menu of the Cockpit app to make sure the correct component has been selected. Otherwise parameters could be changed for the wrong component.
- It is not possible to change prosthesis settings nor to switch to a different mode while the prosthesis battery is being charged. Only the status of the prosthesis can be called up. Instead of the  symbol, the  symbol appears in the bottom row of the screen in the cockpit app.
- The O&P professional's setting is located in the middle of the scale. After making adjustments, this setting can be restored by tapping the **"Standard"** button in the Cockpit app.
- Prosthesis settings should be optimised using the adjustment software. The Cockpit app is not intended for use by the O&P professional to set up the prosthesis. The everyday behaviour of the prosthesis can be changed to a certain extent using the app (e.g. while becoming accustomed to the prosthesis). The O&P professional can use the adjustment software to trace these changes at the next appointment.
- If the settings of a MyMode are to be modified, one must first switch to this MyMode.

9.2.1 Changing the prosthesis setting using the cockpit app

- 1) Once the component is connected and in the desired mode, tap the  symbol in the main menu.
→ The navigation menu opens.
- 2) Tap the menu entry **"Settings"**.
→ A list appears with the parameters for the currently selected mode.
- 3) Change the setting of the desired parameter by tapping the "<", ">" symbols.

INFORMATION: The prosthetist's setting is marked and, after the setting has been changed, can be restored by tapping the "Standard" button.

9.2.2 Overview of adjustment parameters in basic mode

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

The parameters in basic mode describe the dynamic behaviour of the prosthesis in a normal gait cycle. These parameters act as basic settings for automatically adjusting the damping behaviour to the current motion situation (e.g. ramps, slow walking speed, etc.).

The stance function, sitting function and/or stairs and obstacles function can also be activated/deactivated. Further information on the stance function (see Page 80), sitting function (see Page 81), stairs and obstacles function (see Page 82).

The following parameters can be modified:

Parameter	Adjustment software range	Setting range, app	Meaning
Resistance	120 – 180	+/- 10	Resistance against flexion motion e.g. when walking down stairs or when sitting down
Angle	55° – 70°	+/- 3°	Maximum flexion angle in the swing phase
Stance function	Deactivated Activated	0 - deactivated 1 - activated	Activation/deactivation of the stance function. For switching with the Cockpit app, this function has to be activated in the adjustment software. Further information (see Page 80).
Sitting function	Deactivated Activated	0 - deactivated 1 - activated	Activation/deactivation of the sitting function. For switching with the Cockpit app, this function has to be activated in the adjustment software. Further information (see Page 81).
Stair Function	Deactivated Activated	0 - deactivated 1 - activated	Activating/deactivating the stair and obstacle function. For switching with the Cockpit app, this function has to be activated in the adjustment software. Further information (see Page 82).
Pitch	1000 Hz – 4000 Hz	1000 Hz – 4000 Hz	Pitch of beep signal for confirmation tones
Volume	0 – 4	0 – 4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals in case of errors are still generated.

9.2.3 Overview of adjustment parameters in MyModes

⚠ CAUTION

Improper use of the setting parameters in the MyModes

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Get the O&P professional and/or therapist to instruct you regarding the functionality and adjustment options for **all parameters** of the MyModes.

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

The parameters in the MyModes describe the static behaviour of the prosthesis for a specific motion pattern such as cross-country skiing. Damping behaviour is not automatically controlled and adjusted in MyModes.

Parameter	Adjustment software range	Setting range, app	Meaning
Basic flex.	0–200	+/- 20	Level of flexion resistance at the start of flexing the knee joint
Gain	0–100	+/- 10	Increase in flexion resistance (starting with the parameter " Basic flex. ") when flexing the knee joint. The knee joint locks at a certain flexion angle, which depends on the setting of the parameters " Basic flex. " and " Gain ".
Basic ext.	0–60	+/- 20	Level of extension resistance
Locking angle	0–90	+/- 10	Angle up to which the knee joint can be extended. Information: If this parameter is >0, the knee joint is locked in a flexed position in the extension direction. To unlock it, unload the prosthesis and tilt it back for at least 1.5 seconds. This makes extension of the joint possible independently of the settings for the parameters " Basic ext. " and " Locking angle ". This may be necessary to switch to basic mode using a movement pattern.
Pitch	1000 Hz–4000 Hz	1000 Hz–4000 Hz	Pitch of beep signal for confirmation tones
Volume	0–4	0–4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals in case of errors are still generated.

9.3 Turning Bluetooth on the prosthesis on/off

INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (function only available in basic mode) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 86).

INFORMATION

Basic mode (mode 1) has to be active to turn off Bluetooth. If a MyMode is activated, one has to switch to basic mode to turn off Bluetooth.

9.3.1 Switching Bluetooth off/on using the cockpit app

Switching off Bluetooth

- 1) When the component is connected, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".
- 3) Tap the entry "**Deactivate Bluetooth**".
- 4) Follow the on-screen instructions.

Switching on Bluetooth

- 1) Turn the component over or connect/disconnect the battery charger.
→ Bluetooth is switched on for approx. 2 minutes. During this time, the app must be started in order to establish a connection to the component.
- 2) Follow the on-screen instructions.
→ If Bluetooth is switched on, the (⊕) symbol appears on the screen.

9.4 Querying the prosthesis status

9.4.1 Query status through cockpit app

- 1) When the component is connected, tap the ☰ symbol in the main menu.
- 2) In the navigation menu, tap the entry "**Status**".

9.4.2 Status display in the cockpit app

Menu option	Description	Possible actions
Trip	Daily step counter (steps with the prosthesis side)	Reset the counter by tapping the " Reset " button.
Step	Total step counter (steps with the prosthesis side)	For informational purposes only
Service	Display of the next maintenance date	For informational purposes only
Batt.	Current prosthesis charge level, as a percentage	For informational purposes only
Stb/Act: 58/29	Estimated remaining operating time of prosthesis in hours. Rest mode (Stb) e.g. 58 hours, active use (Act) e.g. 29 hours	For informational purposes only

9.5 Mute mode (silent mode)

Activating mute mode (silent mode) turns off the audible feedback signals and the vibration signals. However, warnings in case of component errors are still generated (see Page 99).

Mute mode can be activated/deactivated using the Cockpit app.

INFORMATION

Connecting the battery charger automatically deactivates mute mode again.

9.5.1 Turning mute mode on/off using the Cockpit app

- 1) When the component is connected, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".
- 3) Tap the entry "**Mute mode**".
- 4) Follow the on-screen instructions.

9.6 OPG function (Optimised Physiological Gait)

INFORMATION

The O&P professional can turn the "PreFlex" function on or off using the adjustment software. All other parameters of the OPG function are always active and cannot be influenced.

The OPG function minimises the prosthesis wearer's prosthetic deviations from a harmonious gait pattern and promotes more biomechanically correct walking. This function enables the following features:

PreFlex

PreFlex ensures the knee is in 4° of flexion at the end of swing phase in preparation for loading response. This makes initiating stance phase flexion easier and forward movement is less restricted.

Adaptive yielding control

The knee joint has auto-adaptive stance and swing extension resistance. The stance flexion resistance experienced by the user is dependent on the slope or incline when walking downhill. When walking on a ramp, adaptive yielding control manages flexion depending on the angle of the ramp. The knee joint flexes slowly if the ramp is flat, and flexes quickly if the ramp is steep.

Dynamic stability control (DSC)

DSC ensures the knee will not release stance resistance during biomechanically unstable static and dynamic conditions. Constantly checking multiple parameters, DSC ensures the optimally timed decision for the knee to safely switch from stance to swing. Because DSC is always monitoring knee function, multi-directional movement and walking backward are also possible without risk of stance resistance releasing.





Adaptive swing phase control

Instantaneous adaptation to varied walking cadences and to changes of the pendular mass (e.g. varying footwear) ensures the knee always achieves the swing flexion target angle within (+/-) one degree. The swing phase extension and flexion resistance experienced by the user are auto-adaptive.

The flexed and partially loaded knee will also disable the stance phase on slopes and ramps to allow for greater knee flexion and more ground clearance in the swing phase.

10 MyModes

With the help of adjustment software, the O&P professional can activate and configure up to 5 MyModes in addition to the basic mode. These can be selected by using the Cockpit app. Only the first 3 MyModes can be selected using movement patterns. Switching by using movement patterns has to be activated in the adjustment software by the O&P professional.


These modes are intended for specific motion patterns or postures (e.g. inline skating, running (jogging) ...). Settings can be adjusted using the Cockpit app (see Page 85).

Only when "DUAL" marking is not present (see Page 60)

For information on switching MyModes with the remote control (optional accessory), see the remote control instructions for use.

10.1 Running function as configured MyMode

For running over longer periods of time, the O&P professional can configure a "Running" MyMode which can be activated using the cockpit app or a movement pattern.

In this mode, every step will be performed as a running step with larger swing phase angle and no preflexion at heel strike (PreFlex) (see Page 88).

INFORMATION

The running function will work with specialised running feet such as the 1E95 Challenger as well as with prosthetic feet with axial compression such as the 1C61 Triton Vertical Shock. For details on assembly and alignment please refer to the instructions for use of the foot. Feet without axial compression are generally not suited for running.

10.2 Switching MyModes with the cockpit app

INFORMATION


Bluetooth on the prosthesis must be turned on in order to use the Cockpit app.

If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (function only available in basic mode) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 86).

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app (see Page 83) or if mute mode (silent mode) is activated, there are no beep signals.

Once a connection to a prosthesis has been established, the cockpit app can be used to switch between the MyModes.

- 1) Tap the symbol of the MyMode (1) you want in the main menu of the app.
→ A security question for changing the MyMode appears.
- 2) If you want to change the mode, tap the "OK" button.
→ A beep signal sounds to confirm the switch.
- 3) After switching, a symbol (2) is displayed to identify the active mode.
→ The current mode is also indicated by the name on the lower edge of the screen (3).

10.3 Switching MyModes using motion patterns

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app (see Page 83) or if mute mode (silent mode) is activated, there are no beep signals.

Information on switching


- Switching and the number of motion patterns must be activated by the prosthetist in the adjustment software.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Switching using movement patterns has to be enabled by the prosthetist.
- Position the prosthetic leg back slightly and bounce on the forefoot with the leg extended while maintaining constant contact with the floor.
- Weight must be placed on the forefoot while bouncing.
- The load may not be taken off the foot fully when relieving the load.

Switching process

- 1) Position the prosthetic leg back slightly.
- 2) While maintaining constant contact with the floor, bounce on the forefoot with the leg extended a number of times in one second depending on the desired MyMode (MyMode 1 = 3 times, MyMode 2 = 4 times).
- 3) Keep the prosthetic leg still in this position with no load.

→ A beep and vibration signal will sound to confirm that the movement pattern has been recognised.

INFORMATION: If this beep and vibration signal does not sound, the requirements were not observed when bouncing the foot or mute mode (silent mode) is activated. For more information about mute mode, see the section "Mute mode" (see Page 87).

- 4) After the beep and vibration signal sounds, keep the prosthetic leg extended and still for 1 second.

→ A confirmation signal will sound to indicate that the prosthesis has successfully switched to the corresponding MyMode (2 times = MyMode 1, 3 times = MyMode 2).

INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly kept still or mute mode (silent mode) is activated. Repeat the process to correctly switch to the required mode. For more information about mute mode, see the section "Mute mode" (see Page 87).

10.4 Switching from a MyMode back to basic mode

Information on switching


- Regardless of the configuration of additional MyModes in the adjustment software, it is always possible to switch back to basic mode (mode 1) with a motion pattern.
- It is always possible to switch back to basic mode (mode 1) by connecting/disconnecting the battery charger.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Position the prosthetic leg back slightly and bounce on the forefoot with the leg extended while maintaining constant contact with the floor.
- Weight must be placed on the forefoot while bouncing.
- The load may not be taken off the foot fully when relieving the load.

Switching process

- 1) Position the prosthetic leg back slightly.
- 2) While maintaining constant contact with the floor and with the leg extended, bounce on the forefoot at 3 times or more.
- 3) Keep the prosthetic leg still in this position with no load.
→ A beep and vibration signal will sound to confirm that the movement pattern has been recognised.

INFORMATION: If this beep and vibration signal does not sound, the requirements were not observed when bouncing the foot or mute mode (silent mode) is activated. For more information about mute mode, see the section "Mute mode" (see Page 87).

- 4) After the beep and vibration signal sounds, keep the prosthetic leg extended and still for 1 second.
→ A confirmation signal will sound to indicate that the prosthesis has successfully switched over to basic mode.

INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly kept still or mute mode (silent mode) is activated. Repeat the process to correctly switch to the required mode. For more information about mute mode, see the section "Mute mode" (see Page 87).

11 Additional operating states (modes)

11.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level is 5% (see Page 99). During this time, damping settings are set to their safety mode values. This may be low or high depending on the setting in the adjustment software. The prosthesis is then switched off. You can switch back to basic mode (mode 1) from empty battery mode by charging the product.

11.2 Mode for charging the prosthesis

The product is non-functional during charging.

The product is set to the flexion resistance of safety mode. This may be low or high depending on the setting configured by the prosthetist.

11.3 Safety mode

The product automatically switches to safety mode if a critical system fault occurs (e.g. failure of a sensor signal). Safety mode remains in effect until the error has been rectified.

The switch to safety mode is indicated by beeps and vibration signals immediately prior to switching (see Page 99).

Safety mode can be disabled by connecting and disconnecting the battery charger. If the product switches into safety mode again, this means a permanent error exists. The product must be inspected by an authorised Ottobock Service Centre.

Depending on the type and severity of the error, different remaining functionality is offered in safety mode. This makes limited walking possible for the user depending on the type of error.

The following remaining functionality is available:

- **Minor error:** Continuous stance phase flexion resistance is set, with the ability to initiate the swing phase.
- **Moderate error:** Continuous stance phase flexion resistance is set, with the ability to initiate the swing phase. The swing phase control and stance phase extension resistance may be available or not depending on the type of error.
- Safety mode flexion resistance is configured. This may be low or high depending on the setting configured by the prosthetist.

The following functions are deactivated in safety mode:

- OPG function
- Stairs and obstacles function
- Stance function
- Sitting function

11.4 Overheating mode

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

When the hydraulic unit overheats due to uninterrupted, increased activity (e.g. extended walking downhill), damping is increased along with the rising temperature in order to counteract the overheating. When the hydraulic unit cools down, the product switches back to the damping settings that existed before the overheating mode.

Overheating mode is not activated in the MyModes.

Overheating mode is indicated by a long vibration every 5 seconds.

The following functions are deactivated in overheating mode:

- Sitting function
- Display of the battery charge level without additional equipment
- Switching to a MyMode

12 Storage and bleeding

Air may accumulate in the hydraulic unit if the product is stored for longer periods and not in an upright position. This is noticeable through sounds and irregular damping behaviour.

The automatic bleeding mechanism ensures that all functions of the product are again intact after approximately 10 - 20 steps.

Storage

- When storing the knee joint, press the top of the knee to the flexion stop.
- Avoid extended disuse of the product (use the product regularly).

13 Cleaning

- 1) Clean the product with a damp cloth (fresh water) when needed.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

14 Maintenance

INFORMATION

This component was tested for three million load cycles in accordance with ISO 10328. Depending on the activity level, this corresponds to a service life of three to five years. The duration of use can be individually extended depending on the intensity of use by performing regular service inspections.

Regular maintenance (service inspections) is mandatory in the interest of one's safety and in order to maintain operating reliability and protect the warranty. This maintenance includes an inspection of the sensors and replacement of wear and tear parts.

If maintenance is due, this is indicated after the service interval has expired by three brief beep and vibration signals after disconnecting the battery charger. For further information, please contact your O&P professional.

To have a service inspection carried out, please submit the product as well as the battery charger and power supply unit to the prosthetist.

15 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

15.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

15.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

15.3 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

This product meets the requirements of the European Directive 1999/5/EC for radio equipment and telecommunications terminal equipment. The conformity assessment was drawn up by the manufacturer in accordance with Annex IV of the directive.

15.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
3820 West Great Lakes Drive
Salt Lake City, Utah 84120-7205 USA
Phone + 1-801-956-2400
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website
<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.

5470 Harvester Road

L7L 5N5 Burlington, Ontario

Canada

Phone + 1-800-665-3327

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

16 Technical data

Ambient conditions	
Storage and transport in original packaging	-10°C/+14°F to +60°C/+140°F
Storage and transport without packaging	-10°C/+14°F to +60°C/+140°F Max. 80% relative humidity, non-condensing
Operation	-10°C/+14°F to +60°C/+140°F Max. 80% relative humidity, non-condensing
Charging the battery	0°C/+32°F to +45°C/+113°F

Product	
Reference number	3B5-X2/3B5-X2=ST
Mobility grade according to MOBIS	3 and 4
Maximum body weight	125 kg
Water resistance	Not waterproof, not corrosion-resistant
Range of Bluetooth connection to mobile device	max. 10 m
Weight of the prosthesis without tube adapter and Protective Cover	approx. 1520 g

Tube adapter				
Reference number	2R20	2R21 (with torsion unit)		
Weight	190–300 g/0.42–0.66 lbs	435–545 g/0.96–1.20 lbs		
Material	Aluminium			
Max. body weight	150kg	125kg		
Protection rating	IP67	IP54		
Approved set screws				
Length	10 mm	12 mm	14 mm	16 mm
Reference number	506G3= M8x10	506G3= M8x12	506G3= M8x14	506G3= M8x16

Prosthesis battery	
Battery type	Li-Ion
Charging cycles (charging and discharging cycles) after which at least 80% of the original battery capacity remains available	500
Charge level after 1 hour charging time	30 %
Charge level after 2 hours charging time	50 %
Charge level after 4 hours charging time	80 %
Charge level after 8 hours charging time	Fully charged
Product behaviour during the charging process	The product is non-functional
Operating time of the prosthesis with new, fully charged battery at room temperature	Approx. 5 days with average use

Power supply	
Reference number	757L16*
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10 % to 93 % relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F max. 90 % relative humidity, non-condensing
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ===

Battery charger	
Reference number	4E60*
Storage and transport in original packaging	-25 °C to 70 °C/-13 °F to 158 °F
Storage and transport without packaging	-25 °C to 70 °C/-13 °F to 158 °F Max. 93% relative humidity, non-condensing
Operation	5 °C to 40 °C/41 °F to 104 °F Max. 93% relative humidity, non-condensing
Protection rating	IP40
Input voltage	12 V ===

Only when "DUAL" marking is not present (see Page 60)

Cockpit app	
Article number	Cockpit 4X441-Andr=V*
Supported operating system	Android 4.0.3 or higher
Website for download	http://www.ottobock.com/cockpitapp

Data transfer	
Wireless technology	Bluetooth 2.0
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK

Data transfer	
Data rate (over the air)	704 kbps
Maximum output power (EIRP)	-1.78 dBm (0.66 mW)

Only when "DUAL" marking is present (see Page 60)

Cockpit app	
Reference number	4X441-IOS=* / 4X441-Andr=V* Cockpit
Supported operating system	iOS 9.3 / Android 4.1 or higher
Website for download	http://www.ottobock.com/cockpitapp

Data transfer	
Wireless technology	Bluetooth Smart Ready
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK, $\pi/4$ DQPSK, 8DPSK
Data rate (over the air)	2178 kbps (asymmetrical)
Maximum output power (EIRP):	+8.5 dBm

17 Appendices

17.1 Symbols Used

Manufacturer

Compliance with the requirements according to "FCC Part 15" (USA)

Non-ionising radiation

IP40

Protection against penetration of solid foreign objects with a diameter greater than 1 mm, no protection against water

IP54

Protected against dust, protected against splashing water

IP67

Dustproof, protection against temporary submersion

DUAL

The product's Bluetooth wireless module can establish a connection to devices with the following operating systems: iOS (iPhone, iPad, iPod,...) and Android

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please

observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Serial number

Lot number

Caution, hot surface

Protect from moisture

17.2 Operating states/error signals

The prosthesis indicates operating states and error messages through beeps and vibration signals.

17.2.1 Signals for operating states

Battery charger connected/disconnected

Beep signal	Vibration signal	Event
—	3x long	Charging mode started (3 sec. after connecting the battery charger)
1x short	1x short	Self-test completed successfully, product is operational

Mode switching

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see Page 83).

Beep signal	Vibration signal	Additional action performed	Event
1x short	1x short	Mode switching using the Cockpit app	Mode switching is performed using the Cockpit app.
1x long	1x long	Bouncing on the forefoot followed by unloading the prosthetic leg	Bouncing pattern recognised.
1x short	1x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to basic mode (mode 1) carried out.
2x short	2x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to MyMode 1 (mode 2) carried out.
3x short	3x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to MyMode 2 (mode 3) carried out.

17.2.2 Warnings/error signals



Error during use



Beep signal	Vibration signal	Event	Required action
—	1 x long at intervals of approx. 5 seconds (if mute mode (silent mode) is activated, this signal is not generated)	Overheated hydraulic unit	Reduce activity.

Beep signal	Vibration signal	Event	Required action
–	3 x long	Charge level under 25%	Charge battery soon. Remaining operating time approx. 24 hours
–	5 x long	Charge level under 10%	Charge battery soon Remaining operating time approx. 6 hours
5 x long	5 x long repeated every 60 seconds	Medium error (see Page 92) e.g. a sensor is not operational	Walking possible with restrictions. Please note the change in flexion resistance. The product must be inspected immediately by an O&P professional.
10 x long	10 x long	Battery charge level 5% After the beep and vibration signals, the product switches to empty battery mode and then switches off.	Charge the battery.
30 x long	1 x long, 1 x short repeated every 3 seconds	Severe error/indication of safety mode activation (see Page 92) e.g. one or more sensors are not operational	Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohibited. The product must be inspected immediately by an O&P professional.
–	Continuous	Total failure Electronic control no longer possible. Safety mode active or undetermined valve state. Undetermined product behaviour.	Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohibited. The product must be inspected immediately by an O&P professional.

Error while charging the product

LED on power supply	Status LED on battery charger	Error	Resolution
○		Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.

LED on power supply	Status LED on battery charger	Error	Resolution
○		Non-functional socket	Check socket with another electric appliance.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●		No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.

	Status LED	Charging status indicator (5 LEDs)	Error	Resolution
	The LED ring is lit in weak purple	No LED is lit	Distance between battery charger and receiver of the charging unit on the prosthesis too great. If the distance is more than 2 mm, the prosthesis cannot be charged.	Reduce distance between battery charger and receiver of the charging unit.
	The LED ring is lit up yellow	LED 2 and 4 are lit up	Charger excess temperature	Check whether the specified ambient conditions for charging the battery are met (see Page 95).
		LED 1, 3 and 5 are lit up	Excessively high or low temperature of the prosthesis	
		LED 3 is lit up	The prosthesis is not being charged Distance between battery charger and receiver of the charging unit too great.	Connection may be improved by reducing the distance between the charger and receiver of the charging unit.
	The LED ring is lit up green		The battery charger is functional but not connected to the receiver, or the distance between the battery charger and receiver of the charging unit is too great.	Connect the battery charger or reduce the distance between the charger and the receiver of the charging unit on the prosthesis.

	Status LED	Charging status indicator (5 LEDs)	Error	Resolution
	The LED ring is flashing red		The prosthesis is not being charged Defective battery charger.	Reset the error by disconnecting and connecting the power supply. If the error persists, the battery charger and power supply must be inspected by an authorised Ottobock Service Center.


17.2.3 Error messages while establishing a connection with the cockpit app

Error message	Cause	Correction
Component was connected to another device. Establish connection?	The component was connected to another device	To disconnect the original connection, tap the "OK" button. If the original connection is not to be disconnected, tap the "Cancel" button.
Mode change failed	An attempt was made to switch to a different MyMode while the component was in motion (e. g. while walking)	For safety reasons, switching MyModes is only permitted when components are at rest, e. g. while standing or sitting.
	A current connection to the component was interrupted	Check the following points: <ul style="list-style-type: none"> • Distance from the component to the device • Charge level of the component's battery • Bluetooth of the component switched on? (Switching Bluetooth of the component on/off) • Hold the component with the sole of the foot facing up to make the component "visible" for 2 minutes. • If multiple components were stored, was the correct component selected?

17.2.4 Status signals

Battery charger connected

LED on power supply	Status LED on battery charger	Event
		Power supply and battery charger operational. Battery charger not yet connected to receiver.

LED on power supply	Status LED on battery charger	Event
		Battery charger connected to receiver, good connection. This display turns off automatically after approximately one minute to avoid bothersome light at night. The charging process is not affected.

Battery charger disconnected

Beep signal	Vibration signal	Event	Resolution
1x short	1x short	Self-test completed successfully. Product is operational.	
3x short	3x short	Maintenance note: e.g. service interval has been exceeded, temporary disruption of a sensor signal	<ul style="list-style-type: none"> • Use the Cockpit app to check the next maintenance date for the prosthesis (see Page 87). Should the date lie within the next month, arrange a maintenance appointment with the O&P professional. For this appointment, the O&P professional has to be provided with the battery charger and power supply in addition to the prosthesis with tube adapter. • Conduct the self-test again by connecting/disconnecting the battery charger. • If the beep signal sounds again and the maintenance date has not yet been reached or exceeded, visit the O&P professional soon. If necessary, the prosthetist will forward the prosthesis to an authorised Ottobock Service Center. • The product can be used without restrictions. However, vibration signals may not be generated.

Battery charge level

During the charging process, the current battery charge level is indicated by the LEDs lit on the side of the charger.

LEDs	0	1	2	3	4	5
Battery charge level	0%-10%	10%-30%	30%-50%	50%-70%	70%-90%	>90%

17.3 Directives and manufacturer's declaration

17.3.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

The customer or user of the product must ensure that it is operated in such an environment.

Observe the safety notices in the section "Information on proximity to certain areas" (see Page 68).

The following tables outline the test levels of the tests to be carried out. The higher value applies if there is a difference between the operating environments listed in the tables below.

Table 1 – Electromagnetic emissions for all devices and systems

Directives and manufacturer's declaration – electromagnetic interference		
The product 'X2 Bionic Prosthetic System 3B5-X2' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'X2 Bionic Prosthetic System 3B5-X2' must ensure that the device is operated in such an environment.		
Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1	The product 'X2 Bionic Prosthetic System 3B5-X2' uses HF energy exclusively for its internal functioning. Therefore its HF emissions are very low and interference with neighbouring electronic devices is unlikely.
HF emissions according to CISPR 11	Class B	The product 'X2 Bionic Prosthetic System 3B5-X2' is intended for use in other than residential facilities and in facilities connected directly to a public low-voltage network that also supplies residential buildings.
Harmonics according to IEC 61000-3-2	N/A	
Voltage fluctuations/flicker according to IEC 61000-3-3	N/A	

Table 2 – Electromagnetic interference immunity for all devices and systems

Directives and manufacturer's declaration – electromagnetic interference immunity			
The product 'X2 Bionic Prosthetic System 3B5-X2' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'X2 Bionic Prosthetic System 3B5-X2' must ensure that the device is operated in such an environment.			
Interference immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 6 kV discharge on contact ± 8 kV discharge through air	± 6 kV discharge on contact ± 8 kV discharge through air	Floors should be wood, concrete or ceramic tile. If the floor covering is made of synthetic material, the relative humidity must be at least 30 %.

Interference immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Electrical fast transient/bursts according to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply.
Surges according to IEC 61000-4-5	± 1 kV phase-to-phase voltage ± 2 kV even voltage	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply.
Voltage drops, short interruptions and fluctuations of the supply voltage according to IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ drop of U_T) for 1/2 period $40\% U_T$ (60% drop of U_T) for 5 periods $70\% U_T$ (30% drop of U_T) for 25 periods $< 5\% U_T$ ($> 95\%$ drop of U_T) for 5 s	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply. If the user of the product 'X2 Bionic Prosthetic System 3B5-X2' demands continuous operation even in case of an interruption in the supply voltage, we recommend supplying the product 'X2 Bionic Prosthetic System 3B5-X2' with an uninterruptible power supply or battery.
Magnetic field at mains frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should be equivalent to the typical levels for commercial or hospital environments.
Note: U_T is the rated voltage before applying the test levels.			

Table 4 – Electromagnetic interference immunity for non-life-sustaining devices and systems

Directives and manufacturer's declaration – electromagnetic interference immunity

The product 'X2 Bionic Prosthetic System 3B5-X2' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'X2 Bionic Prosthetic System 3B5-X2' must ensure that the device is operated in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Conducted HF interference according to IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	Not applicable	Portable and mobile wireless equipment should not be used at a lesser distance from the product 'X2 Bionic Prosthetic System 3B5-X2' including the lines than the recommended safety distance calculated using the equation applicable for the transmission frequency. Recommended safety distance: $d = 0.3 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ for 800 MHz to 2.5 GHz where P is the nominal output of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended safety distance in metres (m). The field strength of stationary radio transmitters at all frequencies according to an on-site investigation ^a should be less than the compliance level. Interference is possible in the vicinity of devices that bear the following symbol. 
Radiated HF interference according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	12 V/m 80 MHz to 1 GHz 3 V/m 1 GHz to 2.5 GHz	

Note 1: The higher frequency range applies at 80 MHz and 800 MHz.

Note 2: These directives may not be applicable in all cases. The propagation of electromagnetic factors is influenced by absorption and reflection by buildings, objects and people.

^a The field strength of stationary transmitters such as base stations of radio telephones and mobile land radio equipment, amateur radio stations, AM and FM radio and television stations cannot be precisely determined theoretically in advance. A site survey is recommended to establish the electromagnetic conditions as a result of stationary HF transmitters. If the field strength determined at the site of the product 'X2 Bionic Prosthetic System 3B5-X2' exceeds the compliance level specified above, the product 'X2 Bionic Prosthetic System 3B5-X2' has to be monitored in regards to normal operation at each place of use. If unusual performance characterist-

ics are noted, taking additional steps may be necessary, for example changing the orientation or moving the product 'X2 Bionic Prosthetic System 3B5-X2'

Table 6 – Recommended safe distances between portable and mobile HF telecommunication equipment and the product 'X2 Bionic Prosthetic System 3B5-X2' – for non-life-sustaining devices and systems

Directives and manufacturer's declaration – Recommended safe distances between portable and mobile HF telecommunication equipment and the product 'X2 Bionic Prosthetic System 3B5-X2'

The product 'X2 Bionic Prosthetic System 3B5-X2' is designed for operation in an electromagnetic environment where HF interference is monitored. The customer or user of the product 'X2 Bionic Prosthetic System 3B5-X2' can help prevent electromagnetic interference by complying with the minimum distances between wearable and mobile HF communication equipment (transmitters) and the product 'X2 Bionic Prosthetic System 3B5-X2', as below according to the maximum output of the communication equipment.

Nominal output of transmitter [W]	Safety distance [m] by transmission frequency		
	150 kHz to 80 MHz	80 MHz to 800 MHz $d = 0.3 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	N/A	0.03	0.23
0.1	N/A	0.09	0.73
1	N/A	0.3	2.3
10	N/A	0.95	7.3
100	N/A	3	23

For transmitters with no maximum nominal output specified in the table above, the recommended safety distance d in metres (m) can be determined using the equation in the respective column, where P stands for the maximum nominal output of the transmitter in watts (W) according to the transmitter manufacturer's information.

Note 1: The higher frequency range applies at 80 MHz and 800 MHz.

Note 2: These directives may not be applicable in all cases. The propagation of electromagnetic factors is influenced by absorption and reflection by buildings, objects and people.

The product "X2 - Bionic Prosthetic System 3B5-X2/3B5-X2=ST" is covered by the following patents:

Canada:	CA 2 651 124; CA 2 714 469; CA 2 780 192; CA 2 780 511; CA 2 779 784; CA 2 704 792; CA 2 678 987; CA 2 626 738
China:	CN 101 453 963; CN 101 909 553; CN 101 938 958; CN 102 711 672; CN 102 647 963; CN 101 346 110; CN 102 740 804; CN 101 663 002; CN 102 740 803; CN 102 724 936; CN 102 762 171; CN 107 530 173
Finland:	FI 110 159
Germany:	DE 10 2008 010 281; DE 10 2015 106 384; DE 10 2009 052 887
Japan:	JP 4 718 635; JP 5 619 910; JP 5 547 091; JP 5 394 579; JP 5 678 079; JP 5 968 591
Russia:	RU 2 404 730; RU 2 484 789; RU 2 533 967; RU 2 488 367; RU 2 508 078; RU 2 572 741
South Korea:	KR 101 286 894; KR 101 509 265
Taiwan:	R.O.C. Invention Patent No. . I386194; I459936; I442912; I517845; I551277; I551278; I530278; I542335; I519292; I494095
USA:	US 7 731 759; US 6 908 488; US 8 083 807, US 8 474 329; US 8 876 912; US 8 814 948; US 9 066 818; US 9 572 690; US 9 278 013; US 9 877 849
European Patent	EP 1237513 in DE, FR, GB EP 1940327 in DE, FR, GB, IS, IT, NL, SE, TR EP 2015712 in DE, ES, FR, GB, IT, NL, SE, TR EP 2129340 in DE, FR, GB, IS, IT, NL, SE, TR EP 2222253 in DE, FR, GB, IS, IT, NL, SE, TR EP 2240124 in DE, FR, GB, IT, NL, SE, TR EP 2254525 in DE, FR, GB, IS, IT, NL, TR EP 2498724 in DE, FR, GB, IS, IT, NL, SE, TR EP 2498725 in DE, FR, GB EP 2498726 in DE, FR, GB, IS, IT, NL, SE, TR EP 2498727 in DE, FR, GB, IS, IT, NL, SE, TR EP 2498728 in DE, FR, GB EP 2498730 in DE, FR, GB EP 2772232 in DE, FR, GB, IS, IT, NL, SE, TR

Patents pending in Canada, China, EPA, Germany, Japan, Russia and USA.

Otto Bock Healthcare Products GmbH
Kaiserstraße 39 · 1070 Wien · Austria
T +43 1 523 37 86 · F +43 1 523 22 64
info.austria@ottobock.com · www.ottobock.com