## 3F1=1/3F1=2


<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	9
<b>EN</b> Instructions for use (qualified personnel) .....	29
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	48
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	69
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	89
<b>PT</b> Manual de utilização (Pessoal técnico) .....	110
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	130
<b>SV</b> Bruksanvisning (Fackpersonal) .....	150
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	170
<b>NO</b> Bruksanvisning (Fagpersonell) .....	190
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	210
<b>PL</b> Instrukcja użytkowania (Personel fachowy) .....	228
<b>HU</b> Használati utasítás (szakszemélyzet) .....	249
<b>CS</b> Návod k použití (Odborný personál) .....	269
<b>HR</b> Upute za uporabu (Stručno osoblje) .....	289
<b>SL</b> Navodila za uporabo (Strokovno osebje) .....	307
<b>SK</b> Návod na používanie (Odborný personál) .....	326
<b>TR</b> Kullanma talimatı (Uzman personel) .....	345
<b>RU</b> Руководство по применению (Квалифицированный персонал) .....	364
<b>JA</b> 取扱説明書 (有資格担当者) .....	385
<b>ZH</b> 使用说明书 (专业人员) .....	404

# Quick Reference Guide

de | Dieser "Quick Reference Guide" ersetzt nicht die Gebrauchsanweisung  
 en | This "Quick Reference Guide" does not replace the instructions for use  
 fr | Ce « Quick Reference Guide » ne remplace pas les instructions d'utilisation  
 it | Questa "Quick Reference Guide" non sostituisce le istruzioni per l'uso  
 es | Esta "Quick Reference Guide" no sustituye a las instrucciones de uso  
 pt | Este "Quick Reference Guide" não substitui o manual de utilização  
 nl | Deze "Quick Reference Guide" is geen vervanging voor de gebruiksaanwijzing  
 sv | Denna 'Quick Reference Guide' ersätter inte bruksanvisningen  
 da | Denne "Quick Reference Guide" erstatter ikke brugsanvisningen  
 no | Denne "Quick Reference Guide" erstatter ikke bruksanvisningen  
 fi | Tämä tuotteen mukana toimitettu "Quick Reference Guide" ei korvaa käyttöohjetta  
 pl | Niniejszy "Quick Reference Guide" nie zastępuje instrukcji użytkowania  
 hu | Ez "Quick Start Guide" nem helyettesíti a használati utasítást  
 cs | Tento "Rychlý návod k použití" nenahrazuje Návod k použití  
 hr | Ovaj letak „Quick Reference Guide“ ne zamjenjuje upute za uporabu  
 sk | Táto "Quick Reference Guide" ne nadomešča navodila za uporabo  
 sl | Tento „Quick Reference Guide“ nenahrádza návod na používanie  
 tr | Bu "Hızlı başvuru kılavuzu", kullanım kılavuzunun yerini almaz  
 ru | Это руководство 'Quick Reference Guide' не заменяет руководство по применению  
 ja | 本 "Quick Reference Guide" (クイックガイド) は取扱説明書の代わりとなるものではありません  
 zh | 该 "Quick Reference Guide" (快速参考指南) 并不能取代使用说明书


**99B120=\***


②	①	<b>40-46 cm</b>	<b>46-56 cm</b>
<b>38-48 cm</b>		S	L <sup>1</sup>
<b>&gt; 48 cm</b>		L	L

\* siehe Kapitel 2.2 / see chapter 2.2 / voir chapitre 2.2 / vedere capitolo 2.2 / véase el capítulo 2.2 / consulte o capítulo 2.2 / zie hoofdstuk 2.2 / se kapitel 2.2 / se kapitel 2.2 / se kapitell 2.2 / katso luku 2.2 / patrz rozdział 2.2 / lásd a 2.2 fejezetet / viz kapitola 2.2 / vidi poglavlje 2.2 / glejte poglavje 2.2 / pozri kapitulu 2.2 / bkz. Bölüm 2.2 / см. раздел 2.2 / 第2.2章を参照してください。 / 参见章节 2.2


 6.2


 6.3


 6.4

 6.5

## 6.6

**Angewandte Symbole / Symbole Used / Symboles utilisés / Simboli utilizzati / Símbolos utilizados / Símbolos utilizados / Gebruikte symbolen / Symboler som används / Anvendte symboler / Benyttede symboler / Käytetyt symbolit / Stosowane symbole / Alkalmazott jelképek / Použité symboly / Rabljeni simboli / Uporabljeni simboli / Použité symboly / Kullananlilan semboller / Применяемые символы / 本取扱説明書で使用している記号 / 使用的图标**

de | Kapitel der Gebrauchsanweisung lesen  
 en | Read the section in the instructions for use  
 fr | Veuillez lire le chapitre des instructions d'utilisation  
 it | Leggere il relativo capitolo nelle istruzioni per l'uso  
 es | Leer el capítulo de las instrucciones de uso  
 pt | Ler o capítulo do manual de utilização  
 nl | Lees het hoofdstuk/de hoofdstukken van de gebruiksaanwijzing  
 sv | Läs kapitlet i bruksanvisningen  
 da | Læs kapitlet i brugsanvisningen  
 no | Les kapitlet i bruksanvisningen  
 fi | Lue käyttöohjeen luku  
 pl | Przeczytać rozdział instrukcji użytkowania  
 hu | Olvassa el a használati utasítás fejezetét  
 cs | Přečtěte si příslušnou kapitulu v Návodu k použití  
 hr | Pročitajte poglavlje uputa za uporabu  
 sl | Preberite poglavje navodil za uporabo  
 sk | Prečítajte si kapitolu návodu na používanie  
 tr | Kullanim kılavuzunun bölümlerini okunmalıdır  
 ru | Читать соответствующую главу руководства по применению  
 ja | 取扱説明書のセクションを参照してください  
 zh | 请阅读使用说明书章节

de   Richtig	sv   Rätt	hr   Pravilno
en   Right	da   Rigtigt	sl   prav
fr   Correct	no   Riktig	sk   Správne
it   Corretto	fi   Oikein	tr   Doğru
es   Correcto	pl   Prawidłowo	ru   Правильно
pt   Correto	hu   Helyes	ja   正しい
nl   Goed	cs   Správně	zh   正确

de		Falsch	sv		Fel	hr		Nepravilno
en		Wrong	da		Forkert	sl		false
fr		Incorrect	no		Feil	sk		Nesprávne
it		Errato	fi		Väärin	tr		Yanlış
es		Incorrecto	pl		Nieprawidłowo	ru		Неправильно
pt		Incorreto	hu		Téves	ja		誤り
nl		Fout	cs		Špatně	zh		错误

de		Schraubendreher mit Fähnchengriff
en		Screwdriver with wing handle
fr		tournevis avec manche drapeau
it		cacciavite con impugnatura a bandiera
es		destornillador con mango de bandera
pt		chave de fenda com cabo tipo bandeira
nl		schroevendraaier met vlaggetjesgreep
sv		skruvmejsel med flaggrepp
da		skruetrækker med håndtag
no		skrutrekker med T-håndtak
fi		ruuviavain siipikahvalla
pl		wkrętak z uchwytem flagowym
hu		hasábmarkolatú csavarhúzó
cs		Šroubovák s praporkovou rukojetí
hr		odvijača s drškom u obliku zastavice
sl		kompleta izvijačev z ročajem z zastavico
sk		skrutkovač s práporkovou rukoväťou
tr		tutamaklı tornavida
ru		отвертка с флагообразной ручкой
ja		ウィングハンドル付きドライバー
zh		旗状手柄螺丝刀

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### INFORMATION

Date of last update: 2019-02-05

- ▶ Please read this document carefully before using the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The product "3F1=1 functional cosmesis C-Leg / 3F1=2 functional cosmesis Genium" is called the product/functional cosmesis in the following.

These instructions for use provide you with important information on the use, adaptation and handling of the product.


Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

The functional cosmesis consists of the following components (example 3F1=1):

1. Functional knee part (referred to as "functional knee part" in the following)
2. Functional shank (referred to as "functional shank" in the following)
3. Functional stocking (not included in the scope of delivery)
4. Sheath (with distal stitched edge)

## 2.2 Function

The product helps create the most natural possible appearance. The product forms the exterior finish of a prosthesis. The product is not required for the proper and safe functioning of the prosthesis.

The functional shank is shaped by the O&P professional to create the most natural possible appearance. Streaks, colour differences and small air pockets (up to 1 cm) in the functional shank are concealed by the functional stocking. These differences, which are due to production reasons, do not affect the function of the product and therefore do not constitute grounds for complaints.

The functional stocking finishes the outside of the functional cosmesis. Various colours and sizes are available.

### INFORMATION

Size S of the functional stocking cannot be used for all users due to the body height (see Page 2, dimension points 1 and 2). Additional socket volume has to be created for size L<sup>1</sup> so the adhesion of the functional stocking on the socket can be ensured. Fabricating an individual thigh part with a functional cosmesis for the user is generally recommended in many cases (note the 647G1294 additional documentation).

The sheath improves the sliding capability of the functional stocking.

## 2.3 Combination possibilities

This product can be combined with the following Ottobock components:

### Knee Joint Combination Possibilities

- 3C98-3, 3C88-3 C-Leg
- Genium: 3B1-2, 3B1-2=ST, 3B1-3, 3B1-3=ST

### Prosthetic Foot Combination Possibilities

- 1C60 Triton
- 1C61 Triton Vertical Shock
- 1C63 Triton Low Profile
- 1C64 Triton Heavy Duty

- 1C30 Trias
- 1D10 Dynamic foot
- 1D11 Dynamic foot
- 1M10 Adjust
- Terion K2: 1C11
- 1A30 Greissinger plus
- 1D35 Dynamic Motion
- 1C40 C-Walk
- 1E56 Axtion
- 1E57 Lo Rider
- 1B1 Meridium
- Meridium: 1B1-2
- 1C66 Triton smart ankle
- 1C10 Terion

**Adapter combination possibilities**

- 4R57, 4R57=ST rotation adapter

**Limited compatibility:**

C-Leg: Collision with functional knee part from a tilt angle of 4.3° in the direction of the 4R57 and 4R57=ST rotation adapter button

Genium: Collision with functional knee part from a tilt angle of 3.2° in the direction of the 4R57 und 4R57=ST rotation adapter button

### 3 Application

#### 3.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

The product serves to finish the prosthesis in the sense of restoring a natural appearance that is as inconspicuous as possible.

#### 3.2 Area of application

The product is intended exclusively for wearers of the 3C98-3/3C88-3 C-Leg, 3B1-2/3B1-2=ST Genium and 3B1-3/3B1-3=ST Genium.

The product can be used with a knee axis-floor measurement of up to 560 mm. The minimum knee axis-floor measurement that can be fitted depends on the foot that is used with the tube adapter at maximum insertion.

#### 3.3 Conditions of use

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

The product was developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

#### 3.4 Indications

- For patients with knee disarticulation or transfemoral amputation
- The patient must be able to understand usage and safety messages and put them into practice.

#### 3.5 Qualification

The fitting of patients with the product may only be carried out by qualified personnel (e.g. O&P professional).

### 4 Safety

#### 4.1 Explanation of warning symbols

 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

## 4.2 Structure of the safety instructions

### CAUTION

#### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

## 4.3 General safety instructions

### CAUTION

#### **Non-observance of safety instructions**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety instructions and the stated precautions in this accompanying document.

### CAUTION

#### **Incorrect protective equipment**

Injury due to defective or lack of protective equipment.

- ▶ Wear protective goggles and a dust mask during grinding.

### CAUTION

#### **Use of unsuitable prosthesis components**

Falling due to unexpected behaviour of the product or breakage of load-bearing components.

- ▶ Use the product only in combination with components listed in the section "Combination possibilities" (see Page 30).

### CAUTION

#### **Contact of the product with embers or fire**

Injuries (e.g. burns) and damage to the product.

- ▶ Keep the product away from embers and open flames.

### CAUTION

#### **Wearing the product without a functional stocking**

- > Injuries due to pinching of body parts.
- > Damage due to pinching of garments.
- > Damage to the product.
- > Loss of functionality due to the penetration of moisture into the foam.

- ▶ The product may only be worn with the functional stocking.

### CAUTION

#### **Improper product care**

Injury due to the use of improper cleaning agents.

- ▶ Clean the product only as described in the section "Cleaning and care" (see Page 46).

**⚠ CAUTION****Skin irritation due to irregular care of the product**

> Skin irritation due to contact with a soiled functional shank.

> Skin irritation due to contact with a soiled sheath.

- ▶ Clean the product only as described in the section "Cleaning and care" (see Page 46).
- ▶ The product may only be worn with the functional stocking.

**⚠ CAUTION****Failure to observe warnings/error signals**

Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ Observe the warnings/error signals and corresponding change in damping settings.

**⚠ CAUTION****Impaired perception of the feedback signals**

Falling due to impaired perception of the feedback signals resulting from changed damping behaviour.

- ▶ The audibility of the feedback signals must be checked with the user and the corresponding volume has to be set.

**⚠ CAUTION****Thermal overload due to increased heat accumulation**

Falling due to altered damping behaviour as a result of earlier switching to overheating mode.

- ▶ Observe the warnings/error signals of the prosthesis (overheating mode).

**NOTICE****Independent user changes or modifications made to the product**

Damage to the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.

**NOTICE****Mechanical damage to the product**

Change in or loss of functionality due to damage.

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**INFORMATION**

The penetration of dirt and humidity into the product leads to a loss of comfort. Note the instructions for use of the knee joint (647G890/647G1148).

**INFORMATION**

The product can limit the range of the Bluetooth signal. Note the instructions for use of the knee joint (647G890/647G1148).

## 5 Scope of Delivery and Accessories

### 5.1 Scope of delivery

- 1 pc. Functional shank
- 1 pc. Functional knee part
- 1 pc. Sheath
- 1 pc. Installation ring for 4X259 inductive charger (only for 3F1=2)
- 1 pc. Screwdriver with wing handle
- 4 pc. Cable ties
- 4 pc. Blind plug (plastic set screw)
- 1 pc. 647G1290 instructions for use (user)
- 1 pc. 647G1289 instructions for use (qualified personnel)
- 1 pc. 647G1294 additional documentation

### 5.2 Accessories

The following components are not included in the scope of delivery and must be ordered separately:

- 1 pc. 99B120=\* functional stocking

#### Accessories for the 3F1=1 functional cosmesis C-Leg

- 1 pc. 4X157 charger extension cable – knee

## 6 Preparation for use

### 6.1 Required tools

- Adhesive tape
- 636N9 contact adhesive
- Measurement equipment: folding rule
- 710D20 torque wrench
- Allen key size 2
- Allen key size 4
- 4H105 knee extender (only for C-Leg)  
Read the 647G1375 instructions for use.
- Double-sided adhesive tape
- Protective goggles
- Scissors
- Scriber
- 4X258 installation tool for inductive charger
- Dust mask
- Pedilen or Plastazote (thickness 5 mm)

### 6.2 Check for collision and adapt functional knee part

#### CAUTION

##### Incorrect protective equipment


Injury due to defective or lack of protective equipment.

- ▶ Wear protective goggles and a dust mask during grinding.

Complete the following steps to put on the functional knee part:


- 1) Loosen the screws on the functional knee part.
- 2) Slightly expand the functional knee part and pull it over the frame of the joint.



- 3) Secure the position of the upper knee part with the Allen key size 2.  
**NOTICE! This secures the upper knee part in the required position.**

- 4) Fully extend the knee joint with the assembled socket.  
**INFORMATION: The 4H105 knee extender has to be used with the C-Leg.**
- 5) Check the upper knee part and the ball top for collisions with the socket.
- 6) Mark the collision points on the upper knee part and ball top.

- 7) If necessary, separate the upper knee part from the lower knee part for processing. Push the two blind plugs out of the clips from the outside to the inside using an Allen key size 2.

8) Press the clips out of the functional knee part using an Allen key size 4 and separate the upper knee part from the lower knee part.

9) Sand the upper knee part and ball top in the required area (only permitted in areas marked grey in the illustration).


**INFORMATION: Sanding the white area is not permitted.**

10) Check the upper knee part and ball top for collisions with the socket again. Note steps 2-4!


→ If the socket and upper knee part continue to collide even though the permitted area has been entirely sanded, the upper knee part has to be removed together with the clips and axles.

→ If securing the upper knee part using the Allen key size 2 is no longer possible after processing the functional knee part, the correct position of the upper knee part has to be checked with the help of the functional shank. To do so, insert the joint into the functional shank. After processing has been completed, the upper knee part has to be flush with the edge of the functional shank without contacting the extended socket.

→ If the ball top and socket still collide, the ball top can be removed. (see fig. 7).


- 11) Carefully deburr the sanded areas.
- 12) Reassemble the functional knee part.  
**INFORMATION: Ensure that the axles and clips engage correctly.**

- 13) Reinsert the two blind plugs flush into the clips from the inside.

### 6.3 Only for C-Leg: install charger extension cable

Complete the following steps to install the charger extension cable for the C-Leg (4X157):


- 1) Install the charger extension cable (see the 647G1146 instructions for use).
- 2) Assemble the charger extension cable and O-rings so that the ribs of the functional knee part do not press on the cables and O-ring.

**INFORMATION: The functional knee part has corresponding recesses to accept the O-ring and charging cable.**

### 6.4 Assembling the functional knee part

Complete the following steps to assemble the functional knee part:

- 1) Check the alignment of the actuator stirrup (the screw was already loosened while checking for collisions).

**INFORMATION: The stirrup points up, the wing points down.**

- 2) Assemble the functional knee part on the joint.

**INFORMATION: Verify the form fit of the functional knee part on the joint.**

- 3) Secure the position of the functional knee part with the Allen key size 2.

→ If securing with the Allen key is no longer possible due to processing, the correct position of the upper knee part has to be verified with the help of the functional shank. To do so, insert the joint into the functional shank. The upper knee part has to be flush with the edge of the functional shank.

- 4) Bring the joint to full extension.

**INFORMATION: The 4H105 knee extender has to be used with the C-Leg.**


- 5) Hand-tighten the cable tie.

→ This causes the actuator stirrup to lie against the installed adapter.

→ The upper knee part is held in its position by the Allen key.

- 6) Verify the proper fit of the cable tie.

- 7) Hand-tighten the two screws on the functional knee part using two fingers.


- 8) Check the position of the end of the actuator stirrup.

**INFORMATION: The end of the actuator stirrup has to be flush with the face end of the circular recess in the functional knee part.**

- 9) Remove the Allen key and knee extender (C-Leg).

- 10) Trim the cable tie.

→ When flexing the knee joint, the upper knee part moves noiselessly with the socket.

## 6.5 Finishing the functional shank

### **⚠ CAUTION**


#### **Incorrect protective equipment**

Injury due to defective or lack of protective equipment.

- ▶ Wear protective goggles and a dust mask during grinding.

#### **Shortening the functional shank**

Complete the following steps to shorten the functional shank:

- 1) Bring the prosthesis to full extension.
- 2) Measure from the bottom edge of the upper knee part to the top edge of the foot connection plate.

- 3) Add approx. 15 mm to the measured length and transfer this to the functional shank.


**INFORMATION: On the functional shank, measure from the top to the bottom.**

- 4) Cut off the functional shank at the marked location.

**INFORMATION: Note the foot position for an optimum fitting result.**

#### **Adaptation for fitting users with feet with a high structural height**

For a fitting with a foot with a high structural height, the inner diameter on the lower side of the functional shank has to be adapted accordingly.

The following steps have to be completed to adapt the functional shank with the 1C61 Triton Vertical Shock:

- 1) Enlarge the opening on the lower side of the functional shank.
- 2) Sand the diameter in the foam to 60 mm.
- 3) Foam material (e.g. Pedilen) can be applied in the ankle area of the 1C61 Triton Vertical Shock (below the pyramid adapter and functional ring) to ensure a better hold of the functional shank in the ankle area.

The following steps have to be completed to adapt the functional shank with the 1C60 Triton:


- ▶ Expose the foam on the inside rear for the carbon spring of the prosthetic foot.

#### **Adaptation for fitting users with a small ankle circumference (optional)**

Material can be applied to the tube adapter for fittings with feet that have a low structural height (such as the 1D35, Triton LP).

Complete the following steps to adapt the tube adapter:

- 1) Measure the length of the tube adapter (from the bottom edge of the joint to the start of the pyramid receiver of the tube adapter).
- 2) Cut a rectangle of the measured length and a width of 122 mm out of the material Pedilen or Plastazote (5 mm thickness).


- 3) Attach the material to the tube adapter with double-sided adhesive tape.

**INFORMATION: Do NOT glue the material to the functional shank.**

### Access to the set screws for feet with a high structural height

Complete the following steps to gain access to the set screws:

- 1) Remove the prosthetic foot.  
**INFORMATION: Using two Allen keys size 4 for the disassembly and assembly of the prosthetic foot is recommended.**
- 2) Secure the two adjacent set screws to prevent changing the foot position.
- 3) Turn the other two set screws completely into the tube adapter.
- 4) Insert the knee joint with tube adapter into the functional shank.

- 5) Measure the distance between the bottom edge of the functional shank and the set screws.

- 6) Mark the measured distance to the 4 set screws on the outside of the functional shank.
- 7) Remove the knee joint and tube adapter.
- 8) With a pointed needle or a metal pin heated with a blow dryer, drill a hole in the functional shank at the marked locations.
- 9) Insert the knee joint with the tube adapter into the functional shank.
- 10) Assemble the prosthetic foot through the previously drilled holes.

## 6.6 Sanding the functional shank

### ⚠ CAUTION

#### Incorrect protective equipment

Injury due to defective or lack of protective equipment.

- ▶ Wear protective goggles and a dust mask during grinding.


### INFORMATION

Use a new abrasive sheet to sand the foam.

Sanding with a rounded cone ensures that no grooves are left in the foam in case of slipping off.

The following steps should be observed for sanding the functional shank:


- 1) Protect the footshell with adhesive tape.

- 2) In the area of the ankle at the front and back, do not sand the functional shank flush with the foot cosmesis.

**INFORMATION: A flush transition to the footshell is possible on the sides.**

**INFORMATION: Leave more material at the front and rear. This is compressed in the ankle area by the functional stocking.**

- 3) Sand the functional shank.

**INFORMATION: Do not sand the transition (area marked in white).**


- 4) **For the 3F1=2 functional cosmesis Genium:** In the calf area, sand the leg contour so that it is flush with the installation ring for the inductive charger (4X259).

## 6.7 Assembling the sheath

Complete the following steps to assemble the sheath:

- 1) Pull the sheath over the functional shank so that the stitched side lies below the distal edge of the functional shank.

**INFORMATION: The sheath has to fit over the full length with no longitudinal tension on the foam.**

- 2) Cut off the sheath at the top with an overhang of 6 cm.

- 3) Secure the sheath on the inside all around on the proximal end of the foam with contact adhesive.

**INFORMATION: Note the setting time of the contact adhesive depending on the application thickness and temperature.**

- 4) Fold over the distal end of the sheath and slide it into the recesses of the tube adapter.

**INFORMATION: The sheath must not be glued or cut off in the ankle area.**

- 5) Punch the 4 holes for the tube adapter screws in the sheath as well.


**INFORMATION: Before punching, apply contact adhesive to the corresponding areas to prevent runs in the sheath.**

### 6.8 Only for Genium: preparing the charging unit

Complete the following steps to assemble the charging unit for the Genium:

- 1) Make a longitudinal cut in the sheath.  
**INFORMATION: The slit has to be above the opening in the foam.**


- 2) Secure the sheath on the inside of the opening with contact adhesive.

- 3) Press the installation ring for the inductive charger (4X259) into the opening so it is flush with the foam and sheath.  
**INFORMATION: Please note the 647G1148 instructions for use.**
- 4) Glue in the charging unit flush on the outside of the foam.  
**INFORMATION: Please note the 647G1148 instructions for use.**

## 6.9 Functional Test

After assembly, the following points have to be checked:

1. Is the cable tie correctly installed?
2. Are the functional knee part and functional shank flush with each other?
3. Are the axles and clips correctly engaged?
4. Are the ends of the actuator stirrup flush with the functional knee part?
5. Is the position of the actuator stirrup correct?
6. Can the feedback signals be heard clearly enough?

### 6.9.1 Corrective actions

#### 1. Is the cable tie correctly installed?

Event	Cause	Correction
Noise development	Cable tie too tight	Loosen the cable tie and then slightly tighten it.
	Incorrect position of the cable tie	Check the position of the cable tie (see fig. 11).

#### 4. Are the functional knee part and functional shank flush with each other?

Event	Cause	Correction
Gap between functional knee part and foam or Deformation of the foam due to the functional knee parts	Knee extender not used	1) Loosen the cable tie. 2) Fully extend the knee joint. 3) Use the knee extender.
	Sheath pulled too tight	Put the sheath on again (see Page 41).
	Incorrect position of the actuator stirrup	1) Remove the functional knee part from the joint. 2) Position the actuator stirrup (see fig. 10). 3) Attach the functional knee part to the joint.

#### 3. Are the axles and clips correctly engaged?

Event	Cause	Correction
Noise development	Axles and clips not correctly installed.	1) Install the axles and clips, and verify that they engage correctly (see Page 34). 2) Insert the two blind plugs flush into the clips from the inside (see fig. 8).

**4. Are the ends of the actuator stirrup flush with the functional knee part?**

Event	Cause	Correction
<b>Noise development</b>	Incorrect position of the actuator stirrup ends	<ol style="list-style-type: none"> <li>1) Loosen a screw of the actuator stirrup on the functional knee part.</li> <li>2) Reposition the actuator stirrup (see fig. 12).</li> <li>3) Tighten the screws of the actuator stirrup.</li> <li>4) Repeat steps 1-3 with the second screw.</li> </ol>

**5. Is the position of the actuator stirrup correct?**

Event	Cause	Correction
<b>Collision between actuator stirrup and joint frame</b>	Incorrect position of the actuator stirrup	<ol style="list-style-type: none"> <li>1) Remove the functional knee part from the joint.</li> <li>2) Position the actuator stirrup (see fig. 10).</li> <li>3) Attach the functional knee part to the joint.</li> </ol>

**6. Can the feedback signals be heard clearly enough?**

Event	Cause	Correction
<b>Feedback signals cannot be heard clearly enough</b>	Improper volume setting on the joint	<p>► Adjust the volume of the joint feedback signals (note the instructions for use:</p> <p style="text-align: right;">Genium</p> <p>3B1-2/3B1-2=ST "647G1148", Genium</p> <p>3B1-3/3B1-3=ST "647G1380" and C-Leg</p> <p>3C98-3/3C88-3 "647G1375").</p>

**6.10 Finishing the functional cosmesis**

Complete the following steps to finish the functional cosmesis:

1) Assemble the functional shank, knee joint and prosthetic foot.

2) Fold over the coated ankle part of the functional stocking (coating faces out).

28

- 3) The functional stocking is secured distally 1 cm above the footshell.
- 4) Pull the functional stocking up.  
**INFORMATION: The label field in the functional stocking identifies the back of the stocking.**  
**INFORMATION: Note the correct positioning of the functional zones.**
- 5) For the best possible adhesion, position the end of the functional stocking directly on the socket surface if possible.

#### **INFORMATION**

Together with the user, vary how far the functional stocking is pulled up. The tension is not permitted to noticeably impair the use of the joint. For example, secure the functional stocking on the socket while the prosthesis is angled.

The adhesion of the functional stocking is best directly on the socket material. This needs to be considered when shaping the socket volume.

#### **INFORMATION**

The inductive charger for Genium can be used after putting on the functional stocking. Exposing the charging unit is not necessary.

## **7 Cleaning and care**

### **7.1 Cleaning the functional cosmesis**

- ▶ Clean the plastic parts with a damp, soft cloth and then let them dry.

### **7.2 Cleaning the functional stocking**

- ▶ Wash the product at 30°C/86°F in a washing machine  
Setting: gentle cycle (for delicates).

#### **INFORMATION**

Note the washing symbols on the label.

### **7.3 Care of the functional stocking**

- Wipe the silicone coatings on the inside of the functional stocking with a damp cloth every three days.

## **8 Legal information**

### **8.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for

damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

### 8.3 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 9 Technical data

<b>Environmental conditions</b>	
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F Max. 93% relative humidity, non-condensing
Operation	-10 °C/+14 °F to +60 °C/+140 °F 15% to 93% relative humidity, non-condensing

<b>General information</b>	
Reference number	3F1=1 / 3F1=2
Service life of the functional cosmesis	2 years The service life diminishes or increases depending on the individual level of use. Only when the instructions for use are observed it is possible to achieve the maximum service life, in dependence of the level of use.
Service life of the functional stocking	6 months Wear part, subject to natural wear and tear. The service life diminishes or increases depending on the individual level of use. Only when the instructions for use are observed it is possible to achieve the maximum service life, in dependence of the level of use.
Weight of the functional cosmesis (without adaptations)	910 g

## 10 Appendix

### 10.1 Symbols Used

Manufacturer

LOT:PPPP YYYY WW

Lot number

Declaration of conformity according to the applicable European directives

Do not dry clean

Iron at medium temperature, max. 150°

Do not tumble dry

Gentle cycle at 30°

Do not bleach

Not made with  
natural rubber latex.

Natural rubber latex was not used as a material

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Otto Bock Healthcare Products GmbH