
CE

17CF1=*

DE	Gebrauchsanweisung	7
EN	Instructions for use	11
FR	Instructions d'utilisation	16
IT	Istruzioni per l'uso	21
ES	Instrucciones de uso	25
PT	Manual de utilização	30
NL	Gebruiksaanwijzing	35
SV	Bruksanvisning	40
DA	Brugsanvisning	44
NO	Bruksanvisning	49
PL	Instrukcja użytkowania	53
HU	Használati utasítás	58
CS	Návod k použití	63
SK	Návod na používanie	67
TR	Kullanma talimatı	72
EL	Οδηγίες χρήσης	76
RU	Руководство по применению	81
JA	取扱説明書	86

1

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INFORMATION

Last update: 2019-04-23

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

These instructions for use provide you with important information on the processing of the 17CF1 Carbon Ankle Seven.

2 Product description

The 17CF1 Carbon Ankle Seven is designed to initiate extension moments of the hip and knee during walking and thus achieve extension and knee stability in the lower limbs. The energy generated during heel strike is stored in the carbon matrix and returned at toe-off. In contrast to conventional ankle joint orthoses, the limits in the plantar and dorsal direction are dynamic and are

not subject to static restriction. Due to its design, the carbon spring has a natural outward foot rotation of 7°.

2.1 Available sizes

Selection of carbon springs according to activity level for AFOs

INFORMATION

Installation of the carbon spring in a KAFO:

Use a weaker carbon spring

		Activity level	
Body weight		Normal activity level	High activity level
up to 100 kg		17CF1=L/R1	17CF1=L/R1
up to 90 kg		17CF1=L/R2	17CF1=L/R1
up to 80 kg		17CF1=L/R3	17CF1=L/R2
up to 70 kg		17CF1=L/R4	17CF1=L/R3
up to 60 kg		17CF1=L/R5	17CF1=L/R4
up to 50 kg		17CF1=L/R6	17CF1=L/R5
up to 40 kg		17CF1=L/R8	17CF1=L/R7
up to 30 kg		17CF1=L/R9	17CF1=L/R8
up to 20 kg		17CF1=L/R11	17CF1=L/R10
up to 10 kg		17CF1=L/R12	17CF1=L/R11

2.2 Components

Scope of delivery (see fig. 1)

Item	Quant- ity	Designation	up to 60 kg	up to 100 kg
1	1	Carbon spring		
2	4	Countersunk head screws	501S86=M4x12	501S86=M5x14
3	4	Countersunk washers	507U9=M4	507U9=M5
4	4	Welding nuts	502E3=M4x6	502E3=M5x7.5
5	1	Template		

3 Intended use

3.1 Indications for use

The product is intended **solely** for orthotic fittings of the lower limbs.

3.2 Indications

- Paralysis or weakness of the dorsiflexor and/or plantar flexor muscles.

Conditions of use

The knee extensor musculature must be functionally intact or exhibit a minimum residual function (muscle strength 3 to 5 according to Janda), or be secured with a locked knee joint in a KAFO. The ankle joint has to exhibit residual mobility of 0-0-15° dorsal extension according to the neutral-0 method.

Indications must be determined by the physician.

3.2.1 Contraindications

3.2.1.1 Absolute Contraindications

Not known.

3.2.2 Relative Contraindications

- Structural joint contractures

3.3 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

△ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 Safety instructions

△ CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

△ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

△ CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

△ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Contact with heat, embers or fire

Damage to the product due to high temperatures

- Keep the product away from open flames, embers and other sources of heat.

NOTICE

Thermal overloading of product components

Loss of function due to improper thermal treatment

- Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- Avoid storage in condensing ambient humidity.
- Avoid contact with abrasive substances (e.g. sand, dust).
- Do not expose the product to temperatures below -10 °C or above +40 °C (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparing product for use

5.1 Processing

NOTICE

Damage due to incorrect processing

Premature wear and tear, loss of functionality due to breakage of the carbon spring

- Smooth the ground areas with an abrasive flap disc
- Centre the bores in the carbon spring
- Use a suitable counterbore to deburr the bore holes

INFORMATION

Orthosis alignment

Ottobock recommends checking the alignment of the orthosis using the L.A.S.A.R. Posture or L.A.S.A.R. Posture 3D and adapting it as needed.

Fabricating a plaster model

Due to the design of the carbon spring, make sure the plaster cast is taken with slight outward rotation. Request the information on the plaster casting technique if required (Fabricating a plaster cast – KAFO (knee-ankle-foot orthosis): **646T5411***). Raise the heel slightly by approx. 15 to 20 mm (see fig. 2) for better positioning of the carbon spring below the foot in the plaster cast. The enclosed template may be used for the sagittal positioning of the lower leg (use small template for sizes 10 to 12).

Installing the carbon spring

Prerequisite: The patient-specific plaster model has been fabricated.

- 1) Pull a layer of Perlon stockinette over the plaster model and fabricate the lining using **5 mm** of Pedilin (see fig. 4). The **5 mm** material thickness can be used for all size models.
- 2) Mark the fixation zone and movement zone for the carbon spring, dividing the model into the zones **1** and **2** (see fig. 3, see fig. 5).

- 3) Mark the carbon spring below the foot, with the length extending up to approx. **20 mm** before the rollover zone of the foot (see fig. 6).
- 4) **Optional:** Shorten the carbon spring, at most up to the mark. Deburr and round the cut edges (see fig. 7).
- 5) Create the structure for the carbon spring on the lining using several layers of Pedilin and by sanding (see fig. 8).
- 6) Align the carbon spring in the foot component in the direction of the second foot ray (see fig. 9).
- 7) Contour the Pedilin structure under the carbon spring, avoiding cavities under the carbon spring (see fig. 10).
- 8) Secure the carbon spring on the model with PVA adhesive tape and pull on a layer of Perlon stockinette (see fig. 11).
- 9) For vacuum forming, use homo-polypropylene with a thickness of **4 mm** for all carbon spring sizes (see fig. 12).
- 10) Mark the course of the edge and remove the orthosis from the cast. Then separate the calf component from the foot component according to the markings (see fig. 13).
- 11) Use a ruler to mark the centre of the bores on the calf component (see fig. 14).
- 12) Mark the bores on the foot component. The proximal bore runs in front of the curve of the heel (see fig. 15).
- 13) For the holes in the orthosis shells, use a **4 mm** drill bit for an M4 thread and a **5 mm** drill bit for an M5 thread (see fig. 16).
- 14) Transfer the bore holes to the carbon spring (see fig. 17).
- 15) For the 502E3 welding nut, drill a **5 mm** hole for an M4 thread and a **6.5 mm** hole for an M5 thread (see fig. 18).
- 16) Deburr the bore with a 726S9=90x11.5 counterbore (see fig. 19).
- 17) Press in the welding nut (see fig. 20). To achieve an optimal fit of the bores, screw on the upper section on a trial basis and check the positions of the remaining bore holes. Repeat steps 14 through 16.
- 18) Assemble the carbon spring with the orthosis using the supplied screws, countersunk washers and welding nuts and secure with Loctite 241 (see fig. 21).
- 19) Ottobock recommends positioning the sagittal load line approx. **15 mm** in front of the knee compromise pivot point according to Nietert (see fig. 22).

6 Cleaning

Clean the joint promptly after contact with water containing salt, chlorine or soap, or if it gets dirty (e.g. sand).

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a cloth or allow it to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Spare parts are listed in the section "Components/design".

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

9.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

9.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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