

## 28U90

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	7
<b>FR</b> Instructions d'utilisation .....	11
<b>IT</b> Istruzioni per l'uso .....	15
<b>ES</b> Instrucciones de uso .....	19
<b>PT</b> Manual de utilização .....	23
<b>NL</b> Gebruiksaanwijzing .....	27
<b>SV</b> Bruksanvisning .....	30
<b>DA</b> Brugsanvisning .....	34
<b>PL</b> Instrukcja użytkowania .....	38
<b>HU</b> Használati utasítás .....	41
<b>CS</b> Návod k použití .....	45
<b>RO</b> Instrucțiuni de utilizare .....	49
<b>BG</b> Инструкция за употреба .....	53
<b>EL</b> Οδηγίες χρήσης .....	57
<b>RU</b> Руководство по применению .....	61
<b>ZH</b> 使用说明书 .....	65

	<b>Artikelnummer / Article Number</b>	<b>Seite / Side</b>	<b>Schuhgröße/ Shoe Size</b>	<b>Farbe/ Colour</b>	<b>Farbe Polster/ Colour Pad</b>
	28U90=L35-37-0	links/left	L35-37	Transparent/transparent	Beige/beige
	28U90=L35-37-7	links/left	L35-37	Schwarz/black	Schwarz/black
	28U90=R35-37-0	rechts/right	R35-37	Transparent/transparent	Beige/beige
	28U90=R35-37-7	rechts/right	R35-37	Schwarz/black	Schwarz/black
	28U90=L37-39-0	links/left	L37-39	Transparent/transparent	Beige/beige
	28U90=L37-39-7	links/left	L37-39	Schwarz/black	Schwarz/black
	28U90=R37-39-0	rechts/right	R37-39	Transparent/transparent	Beige/beige
	28U90=R37-39-7	rechts/right	R37-39	Schwarz/black	Schwarz/black
	28U90=L39-41-0	links/left	L39-41	Transparent/transparent	Beige/beige
	28U90=L39-41-7	links/left	L39-41	Schwarz/black	Schwarz/black
	28U90=R39-41-0	rechts/right	R39-41	Transparent/transparent	Beige/beige
	28U90=R39-41-7	rechts/right	R39-41	Schwarz/black	Schwarz/black
	28U90=L41-44-0	links/left	L41-44	Transparent/transparent	Beige/beige
	28U90=L41-44-7	links/left	L41-44	Schwarz/black	Schwarz/black
	28U90=R41-44-0	rechts/right	R41-44	Transparent/transparent	Beige/beige
	28U90=R41-44-7	rechts/right	R41-44	Schwarz/black	Schwarz/black
	28U90=L44-47-0	links/left	L44-47	Transparent/transparent	Beige/beige
	28U90=L44-47-7	links/left	L44-47	Schwarz/black	Schwarz/black
	28U90=R44-47-0	rechts/right	R44-47	Transparent/transparent	Beige/beige
28U90=R44-47-7	rechts/right	R44-47	Schwarz/black	Schwarz/black	

<b>Material</b>	PA, PP
-----------------	--------

## INFORMATION

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### INFORMATION

Last update: 2019-05-20

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use contain important information for fitting and setting up your 28U90 ankle foot orthosis.

## 2 Intended use

### 2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

### 2.2 Indications

- Dorsal flexor weakness, (e.g. after stroke, traumatic brain injury, in multiple sclerosis, neuromuscular atrophy, or peroneal paralysis)

Indications must be determined by the physician.

### 2.3 Contraindications

#### 2.3.1 Absolute Contraindications

- Moderate to severe spasticity in the lower leg
- Moderate to severe oedema
- Leg ulcers
- Moderate to severe foot deformities

#### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

### 2.4 Further Usage Restrictions

The orthosis must not be used for the prevention of contractures (shortening) etc.

The orthosis must not be used for sports activities that involve jumping, sudden movements or rapid step sequences (e.g. basketball, badminton, riding sports). Sports activities should generally be discussed with the patient.

### 2.5 Effects

The ankle foot orthosis provides lifting of the foot and limits plantar flexion.

## 3 Safety

### 3.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 3.2 General safety instructions

 <b>CAUTION</b>
<b>Reuse on other persons and improper cleaning</b>
Skin irritation, formation of eczema or infections due to contamination with germs
▶ The product may be used by one person only.
▶ Clean the product regularly.

**⚠ CAUTION**

**Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

**⚠ CAUTION**

**Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

**NOTICE**

**Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 4 Handling

**⚠ WARNING**

**Operation of motor vehicles**

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.

**INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

### 4.2 Adaptation

**⚠ CAUTION**

**Improper moulding or application**

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

## **⚠ CAUTION**

### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

## **NOTICE**

### **Use without suitable footwear**

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.

### **Adapting the orthosis size:**

- 1) **If necessary:** Cut or sand down the orthosis to size according to shoe size.
- 2) Then, deburr the worked edges and polish out any notches.

### **Adapting the orthosis:**

- 1) **If necessary:** Thermoplastically re-form the orthosis by applying gentle heat to a limited extent.

**NOTICE! Re-forming of the orthosis in cold conditions or at temperatures above 80°C is not permitted.**

- 2) Attach the calf pad and the hook-and-loop strap (see fig. 2).

**INFORMATION: The hook-and-loop strap can be fastened medially or laterally.**

- 3) Once fitting is complete, reattach the hook-and-loop strap to the calf pad (see fig. 3).

### **Final inspection:**

- ▶ Conduct trial walking with the patient. Testing on inclines, ramps and stairs is mandatory. Fine tune as required.

## **4.3 Application**

## **NOTICE**

### **Mechanical overload**

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

1. Open the hook-and-loop strap.
2. Insert the orthosis into the shoe.
3. Put the orthosis on together with the shoe.
4. Tighten and close the hook-and-loop strap.

## **5 Cleaning**

Clean the orthosis regularly:

### **Textile pads**

- 1) Remove the textile pads from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the textile pads in warm water at **30 °C** with standard mild detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## Orthosis frame

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

### 7.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 7.3 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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