12K41, 12K42, 12K44, 12K50

[DE] Gebrauchsanweisung (Benutzer)	2
[EN] Instructions for use (user)	13
[FR] Instructions d'utilisation (Utilisateur)	23
[IT] Istruzioni per l'uso (Utilizzatore)	34
[ES] Instrucciones de uso (Usuario)	45
[NL] Gebruiksaanwijzing (Gebruiker)	56
[SV] Bruksanvisning (Användare)	67
[RU] Руководство по применению (Пользователь)	77
[JA] 取扱説明書(ユーザー用)	88

INFORMATION

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INFORMATION

Date of last update: 2019-07-18

- Please read this document carefully before using the product.
- Have yourself instructed by qualified personnel in the proper and safe use of the product.
- Please contact the qualified personnel if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). The contact information for the manufacturer is found on the back page.
- Please keep this document in a safe place.

The "12K41=*" ErgoArm, 12K42=*" ErgoArm plus, 12K44=*" ErgoArm Hybrid plus, 12K50=*" ErgoArm Electronic plus" product is referred to simply as the product/elbow joint below.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Construction and Function

2.1.1 Design

1. Easy Plug
2. Strap clamp
3. Bowden cable (mechanical lock)
4. Dial
5. Humeral rotation feature adjusting screw

2.1.2 Function

The product is a mechanically/electronically driven elbow joint. It helps the wearer perform everyday tasks in combination with other prosthetic components (Combination possibilities).

Various control programmes make optimal adaptation to the individual needs and abilities possible.

The ErgoArm models differ in the following characteristics:

Models	AFB	Easy Plug	Lock	Slip-stop
12K41=* ErgoArm	-	-	Mechanical	Mechanical
12K42=* ErgoArm plus	x	-	Mechanical	Mechanical

Models	AFB	Easy Plug	Lock	Slip-stop
12K44=* ErgoArm Hybrid plus	x	x	Mechanical	Mechanical
12K50=* ErgoArm Electronic plus	x	x	Mechanical or Electronic	Mechanical or Electronic

2.1.2.1 Definition of terms

AFB (Automatic Forearm Balance)

The AFB flexion aid stores the energy that is released when the prosthetic arm is extended and uses it to assist with flexion. The complex mechanism of the flexion aid controls the dynamics of forearm movement and equalises the leverage forces that develop via a compensation gear. As a result, the patient only has to exert minimal force on the flexion cable or through body momentum. The compensation coefficient can be easily adjusted to the individual weight of the prosthetic forearm and different clothing.

Lock

The product is either locked mechanically or electronically (see Page 13). The internal ratchetless lock can be unlocked or locked in any position, even under a load. The product has a load capacity of up to 230 N in the locked state (this value may decline over the course of several years due to wear and tear). The lock will slip under higher loads.

Mechanical lock

A strong pull on the cable (audible clicking sound) locks or unlocks the joint.

Thanks to the integrated compensation mechanism, the lock can be released with little effort, even under load. The integrated strain relief reduces the risk of damage when the pull cable is pulled too hard.

Electronic lock

The lock can be released and engaged using electrode signals or switch signals (depending on the chosen switching version).

Slip-stop function

The slip-stop function makes it possible to lower the forearm in a controlled manner without having to fully release the lock.

Slip-stop function, mechanical

A gentle pull on the cable (no audible click) unlocks the product so the forearm can be moved continuously. When the pull on the cable stops, the product is locked again.

Slip-stop function, electronic

The slip-stop function is activated by the electrode signals (depending on the chosen control version). When the electrode signal stops, the product is locked again.

3 Application

3.1 Indications for use

The product is intended **exclusively** for upper limb exoprosthetic fittings.

3.2 Conditions of use

The product is designed for use by unilateral or bilateral amputees.

The product is intended **exclusively** for use on adults.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

4 Safety

4.1 Explanation of warning symbols

⚠ WARNING Warning regarding possible serious risks of accident or injury.

⚠ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 Structure of the safety instructions

⚠ CAUTION

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

⚠ WARNING

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

⚠ WARNING

Operating the product near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.

- ▶ When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ WARNING

Use of the product while operating a vehicle or machinery

- > Risk of accident due to unexpected behaviour of the product.
- > Risk of injury due to faulty control or malfunction of the product.
- ▶ Applicable national legal regulations must be complied with when operating vehicles of any kind or when operating machinery with the product.
- ▶ Have your driving ability examined and verified by an authorised test centre.
- ▶ Switch off the product before driving.

INFORMATION

The ability of prosthesis wearers to operate a vehicle is determined on a case-by-case basis. Criteria include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, prosthesis design) and the amputee's individual abilities.

⚠ CAUTION

Changing prosthetic components when switched on

Risk of injury due to faulty control or malfunction of the prosthesis solution.

- ▶ Remove the rechargeable battery from the battery receptacle or switch off the prosthesis solution before changing prosthetic components (such as the terminal device).
- ▶ The 12K50= ErgoArm can be locked and unlocked by means of the mechanical pull cable even when it is switched off.

⚠ CAUTION

Overloading due to unusual activities

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional. If necessary, he or she will forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Manual unlocking of elbow lock under load

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

⚠ CAUTION

Penetration of dirt and humidity into the product

Risk of injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.
- ▶ Do not expose the product, and especially the elbow joint, to dripping or splashing water.
- ▶ Wear the product, and especially the elbow joint, under suitably resistant clothing in the rain.

⚠ CAUTION

Signs of wear and tear on the product components

Injury due to faulty control or malfunction of the product

- ▶ In the interest of preventing injuries and to maintain product quality, we recommend performing regular service.
- ▶ Please contact your O&P professional for more information.

⚠ CAUTION

Unintentional unlocking of the terminal device

Risk of injury due to releasing the terminal device from the forearm (e.g. while carrying objects).

- ▶ When using a quick-disconnect wrist unit, position the terminal device so that slight twisting does not release the terminal device from the forearm.

⚠ CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the product components.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Independent manipulation of the product

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

⚠ CAUTION

Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.99 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
 - DECT cordless phones incl. base station: 0.35 m
 - WiFi (routers, access points,...): 0.22 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

⚠ CAUTION

Risk of pinching in the joint flexion area

Injuries due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the joint.

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product and the contacts with a damp cloth and mild soap only.

NOTICE

Damage to the battery

Damage to the battery due to using both connection options.

- Only use one of the two connections for the battery (connection on the forearm cable or Easy Plug connection).

5 Scope of Delivery and Accessories

12K41 ErgoArm plus scope of delivery

- 1 pc. 12K41= * ErgoArm
- 1 pc. 647G1475 Instructions for use (user)

12K42 ErgoArm plus scope of delivery

- 1 pc. 12K42= * ErgoArm plus
- 1 pc. 647G1475 Instructions for use (user)

12K44 ErgoArm Hybrid plus scope of delivery

- 1 pc. 12K44= * ErgoArm Hybrid plus
- 1 pc. 647G1475 Instructions for use (user)

12K50 ErgoArm Electronic plus scope of delivery

- 1 pc. 12K50= * ErgoArm Electronic plus
- 1 pc. 647G1475 Instructions for use (user)

6 Use

6.1 Adjusting the level of compensation

⚠ CAUTION

Incorrect donning and doffing of the prosthesis

Risk of injury due to sudden flexion of the prosthesis.

- Switch off the prosthesis before donning or doffing.
- The prosthesis must always be flexed for donning and doffing.

⚠ CAUTION

Sudden flexion of the product due to high level of compensation

Risk of injury due to sudden flexion of the prosthesis.

- Note the level of compensation when taking off heavy clothing.
- Note the setting for the level of compensation when the prosthesis is not being worn.
- The prosthesis must always be flexed for donning and doffing.


Complete the following steps to set the level of compensation:

- 1) Flex the arm.
- 2) Use the dial to adjust the compensation.
INFORMATION: The adjustment range is limited by stops.
INFORMATION: The adjustment mechanism is protected by a slip clutch.

There are three different compensation ranges:

2

**Range 1 = low compensation level**

Permits the arm to swing freely while walking.

Range 2 = moderate compensation level

Balance is positively supported when flexing the arm. Compensation is reduced again during extension.

Range 3 = constant compensation level

If set correctly, the weight of the forearm is balanced by the compensation level.

Range 4 = low compensation level

Low compensation until flexion stop.

6.2 Adjusting the upper arm rotation joint

Complete the following steps to rotate and adjust the upper arm rotation joint:

3

- 1) Loosen the adjusting screw.
- 2) Move the upper arm rotation joint to the desired position.
- 3) Re-tighten the adjusting screw.

7 Charging the battery

The following points must be observed when charging the battery:


- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product once a day when used by the patient on a daily basis.
- When used daily, the complete charging unit (power supply – battery charger) may remain plugged into the wall socket.

7.1 757L20 battery charger and 757B2* EnergyPack

NOTICE
Deep discharge of the battery

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes ☀.
- 5) Remove the rechargeable battery after the charging process is complete.

7.2 757L35 battery charger and 757B35-* rechargeable battery

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (Status signals).
- 2) The charging process starts.
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

8 Cleaning

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed.
Ensure that no liquid penetrates into the system component(s).
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

9 Maintenance

We recommend performing regular maintenance (service inspections) to prevent injuries and maintain the quality of the product. This should be carried out when the maintenance interval expires, i.e. after 24 months. The manufacturer grants a grace period of no more than one month before or three months after the due date. The maintenance includes inspecting the sensors and replacing wear and tear parts.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

10 Legal information

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE conformity

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

11 Technical data

Environmental conditions	
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-20 °C/-4 °F to +40 °C/+104 °F max. 93 % relative humidity, non-condensing
Operation	+5 °C/+41 °F to +40 °C/+104 °F max. 93 % relative humidity, non-condensing

Reference number	12K41=*	12K42=*	12K44=*	12K50=*
Current draw (in mA)	-	-	-	3.5
Power supply	-	-	-	757B20 757B21 757B25=* 757B35=*
Max. load in locked state (N)			230	
Flexion angle (in °)			15°–145°	
Service life			5 years	

12 Appendices

12.1 Symbols on the product

Declaration of conformity according to the applicable European directives

SN YYYY WW NNN Serial number

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Manufacturer

12.2 Display of the current charge level

12.2.1 757B2* EnergyPack

LED indicator	Event
	Charge level above 50% (illuminated in green)
	Charge level below 50% (flashing alternating green and orange)
	Charge level below 5% (illuminated in orange)

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

12.2.2 757B35-* Battery

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

Charging receptacle	Event
	Battery fully charged (illuminated in green)
	Battery 50% charged (illuminated in yellow)
	Battery drained (illuminated in orange)

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Otto Bock Healthcare Products GmbH