






## 50P50 Manu Junior 50P51 Manu Junior Pollex

<b>DE</b> Gebrauchsanweisung .....	5	<b>FI</b> Käyttöohje.....	33
<b>EN</b> Instructions for use .....	8	<b>PL</b> Instrukcja użytkowania.....	36
<b>FR</b> Instructions d'utilisation.....	10	<b>HU</b> Használati utasítás .....	39
<b>IT</b> Istruzioni per l'uso .....	13	<b>CS</b> Návod k použití.....	41
<b>ES</b> Instrucciones de uso .....	16	<b>SK</b> Návod na používanie .....	44
<b>PT</b> Manual de utilização.....	19	<b>ET</b> Kasutusjuhis .....	47
<b>NL</b> Gebruiksaanwijzing.....	22	<b>RU</b> Руководство по применению .....	50
<b>SV</b> Bruksanvisning .....	25	<b>JA</b> 取扱説明書.....	53
<b>DA</b> Brugsanvisning .....	28	<b>ZH</b> 使用说明书.....	55
<b>NO</b> Bruksanvisning.....	30		

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### INFORMATION

Date of last update: 2020-03-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50P50 Manu Junior and 50P51 Manu Junior Pollex wrist supports for children.

## 2 Intended use

### 2.1 Indications for use

The support is intended **exclusively** for treating the upper limbs and **exclusively** for contact with intact skin.

The support must be used according to the indication.

### 2.2 Indications

- Congenital malformations of the wrist/carpus
- Tendinopathy
- Ligament injury, ligament strain in the wrist
- Wrist instability
- Impingement syndrome
- Scaphoid pseudoarthrosis
- Juvenile rheumatoid arthritis
- Irritation in the wrist area
- Wrist drop due to paresis or paralysis of the radial nerve

- Carpal tunnel syndrome

### 50P51 Manu Arexa Pollex additional indications

- Irritation and instability in the area of the thumb saddle joint

The indication must be determined by the physician.

## 2.3 Contraindications

### 2.3.1 Absolute Contraindications

None known.

### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia in the area treated; lymphatic flow disorders, including unclear soft tissue swelling distal to the device; sensory and circulatory disorders in the area of the forearm and the hand.

## 2.4 Mechanism of Action

The supports stabilize the wrist joint, and the 50P51 Manu Junior Pollex also stabilizes the thumb joint.

## 3 Safety

### 3.1 Explanation of warning symbols

 **CAUTION**

Warning regarding possible risks of accident or injury.

 **NOTICE**

Warning regarding possible technical damage.

### 3.2 General safety instructions

 **CAUTION**

#### Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.

- ▶ Clean the product regularly.

**CAUTION**

**Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

**NOTICE**

**Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

**NOTICE**

**Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

**NOTICE**

**Improper molding or application**

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ The product may only be molded by qualified personnel.
- ▶ Do not make any improper changes to the product.
- ▶ Always apply the product according to the information in the instructions.

## 4 Handling

**INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Size selection

- 1) Measure the wrist circumference.
- 2) Determine the size of the support (see size chart).

### 4.2 Fitting and Application

**CAUTION**

**Incorrect or excessively tight application**

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

An anatomically pre-shaped metal splint is located in the outer pocket of the support. The shape of the splint can be customized to fit the user.

- 1) Open the support; the hook-and-loop closure at the thumb remains fastened.
- 2) Remove the splint from the pocket (see fig. 3).
- 3) Mould the splint by hand so that it fits the anatomical shape of the patient's palm (see fig. 4).
- 4) Insert the splint into the pocket with the outer spoon-shaped curve facing the palm.

- 5) Guide the thumb through the thumb opening and fasten the hook-and-loop closures (see fig. 5).
  - 6) Check the fit of the support and adjust the hook-and-loop closure at the thumb if necessary (see fig. 6, see fig. 7).
  - 7) **Optional:** Trim the hook-and-loop closures with scissors (see fig. 8).
- The support fits snugly without pinching.

### 4.3 Removal

- 1) Open the hook-and-loop closures on the support one after the other; the hook-and-loop closure for the thumb remains fastened.
- 2) Remove the support.

## 5 Cleaning

### NOTICE

#### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

- 1) Remove the metal splint from the pocket of the support.
- 2) Fasten all hook-and-loop closures.
- 3) Recommendation: Use a laundry bag or net.
- 4) Wash the support in warm water at **40 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 6) Insert the metal splint into the pocket of the support.

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

## 7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

## 7.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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