

## 17LA3N


|  |     |
|--|-----|
| <b>DE</b> Gebrauchsanweisung .....         | 5   |
| <b>EN</b> Instructions for use .....       | 12  |
| <b>FR</b> Instructions d'utilisation ..... | 20  |
| <b>IT</b> Istruzioni per l'uso .....       | 28  |
| <b>ES</b> Instrucciones de uso .....       | 36  |
| <b>PT</b> Manual de utilização .....       | 44  |
| <b>NL</b> Gebruiksaanwijzing .....         | 52  |
| <b>SV</b> Bruksanvisning .....             | 60  |
| <b>DA</b> Brugsanvisning .....             | 67  |
| <b>NO</b> Bruksanvisning .....             | 75  |
| <b>FI</b> Käyttöohje .....                 | 82  |
| <b>CS</b> Návod k použití .....            | 89  |
| <b>SK</b> Návod na používanie .....        | 97  |
| <b>RU</b> Руководство по применению .....  | 104 |
| <b>ZH</b> 使用说明书 .....                      | 113 |

1

2

## INFORMATION

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### INFORMATION

Date of last update: 2020-03-25

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the processing of the 17LA3N ankle joint.

## 2 Product description

### 2.1 Available sizes

The size of the product is selected according to the planned use (unilateral or bilateral) and the body weight.

|            | Maximum body weight in kg (lbs)    |                                    |
|------------|------------------------------------|------------------------------------|
|            | Unilateral use (lateral or medial) | Bilateral use (lateral and medial) |
| 17LA3N=10* | 15 (33)                            | 25 (55)                            |
| 17LA3N=12* | 20 (44)                            | 40 (88)                            |
| 17LA3N=14* | 50 (110)                           | 80 (176)                           |
| 17LA3N=16* | 85 (187)                           | 120 (264.5)                        |
| 17LA3N=20* | 110 (242.5)                        | 160 (352.5)                        |

#### INFORMATION

- ▶ Use the next larger size if the product is being worn unilaterally and in case of flexion contractures in the knee or hip  $>10^\circ$ , distinct torsion or varus/valgus instabilities or increased physical activity.
- ▶ Use the product bilaterally if there is an ischial support.

### 2.2 Components/design (see Fig. 1)

| Item | Component   | Reference number  |
|------|-------------|---|
| 1    | Set screw   | 30Y349=12 (size 10, 12)<br>30Y309=12 (size 14)<br>30Y309=20 (size 16, 20)                       |
| 2    | Stop bumper | <b>Stop set</b><br>17LA35N=12 (size 10, 12)<br>17LA35N=14 (size 14)<br>17LA35N=20 (size 16, 20) |
| 3    | Stop pin    |   |
| 4    | Stop sleeve |   |
| 5    | Bearing nut |   |
| 8    | Joint screw | <b>Screw connection</b><br>17LA32N=10 (size 10)   |

| Item | Component               | Reference number   |
|------|-------------------------|--|
|      |                         | 17LA32N=12 (size 12)<br>17LA32N=14 (size 14)<br>17LA32N=16 (size 16)<br>17LA32N=20 (size 20)   |
| 6    | Base body               |  |
| 7    | Axial washers           | <b>Axial washers set</b><br>17LA33N=10 (size 10)<br>17LA33N=12 (size 12)<br>17LA33N=14 (size 14)<br>17LA33N=16 (size 16)<br>17LA33N=20 (size 20) |
| 9    | Countersunk head screws | 501S137=M4x8 (size 10, 12)<br>501S137=M4x10-2 (size 14)<br>501S137=M5x10-1 (size 16, 20)   |
| 10   | Ball                    | <b>Compression spring set</b><br>17LA34N=12 (size 10, 12)<br>17LA34N=14 (size 14)<br>17LA34N=16 (size 16)<br>17LA34N=20 (size 20)                |
| 11   | Spring guide housing    |  |
| 12   | Compression spring      |  |
|      |                         |  |

### Not included in scope of delivery

| Item            | Component   | Reference number |
|-----------------|---|------------------|
| not illustrated | Lamination bar  | 17LS3=*          |
|                 | or side bar   | 17LV3=*          |
|                 | Foot stirrup  | 17LF3N=*         |
|                 | or foot stirrup   | 17LF31N=*        |
| see fig. 2      | <b>Lamination dummy set</b><br>Lamination dummy – lower part (item 1)<br>Shoulder screw (item 3)<br>The cover for the lamination dummy (item 2) is not included in the set. | 17LD1N=*         |

## 3 Intended use

### 3.1 Indications for use

System ankle joints are intended to be used **exclusively** for providing a patient with a lower limb orthosis.

In case of unilateral use, the manufacturer recommends that the system ankle joint be processed by means of the lamination resin technique or carbon prepreg technology.

## 3.2 Indications

Partial or total paralysis of the leg muscles.

Indications must be determined by the physician.

## 4 Safety

### 4.1 Explanation of warning symbols

**⚠ CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 4.2 Safety instructions

**⚠ CAUTION**

#### **Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**⚠ CAUTION**

#### **Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

**⚠ CAUTION**

#### **Excessive strain due to use on more than one patient**

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

**NOTICE**

#### **Thermal overloading of the product**

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

#### NOTICE

##### **Thermal overloading of product components**

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

#### NOTICE

##### **Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

## 5 Preparation for use

#### CAUTION

##### **Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

### 5.1 Processing

#### INFORMATION

The instructions for fabricating an orthosis with unilateral joints are included in the 646T5441 technical information. The laminate structure described here is a recommendation from Ottobock and must be adjusted to the needs of the patient.

- > **Required materials:** 17LD1N=\* lamination dummy with shoulder screw, 17LD1N=-P dummy cover, 17LF3N=\* or 17LF31N foot stirrup, 17LS3=\* lamination bar, 501S137=\* countersunk head screws, 636K8 Plastaband, 743Y56=4 alignment adapter
- > The plaster model has been produced in the neutral position.
- 1) Install the carbon fibre-reinforced lamination bar with the two countersunk head screws on the lower part of the lamination dummy.
- 2) The bearing bushing must be removed from the foot stirrup before the foot stirrup can be mounted on the lamination dummy. **INFORMATION: The bearing bushing can be easily pushed out of the foot stirrup by means of the lamination dummy cover. Important! Set the bearing bushing aside in a safe place and reinsert it in the foot stirrup prior to final assembly of the orthosis.**
- 3) Place the reinforced foot stirrup in the lamination dummy cover.
- 4) Set the cover onto the lower section of the lamination dummy and secure it with the shoulder screw. The design of the cover causes it to be moved into the neutral position.
- 5) Attach the lamination dummy to the plaster model by inserting the shoulder screw through the cover, foot stirrup and lower section and fastening it to the alignment adapter.
- 6) Seal the screws of the bar and dummy with Plastaband.
- 7) Laminate the orthosis according to the patient's needs.
- 8) After the laminate hardens, use a knife to cut free the lamination dummy and remove it.
- 9) Remove the Plastaband from the screws.
- 10) Remove the two countersunk head screws of the lamination bar.
- 11) Loosen the lamination dummy screws and remove the cover.
- 12) Use the oscillating saw to expose the edges for removal.
- 13) Remove the components from the plaster.

## 5.2 Assembling the joint

- ▶ Install the ankle joint for the first fitting.

### Mounting the axial washers

- 1) Select and insert the axial washers (see the following table for the thickness of the axial washers).
- 2) Optimise the movement of the orthosis joint by replacing the axial washers step-by-step: Use thicker washers if there is play in the medial/lateral direction, thinner washers if there is resistance in plantar flexion/dorsal extension.

| Orthosis joint size   | Thickness of axial washers [mm]           |
|---|---|
| 17LA3N=10*, 17LA3N=12*  | 0.70, 0.75, <b>0.80</b> , 0.85, 0.90      |
| 17LA3N=14*, 17LA3N=16*, 17LA3N=20*  | 0.85, 0.90, <b>0.95</b> , 1.0, 1.05, 1.10 |
| <b>Start by inserting 2 axial washers of the thickness indicated in bold.</b> |   |

### Mounting the foot stirrup

- 1) Place the bearing bushing (item 1) in the foot stirrup (item 2) (see fig. 3).
- 2) Insert the foot stirrup together with the axial washers (item 3) into the base body (see fig. 4). The stops of the axial washers must be positioned in the receivers on the underside of the joint for this purpose (see fig. 5).
- 3) Slide the bearing nut from the back side through the base body and the foot stirrup and use the joint screw to screw it from the front.
- 4) Check whether the foot component can be moved smoothly and without play. Use axial washers of varying thicknesses if necessary.
- 5) Make sure that the axial washers prevent direct contact between the metal base body and the foot stirrup.
- 6) For the final assembly, secure the screw connection with Loctite 241 and the corresponding torque.

|                    | Torque in Nm (lbf in.) |             |             |             |             |
|--------------------|------------------------|-------------|-------------|-------------|-------------|
|                    | 17LA3N=1-0*            | 17LA3N=1-2* | 17LA3N=1-4* | 17LA3N=1-6* | 17LA3N=2-0* |
| <b>Joint screw</b> | 3 (27)                 | 3 (27)      | 5 (44.5)    | 8 (71)      | 8 (71)      |

### Mounting the bar

- Fasten the bar to the basic piece using the countersunk head screws and secure the screw connections with Loctite 241 and the corresponding torque.

|                                | Torque in Nm (lbf in.) |             |             |             |             |
|--------------------------------|------------------------|-------------|-------------|-------------|-------------|
|                                | 17LA3N=1-0*            | 17LA3N=1-2* | 17LA3N=1-4* | 17LA3N=1-6* | 17LA3N=2-0* |
| <b>Counter-sunk head screw</b> | 3,5 (31)               | 3,5 (31)    | 3,5 (31)    | 7 (62)      | 7 (62)      |

## 5.3 Installing the modules

The individual combination options of the modules, the stop pin and the compression spring are used to adapt the system ankle joint to the patient.

### Mounting the compression spring

- 1) Place the ball together with the spring guide housing in the base body; the ball is positioned in the receiver of the spring guide housing.
- 2) Mount the compression spring.
- 3) For the final assembly, apply Loctite 241 to secure the thread of the set screw.
- 4) Close the channel on the base body with the set screw, observing the angle setting of the ankle joint as you do so.

### Mounting the stop pin

- 1) Slide the stop pin into the stop sleeve and insert them into the base body together. **INFORMATION: When installing the stop, ensure that the tapered side of the stop is resting on the foot stirrup.**
- 2) Insert the stop bumper.
- 3) For the final assembly, apply Loctite 241 to secure the thread of the set screw.
- 4) Close the channel on the base body with the set screw, observing the angle setting of the ankle joint as you do so.

### Adjusting the stop and spring pre-tension

- ▶ To adjust the pre-tension, use a Torx wrench to screw the set screw deeper into the channel or unscrew it.

## 6 Cleaning

Clean the orthosis joint promptly after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the orthosis joint with clean, fresh water.
- 2) Dry the orthosis joint with a cloth or allow it to air-dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

## 7 Maintenance

### INFORMATION

**The product may be exposed to increased loads by the patient.**

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Only use 633F7 special lubricant.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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