


17AD1*, 17AF2*

DE Gebrauchsanweisung.....	4
EN Instructions for use.....	8
FR Instructions d'utilisation	11
IT Istruzioni per l'uso.....	15
ES Instrucciones de uso	19
PT Manual de utilização	23
NL Gebruiksaanwijzing	27
SV Bruksanvisning	30
DA Brugsanvisning	34
NO Bruksanvisning	37
PL Instrukcja użytkowania	41
HU Használati utasítás.....	45
CS Návod k použití	49
TR Kullanma talimatı	52
EL Οδηγίες χρήσης.....	56
RU Руководство по применению.....	60

INFORMATION

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INFORMATION

Last update: 2020-04-01

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

Explanation of symbols

 **CAUTION** Warnings regarding possible risks of accident or injury.

 **NOTICE** Warnings regarding possible technical damage.

 **INFORMATION** Additional information on the fitting / use.

1 Application

1.1 Indications for use

The Ankle Joints are to be used **exclusively** for the fitting with ankle foot orthoses (AFOs) made of thermoplastic material and must be used in pairs.

Size	Weight	Lower Leg length (MTP/floor)
17AD1=120 / 17AF2=77	up to 100 kg / 220 lbs	up to 55 cm
17AD1=93 / 17AF2=63	up to 35 kg / 77 lbs	up to 35 cm

1.2 Indication

Paralysis or partial paralysis of the lower leg muscles, for guidance of the ankle joints. Specific indications must be determined by the physician.

1.3 Safety instructions

CAUTION

Risk of injury due to excessive load on the ankle joint. A deviation from specified values may lead to breaks caused by excessive flexion loads. This can cause the patient to fall. Observe the alignment instructions and adhere to the specified structural limits.

CAUTION

Risk of injury due to misaligned joint installation. Deviations in the frontal and / or sagittal plane may cause joint stiffness. Ensure the orthosis ankle joints are always installed parallel.

CAUTION

Results of product deterioration. Wear and tear on system components can lead to malfunction of the orthotic joint. This can cause the patient to fall. In the interest of the patient's safety, we recommend inspecting the orthotic joints every six months (refer to Section 3.2 Maintenance).

NOTICE

Risk of breakage due to use of improper materials. The ankle joints have been developed for the fabrication of AFOs made of thermoplastic material. It is not recommended to use these ankle joints in lower leg orthoses made in carbon lamination resin technique. Static loads would be transferred directly into the joints and the stability would be jeopardised.

1.4 Design/Replacement parts

The 17AD1*, 17AF2* Ankle Joints are made of stainless steel. The joints are provided with bearing pins (17AD1*) or bearing washers (17AF2*).

The 17AD1* Ankle Joint is to be equipped either with stop screws for limitation of movement or with compression springs for influencing the movement, e.g. for dorsal flexion.

The 17AF2* is a free-motion ankle joint without limitation of movement.

Replacement parts (fig 1 + 2):

1.	501S89=M5×6	Fixing screws	17AD1*, 17AF2*
2.	501S89=M5×8	Fixing screws for thicker plastic material	17AD1*, 17AF2*
3.	513D83=1.1×3.7	Compression spring	17AD1*
4.	30M8	Bearing pin	17AD1*
5.	501S32=M4X7.5X10	Bearing screw	17AD1*
	507U5=9.8X7.4X1.5	Red washer	17AD1*
	17Y93=M4X6X6.1	Bearing nut	17AD1*
6.	509Y1=5.0	Pressure ball	17AD1*
7.	506G21=M6×6	Set screw	17AD1*
8.	506A8=5x18	Stop pin	17AD1*
9.	29PK5	Alignment tubing	17AD1*, 17AF2*

2 Technical data

Size	17AD1=120	17AD1=93	17AF2=77	17AF2=63
Weight (~ g; pair)	160 g	130 g	60 g	30 g
Max. body weight	100 kg / 220 lbs	35 kg / 77 lbs	100 kg / 220 lbs	35 kg / 77 lbs
Lower leg length (MTP/floor)	55 cm	35 cm	55 cm	35 cm

3 Handling

3.1 Processing

⚠ CAUTION

Risk of injury due to use of improper materials. The ankle joints have been developed for the fabrication of AFOs made of thermoplastic material. It is not recommended to use these ankle joints in lower leg orthoses made in carbon lamination resin technique. Static loads would be transferred directly into the joints and the stability would be jeopardised.

Please note that the Ottobock Ankle Joints, due to the different alignment parameters such as the created forefoot lever load, are subjected to different loads.

For mounting the Ankle Joints please observe the following instructions.

- Determine the compromise pivot axis of the ankle joint, using the 743A7 Pivot Point Adjustment Aid for example, and mark it on the plaster negative (fig. 3).
- Pierce through the marked pivot points and position the Alignment Tube (fig. 4 + 5).
- The reinforcement bar used to fix the plaster positive must not be in contact with the tube to prevent it from distorting (fig. 6).
- Cut the tube ends on the plaster positive flush with the plaster's contours (fig. 7).
- Clamp the plaster positive, for example in the 743A7 Pivot Point Adjustment Aid, and determine the joints' alignment lines (fig. 8).
- Fully screw the long alignment screws through the threaded pivot holes of the joints and insert them into the tubes. Lay up the joints to the alignment lines on the plaster positive and using Bending Irons adapt the joint bars to the shape of the cast (fig. 9 + 10).
- Attach the adapted joint bars with Moulding Clay or similar and small pins to the plaster positive.
- The Moulding Clay must be trimmed flush with the joint bars' contours (fig. 11).
- Prepare for vacuum forming with 99B25 Nylon Stockinette (2 layers) or similar.
- Vacuum form as usual with the desired thermoplastic material (fig 12).
- Mark and cut to the trim lines and grind the cut edges smooth (fig. 13).

INFORMATION

At the foot stirrup, the cutting edge must be directly below the joint to achieve as much contact of the foot stirrup as possible.

- Pierce through the screw holes for fixing the joints using a pre-heated bradawl, for example, and rework the holes to fit the screws. (fig. 14)
- Attach the joints with screws. Do not forget to apply 636K13 Loctite before final delivery of the orthosis. Shorten any protruding screws (fig. 15).

INFORMATION

Make sure to fully screw the attachment screws into the threads!

Make sure to secure the joints' fixing screws with 636K13 Loctite before delivering the orthosis.

INFORMATION

Make sure that with the 17AD1* the compression spring is not tightened up to the stop, as this can lead to increased wear.

3.2 Maintenance

The free-motion 17AF2=77 and 17AF2=63 Ankle Joints need no maintenance.

For lubrication of the 17AD1=120 and 17AD1=93 double-function joints, please use 633F7 Ot-tobock Lubricant.

We recommend to check orthoses with these functional components every 6 months for proper function.

3.3 Wearing Parts for maintenance

Only the 17AF2=77 and 17AF2=63 double-function ankle joints require replacement parts.

Please refer to the list of replacement parts in Section 1.4 "Design/Replacement Parts".

4 Disposal

Dispose of the product in accordance with national regulations.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website..

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