
CE

## 17CF1=\*

DE	Gebrauchsanweisung .....	7
EN	Instructions for use .....	11
FR	Instructions d'utilisation .....	16
IT	Istruzioni per l'uso .....	21
ES	Instrucciones de uso .....	26
PT	Manual de utilização .....	31
NL	Gebruiksaanwijzing .....	36
SV	Bruksanvisning .....	41
DA	Brugsanvisning .....	45
NO	Bruksanvisning .....	50
PL	Instrukcja użytkowania .....	54
HU	Használati utasítás .....	59
CS	Návod k použití .....	64
SK	Návod na používanie .....	69
TR	Kullanma talimatı .....	73
EL	Οδηγίες χρήσης .....	78
RU	Руководство по применению .....	83
JA	取扱説明書 .....	88

1

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### INFORMATION

Date of last update: 2020-04-02

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17CF1 Carbon Ankle Seven.

## 2 Product description

The 17CF1 Carbon Ankle Seven is designed to initiate extension moments of the hip and knee during walking and thus achieve extension and knee stability in the lower limbs. The energy generated during heel strike is stored in the carbon matrix and returned at toe-off. In contrast to conventional ankle joint orthoses, the limits in the plantar and dorsal direction are dynamic and are not subject to static restriction. Due to its design, the carbon spring has a natural outward foot rotation of 7°.

### 2.1 Available sizes

Selection of carbon springs according to activity level for AFOs

#### INFORMATION

##### Installation of the carbon spring in a KAFO:

Use a weaker carbon spring

		Activity level	
Body weight	Normal activity level	High activity level	Spring width
up to 100 kg	17CF1=L/R1	17CF1=L/R1	30 mm
up to 90 kg	17CF1=L/R2	17CF1=L/R1	30 mm
up to 80 kg	17CF1=L/R3	17CF1=L/R2	30 mm
up to 70 kg	17CF1=L/R4	17CF1=L/R3	30 mm
up to 60 kg	17CF1=L/R5	17CF1=L/R4	30 mm
up to 50 kg	17CF1=L/R6	17CF1=L/R5	25 mm/30 mm
up to 40 kg	17CF1=L/R8	17CF1=L/R7	25 mm
up to 30 kg	17CF1=L/R9	17CF1=L/R8	25 mm
up to 20 kg	17CF1=L/R11	17CF1=L/R10	22 mm
up to 10 kg	17CF1=L/R12	17CF1=L/R11	22 mm

### 2.2 Components

#### Scope of delivery (see fig. 1)

Item	Quant- ity	Designation	up to 60 kg	up to 100 kg
1	1	Carbon spring		
2	4	Countersunk head screws	501S86=M4x12	501S86=M5x14
3	4	Countersunk washers	507U9=M4	507U9=M5
4	4	Welding nuts	502E3=M4x6	502E3=M5x7.5

<b>Scope of delivery (see fig. 1)</b>				
<b>Item</b>	<b>Quant- ity</b>	<b>Designation</b>	<b>up to 60 kg</b>	<b>up to 100 kg</b>
5	1	Template		

## **3 Intended use**

### **3.1 Indications for use**

The product is intended **solely** for orthotic fittings of the lower limbs.

### **3.2 Indications**

- Paralysis or weakness of the dorsiflexor and/or plantar flexor muscles.

### **Conditions of use**

The knee extensor musculature must be functionally intact or exhibit a minimum residual function (muscle strength 3 to 5 according to Janda), or be secured with a locked knee joint in a KAFO. The ankle joint has to exhibit residual mobility of 0-0-15° dorsal extension according to the neutral-0 method.

Indications must be determined by the physician.

### **3.2.1 Contraindications**

#### **3.2.1.1 Absolute Contraindications**

Not known.

#### **3.2.2 Relative Contraindications**

- Structural joint contractures

## **3.3 Qualification**

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

## **4 Safety**

### **4.1 Explanation of warning symbols**

**△ CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### **4.2 Safety instructions**

**△ CAUTION**

#### **Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- Only use the product for the defined area of application.
- If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**△ CAUTION**

#### **Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- Use caution when working with the product.
- Check the product for proper function and readiness for use.

- In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

### **⚠ CAUTION**

#### **Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- Observe the alignment, assembly and adjustment instructions.

### **⚠ CAUTION**

#### **Excessive strain due to use on more than one patient**

Risk of injury and loss of functionality as well as damage to the product

- Use the product on only one patient.
- Observe the maintenance recommendations.

### **NOTICE**

#### **Contact with heat, embers or fire**

Damage to the product due to high temperatures

- Keep the product away from open flames, embers and other sources of heat.

### **NOTICE**

#### **Thermal overloading of product components**

Loss of function due to improper thermal treatment

- Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

### **NOTICE**

#### **Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- Avoid storage in condensing ambient humidity.
- Avoid contact with abrasive substances (e.g. sand, dust).
- Do not expose the product to temperatures below -10 °C or above +40 °C (e.g. sauna, excessive sunlight, drying on a radiator).

## **5 Preparing product for use**

### **5.1 Processing**

### **NOTICE**

#### **Damage due to incorrect processing**

Premature wear and tear, loss of functionality due to breakage of the carbon spring

- Smooth the ground areas with an abrasive flap disc
- Centre the bores in the carbon spring
- Use a suitable counterbore to deburr the bore holes

## INFORMATION

### Orthosis alignment

Ottobock recommends checking the alignment of the orthosis using the L.A.S.A.R. Posture or L.A.S.A.R. Posture 3D and adapting it as needed.

### Fabricating a plaster model

Due to the design of the carbon spring, make sure the plaster cast is taken with slight outward rotation. Request the information on the plaster casting technique if required (Fabricating a plaster cast – KAFO (knee-ankle-foot orthosis): **646T5411\***). Raise the heel slightly by approx. 15 to 20 mm (see fig. 2) for better positioning of the carbon spring below the foot in the plaster cast. The enclosed template may be used for the sagittal positioning of the lower leg (use small template for sizes 10 to 12).

### Installing the carbon spring

**Prerequisite:** The patient-specific plaster model has been fabricated.

- 1) Pull a layer of Perlon stockinette over the plaster model and fabricate the lining using **5 mm** of Pedilin (see fig. 4). The **5 mm** material thickness can be used for all size models.
- 2) Mark the fixation zone and movement zone for the carbon spring, dividing the model into the zones **1** and **2** (see fig. 3, see fig. 5).
- 3) Mark the carbon spring below the foot, with the length extending up to approx. **20 mm** before the rollover zone of the foot (see fig. 6).
- 4) **Optional:** Shorten the carbon spring, at most up to the mark. Deburr and round the cut edges (see fig. 7).
- 5) Create the structure for the carbon spring on the lining using several layers of Pedilin and by sanding (see fig. 8).
- 6) Align the carbon spring in the foot component in the direction of the second foot ray (see fig. 9).
- 7) Contour the Pedilin structure under the carbon spring, avoiding cavities under the carbon spring (see fig. 10).
- 8) Secure the carbon spring on the model with PVA adhesive tape and pull on a layer of Perlon stockinette (see fig. 11).
- 9) For vacuum forming, use homo-polypropylene with a thickness of **4 mm** for all carbon spring sizes (see fig. 12).
- 10) Mark the course of the edge and remove the orthosis from the cast. Then separate the calf component from the foot component according to the markings (see fig. 13).
- 11) Use a ruler to mark the centre of the bores on the calf component (see fig. 14).
- 12) Mark the bores on the foot component. The proximal bore runs in front of the curve of the heel (see fig. 15).
- 13) For the holes in the orthosis shells, use a **4 mm** drill bit for an M4 thread and a **5 mm** drill bit for an M5 thread (see fig. 16).
- 14) Transfer the bore holes to the carbon spring (see fig. 17).
- 15) For the 502E3 welding nut, drill a **5 mm** hole for an M4 thread and a **6.5 mm** hole for an M5 thread (see fig. 18).
- 16) Deburr the bore with a 726S9=90x11.5 counterbore (see fig. 19).
- 17) Press in the welding nut (see fig. 20). To achieve an optimal fit of the bores, screw on the upper section on a trial basis and check the positions of the remaining bore holes. Repeat steps 14 through 16.
- 18) Assemble the carbon spring with the orthosis using the supplied screws, countersunk washers and welding nuts and secure with Loctite 241 (see fig. 21).
- 19) Ottobock recommends positioning the sagittal load line approx. **15 mm** in front of the knee compromise pivot point according to Nieter (see fig. 22).

## **6 Cleaning**

Clean the joint promptly after contact with water containing salt, chlorine or soap, or if it gets dirty (e.g. sand).

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a cloth or allow it to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

## **7 Maintenance**

### **INFORMATION**

**The product may be exposed to increased loads by the patient.**

- Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Spare parts are listed in the section "Components/design".

## **8 Disposal**

Dispose of the product in accordance with national regulations.

## **9 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **9.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### **9.2 CE conformity**

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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