17B203=* E-MAG Active

DE Gebrauchsanweisung (Benutzer)	3
EN Instructions for use (user)	9
FR Instructions d'utilisation (Utilisateur)	16
IT Istruzioni per l'uso (Utilizzatore)	23
ES Instrucciones de uso (Usuario)	30
PT Manual de utilização (Usuário)	36
NL Gebruiksaanwijzing (Gebruiker)	43
SV Bruksanvisning (Användare)	50
DA Brugsanvisning (Bruger)	56
NO Bruksanvisning (Bruker)	62
FI Käyttöohje (Käyttäjä)	69
CS Návod k použití (Uživatel)	75
SK Návod na používanie (Uživateľ)	82
RU Руководство по применению (Пользователь)	88
ZH 使用说明书 (用户)	96

1 Vorwort

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2020-04-14

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the use of the 17B203 E-MAG Active orthotic knee joint.

2 Product description

2.1 Construction and Function

The E-MAG Active is an orthotic knee joint with stance phase control which offers you a free swing phase and locks the knee joint for the stance phase prior to heel strike. The orthotic knee joint determines the position of your leg with the help of the electronics and enables the swing phase. Read the 646D1043 "Quick start guide" regarding this.


The E-MAG Active is an orthotic knee joint with a concealed locking mechanism, which means that the joint is less susceptible to dirt and also prevents clothing from becoming trapped. Not only is this gentle on your clothing but, above all, offers you enhanced safety.

The E-MAG Active has an additional locking function (PreLock) that is activated at 15° knee flexion. This secures the orthotic knee joint before the swing phase ends. The orthotic knee joint is then fully locked when it reaches the extension stop.

Two buttons allow you to electronically unlock the orthotic knee joint (once) or engage the lock (permanently). In addition, the mechanical (temporary) unlocking feature supports a larger area of application, for example when cycling. These instructions for use describe the functions and operation of the E-MAG Active as well as how to charge it and care for it.

Your O&P professional will be happy to answer any technical questions you may have.

2.2 Product overview

(see fig. 1), item	Component
1	Orthotic knee joint
2	Unlocking mechanism
3	Battery
4	Electronics

3 Safety

3.1 Explanation of warning symbols

CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Exposure of the product to unsuitable environmental conditions

Patient injury, damage, brittleness or destruction due to improper handling

- ▶ Do not expose the product to condensing ambient humidity or liquids.
- ▶ Do not expose the product to abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

⚠ CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

⚠ CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION

Improper handling

Risk of falling due to accidentally opening or engaging the lock

- ▶ Before placing weight on the orthosis, verify that the lock is in the desired position.
- ▶ Avoid accidentally opening or engaging the lock.
- ▶ Instruct the patient in the correct handling of the orthotic joints.

⚠ CAUTION

Knee not fully extended at heel strike

Risk of injury due to unlocked orthosis joint

- ▶ Fully extend the knee prior to heel strike.

⚠ CAUTION

Risk of injury due to electromagnetic interference

Malfunctions may occur in the vicinity of high-voltage power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores). This can cause the patient to fall.

Avoid proximity to strong sources of magnetic and electric interference (e.g. transformer stations, transmitters).

⚠ CAUTION

Risk of falling due to improper operation of the service switches

The A and T service switches may only be adjusted by qualified personnel.

⚠ CAUTION

Risk of falling while using the control unit

To avoid the risk of falling, you may only activate the locking and sitting function while standing and not while walking.

⚠ CAUTION

Risk of falling due to use of grease (or other lubricants)

This could contaminate the joint and jam the lock. The bearing washers and joint locking mechanism should never be lubricated.

NOTICE

Risk of damage due to incorrect environmental conditions

Solid particles and fluid in the joint, lock or electronics of the E-MAG Active may lead to malfunctioning. Do not expose the components to excessive smoke or dust in the environment, mechanical vibrations or shocks.

4 Preparation for use

INFORMATION

The E-MAG Active can only be operated with the intended 317B20 battery. It was developed especially for the E-MAG Active. Both the receptacle and lock unit are designed to fit with the form of this battery.

INFORMATION

Before putting the E-MAG Active into operation, read the instructions for the 317B20 battery and the 317L20 battery charger.

This section will describe how the E-MAG Active can be activated/deactivated and how it can be reactivated from the energy-saving mode. The section also explains how the electronic system starts and which signals the system emits during activation.

4.1 Charging the battery

INFORMATION

After charging, the E-MAG Active can be used for around 5,000 steps with one battery charge. Nevertheless, charging the battery daily (ideally at night) is recommended.

INFORMATION

In daily use, the charging time depends on the remaining capacity and maintenance of the battery. The charger can charge two batteries at the same time.

INFORMATION

The battery should always be removed when the orthosis is not in use.

INFORMATION

The LEDs on the charger must not be illuminated prior to inserting the battery. If they are illuminated, please submit the charger to an authorised Ottobock Service Centre.

First plug the battery charger into a power outlet and then insert the battery into the charger (see fig. 2). One full charging cycle takes **8 to 10 hours**. Place the charged battery into its receptacle; it should snap easily into the holder. The contacts have to face up and in. In order to do so, insert the battery into the receptacle at an angle and press it against the orthosis. You can hear the battery lock into place when it is properly inserted.

4.2 Activating/deactivating the E-MAG Active

⚠ CAUTION

Risk of falling

Missing signals during activation indicate a failure of the E-MAG Active. Remove the battery immediately and then reinsert. Please contact your specialist workshop immediately if the error occurs again.

ℹ INFORMATION

Failure to unlock the knee joint for more than one hour automatically causes the control electronics to switch to stand-by mode. To reactivate the system, simply push the top button on the control electronics.

Activation is performed by inserting the battery. Every time the E-MAG Active is activated, the electronic control unit runs a short self-test that emits the following signals:

Signal	Visual feedback (LEDs)	Audible feedback	Tone/length
System ready > knee joint locked	All LEDs flash briefly (white)		Rising
System not ready > knee joint flexed (open)	Both LEDs flash continuously (white and green)		Tone-pause-tone series

When the E-MAG Active is flexed, the tone-pause-tone series is emitted as a warning signal until the E-MAG Active is extended/locked. Deactivation is performed by removing the battery.

5 Use

This section explains how to operate the knee joint control unit by button and tells you about the warnings used by the feedback system of the E-MAG Active to inform you of malfunctions.


5.1 Control and use of 17B203 electronics

⚠ CAUTION

Risk of falling due to unlocking / locking

Both permanent and one-time unlocking as well as permanent locking of the system can result in an elevated risk of falling. The patient should not use these functions whilst walking.

The electronic control unit measures the position of the orthosis during walking. It unlocks the orthotic knee joint before toe-off and locks it prior to heel strike at the end of the swing phase so the user can step down safely. Your O&P professional adapted the electronic control unit of the E-MAG Active to your needs during the initial fitting. Two manual switching possibilities are available to you to give you full control over the E-MAG Active (see fig. 3):


5.1.1 Electronic locking (permanent)

If you feel unsafe, you can lock the joint at any time using the top button for locking the knee joint. This may be the case when walking on a slope or if you want to stand for an extended period of time without focusing on the electronic control unit. In order to do so, briefly press the upper button twice in quick succession (like double-clicking a computer mouse). This function is confirmed by an audible signal. The locking function is deactivated by pressing the upper button twice in quick succession until you hear an audible signal.

5.1.2 Electronic unlocking (once)

The patient can use the lower button to unlock the orthotic knee joint once (e.g. in order to sit down). In order to do so, briefly press the lower button twice in quick succession (like double-clicking a computer mouse). This function is **not** confirmed by an audible signal.

5.1.3 Mechanical unlocking (temporary)

For certain activities (such as cycling), it may be beneficial to temporarily unlock the orthotic knee joint. Mechanical unlocking of the orthotic knee joint is indicated by the symbol of a cyclist (switch to "🚲").

The unlocking feature needs to be cancelled in order to restore the functionality of the orthotic knee joint. In order to secure the orthotic knee joint again, the switch is operated once more in the opposite direction. The normal state of the joint is indicated by the symbol of a walking person (switch to "🚶").

5.2 Warning signals, use of 17B203 electronics

5.2.1 Battery charge level

INFORMATION

The electronics will turn off if the battery is discharged completely. The lock can be released manually.

Should the battery charge level fall below the preset charge capacity, the electronics will emit the warning signals listed here:

Signal	Visual feedback (LEDs)	Acoustic feedback	Tone/length
Battery low	Alternate flashing (white and red)		Series of high-low tones/individual signals
Battery empty			Falling/long

The warning signal is repeated every 15 minutes. Replacing the battery immediately is highly recommended. How quickly the battery is discharged depends on the ambient temperature: The colder it is, the faster the battery will be drained. The electronics will turn off if the battery power is insufficient. If this occurs, the orthotic knee joint can only be unlocked manually (see fig. 4).

5.2.2 Failure of the electronics

The E-MAG Active is not vulnerable to failure. However, should disruptions occur, the electronics recognise the failure and send a warning signal.

Signal	Visual feedback (LEDs)	Audible feedback	Tone/length
Electronics malfunction	Continuously lit (all LEDs lit, except green)		High/continuous tone

6 Cleaning

- 1) Keep your orthosis away from damp environments. Use a slightly moistened cloth to wipe the orthosis, E-MAG Active or any of the components clean from outside.
- 2) Due to the construction of the orthotic knee joint, clothing fibres cannot get into the mechanics. However, please observe the six-month maintenance intervals, because fine dust can settle in the electronics and mechanism, which can lead to malfunctions of the E-MAG Active.

7 Maintenance

Have your O&P professional check your E-MAG Active for wear and tear as well as proper functioning of the electronic components every six months.

8 Technical data & accessories

Storage temperature range	-20 °C to +70 °C (-4 °F to 176 °F)
Operating temperature range	-15 °C to +50 °C (-5 °F to 122 °F)
Relative humidity for both ranges	15% to 93%
Power supply for joint	NiMh battery, 4.8 V nominal voltage
Range of a battery	About 5,000 steps

9 Disposal

Dispose of the product in accordance with national regulations.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

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