

21Y21

DE	Gebrauchsanweisung	5
EN	Instructions for use	10
FR	Instructions d'utilisation	16
IT	Istruzioni per l'uso	21
ES	Instrucciones de uso	26
PT	Manual de utilização.....	31
NL	Gebruiksaanwijzing.....	37
SV	Bruksanvisning	42
DA	Brugsanvisning.....	47
NO	Bruksanvisning.....	52
PL	Instrukcja użytkowania.....	57
HU	Használati utasítás	63
CS	Návod k použití.....	68
RO	Instrucțiuni de utilizare.....	73
TR	Kullanma talimatı.....	79
EL	Οδηγίες χρήσης	84
RU	Руководство по применению	90
JA	取扱説明書.....	96
ZH	使用说明书.....	100

INFORMATION

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INFORMATION

Last update: 2020-04-24

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

Explanation of Symbols

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting/use.

1 Description

1.1 Indications for use

The Ottobock ClickValve is to be used **exclusively** for the prosthetic fitting of lower limb amputees.

1.2 Field of application

The ClickValve is suitable for amputees with the following requirements:

- Easy handling (threadless).
- Lateral gripping surfaces facilitating valve opening in case of restricted finger functionality.
- The valve opening is large enough to allow easy removal of donning assists.
- Easier position control of the residual limb.
- Audible feedback when the valve has been positioned correctly.

1.3 Safety instructions

Please forward the following safety instructions to your patients:

NOTICE

Risk of damage if using powder. Any kind of powder on the inner side of the socket or on the skin can restrict or disable the functionality of the valve during use. Avoid using powder in areas close to the valve.

1.4 Function

The ClickValve is a threadless one-way valve. It is opened and closed by pressing together two lateral gripping surfaces. It is used to evacuate air from the prosthetic socket. The included safety shackle (4) enables the patient to have the upper part of the valve close to the prosthesis at any time.

2 Components

21Y21 ClickValve (Fig. 1)
(1) Lower screw ring
(2) Upper screw ring
(3) Valve upper part
(4) Safety shackle

21Y21 ClickValve (Fig. 1)

(5) Air release button

(6) Vacuum forming dummy

(7) Lamination dummy

(8) 501S86=M5x12 Countersunk Head Screw

(9) 503S3 Screw

(10) Pin wrench

(11) O-ring

3 Handling

The dummy provided is used during vacuum forming or lamination and allows for precise fitting into the socket made of different materials. Please always make sure that the full dummy surface rests on the plaster positive.

3.1 Setup

NOTICE

Risk of damage to the valve. The valve is not intended for the fixation of inner socket and frame socket. It has no holding function. Only attach the valve to the test socket / to the flexible inner socket.

3.1.1 Installation preparation for sockets made of ThermoLyn

The enclosed vacuum forming dummy (6) provides for that the lower screw ring is flush on the inner side of the socket.

- Soak 99B25 Nylon Stockinette with 85F1 Plaster Parting Agent and pull it over the plaster positive.
- Secure the vacuum forming dummy to the plaster positive with the enclosed 503S3 Screw (Fig. 2).

3.1.1.1 Completion of a test socket

- After vacuum forming of the test socket of ThermoLyn, grind away the material covering the vacuum forming dummy until it is flush (Fig. 3).
- Remove the 503S3 Screw and vacuum forming dummy.

- Disassemble the valve and either apply 617H46 Bonding Agent to the contact surface of the lower screw ring (Fig. 4) or place the included O-ring (11) in this position. Position the lower screw ring in the socket (Fig. 5). Place the ring of the safety shackle on the lower screw ring.
- Tighten the upper screw ring. The pin wrench can be used for tightening the upper screw ring (Fig. 6).
- Affix the upper part of the safety shackle to the upper part of the valve. This will prevent that the upper part of the valve can get lost.

3.1.1.2 Completion of a flexible inner socket with carbon frame

The enclosed lamination dummy (7) provides for the space required for the upper screw ring in the laminated socket.

- After vacuum forming of the definitive inner socket of ThermoLyn, grind away the material covering the vacuum forming dummy until it is flush (Fig.3).
- Remove the 503S3 Screw.
- Put the lamination dummy onto the vacuum forming dummy and secure with the enclosed countersunk head screw (Fig. 7). The 5Z1 Lamination Kit with the aid of the Technical Information 2.3.5 or the 5Z4 Lamination Kit for CAT-CAM can be used for reinforcement of the frame.
- After laminating, grind away the material covering the lamination dummy until it is flush and remove the countersunk head screw.
- Remove the carbon frame with inner socket from the plaster positive.
- Disassemble the valve and either apply 617H46 Bonding Agent to the contact surface of the lower screw ring (Fig. 4) or place the included O-ring (11) in this position. Position the lower screw ring in the socket (Fig. 5). Place the ring of the safety shackle on the lower screw ring.
- Tighten the upper screw ring. The pin wrench can be used for tightening the upper screw ring (Fig. 6).
- Position the flexible inner socket in the carbon frame. Affix the upper part of the safety shackle to the upper part of the valve. This will prevent that the upper part of the valve can get lost.
- Make sure that the valve upper part can be easily installed into and removed from the valve base part. Remove more excessive material if necessary.

3.1.2 Installation preparation for laminated sockets

The enclosed vacuum forming dummy provides for that the lower screw ring is flush with the inner side of the socket.

- Pull 99B25 Nylon Stockinette over the plaster positive.
- Secure the vacuum forming dummy to the plaster positive with the enclosed 503S3 Screw (Fig.2).
- Pull 99B81 PVA Bag over the plaster positive and reinforce the laminated socket adequately.

3.1.2.1 Completion of a laminated socket without flexible inner socket

- After laminating the lamination resin socket, grind away the material covering the vacuum forming dummy until it is flush.
- Remove the 503S3 Screw and vacuum forming dummy.
- Disassemble the valve and either apply 617H46 Bonding Agent to the contact surface of the lower screw ring (Fig. 4) or place the included O-ring (11) in this position. Position the lower screw ring in the socket (Fig. 5). Place the ring of the safety shackle on the lower screw ring.
- Tighten the upper screw ring. The pin wrench can be used for tightening the upper screw ring (Fig. 6).
- Affix the upper part of the safety shackle to the upper part of the valve. This will prevent that the upper part of the valve can get lost.
- The lamination dummy is not needed here.

3.2 Instructions for use

INFORMATION

Silicone spray can be applied to the rubber sealing ring to make it easier to insert the upper part of the valve into the base part.

Place the ClickValve upper part into the valve's base part by simply pressing. An audible feedback is emitted when the valve upper part has been pressed in properly. To remove the upper part of the valve, press the air release button, then press together the two lateral gripping surfaces with the thumb and index finger and pull the upper part of the valve out of the base part.

3.3 Maintenance instructions

You may choose between two different possibilities:

- The valve should be cleaned regularly with clear, warm water.
- There are two lateral air exhaust openings in the upper part of the valve. Cover one of it airtight and blow compressed air through the other one (Fig. 8).

4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

4.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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