
8066 Genu Immobil Vario T

DE Gebrauchsanweisung	4
EN Instructions for use	8
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	15
ES Instrucciones de uso	19
PT Manual de utilização	23
NL Gebruiksaanwijzing	27
SV Bruksanvisning	31
DA Brugsanvisning	35
NO Bruksanvisning	39
FI Käyttöohje	42
PL Instrukcja użytkowania	46
HU Használati utasítás	50
CS Návod k použití	54
RO Instrucțiuni de utilizare	58
HR Upute za uporabu	62
SK Návod na používanie	65
BG Инструкция за употреба	69
TR Kullanma talimatı	73
EL Οδηγίες χρήσης	77
RU Руководство по применению	81
JA 取扱説明書	85
ZH 使用说明书	89

**Material**

Al/aluminium; CO/cotton; PA/polyamide; PES/polyester; PP/polypropylene; PU/polyurethane; Schmelzkleber/hotmelt adhesive; Stahl/steel (Nickel Legierung/nickel alloy)

Beidseitig anwendbar / Suitable for either side	Art.-Nr./Art. No. Größe / Size	Länge / Length		A Unterschenkelumfang / Lower leg circumference		B Oberschenkelumfang / Thigh circumference	
		[cm]	[inch]	[cm]	[inch]	[cm]	[inch]
	8065 Universal	55	21.7	max. 52	max. 20.5	max. 70	max. 27.6
	8066 Universal	55 - 64	21.7 - 25.2	max. 52	max. 20.5	max. 70	max. 27.6

INFORMATION

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INFORMATION

Date of last update: 2020-05-05

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the fitting and application of the 8066 Genu Immobil Vario T, a ROM knee orthosis with length-adjustable aluminium splints for immobilisation and gradual buildup of movement with flexion and extension limits at various degrees.

2 Product description

8066 Genu Immobil Vario T (see fig. 1)	
Item	Description
1	Thigh pad
2	Hook-and-loop strap
3	Strap guide loop
4	Release button
5	Orthotic joint
6	Hook strap on the back of the splint
7	Adjustment wheel for flexion and extension
8	Splint, length-adjustable
9	Calf pad

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Postoperative and posttraumatic immobilisation of the knee joint, with the option of limiting the range of motion
- Acute and subacute knee injuries
- Knee ligament injuries and grafts
- Patella dislocation

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

None.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

3.4 Effects

- Immobilises the knee joint in adjustable positions (ROM – range of motion)
 - Flexion: 0° to 120°
 - Extension: 0° to 90°
- Allows early mobilisation of the patient through gradual increases in range of motion in 10° increments
- Stabilises the capsular ligaments in the knee joint

4 Safety

4.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 General safety instructions

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.

- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

5.1 Size selection

The product is available in a universal size.

The length can be adapted to the leg using the telescoping bars (see fig. 5).

5.2 Adaptation

The initial fitting must be carried out by trained, qualified personnel.

Adapting the telescoping bars to the leg length

- > The product has been applied to the patient.
- ▶ To adapt the telescoping bars to the patient's leg length, turn the gear wheel between the thigh and lower leg sleeves.

Adjusting the flexion angle

- 1) Slide the release button up and hold it in position.
- 2) Turn the adjustment wheel to set the desired flexion angle, 1. flexion or 2. extension, on the sliders.
- 3) Let go of the release button to lock the chosen position. **CAUTION! Make sure the setting is the same for both orthotic joints.**

5.3 Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

> All straps on the orthosis are opened.

- 1) Pass the thigh and calf pads around the leg and fasten them with the hook-and-loop closures at the front of the leg.
- 2) Position the lateral orthotic joint at the height of the middle of the patella and secure it to the pad with the hook strap on the back of the splint.
- 3) Repeat this step for the medial orthotic joint. **CAUTION! Adjust the orthosis so the splints are parallel to each other and the medial and lateral orthotic joints are at the same height.**
- 4) Pass the straps through the strap guide loops. Close the straps first below and then above the patella. Make sure that the hook-and-loop sections are placed over one another.
 - The orthosis should fit firmly but is still comfortable to wear.

5.4 Trial fitting

- ▶ Check that the orthosis is fitted properly on the leg before handing over to the patient.

- ▶ Position the orthosis on the leg in such a way that the orthotic joints are aligned at the height of the middle of the patella; alternatively, adjust them in the compromise pivot point according to Nietert.

5.5 Removal

- ▶ To take off the orthosis, open the straps and pads without taking off the splints. This prevents ongoing readjustment.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

Pads

- 1) Remove the pads from the orthosis.
- 2) Wash the pads by hand in warm water at **30 °C** with standard mild detergent. Rinse thoroughly. Do not use fabric softener.
- 3) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 4) Reattach the pads to the orthosis.

7 Disposal

Dispose of the product in accordance with national regulations.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

8.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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