
50K13 Genu Arexa






DE Gebrauchsanweisung	5
EN Instructions for use	10
FR Instructions d'utilisation	15
IT Istruzioni per l'uso	20
ES Instrucciones de uso	25
PT Manual de utilização	30
NL Gebruiksaanwijzing	35
SV Bruksanvisning	41
DA Brugsanvisning	46
FI Käyttöohje	51
PL Instrukcja użytkowania	55
HU Használati utasítás	61
CS Návod k použití	66
SK Návod na používanie	71
BG Инструкция за употреба	75
TR Kullanma talimatı	81
EL Οδηγίες χρήσης	86
RU Руководство по применению	91
JA 取扱説明書	97
ZH 使用说明书	102



	Seite/ Side	Größe/ Size	Umfang/Circumference			
			A (cm)	A (inch)	B (cm)	B (inch)
	L/R	S	32-37	12,6-14,6	39-46	15,4-18,1
	L/R	M	37-41	14,6-16,1	46-53	18,1-20,9
	L/R	L	41-46	16,1-18,1	53-61	20,9-24,0
	L/R	XL	46-50	18,1-19,7	61-69	24,0-27,2
L/R	XXL	50-56	19,7-22,1	69-79	27,2-31,1	

Material	AL, SS, klettbarer Velours aus PA/velastic velour of PA, PU-Schaum/PU-foam, EVA-Schaum/EVA-foam, PA, POM, TPU, PP
-----------------	---

INFORMATION

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INFORMATION

Date of last update: 2020-05-05

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K13 Genu Arexa knee joint orthosis.

2 Product description

see fig. 1

Item	Component	Item	Component
1	Strap guide loop	6	Pad retainer
2	Pad button	7	Extension stop (square grip)
3	Clip rivet	8	Flexion stop (round grip)
4	Tibia pad	-	Straps 1 through 5
5	Condyle pad	-	Strap 6 (CI version, optional ACL/PCL version)

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Conservative treatment of injuries (e.g. after distension, partial ruptures, ruptures) of the anterior and/or posterior cruciate ligament and/or the medial and/or lateral collateral ligament of the knee joint
- Functional follow-up treatment after ligament surgery of the anterior and/or posterior cruciate ligament and/or the medial and/or lateral collateral ligament of the knee joint
- Conservative therapy of chronic complex knee joint instabilities

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

Not known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the area of the legs, e.g. in case of diabetic neuropathy.

3.4 Effects

The orthosis relieves and stabilises the knee joint.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 General safety instructions

 **CAUTION**

Use of the product on damaged skin

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product on healthy skin.

 **CAUTION**

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

 **CAUTION**

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

 **CAUTION**

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

 **NOTICE**

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Incorrect setting

Damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ The orthosis may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

5.1 Size Selection

1. Measure the circumference of the knee joint (circumference A).
2. Measure the circumference of the thigh (circumference B) **15 cm** above the centre of the patella.
3. Determine the size of the orthosis (see size chart).

5.2 Fitting and Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

5.2.1 Adjusting the Range of Motion (Flexion/Extension)

Available extension stops (square grip):	0°, 10°, 20°, 30°, 45°
Available flexion stops (round grip):	0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

INFORMATION

The joint stops (medial and lateral) must have the same degree settings.

- 1) Remove the pad retainer (and the spacer rings if applicable) (see fig. 5).
- 2) Insert the desired extension and flexion stops into the joint (see fig. 6).
→ The stops will engage with slight resistance.
- 3) **Optional:** Snap the required number of spacer rings onto the joint (see fig. 7).
- 4) Attach the pad retainer.

5.2.2 Adjusting the Strap Height

The height of straps 1 and 5 can be adjusted for knee joint relief in case of post-operative therapy and/or to adapt the position of the straps to the anatomical shape of the leg.

- 1) Remove the clip rivet of the strap holder for straps 1 and 5 on both sides.
- 2) Reposition the strap holder.
- 3) Attach the strap holder to the orthosis frame with the supplied new clip rivets.

5.2.3 Applying the ACL Version

> **Prerequisite:** The orthosis is in the delivery condition (ACL version) (see fig. 2).

> **Prerequisite:** The patient is sitting on the front edge of a chair and bends the knee at **45°**.

- 1) Position the orthosis with the middle of the joints approximately even with the top of the patella.
 - Although this puts the orthosis higher than its final position, it settles into the correct position on its own.
- 2) Fasten strap 1 (see fig. 2).
- 3) Fasten strap 2.
- 4) Engage strap 3 on the pad button with the strap guide loop and tighten the strap.
 - The strap guide loop engages audibly.
- 5) Tighten strap 4.
- 6) Push the orthosis towards the thigh.
- 7) Fasten strap 5.
- 8) Tighten the orthosis straps so that the orthosis fits firmly but is comfortable to wear.
- 9) After **15 min** of active movement (e.g. slight knee bends): if necessary, slightly retighten the straps in the specified sequence (see fig. 2, straps 1-5).

Optional: Use the included strap 6 if the knee joint requires more stabilisation. Open the hook-and-loop closure of the strap holder for strap 6 on both sides and mount the strap on the front.

5.2.4 Conversion to and application of the PCL version

> **Prerequisite:** The orthosis is in the delivery condition (ACL version) and has to be converted to the PCL version (see fig. 3).

- 1) Turn the strap holder for strap 2 from the rear to the front of the orthosis and align the strap (see fig. 3).
- 2) **Optional:** Use the included strap 6 if the knee joint requires more stabilisation. In doing so, strap 2 remains on the front of the orthosis and strap 6 is attached with hook-and-loop to the strap holder on both sides on the back of the orthosis.
- 3) Put on the orthosis and fasten the numbered straps in sequence (see fig. 3, straps 1-5, **optional** strap 6).

5.2.5 Conversion to and application of the CI version

> **Prerequisite:** The orthosis is in the delivery condition (ACL version) and has to be converted to the CI version (see fig. 4).

- 1) Attach the supplied strap 6 to the strap holder on both sides on the front of the orthosis with hook-and-loop (see fig. 4).
- 2) Put on the orthosis and fasten the numbered straps 1-6 in sequence (see fig. 4).

5.3 Delivery

NOTICE

Incorrect setting

Damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ The orthosis may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.

- ▶ Check for proper fit of the orthosis on the leg at the time of delivery.

→ The centre of the patella and the centre of the joints are approximately lined up with each other.

5.4 Cleaning

INFORMATION

- ▶ Ensure proper alignment when inserting the tibia pad. The wide side of the pad has to face up.

- 1) Remove the straps and pads from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 5) If the orthosis comes into contact with salt water or dirt, rinse with clear water and allow to air dry.

5.5 Information on Re-Use (Applies for Germany Only)

The product is approved for re-use **exclusively** in Germany. Products that are being re-used are subject to increased strain. Re-use is only permitted up to a maximum total period of use of **2 years**.

Reconditioning

Hygienic reconditioning of the product (according to the Robert-Koch Institute guidelines) is required prior to re-use.

All components subject to pronounced wear or that come into direct contact with skin are not suitable for re-use. These components are contained in the 29K110 replacement kit.

All components suitable for re-use must be checked for proper functionality and damage prior to re-use.

5.5.1 Replacing the Straps and Pads

Straps:

- 1) Release straps 3 and 4 from the frame with the quick-release closures and replace them.
- 2) Loosen straps 1, 2, 5 and 6 from the strap guide loops and replace them. Note section 5.2.2 for the adjustment of straps 1 and 5.

Strap and Frame Pads

The strap and frame pads are attached with hook-and-loop closures and are easy to replace.

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

7.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

7.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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