## 88FO\* SAFO

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	7

1

2

3

## INFORMATION

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### INFORMATION

Date of last update: 2020-05-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

### 1.1 Construction and Function

The SAFO ONE and SAFO Walk silicone ankle foot orthoses are made of silicone. They are custom made according to a plaster cast.

#### How it works

Reinforcement at the back of the foot keeps the foot from lowering during the swing phase and prevents stumbling. Compensating movements to increase the ground clearance are no longer necessary. The thin, flexible sole makes it possible to feel the ground conditions and perceive even small irregularities.

The orthosis enables a more natural gait pattern and optimum walking comfort for the user.

## 2 Intended use

### 2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

## 2.2 Indications

- Dorsiflexor weakness, with or without mild spasticity (e.g. post CVA, traumatic brain injury, multiple sclerosis, neuromuscular atrophy, peroneal paralysis or infantile cerebral palsy (ICP))
- Minor impairment of plantarflexors

Indications must be determined by the physician.

## 2.3 Contraindications

### 2.3.1 Absolute Contraindications

Not known.

### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the area of the lower leg and the foot, e.g. in case of diabetic neuropathy.

## 3 Safety

### 3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

### 3.2 General safety instructions

 **CAUTION**

#### **Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

 **CAUTION**

#### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

### NOTICE

#### Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

## 4 Handling

### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

### 4.1 Application

### ⚠ CAUTION

#### Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

### NOTICE

#### Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

### INFORMATION

The orthosis is worn directly on the skin. To reduce friction while putting on the orthosis, you can apply body powder or Ottobock Pro Comfort Gel to the inside.

1) **NOTICE! Pay special attention to the area of the sole during inspection.**

Inspect the orthosis for damage. Do not use the orthosis if it is damaged.

- 2) Open all hook-and-loop closures on the orthosis.
- 3) Pull the orthosis over the foot so it fits snugly (see fig. 1).
- 4) Fasten all hook-and-loop closures on the orthosis (see fig. 2).
- 5) Check the fit of the orthosis (see fig. 3). If the orthosis does not fit correctly, contact your orthopaedic technician.

## 4.2 Cleaning

### NOTICE

#### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

> **Recommended cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)

- 1) Clean the product with clear fresh water and a pH neutral soap.
- 2) Rinse the soap away with clear fresh water.
- 3) Dry the product with a soft cloth.
- 4) Allow to air dry in order to remove residual moisture.

## 5 Disposal

Dispose of the product in accordance with national regulations.

## 6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

Ottobock SE & Co. KGaA  
Max-Näder-Straße 15 · 37115 Duderstadt · Germany  
T +49 5527 848-0 · F +49 5527 848-3360  
healthcare@ottobock.de · www.ottobock.com