

8060 Genu Immobil 0°

DE Gebrauchsanweisung	4
EN Instructions for use	7
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	14
ES Instrucciones de uso	18
PT Manual de utilização	21
NL Gebruiksaanwijzing	25
SV Bruksanvisning	29
DA Brugsanvisning	32
NO Bruksanvisning	36
PL Instrukcja użytkowania	39
HU Használati utasítás	43
CS Návod k použití	46
HR Upute za uporabu	50
SK Návod na používanie	53
BG Инструкция за употреба	56
TR Kullanma talimatı	60
EL Οδηγίες χρήσης	64
RU Руководство по применению	67
ZH 使用说明书	71

Material Al/aluminium; CO/cotton; PA/polyamide; PP/polypropylene; PU/polyurethane

Beidseitig anwendbar / Suitable for either side	Größe / Size	Länge / Length		A Unterschenkelumfang / Lower leg circumference		B Oberschenkelumfang / Thigh circumference	
		[cm]	[inch]	[cm]	[inch]	[cm]	[inch]
	S	40	15.7	35 - 55	13.8 - 21.6	42 - 62	16.5 - 24.4
	M	50	19.7	40 - 60	15.7 - 23.6	48 - 68	18.9 - 26.7
	L	60	23.6	42 - 62	16.5 - 24.4	42 - 62	20.4 - 28.3
	XL	70	27.6	45 - 65	17.7 - 25.5	45 - 65	21.6 - 29.5

INFORMATION

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INFORMATION

Date of last update: 2020-05-18

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 8060 Genu Immobil 0° immobilisation orthosis.

2 Product description

8060 Genu Immobil 0° (see fig. 1)	
Item	Description
1	Strap guide loop
2	Orthosis shaft with padding

8060 Genu Immobil 0° (see fig. 1)	
Item	Description
3	Side parts with hook-and-loop straps, lateral and medial
4	Aluminium splints
5	Loop strap, narrow

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Postoperative and posttraumatic immobilisation of the knee joint
- Acute knee joint injuries
- Knee ligament injuries and grafts

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

None.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

3.4 Effects

- Immobilises the knee joint at 0° flexion
- Reduces collateral instability

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 General safety instructions

 **CAUTION**

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

 **CAUTION**

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

5.1 Size selection

- 1) Measure the circumference of the thigh 15 cm above the centre of the patella.
- 2) Measure the circumference of the lower leg 15 cm below the centre of the patella.
- 3) Determine the size of the orthosis (see size chart).

5.2 Adaptation

The initial fitting must be carried out by trained, qualified personnel.

Adapting the splints

The lateral and posterior aluminium splints are removable and individually adaptable to the desired flexion angle.

- 1) Remove the splint from the pocket of the orthosis.
- 2) Shape the splint by hand.
- 3) Slide the adapted splint into the pocket.
- 4) Put the orthosis on the patient for trial fitting.
- 5) Repeat the steps until the desired result is obtained.

5.3 Application

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.

▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- > The patient sits on the front edge of a chair and the knee is extended.
- > All straps on the orthosis are opened.
- 1) Pass the orthosis around the knee from behind and apply it to the leg. Centre the patella cutout over the patella.
- 2) **Optional:** For legs with a larger circumference, the side parts with straps and strap guide loops can be positioned on the orthosis shaft.
- 3) Pass the straps through the strap guide loops. Close the straps first below and then above the knee. Make sure that the hook-and-loop sections are placed over one another (see fig. 2).
- 4) **Optional:** Attach the narrow strap to the orthosis so it is centred and passes over the patella, for example to secure an additional cooling element (see fig. 3).
 - The orthosis fits firmly but is still comfortable to wear (see fig. 4).
- ▶ Check that the orthosis is fitted properly on the leg before handing over to the patient.

5.4 Removal

- 1) Open all straps one after the other.
- 2) Remove the orthosis from the leg.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

Pads

- 1) Remove the splint from the pocket of the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the orthosis in warm water at **30 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Avoid exposure to direct heat (e.g. sunlight, oven or radiator).
- 5) Reinsert the splint in the pocket of the orthosis.

7 Disposal

Dispose of the product in accordance with national regulations.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

8.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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