50S20-1 Malleo Immobil Night Splint


DE	Gebrauchsanweisung	4
EN	Instructions for use	8
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	17
ES	Instrucciones de uso	21
PT	Manual de utilização	26
NL	Gebruiksaanwijzing	30
SV	Bruksanvisning	35
DA	Brugsanvisning	39
NO	Bruksanvisning	43
FI	Käyttöohje	47
PL	Instrukcja użytkowania	51
HU	Használati utasítás	55
CS	Návod k použití	60
RO	Instrucțiuni de utilizare	64
HR	Upute za uporabu	68
TR	Kullanma talimatı	72
EL	Οδηγίες χρήσης	77
RU	Руководство по применению	81
JA	取扱説明書	86
ZH	使用说明书	90

5

	Artikelnummer/ Article Number	Größe/ Size	Schuhgröße/ Shoe Size (EU)	Schuhgröße/ Shoe Size (US)	
				♂	♀
	50S20-1	S	bis/ up to 37	5 ½	6 ½
	50S20-1	M	38–42	6–9	7–10
	50S20-1	L	ab/ over 43	9 ½	10 ½

Material	PES, PP, PA, EVA, PU Schaum/PU foam
-----------------	-------------------------------------

INFORMATION

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INFORMATION

Date of last update: 2020-05-20

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50S20-1 Malleo Immobil Night Splint immobilising orthosis.

2 Intended use

2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Tendon extensions in the calf musculature
- Functional contracted pes equinus
- Plantar fasciitis
- Tendinopathy
- Achilles tendon irritation

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

2.4 Effects

The orthosis serves to immobilise the foot/ankle joint area.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

⚠ CAUTION**Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE**Improper use and changes**

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

NOTICE**Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

4.2 Fitting and Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- > The patient is sitting on the front edge of a chair.
 - > Open all closures on the orthosis.
- 1) Position the foot in the middle of the orthosis with the knee flexed (see fig. 1).
 - 2) Fasten the closure above the ankle (see fig. 2).
 - 3) Fasten the closure across the dorsum of the foot (see fig. 3).
 - 4) Fasten the closure at the shin (see fig. 4).
 - 5) Extend the knee to check the fit of the closures.
 - 6) **Optional:** Incrementally adjust the hook-and-loop closures that support dorsal flexion to increase or reduce the angle of the foot (see fig. 5).

4.3 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

Textile component

- 1) Remove the textile pads from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the textile pads in warm water at **30 °C** with standard mild detergent. Rinse thoroughly.

- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

Orthosis frame

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

6.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

6.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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