

## 50K6-1 Agilium Patella

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	8
<b>FR</b> Instructions d'utilisation .....	12
<b>IT</b> Istruzioni per l'uso .....	17
<b>ES</b> Instrucciones de uso .....	21
<b>PT</b> Manual de utilização .....	26
<b>NL</b> Gebruiksaanwijzing .....	30
<b>SV</b> Bruksanvisning .....	35
<b>DA</b> Brugsanvisning .....	39
<b>NO</b> Bruksanvisning .....	43
<b>FI</b> Käyttöohje .....	47
<b>PL</b> Instrukcja użytkowania .....	51
<b>HU</b> Használati utasítás .....	56
<b>SK</b> Návod na používanie .....	60
<b>TR</b> Kullanma talimatı .....	64
<b>JA</b> 取扱説明書 .....	69
<b>ZH</b> 使用说明书 .....	73

	Seite/ Side*	Größe/ Size	Umfang/Circumference**			
			A (cm)	A (inch)	B (cm)	B (inch)
	L/R	XS	32-35	12.6-13.8	40-44	15.8-17.3
	L/R	S	35-38	13.8-15.0	44-48	17.3-18.9
	L/R	M	38-41	15.0-16.1	48-52	18.9-20.5
	L/R	L	41-44	16.1-17.3	52-56	20.5-22.0
	L/R	XL	44-48	17.3-18.9	56-61	22.0-24.0
	L/R	XXL	48-52	18.9-20.5	61-65	24.0-25.6

\* Links/Rechts; left/right

\*\* unterhalb/oberhalb Patellamitte / below/above mid patella

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### INFORMATION

Date of last update: 2020-05-06

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K6 Agilium Patella Pro osteoarthritis brace.

## 2 Product description

Item	Component	Item	Component
1	Donning aid	7	Dynamic tension element
2	Patella support	7.1	High-strength alignment strap
3.1	Upper ratchet hook	7.2	Elastic rubber band
3.2	Lower ratchet hook	7.3	Hook and loop closures for the dynamic tension element
4	Spring elements	8	Guide channel
5.1	Upper ratchet strap	9	Hook and loop guides
5.2	Lower ratchet strap	10	Splint
6	Ratchet strap buckle		

## 3 Intended use

### 3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 3.2 Indications

- Retropatellar osteoarthritis
- Patellofemoral misalignment
- Patellofemoral pain syndrome
- Chondropathy/chondromalacia patellae
- Condition after patella luxations and subluxations
- Condition after injury to the patellar tendon
- Anterior knee pain after surgical procedures on the knee joint (e.g. total knee replacement, lateral release surgery)

The indication must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

None known.

#### 3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries; inflammation; prominent, swollen scars; reddening and hyperthermia of the treated limb; pronounced varicose veins, especially with impaired return flow; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area where the device will be applied; sensory and circulatory disorders in the legs, e. g. associated with diabetic neuropathy.

### 3.4 Further usage restrictions

The brace is suitable for sports that may involve bodily contact with other persons (contact sports) **exclusively** in combination with the 29K129 sports cover.

### 3.5 Mechanism of Action

The dynamic re-alignment mechanism provided by the brace guides the patella within the physiologically correct patellofemoral groove. It does this without increasing the pressure caused by the tracking system in relation to the flexion angle.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
---	---

 <b>NOTICE</b>	Warning regarding possible technical damage.
--	--

### 4.2 General safety instructions

 <b>CAUTION</b>
---

#### Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

 <b>CAUTION</b>
---

#### Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

## NOTICE

### Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## NOTICE

### Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

## 5 Handling

### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

### 5.1 Selecting the Size

- 1) Measure the circumference of the thigh 15 cm above the center of the patella.
- 2) Measure the circumference of the lower leg 15 cm below the center of the patella.
- 3) Determine the size of the brace (see size chart).

### 5.2 Adaptation

#### CAUTION

#### Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

- 1) **Donning:** Hold the brace on both sides with the donning loops (1) and pull it over the knee until the patella cutout is centered on the patella.
- 2) **Positioning** (see fig. 2): Place the patient's knee in a 20°-30° flexion position. Grasp the patella carefully and guide the patella support (2) to the outer edge of the patella.
- 3) **Clasping** (see fig. 3): Hook the upper ratchet hook (3.1), with slight tension applied to the spring elements (4), onto the ratchet strap (5.1). Hook the lower ratchet hook (3.2), with slight tension applied to the spring elements (4), onto the ratchet strap (5.2).

**INFORMATION: The maximum pre-tension is specified by a mark on the ratchet hook. Apply the maximum amount of pre-tension by aligning the mark with the splint (10).**

→ The patella support (2) lies against the outer edge of the patella.

- 4) **Fastening** (see fig. 4): The patella support is held in the upper and lower guide channels (8) by the hook and loop guides (9). The hook and loop guides can be readjusted if necessary.
- 5) **Functional check** (see fig. 5): Carefully flex the knee to the degree approved by the attending physician and extend it to check that the brace is working optimally and is positioned properly.
  - **Does it require adjustment?** See 6 "Readjustment"
  - **Is the brace positioned properly?** See 7 "Shorten"

- 6) **Readjustment (optional)** (see fig. 6):  
**Readjusting at the front:** The ratchet hooks can be unhooked from the ratchet straps and readjusted.  
**Readjusting at the back:** The hook and loop closures on the dynamic tension element (7.3) can be moved medially or laterally. To do so, first move the two hook and loop straps that are positioned one over the other as necessary. Then, readjust only the upper hook and loop strap that is under tension (7.1) and limits realignment.
- 7) **Shorten** (see fig. 7): When shortening, note that every ratchet hook requires 2 ratchet straps. Unhook the ratchet hooks from the ratchet straps and shorten the ratchet straps as necessary to adapt them. Then round and de-burr the ends.  
**Tip:** Count the number of ratchet positions while the hooks are hooked in and mark the point to cut. Note that 2 ratchet straps are always required for each hook.
- 8) **Final check:** Repeat the functional check through careful flexion and extension within the range approved by the physician. The brace can be readjusted at any time using the additional 2 pairs of holes in the patella support.
- 9) **If necessary:** Use the hook and loop points to attach the additional straps around the top and bottom edges of the brace.

**5.3 Application**

**NOTICE**

**Use of a worn or damaged product**  
 Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

**NOTICE**

**Use for contact sports**  
 Impaired function of the brace components due to overloading

- ▶ Inform the patient that the brace should not be used without the sports cover for sports where there may be bodily contact with other persons.

- 1) Open the ratchet strap buckle (6).
- 2) Pull the brace on using the upper donning loops.  
 → The patella cutout is centered over the patella.
- 3) Close the ratchet strap buckle (6).

**5.4 Removal**

- 1) Open the ratchet strap buckle (6).
- 2) Take the brace off.


**5.5 Cleaning**

**NOTICE**

**Use of improper cleaning agents**  
 Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Close the ratchet strap buckle.  
**NOTICE! Do not change the setting of the dynamic tension element (7)!**
- 2) Hand wash the brace in warm water at 30 °C with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 3) Allow to air dry. Do not expose to direct heat (e.g. sunshine, stove or radiator).

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### 7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

### 7.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

### 7.3 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

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