# Lisa

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	41
<b>FR</b> Instructions d'utilisation .....	79
<b>IT</b> Istruzioni per l'uso .....	117
<b>ES</b> Instrucciones de uso .....	155
<b>PT</b> Manual de utilização .....	193
<b>NL</b> Gebruiksaanwijzing .....	231
<b>PL</b> Instrukcja użytkowania .....	269
<b>HR</b> Upute za uporabu .....	307
<b>SK</b> Návod na používanie .....	345
<b>TR</b> Kullanma talimatı .....	383
<b>RU</b> Руководство по применению .....	419
<b>KO</b> 사용 설명서 .....	459

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DE

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# Table of contents

EN

<b>1</b>	<b>Foreword</b> .....	<b>43</b>
<b>2</b>	<b>Product description</b> .....	<b>43</b>
2.1	Function .....	43
<b>3</b>	<b>Intended use</b> .....	<b>43</b>
3.1	Indications for use .....	43
3.2	Indications .....	43
3.3	Contraindications .....	44
3.3.1	Absolute Contraindications .....	44
3.3.2	Relative Contraindications .....	44
3.4	Qualification .....	44
<b>4</b>	<b>Safety</b> .....	<b>44</b>
4.1	Explanation of warning symbols .....	44
4.2	Safety instructions for assembly and adjustments .....	44
4.3	Safety instructions for use .....	44
4.4	Side effects .....	47
4.5	Further information .....	48
4.6	Nameplate and warning labels .....	48
4.6.1	Nameplate .....	48
4.6.2	Warning labels .....	48
<b>5</b>	<b>Delivery</b> .....	<b>49</b>
5.1	Scope of delivery .....	49
5.2	Options .....	49
5.3	Storage .....	49
5.3.1	Storage during daily use .....	49
5.3.2	Storage during extended disuse .....	49
<b>6</b>	<b>Preparing the product for use</b> .....	<b>49</b>
6.1	Unfolding .....	49
6.2	Installing the foot support .....	50
<b>7</b>	<b>Delivery</b> .....	<b>51</b>
7.1	Transport to the customer .....	51
7.2	Handing over the product .....	51
<b>8</b>	<b>Use</b> .....	<b>52</b>
8.1	Instructions for use .....	52
8.2	Basic functions .....	52
8.2.1	Adjusting the back support angle .....	52
8.2.2	Adjusting the foot support .....	53
8.2.3	Adjusting the seat depth .....	53
8.2.4	Operating the wheel lock .....	54
8.2.5	Lifting the rehab buggy .....	54
8.2.6	Using the cover for the back support folding scissor mechanism .....	54
8.3	Positioning belts .....	55
8.3.1	Safety instructions for qualified personnel .....	55
8.3.2	Safety instructions for the user .....	55
8.3.2.1	Warning labels .....	56
8.3.3	Lap belt .....	56
8.3.4	Five-point belt .....	57
8.4	Additional options .....	59
8.4.1	Abductor .....	59
8.4.2	Storage bag .....	60
8.4.3	Buggy board .....	60
8.4.4	Canopy .....	61
8.4.5	Grab rail .....	62
8.4.6	Hook-and-loop foot straps .....	62
8.4.7	Headrests .....	63

---

8.4.8	"Swivelling" front wheels with swivel lock.....	63
8.4.9	Padding.....	64
8.4.10	Slip sack.....	65
8.4.11	Rain cape.....	66
8.4.12	Groin strap.....	66
8.4.13	Lateral guides.....	67
8.4.14	Armrests.....	68
8.4.15	Backpack.....	69
8.5	Disassembly and transport.....	69
8.6	Use in vehicles for transporting persons with reduced mobility.....	70
8.6.1	Required accessories.....	71
8.6.2	Using the product in a vehicle.....	71
8.6.3	Restrictions for use.....	73
8.7	Care.....	73
8.7.1	Cleaning.....	73
8.7.1.1	Cleaning by hand.....	73
8.7.1.2	Cleaning belts/straps.....	74
8.7.2	Disinfection.....	74
<b>9</b>	<b>Maintenance and repair.....</b>	<b>74</b>
9.1	Maintenance.....	74
9.1.1	Maintenance intervals.....	75
9.1.2	Maintenance tasks.....	75
9.2	Repair.....	75
<b>10</b>	<b>Disposal.....</b>	<b>75</b>
10.1	Disposal information.....	75
10.2	Information on re-use.....	75
<b>11</b>	<b>Legal information.....</b>	<b>76</b>
11.1	Liability.....	76
11.2	Warranty.....	76
11.3	Service life.....	76
<b>12</b>	<b>Technical data.....</b>	<b>76</b>
<b>13</b>	<b>Appendices.....</b>	<b>77</b>
13.1	Threshold values for rehab buggies transportable by train.....	77

# 1 Foreword

## INFORMATION

Date of last update: 2020-04-01

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

## INFORMATION

- ▶ New information regarding product safety and product recalls as well as the declaration of conformity can be obtained at [oa@ottobock.com](mailto:oa@ottobock.com) or from the manufacturer's service department (see inside or outside of back cover for addresses).
- ▶ You can request this document as a PDF file at [oa@ottobock.com](mailto:oa@ottobock.com) or from the manufacturer's service department (see inside or outside of back cover for addresses). The PDF file can also be displayed in a larger size.

You have received a product that is very versatile for everyday use at home and outdoors.

In order to exclude injuries of any type, familiarise yourself with the handling, functions and intended use of the product before using it. These instructions for use provide you with the related necessary information.

### Please note the following in particular:

- All users and/or their attendants must be trained by qualified personnel in the use of the product. In particular, users and/or attendants must be informed of the residual risks with the aid of the safety notices in the instructions for use (user).
- The product was adapted to the needs of the user. Subsequent changes may be made only by qualified personnel. We recommend checking the product settings **once per year** to ensure optimal treatment over the long term. Especially for users with a changing anatomy (for example body dimensions, weight), an adjustment at least **once every six months** is recommended.
- Your product may differ from the models shown. In particular, not all the options described in these instructions for use will be installed on your product.
- The manufacturer reserves the right to make technical changes to the model described in these instructions for use.

## 2 Product description

### 2.1 Function

The product is intended exclusively for transporting one child on the seat.

The product can be used on solid ground both indoors and outdoors.

## 3 Intended use

The safe use of the product can only be ensured in case of intended use in accordance with the information contained in these instructions for use. The user is ultimately responsible for accident-free operation.

### 3.1 Indications for use

The product is intended for everyday indoor and outdoor use, by an attendant, of children with temporary or permanent limitations of the ability to walk, inability to walk or difficulty standing up.

The product is suitable for users with intact skin whose anatomy (such as body dimensions and weight) permits the intended use of the product. It is operated by an attendant.

The product may only be used with the options that are listed on the product order form. Ottobock assumes no liability for combinations with medical devices and/or accessories from other manufacturers outside the modular system.

Combinations based on a combination agreement that have been evaluated for effectiveness and safety are an exception to this.

### 3.2 Indications

- Minor to pronounced or complete restrictions of mobility

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

- None known

#### 3.3.2 Relative Contraindications

- Failure to meet physical or mental requirements

### 3.4 Qualification

Installation and adjustment tasks may only be carried out by qualified personnel. Compliance with all manufacturer specifications and all applicable legal provisions is required. Please contact the manufacturer's service department for further information (see inside or outside of rear cover for addresses).

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Safety instructions for assembly and adjustments

 **WARNING**

**Making incorrect adjustments**  
 Risk of falling, tipping over or malposition of the user due to incorrect settings

- ▶ Settings may only be changed by qualified personnel unless otherwise prescribed.
- ▶ Before testing setting changes with the user seated, all screw connections must be firmly tightened.
- ▶ Check the product for safety and functionality after changing the settings.

 **CAUTION**

**Unsecured screw connections**  
 Pinching, crushing, tipping over, falling of user due to assembly errors

- ▶ After all adjusting/readjusting work, retighten the attachment screws/nuts firmly.
- ▶ Please note that all adjustment levers are tightened by hand without the help of tools.

### 4.3 Safety instructions for use

#### Hazards during preparation for use

 **WARNING**

**Independent modification of settings**  
 Serious injuries to the user due to improper changes to the product

- ▶ Do not modify the settings established by the qualified personnel. Only the settings described in the section "Use" in these instructions for use may be adjusted independently.
- ▶ In case of problems with the settings, please contact the qualified personnel who adjusted your product.

 **WARNING**

**Improper handling of packaging materials**  
 Risk of suffocation due to neglect of the duty to supervise

- ▶ Packaging materials must be kept out of the reach of children.

## Hazards when getting in

### CAUTION

#### Wheel lock not engaged when getting in or out

Falling, tipping over of the user due to user error

- ▶ Before the user gets in or out, always engage and check the wheel lock.
- ▶ Engage the wheel lock to prevent the rehab buggy from moving on uneven ground or during transfers (e.g. into a car).

### CAUTION

#### Improper use of leg support while getting in

Tipping over, falling due to user error

- ▶ The leg support must not be used for getting in or out.
- ▶ In necessary, fold up the footplates before getting in.

## Hazards while driving

### WARNING

#### Pushing too fast

Tipping over, falling of the user due to caster wheel wobble

- ▶ The swivelling front wheels may begin to wobble at higher speeds and lead to an abrupt stop, causing the rehab buggy to tip over forwards.
- ▶ Therefore, push the rehab buggy only at normal walking speed or engage the swivel lock (see section "Use"). It is not permissible to let go of the push bar while pushing the buggy or to thrust it away from you.

### WARNING

#### Unsupervised parking

Falling out, falling of the user due to neglect of the duty to supervise

- ▶ Never leave the user unattended, not even when the positioning belts or safety belts are secured and the wheel locks are activated.

### WARNING

#### Impermissible use

Tipping over, rolling over of the user due to failure to observe specifications

- ▶ Using a rehab buggy beyond normal conditions may be dangerous.
- ▶ Please note that this product is not suitable for jogging, running, skating, etc.

### CAUTION

#### Lack of driving experience

Tipping over, falling due to errors in handling the product

- ▶ Practice on level, open ground first.
- ▶ Learn with the support of an assistant how the rehab buggy reacts to changes in the centre of gravity, e.g. downward or upward slopes, inclines or when overcoming obstacles.
- ▶ Always activate the anti-tipper (if any).

### CAUTION

#### Incorrect centre of gravity

Tipping over, overturning due to incorrect centre of gravity

- ▶ When reaching for objects, the user must not lean too far out of the seat.
- ▶ Put the back support into neutral position (0–5°) before negotiating slopes and obstacles on slopes or ramps.
- ▶ Always provide support to the rear when negotiating slopes and obstacles on slopes or ramps.
- ▶ Hold the product tightly by the push handles with both hands when driving on inclines and downgrades. Otherwise, the rehab buggy could tip over sideways when driving diagonally on a slope or roll away when driving on downgrades.
- ▶ Never park the product on slopes with the user in the seat.

**⚠ CAUTION**

**Improper loading**

Tipping over, falling of the user due to shift in centre of gravity

- ▶ Please note that any additional load may have a negative impact on the stability of the rehab buggy.
- ▶ Do not hang heavy bags or the like on the push bar/push handles.

**⚠ CAUTION**

**Incorrect centre of gravity setting**

Falling of the user due to extreme settings

- ▶ Check the standard adjustments of the rehab buggy for stability against tipping and the ergonomically correct sitting position of the user. Avoid any extreme settings.

**⚠ CAUTION**

**Risky operation**

Falling, tipping over backwards due to approaching obstacles incorrectly

- ▶ Push slowly when crossing obstacles (e.g. steps, curbs) and negotiating uphill or downhill slopes and inclines.
- ▶ Never cross obstacles at an angle. Always approach obstacles head on (at an angle of 90°).
- ▶ Raise the front wheels before crossing obstacles.
- ▶ Avoid collisions with obstacles and dropping off curbs/ledges.
- ▶ Avoid riding cross-country.

**Hazards when overcoming obstacles**

**⚠ WARNING**

**Overcoming steps and obstacles without assistance**

Tipping over, falling of the user due to failure to observe transportation instructions

- ▶ Always have accompanying persons help you negotiate steps and other obstacles.
- ▶ Use available facilities (e.g. access ramps or lifts).

**⚠ CAUTION**

**Improper lifting by attendants**

Tipping over, falling of the user due to lifting by removable components

- ▶ Attendants may lift the buggy only by welded or permanently attached components (e.g. front frame tube above the front wheels, rear axle).
- ▶ It is particularly dangerous to lift the buggy by the unfixed push handle/push bar.

**Hazard in case of broken skin**

**⚠ CAUTION**

**Skin damage**

Skin damage or pressure points due to overloading

- ▶ Check your skin for intactness before and during use of the product.
- ▶ Pay attention to diligent skin care and pressure redistribution during interruptions in using the product.
- ▶ If skin damage or other problems occur during use, stop using the product. Consult the qualified personnel.

**Hazards due to fire, heat and cold**

**⚠ WARNING**

**Seat cushions and back support pads may ignite**

Burns due to user error

- ▶ The seat and back support upholstery as well as seat cushions, padding and covers fulfil the normative requirements for flame resistance. However, they may still ignite if fire is handled improperly or negligently.
- ▶ Keep away from all ignition sources, especially lit cigarettes.

**⚠ CAUTION****Extreme temperatures**

Hypothermia or burns due to contact with components, failure of components

- ▶ Do not expose the product to any extreme temperatures (e.g. direct sunlight, sauna, extreme cold).
- ▶ Do not leave the product in the immediate vicinity of heaters.

**Hazards due to improper use of the product****⚠ WARNING****Overloading**

Severe injuries if the product tips over due to overloading, damage to the product

- ▶ Do not exceed the maximum load capacity (see the nameplate and section "Technical data").
- ▶ Please note that certain accessories and add-on components will reduce the remaining load capacity.

**⚠ WARNING****Loss of stability due to overloading**

Falling, tipping over of the user due to changed centre of gravity

- ▶ Please note that any additional load may have a negative impact on the stability of the rehab buggy.
- ▶ Never hang heavy bags or the like on the push bar.

**⚠ WARNING****Exceeding the service life**

Serious injuries due to failure to observe the manufacturer's requirements

- ▶ Using the product beyond the specified expected service life leads to increased residual risk.
- ▶ Observe the specified service life.

**⚠ CAUTION****Use of the product during diagnostic examinations and therapeutic treatment**

Impairment of the examination results or the effectiveness of treatment due to interactions of the product with devices that are used

- ▶ Make sure that examinations and treatments are carried out exclusively under the prescribed conditions.

**⚠ CAUTION****Uncontrolled driving behaviour, unexpected sounds or odours**

Falling, tipping, collision with persons or nearby objects due to defects

- ▶ If any faults, defects or other hazards that can lead to personal injury are detected, the product must be taken out of service immediately. This includes uncontrolled movements as well as sounds that are unexpected or previously not noted or odours that deviate significantly from the state of the product at the time of delivery.
- ▶ Contact the qualified personnel.

**⚠ CAUTION****Lack of maintenance**

User falling or tipping over due to lack of maintenance

- ▶ Check that the product is working correctly before each use.
- ▶ Observe the specified maintenance intervals (see Page 74, "Maintenance" section).

**NOTICE****Use under incorrect environmental conditions**

Damage to product due to corrosion or abrasion

- ▶ Do not use the product in salt water.
- ▶ Make sure that the wheel bearings are not damaged by sand or other particles.

**4.4 Side effects**

The following side effects may occur during use of the product:

- Neck, muscle and joint pain
  - Circulatory disorders, risk of pressure sores
- Contact a doctor or therapist in case of problems.

### 4.5 Further information

**INFORMATION**  
 Even in the event of compliance with all applicable guidelines and standards, alarm systems (e.g. in department stores) may respond to your product. Should this happen, remove your product from the area where the alarm was triggered.

**INFORMATION**  
 Please note that all options and add-on components will reduce the remaining load capacity.

**INFORMATION**  
 The serial number required for enquiries and orders can be found on the nameplate. For explanations of the nameplate, see the section "Nameplate" (see Page 48).

### 4.6 Nameplate and warning labels

#### 4.6.1 Nameplate

The nameplates are found on the frame.

Label	Meaning
	<b>A</b> Manufacturer's product name
	<b>B</b> CE marking
	<b>C</b> Maximum user weight (see section "Technical data")
	<b>D</b> Manufacturer information/address
	<b>E</b> Serial number <sup>1)</sup>
	<b>F</b> Manufacturing date <sup>2)</sup>
	<b>G</b> Symbol for medical device
	<b>H</b> <b>WARNING!</b> Read the instructions for use before using the product. Observe important safety-related information (e.g. warnings, precautions).
	<b>I</b> Manufacturer's reference number for the product variant
	<b>J</b> Serial number (PI) <sup>3),1)</sup>
	<b>K</b> Global Trade Item Number (DI) <sup>4)</sup>

1) YYYY = year of manufacture; WW = week of manufacture; PP = production site; XXXX = sequential production number  
 2) YYYY = year of manufacture; MM = month of manufacture; DD = day of manufacture  
 3) UDI-PI to GS1 standard; UDI = Unique Device Identifier, PI = Product Identifier  
 4) UDI-DI to GS1 standard; UDI = Unique Device Identifier, DI = Device Identifier

#### 4.6.2 Warning labels

Label	Meaning
	Fixation point to attach the product in vehicles for transporting persons with reduced mobility

## 5 Delivery

### 5.1 Scope of delivery

- Folded rehab buggy with installed foot support holders
- Individual foot support
- Options according to the order
- Instructions for use and required tools
- Instructions for use for accessories (depending on equipment)

### 5.2 Options

The standard model can be fitted to the user's personal requirements thanks to a large range of options. For use of these options: see Page 52 ff.

### 5.3 Storage

#### 5.3.1 Storage during daily use

The product should always be stored so it is protected against external influences.

#### 5.3.2 Storage during extended disuse

The product must be stored in a dry place. Specific information for extended storage: see Page 76.

It is not necessary to disassemble the product.

Maintain sufficient clearance from sources of heat. If the product is parked for an extended period of time or the tyres overheat (e.g. in the vicinity of radiators or in case of exposure to strong sunlight behind glass), the tyres may become permanently deformed.

Direct exposure to sunlight or UV radiation causes the tyres to age more quickly. As a result, the tread surface hardens and corner pieces break out of the tread.

The tyres should be replaced every **2 years** regardless of wear and tear.

## 6 Preparing the product for use

### 6.1 Unfolding

#### ⚠ CAUTION

#### **Holding on to the product incorrectly while unfolding and folding**

Pinching, crushing due to lack of caution in danger areas

- ▶ Hold the product only by the described parts.
- ▶ Never reach into the folding scissor mechanism.
- ▶ Avoid pinching yourself on movable components.

#### INFORMATION

On the topic of disassembly/transport: see Page 69.

- 1) Remove the packaging material and set down the rehab buggy (see fig. 1).
- 2) Open the transport safety device (see fig. 2).
- 3) Lift the rehab buggy by the push handles. Place one foot behind a rear wheel and tip the rehab buggy up and over the rear wheel (see fig. 3).
- 4) Move the folding lever under the seat surface (ball grip) down in a semi-circular motion and unfold the rehab buggy (see fig. 4).

**NOTICE! Make sure that the slides are resting completely on the stops.**

- 5) Push the folding lever on the back support (ball grip) downwards to the stop (see fig. 5).

## 6.2 Installing the foot support

### Initial installation

The foot support is not assembled with the two foot plate bars on delivery:

- 1) Pull the pivot bearings up and turn them to the front (see fig. 6, item 1). Allow the pivot bearings to engage downward into place.
- 2) Loosen the clamping screws (see fig. 6, item 2).
- 3) From below, slide the foot plate bars to the same height in the clamping parts of the pivot bearings (see fig. 6, item 3) and tighten the clamping screws slightly.
- 4) Align the foot plate bars so that the foot plate can be joined together with the remaining foot plate bar (see fig. 7).
- 5) Fold down the foot plate.
- 6) Slide the notch in the foot plate over the holder on the foot plate bar (see fig. 8).
- 7) Turn the twistlock a quarter of a turn so that the bolt can engage in the foot plate (see fig. 9).

**NOTICE! Make sure the foot plate has securely locked by giving it a short tug.**


8) Firmly tighten the clamping screws (see fig. 6, item 2).

**INFORMATION: Before tightening the screw connections, make sure that the foot plate bars are at the same height.**

#### Intended use in everyday life

- 1) **If necessary:** Fold the pivot bearings to the front. To do so, pull the pivot bearings up (see fig. 6, item 1) and turn them to the front. Allow the pivot bearings to engage downward into place.
- 2) Fold down the foot plate.
- 3) Slide the notch in the foot plate over the holder on the remaining foot plate bar (see fig. 8).
- 4) Turn the twistlock a quarter of a turn so that the bolt can engage in the foot plate (see fig. 9).

**NOTICE! Make sure the foot plate has securely locked by giving it a short tug.**

## 7 Delivery

### 7.1 Transport to the customer

#### NOTICE

##### Use of unsuitable packaging

Damage to the product caused by transportation using incorrect packaging

- ▶ Use only the original packaging for delivery of the product.

### 7.2 Handing over the product

#### ⚠ WARNING

##### Lack of instruction

Tipping over, falling of the user due to lack of knowledge

- ▶ Instruct the user or the attendant in the proper use of the product when handing it over.

The following steps must be performed for the safe delivery of the product:

- A sitting test has to be conducted with the user of the product. Pay special attention to proper positioning according to medical considerations.
- The user and any attendants must be instructed in the safe use of the product. In particular, the enclosed instructions for use are to be used.
- The instructions for use must be given to the user/attendant upon delivery of the wheelchair.
- **Depending on equipment:** The supplied instructions for use for accessories have to be handed over in addition.

## 8 Use

### 8.1 Instructions for use

- The products in this series generally meet the minimum technical requirements for wheelchairs transportable by train. However, please note that it is possible that not every individual product will fulfil the minimum requirements due to the variation in designs (see Page 77).

### 8.2 Basic functions

#### 8.2.1 Adjusting the back support angle

##### **CAUTION**


##### **Change in centre of gravity following back angle adjustment**

Falling over, falling out of the user due to adjustment errors

- ▶ After each modification of the back angle, check the tilt stability. To do so, hold the product with both hands by the push handle/push bar.

- 1) Engage the wheel lock (see Page 54).
- 2) Secure the back support with one hand (see fig. 10).
- 3) Pull the release cord a few centimetres upwards (see fig. 10) to release the angle adjustment of the back support (see fig. 11, item 1).
- 4) Move the back support to the desired angle while keeping the release cord pulled (see fig. 12, see fig. 13).
- 5) Let go of the release cord again so that the bars on the left and right side lock into place (see fig. 11).

**CAUTION! Check that the lock is firmly engaged by pulling on the back support.**

## 8.2.2 Adjusting the foot support

### Adjusting the foot support height

- 1) Loosen the clamping screws in the pivot bearings on the left and right side (see fig. 14, item 1).
- 2) Slide the foot plate bars with the foot plate up or down (see fig. 14, item 2).
- 3) Firmly retighten the clamping screws (see fig. 14, item 1).



**INFORMATION: The pivot bearings on the front frame must be lowered for lower leg lengths of 340 mm or greater.**

- 4) To do so, loosen the mounting screws on the front frame on both sides and screw them into a lower hole (see fig. 14, item 3).

→ The pivot bearings can now engage into a lower position.


### Adjusting the foot support depth

- 1) Fold up the foot plate.
- 2) Loosen the circlip on the twistlock (see fig. 15, item 1) and pull the twistlock out of the foot plate bar.
- 3) Insert the twistlock into the second hole and secure it with the circlip.
- 4) Remove the screw from the foot plate bar on the other side of the foot support and slide the foot plate to the second hole.
- 5) Tighten the screw firmly in the hole (not illustrated).

## 8.2.3 Adjusting the seat depth

- 1) Loosen and remove the screw connections on both sides of the seat surface (see fig. 16, item 1).
- 2) Adjust the desired seat depth, reinsert the screws and firmly tighten them (see fig. 17).

### 8.2.4 Operating the wheel lock

#### **⚠ WARNING**

#### **Insufficient braking**

Accident, falling of the user due to loss of braking function

► Verify the functionality of the wheel lock prior to every use.

- 1) **Engage the wheel lock:** Press the release lever down with your foot (see fig. 18).
- 2) **Release the wheel lock:** Pull the release lever up with the tip of your foot (see fig. 19).

### 8.2.5 Lifting the rehab buggy

- 1) Screw the tip-assist to the rear frame on the left or right side.
- 2) To tilt the rehab buggy, place one foot on the tip-assist and pull the buggy towards you (see fig. 20).

### 8.2.6 Using the cover for the back support folding scissor mechanism

The folding scissor mechanism of the back support is covered with a pocket-type cover (see fig. 21). The pocket can only carry soft or small objects as these come into close proximity with the frame tubes when the back support is at maximum inclination.

## 8.3 Positioning belts

### 8.3.1 Safety instructions for qualified personnel

#### **⚠ CAUTION**

##### **Incorrect approach to the adjustment process**

Injuries, malpositions, user discomfort due to adjustment errors

- ▶ The belt system may only be installed by qualified personnel.
- ▶ The qualified personnel is responsible for the individual positioning and fitting of the belt system.
- ▶ Adjusting the belt system too tightly may lead to unnecessary pain or user discomfort.
- ▶ Adjusting the belt system too loosely can cause the user to slide into a dangerous position. In addition, the fastening snaps could open unintentionally if they slide against hard parts of clothing (e.g. buttons).
- ▶ For proper installation and adjustment, follow the instructions for use included with the product.

#### **⚠ CAUTION**

##### **Lack of instruction**

Injuries, malpositions, illness of the user due to incorrect information

- ▶ The qualified personnel is responsible for making sure that the user and/or attendant/nursing staff has understood the proper adjustment, use, maintenance and care of the belt system.
- ▶ In particular, ensure that the user and/or attendant/care staff knows how to quickly loosen and open the product to avoid delays in case of emergency.

### 8.3.2 Safety instructions for the user

#### **⚠ WARNING**

##### **Use of the belt system as a passenger restraint system in vehicles for transporting persons with reduced mobility is forbidden**

Serious injuries due to improper handling of the product

- ▶ Under no circumstances may the belts and positioning aids offered with the product be used as part of a passenger restraint system in vehicles for transporting persons with reduced mobility.
- ▶ Note that the belts and positioning aids offered with the product are only intended to help support the user sitting in the product.

#### **⚠ WARNING**

##### **Incorrect application of the belt system/positioning aid**

Strangulation, suffocation due to incorrectly placed positioning belts/positioning aids

- ▶ Follow the instructions for proper application so that the user cannot slip into a position that would endanger him or her.
- ▶ Note that the user's back must be up against the back support pad when securing the belt system/positioning aid and during use of the product.
- ▶ Do not use positioning belts or positioning aids if the user is able to sit without them.

**⚠ CAUTION**

**Improper adjustments**

Injuries, malpositions, illness of the user due to adjustment changes

- ▶ The belt system is an important part of an individual seating unit/seating solution. Do not modify the installation position and basic settings established by the qualified personnel.
- ▶ In case of problems with these adjustments (such as an unsatisfactory sitting position), promptly contact the qualified personnel who fitted the product.
- ▶ **Immediately** consult the qualified personnel if you detect signs of discomfort or fear when using the product.
- ▶ Have the basic settings of the belt system checked regularly. Adjustments may be required due to the growth of the user or because of changes in the course of the disease.

**⚠ CAUTION**

**Improper use**

Falls, user falling out due to improper use

- ▶ The positioning system for the upper body must be put on when getting into the product and used in combination with a positioning lap belt at all times while using the product.
- ▶ Only open the positioning system for the upper body and the lap belt when the user is ready to get out of the product.
- ▶ Do not leave the user unsupervised if the cognitive abilities of the user could lead to unintentional opening of the positioning system.

**8.3.2.1 Warning labels**

Label	Meaning
Posture Support Device Not suitable for use as a passenger restraint in transport.	Positioning system Not suitable for use as a personal restraint system in vehicles for transporting persons with reduced mobility.

**8.3.3 Lap belt**

The standard equipment of the rehab buggy includes a simple lap belt. It provides additional support for the pelvis and keeps the user from sliding out of the seat.


**Operating instructions**

- Place the user in an upright, **90°** seated position (if physiologically possible).
- Ensure that the back is up against the back support pad (if physiologically possible).
- The lap belt should run in the direction of the thighs in front of the pelvic bones if possible (not illustrated).
- The lap belt should be at an angle of about **60°** to the seat surface. The belt strap should run across the thighs in front of the pelvic bones (see fig. 22).

**Possible errors**

- The lap belt is positioned above the user’s pelvis in the area of the soft part of the stomach.
- The user does not sit upright in the seat.
- If the lap belt is too loose, the user can shift/slide out to the front.
- During the installation/adjustment, the lap belt is routed over parts of the seating system (e.g. over arm supports or seat pads). This causes the lap belt to lose its positioning and retaining function.


### Installing and adjusting the belt

- 1) Thread the belt in front of the back support upholstery through the two eyebolts at the lower end of the back support (see fig. 23, item 1).
- 2) Adjust the belt ends at the belt closure. Do so by loosening the belt ends in the buckle, adjusting the length and pulling them tight.

**INFORMATION: In the case of a two-piece lap belt, the belt length is adjusted with the hook-and-loop fasteners on the belt ends. The belt ends must overlap by at least 150 mm.**

### Using the belt

- 1) To unfasten the belt, press the quick-release buckle from above and below. Pull the quick-release buckle out of the closure (see fig. 23).
- 2) Seat the user in the seat.
- 3) Fasten the belt by means of the quick-release buckle until the closure audibly engages.

### 8.3.4 Five-point belt

The rehab buggy can be equipped with a five-point belt. It additionally stabilises the pelvis and shoulders and prevents the user from sliding out of the seat. The buckle makes it easy to get in and out.

The optional five-point belt consists of shoulder harnesses with shoulder pads, the lap belt piece and the centre strap.

#### Operating instructions

- Place the user in an upright, 90° seated position (if physiologically possible).
- Ensure that the back is up against the back support pad (if physiologically possible).

#### Possible errors



- The user does not sit upright in the seat.
- If the positioning belt is too loose, the user can shift/slide out to the front.
- During the installation/adjustment, the positioning belt is routed over parts of the seating system (e.g. over arm supports or seat pads). This causes the lap belt to lose its retaining function.

**Installing and adjusting the shoulder harness without padding (see fig. 24)**

- 1) Pass the shoulder harness behind the padding around the back support and fasten the closure buckle (see fig. 25).
- 2) Adjust the length on the belt closures. Do so by loosening the belt ends in the buckles, adjusting the length and pulling them tight.

**Installing and adjusting the shoulder harness with padding (see fig. 26)**

- 1) At the upper edge of the back support, open the hook-and-loop closure on the padding (see fig. 27, item 1). Push the padding slightly forwards.
- 2) Open the zips on the top and bottom of the padding (see fig. 27, item 2) and pass the shoulder harnesses through the gap between the zips from the front.
- 3) Adjust the zip position to set the desired shoulder height.
- 4) Pass the shoulder harness behind the padding around the back support and fasten the closure buckle (see fig. 25).
- 5) At the top edge of the back support, fasten the hook-and-loop closure on the padding (see fig. 27, item 1).
- 6) Adjust the length on the belt closures. Do so by loosening the belt ends in the buckles, adjusting the length and pulling them tight.

**Installing and adjusting the lap belt**

- 1) Thread the lap belt of the five-point belt in front of the back support upholstery through the two eyebolts at the lower end of the back support (see fig. 28).
- 2) Adjust the length of the lap belt at the buckle.

**Installing and adjusting the centre belt**

- 1) Attach the centre belt of the five-point belt to the loop on the seat surface (see fig. 29).
- 2) Adjust the length of the centre belt using the triglide.

**INFORMATION: The centre belt and lap belt must not be removed before folding up the rehab folding buggy.**

### Using the five-point belt

- 1) To open the belt, press the release button and undo the buckle (see fig. 30).
- 2) Seat the user in the seat.
- 3) Bring all of the belts together in the buckle. The buckle can only be closed if the strap tabs are inserted in the correct order (see fig. 30).

## 8.4 Additional options

### ⚠ CAUTION

#### Improper installation of options/accessories

Falling of the user due to incorrect installation, damage of option components

- ▶ All option/accessory components must be installed by the qualified personnel unless they were pre-installed by the manufacturer. A few notes on installation are provided here for the benefit of the qualified personnel.
- ▶ Always mount the option/accessory components such that the stability of the product is maintained.
- ▶ Inform the attendant that installing and using some of the option components reduces the maximum load.
- ▶ Inform the attendant that installing some of the option components (e.g. storage bag, canopy) may restrict the ability to fold the product. These components should be removed prior to folding if necessary.


### 8.4.1 Abductor

The seating unit can be equipped with a pommel that ensures that the legs are abducted (spread to the sides) from the centre of the body for certain indications.

### ℹ INFORMATION

For use in a vehicle for transporting persons with reduced mobility, the pommel has to be folded away to the front and secured in this position.

- 1) Screw the pommel to the folding scissor mechanism under the seat surface (see fig. 31).
- 2) To make getting in and out easier, you can fold the pommel forward (see fig. 32). Fold the pommel towards the seat surface, press it down and keep it pressed down to fold it forward.

### 8.4.2 Storage bag

#### NOTICE

##### Incorrect attachment

Damage to product due to incorrect installation

- ▶ During the attachment, make sure that the product cannot grind against the wheels.

#### INFORMATION

For use in a vehicle for transporting persons with reduced mobility, the storage bag can remain on the product provided it is emptied and the contents are securely stowed.

- 1) Fasten the storage bag to the frame under the seat (see fig. 33, see fig. 34).
- 2) Fold the storage bag together with the rehab buggy.

**NOTICE! The storage bag may not be loaded with more than 5 kg. Empty the storage bag before folding up the rehab buggy!**

### 8.4.3 Buggy board

#### CAUTION

##### Improper use of the buggy board

Falling, tipping, injuries due to user error

- ▶ Do not exceed the max. load of the buggy board. This is **20 kg**.
- ▶ Do not use the buggy board as a scooter.

#### NOTICE

##### Damage during folding

Damage to the product due to user error

- ▶ Do not fold the rehab buggy up when the buggy board is installed.

**INFORMATION**

All information regarding the attachment, adjustment and use of the buggy board are contained in the manufacturer's instructions for use enclosed with the product.

The buggy board (see fig. 35) allows children to ride safely while standing on their own feet, has spring loaded casters and offers sufficient standing room. It is attached to the back frame and can be folded up or detached using the quick-disconnect units.

**8.4.4 Canopy****⚠ CAUTION****Incorrect use of the canopy**

Falling, tipping of the user due to uncontrolled rolling away

- ▶ An unfolded canopy is a large surface that can be caught by the wind. Close the canopy every time you park the rehab buggy. Never leave the user unattended.

**NOTICE****Improper use of the canopy**

Damage to the fabric due to user error

- ▶ Always remove the canopy before folding the rehab buggy so that it does not get folded up as well.
- ▶ Remove the canopy from the rehab buggy if it is stored for a long period to protect it from damage.

The rehab buggy can be equipped with a canopy to protect the user from sun or rain.

**Mounting the canopy**

- 1) Move the back support to an upright position.
- 2) Attach the canopy to the push handles with the provided snap-in clips (see fig. 36).
- 3) Pull the canopy over the back support and close the hook-and-loop closures around the sides of the push handles.
- 4) In case of rain, pull the supplied transparent rain cover over the canopy and down over the foot support (see fig. 37). The push handles remain uncovered.

**Removing the canopy**

- 1) Move the back support to an upright position.
- 2) Open the hook-and-loop closures of the canopy on the sides around the push handles.
- 3) Remove the snap-in clips from the push handles.
- 4) Remove the canopy.

### 8.4.5 Grab rail

For certain indications, the rehab buggy can be equipped with a grab rail that prevents the user from falling forwards out of the seat.

The grab rail can be removed to make getting into the product easier.

#### Mounting the grab rail

- 1) Screw the receivers for the grab rail to the seat tubes (see fig. 38).
  - 2) Slide the pad onto the grab rail.
  - 3) Press in the tripod springs on the left and right sides of the grab rail.
  - 4) Insert the two ends of the grab rail into the receivers to the stop.
- The arch of the grab rail should point towards the seat.

#### Removing the grab rail

- 1) Press and hold the tripod springs of the grab rail.
- 2) Pull the grab rail out of the receivers (see fig. 39).

### 8.4.6 Hook-and-loop foot straps

#### INFORMATION

For use in a vehicle for transporting persons with reduced mobility, the feet have to be taken out of the hook-and-loop foot straps. In addition, the opened hook-and-loop foot straps should be secured on the foot plate next to the feet.

The seating unit can be equipped with hook-and-loop foot straps that secure the child's feet firmly to the foot plate for certain indications.

The foot straps are attached to the heel part of the foot plate with two rivets or screws positioned beside one another (see fig. 40). Corresponding holes must be made in the foot plate.

In addition, foot straps can also be used for the forefoot area.

#### Using the hook-and-loop foot strap

- 1) Open the hook-and-loop closure.

- 2) Place the foot on the foot plate.
- 3) Pass the strap tightly around the shoe and fully fasten the hook-and-loop closure.

### 8.4.7 Headrests

The rehab buggy has one head support model for use with or without padding respectively.

#### Using the head support without padding

- 1) Pass the two attachment straps of the head support around the back support on the back side (see fig. 41).
- 2) Connect the two closures on the back side and pull the straps tight (not illustrated).

#### Using the head support with padding

- 1) Find a position that is appropriate for the padding (see fig. 42) and screw the head support to the middle section of the back support pad at the desired height (see fig. 50, item 1).
- 2) Remove the head support and the padding prior to folding up the rehab buggy.

### 8.4.8 "Swivelling" front wheels with swivel lock

The rehab buggy can be equipped with swivelling front wheels to make steering easier. If needed, the caster wheels can be locked with the integrated swivel lock. This function increases directional stability on uneven ground.

The swivel lock is installed by the specialist dealer or manufacturer prior to delivery.

- 1) **Activate swivel lock:** Turn the locking pin on the two front wheels by **90°** and allow it to engage in the hole of the caster fork (see fig. 43).
- 2) **Deactivate swivel lock:** Pull the locking pin out of the caster forks on the two front wheels and turn the locking pin by **90°** so that it remains in the open position (see fig. 44).

### 8.4.9 Padding

The padding enhances the user's sitting comfort and also reduces the seat surface in the first months.

#### Attaching simplified padding

- 1) Put the seat pad onto the seat upholstery.
- 2) Guide the flap on the back end of the seat pad to the back, insert it between the seat and back support upholstery and fasten it under the seat pad (see fig. 45).
- 3) Move the lateral pads of the seat pad to an upright position on the right and left sides.
- 4) Put on the back support pad and close the hook-and-loop connection to the seat pad on the lower end across the entire width (see fig. 46).


#### Attaching seat pad with multifunctional back

- 1) Put the seat pad onto the seat upholstery.
- 2) Guide the flap on the back end of the seat pad to the back, insert it between the seat and back support upholstery and fasten it under the seat pad (see fig. 45).
- 3) Put on the back support pad and lean it against the back support upholstery. Close the toggle fasteners on the right and left sides between the back support pad and seat pad (see fig. 47).
- 4) **For first-time installation:** Pull off the adhesive pad on the outside of the lateral pad and attach it to the outside of the frame (see fig. 48).
- 5) Pull the lateral pads to the outside around the front frame on the right and left sides (and around the arm supports as well if applicable) and fasten them.
- 6) Properly adjust the back support pad to the back support upholstery and close the hook-and-loop connections to the lateral pads (see fig. 49).

**INFORMATION:** The back support pad is equipped with zip closures with double sliders to create openings for the head support, the shoulder harnesses of the five-point belt or the groin strap (see fig. 50, item 1).

#### Optional: using a compensation pad

Compensation pads can be used to adjust the seat width if necessary.

- 1) Pull the lateral pads on the left and right sides of the front frame off the hook-and-loop closure (see fig. 48).
- 2) Fasten the compensation pads to the lateral pads (see fig. 51).
- 3) Pull the lateral pads to the outside around the front frame on the right and left sides again (and around the arm supports as well if applicable). Close the hook-and-loop closures (see fig. 48).

#### 8.4.10 Slip sack


##### NOTICE

##### Incorrect attachment

Damage to the product

- ▶ During the attachment, make sure that the product cannot grind against the wheels.

- 1) When needed, make a cut in the slip sack for putting through the lap belt, the pommel and the five-point belt or the groin strap. Use the provided positions with the applied leather applications for this purpose.  
**INFORMATION: We recommend that you first make cutting marks before cutting the leather applications. You should then pre-punch the start and end of the cutting mark with a hole punch and finally make a cut that connects the two holes.**
- 2) Place the slip sack in the rehab buggy and tie the upper fastening cords to the push handles (see fig. 52, item 1).
- 3) Attach the lower fastening cords to the lower folding levers of the rear frame (see fig. 53, item 1).

#### 8.4.11 Rain cape

- 1) Pull the hood over the child's head with the zip open (see fig. 54).
- 2) Unfold the upper edge of the rain cape over the back support and the lower edge around the foot support.

#### 8.4.12 Groin strap

The seating unit can be equipped with a groin strap that positions the pelvis in the seat for certain indications.

##### Putting on the groin strap

- 1) Position the groin strap with the quick-release buckles on the thigh belts pointing downwards on the seat surface (see fig. 55).
- 2) Insert the attachment tabs of the groin strap to the back under the back support, wrap them around the seat edge and fasten them (see fig. 56).  
**NOTICE! If padding is being used, fasten the tabs to the padding.**
- 3) Pass a tensioning strap for the leg straps through the eyebolt and the end of the loop of the tensioning strap on the right and left sides respectively (see fig. 57).
- 4) Tighten the tensioning belts and thread them into the clips (see fig. 58).
- 5) Place the user on the groin strap and thread the leg straps between the legs and over the thighs to the outside (see fig. 59).  
**NOTICE! Make sure that the part of the leg strap on which the buckle is fastened is lying on the pelvic bone.**
- 6) Close the two clips at the leg straps (see fig. 60).
- 7) The length of the leg straps can be adjusted by pulling on the free belt ends in the belt adjustment if needed.

#### 8.4.13 Lateral guides

The lateral guide pads are filled with 10 mm-thick foam inserts which can be variably inserted to enable seat width reductions of **40 mm**, **60 mm** or **80 mm** (see fig. 61).

##### Using the lateral guides

- 1) Push the lateral guide pads with the insert onto the connection point between the back support and frame on both sides (see fig. 62).
- 2) Press the insert from above onto the connection point between the seat tube and front frame on both sides.

#### 8.4.14 Armrests

##### **⚠ CAUTION**

##### **Improper lifting by attendants**

User tipping over or falling due to lifting by removable components

- ▶ The arm supports must not be used to lift the product.

The seating unit can be equipped with arm supports that support the arms for certain indications.


The arm supports are screwed on to the left and right side of the front frame (see fig. 63). The height of the arm supports can be set in two different positions.

##### **Installing the arm supports**

- 1) Hold the arm support adapter against the frame tube.
- 2) Insert each of the included bushings with collars through the arm support adapters and selected frame tube holes from the inside (see fig. 64).
- 3) Insert each of the included bushings without collars through the arm support adapters and selected frame tube holes from the outside (see fig. 65).
- 4) Screw on the arm supports (see fig. 66). Make sure that the screw head is on the outside.

##### **Using the arm supports**

- 1) Fold up the arm support before folding the rehab buggy.
- 2) Fold down the arm support before using the rehab buggy.

### 8.4.15 Backpack

#### ⚠ CAUTION

#### Improper installation/use of the backpack

Falling of the user due to incorrect installation, damage to the product

- ▶ Always mount the backpack such that the stability of the product is maintained.
- ▶ Please note that the maximum load of the overall product must not be exceeded after attaching the backpack as well (see Page 76).

- 1) Hang the top attachment loops of the backpack over the back tubes (see fig. 67, item 1).
- 2) Open the bottom attachment loops of the backpack and guide them around the tube of the folding scissor mechanism at the indicated position (see fig. 67, item 2).
- 3) Close the attachment loops.
- 4) Make sure that the backpack is up against the back support without hanging down. Set the attachment loops to a correspondingly short length to ensure this if necessary.

### 8.5 Disassembly and transport

#### ⚠ CAUTION

#### Holding on to the product incorrectly while unfolding and folding

Pinching, crushing due to lack of caution in danger areas

- ▶ Hold the product only by the described parts.
- ▶ Never reach into the folding scissor mechanism.
- ▶ Avoid pinching yourself on movable components.

**NOTICE****Deformation while folded**

Damage to the product, problems with unfolding due to improper loading

- ▶ Never place heavy objects onto the folded product.
- ▶ Make sure that objects placed in the storage bag cannot be crushed while folding.

- 1) Remove the following accessory components from the rehab buggy:  
Buggy board (see Page 60), canopy (see Page 61), grab rail (see Page 62), head support (see Page 63), padding (see Page 64), rain cape (see Page 66), slip sack (see Page 65), groin strap (see Page 66) and lateral guides (see Page 67).  
If padding is being used, the head support, the shoulder harnesses of the five-point belt and the groin strap may remain attached to the padding when removing it.
- 2) Turn the twistlock on the foot support a quarter of a turn so that the bolt disengages from the foot plate (see fig. 8).
- 3) Fold the foot plate up and away to the side along with the fixed foot plate bar.
- 4) Fold the remaining foot plate bar away to the side.
- 5) **If present:** Fold up the arm supports (see Page 68).
- 6) Release the wheel lock (see Page 54).
- 7) Move the back support to an upright position and pull the lock open (see Page 52).
- 8) Pull the folding lever on the back support (ball grip) upwards (see fig. 5).
- 9) Move the folding lever on the seat surface (ball grip) up in a semi-circular motion and fold the rehab buggy (see fig. 4).
- 10) Set down the rehab buggy and engage the transport safety device (see fig. 2).

**8.6 Use in vehicles for transporting persons with reduced mobility****⚠ WARNING****Use in vehicles for transporting persons with reduced mobility**

Serious injuries in case of accidents due to user error

- ▶ Always use the seats and personal restraint systems in the vehicle for transporting persons with reduced mobility first. This is the only way to ensure optimum protection of passengers in the event of an accident.
- ▶ If the product is to be used as a seat in a vehicle for transporting persons with reduced mobility, the safety elements offered by the manufacturer and appropriate fastening and personal restraint systems must be used. For more information, please refer to our brochure with the order number 646D158.
- ▶ Never transport more than one person in the product.
- ▶ Only use the product in a vehicle for transporting persons with reduced mobility with the back support in a vertical position.
- ▶ Observe the limitations regarding installed options (see Page 73).

**⚠ WARNING****Use of the belt system as a passenger restraint system in vehicles for transporting persons with reduced mobility is forbidden**

Serious injuries due to improper handling of the product

- ▶ Under no circumstances may the belts and positioning aids offered with the product be used as part of a passenger restraint system in vehicles for transporting persons with reduced mobility.
- ▶ Note that the belts and positioning aids offered with the product are only intended to help support the user sitting in the product.

**⚠ WARNING****Prohibited transportation of the passenger with activated back angle adjustment**

Loss of safe restraint in the product due to user error

- ▶ Ensure the passenger is seated in a nearly upright position during transport.
- ▶ If a back angle adjustment mechanism is installed, move the backrest to a nearly upright position prior to travel.
- ▶ Check the locking mechanism on both sides.

**INFORMATION**

For transporting children weighing up to **36 kg**, Ottobock strongly recommends using a child restraint system (child car seat) tested in accordance with ECE R44/04 or higher, or a restraint system developed specially for children with disabilities and which conforms to ECE R44/04. Always observe the regulations for your country or state.

The product has been tested by the manufacturer according to ISO 7176-19 and may be used as a seat in vehicles for transporting persons with reduced mobility subject to the conditions defined below.

The product must be sufficiently secured during transport in vehicles for transporting persons with reduced mobility. The illustrations that follow show an example for anchoring in a motor vehicle.

The manufacturer is not responsible for the fastening systems that are used. Ensure that only fastening systems that meet the applicable legal requirements and are designed for the overall weight of the product including the user are used.

The transport weight of the person to be transported in a vehicle for transporting persons with reduced mobility corresponds to the maximum permissible user weight (see Page 76).

**8.6.1 Required accessories**

The eyebolts on the front and rear frames are already mounted ex works.

**8.6.2 Using the product in a vehicle****⚠ WARNING****Positioning in vehicles for transporting persons with reduced mobility**

Serious injuries in case of accidents due to user error

- ▶ Positioning of the product in vehicles for transporting persons with reduced mobility may only be performed by the qualified personnel.
- ▶ The product must always face forwards when it is used as a seat in a vehicle for transporting persons with reduced mobility.
- ▶ Instruct the qualified personnel regarding the mounting points on your product described below.

**⚠ WARNING****Inadequate transportation safety**

Loss of safe restraint due to failure to observe transportation instructions


- ▶ Observe the following instructions for correct transport safety in the vehicle for transporting persons with reduced mobility.
- ▶ If necessary, instruct the qualified personnel on the following information.

**Securing the product in the vehicle for transporting persons with reduced mobility**

The product is attached to the eyebolts (fixation points) using four vehicle side restraint belts.

The fixation points are marked with stickers. The stickers show where the vehicle side restraint belts have to be engaged.

- 1) Position the product in the vehicle for transporting persons with reduced mobility. For more information, please refer to section 5 in the brochure "Transporting persons with reduced mobility", order number 646D158.
- 2) Engage the hook of the vehicle side restraint belts at the front and rear respectively in the eyebolts and tighten them (front/rear: see fig. 68; rear: see fig. 69).
- 3) Tighten the vehicle side restraint belts at the front and rear as firmly as possible.
- 4) Remove all prohibited components from the product (see Page 73).
  - The product showing the correct positioning of the attachment straps (see fig. 68).

### Information on correct transport safety of the user in the vehicle for transporting persons with reduced mobility

- Using the personal restraint system of the vehicle for transporting persons with reduced mobility is required. The three-point personal restraint system is attached in the vehicle for transporting persons with reduced mobility as follows:
  - The lap belt of the personal restraint system is usually attached by the qualified personnel to the left and right sides respectively on the mounting point/pin of the product.
  - The shoulder harness of the personal restraint system is usually mounted on the vehicle pillar and is attached by the qualified personnel to the corresponding mounting point/pin provided on the lap belt.
- The straps of the personal restraint system must always be routed close to the user's body. The straps must not be routed over the wheels.
- The shoulder harness must always be routed over the user's shoulder.
- The belt strap must not be twisted on the user's body.

### Placement of the personal restraint system integrated in the vehicle for transporting persons with reduced mobility

- 1) Pull each end of the restraint lap belt from the inner side of the seat between the upholstery and frame through to the outside (see fig. 70).
- 2) Engage the end of the restraint lap belt on the mounting point/pin (see fig. 69).
- 3) Secure the shoulder harness on the mounting point/pin provided on the lap belt (not illustrated).
  - The restraint lap belt is pulled through and fastened.

### 8.6.3 Restrictions for use

#### **⚠ WARNING**

#### **Risk of accidents and injury due to using the product with certain settings and/or installed options**

Severe injury in case of accidents due to options coming loose

- ▶ Before using the product as a seat in a vehicle for transporting persons with reduced mobility, remove options that need to be taken off for safe transportation. Please observe the following table.
- ▶ Stow all dismantled components securely in the vehicle for transporting persons with reduced mobility.
- ▶ Please note that certain settings on the product exclude the use of the product in a vehicle for transporting persons with reduced mobility.

Option <sup>1)</sup>	Transportation in a vehicle for transporting persons with reduced mobility not possible	Remove option	Secure option on product
Pommel			X <sup>2)</sup>
Hook-and-loop foot straps			X <sup>3)</sup>
Storage bag		X <sup>4)</sup>	
Buggy board		X	
Canopy		X	
Grab rail		X	
Head support		X	
Slip sack		X	
Rain cape		X	
Lateral guides		X	
Backpack		X	

<sup>1)</sup> Not all of the options named are installed on all products.

<sup>2)</sup> For use in a vehicle for transporting persons with reduced mobility, the pommel has to be folded away to the front and secured in this position (see Page 59).

<sup>3)</sup> The feet are not permitted to be in the straps, or to be secured with them, during transport in a vehicle for transporting persons with reduced mobility (see Page 62).

<sup>4)</sup> For use in a vehicle for transporting persons with reduced mobility, the storage bag can remain on the product provided it is emptied and the contents are securely stowed.

### 8.7 Care

#### **⚠ CAUTION**

#### **Lack of or improper cleaning**

Health hazard due to infections, damage to the product due to user error

- ▶ Clean the product at regular intervals.
- ▶ Do not clean the product with a jet of water or a pressure washer. The penetration of water can cause corrosion and subsequent malfunctions.
- ▶ Check the driving behaviour of the product after cleaning it.

#### 8.7.1 Cleaning

##### 8.7.1.1 Cleaning by hand

- 1) Remove the pads and covers (Detaching and attaching pads, Detaching/attaching the covers) and separate the covers from the pads.
- 2) Wipe the covers and pads with a damp cloth. Allow to dry.
- 3) **If needed:** Wash covers made from spacer fabrics at max. **40 °C**.

**INFORMATION: The covers may also be washed in a washing machine if they are put in a washing bag or pillowcase.**

#### **Important information on cleaning**

- Do not spray the product with a pressure washer.

- Do not use any aggressive cleaners, solvents or hard brushes etc.
- Clean plastic parts, frame parts and the chassis and wheels with a mild cleaner and a damp cloth. Dry thoroughly afterwards.

### 8.7.1.2 Cleaning belts/straps

#### Cleaning a belt system with metal closure

##### INFORMATION

Observe the washing recommendations on the product and the information in the corresponding instructions for use provided for the product.

- Straps with metal closures **may not be washed in the washing machine** as the penetration of water could cause corrosion and subsequent malfunctions.
- Clean the belt straps by gently dabbing them with warm soapy water (with some disinfectant) or carefully wipe with a dry, clean, absorbent cloth.

#### Cleaning a belt system with plastic closure

- Depending on the model, belts/straps with plastic buckles can be washed in the washing machine between **40 °C** and **60 °C**.
- Recommendation: Use a laundry bag or net and mild detergent.

##### INFORMATION

Alternatively, the belt straps can be cleaned by gently dabbing them with warm soapy water (with some disinfectant) or carefully wiped with a dry, clean, absorbent cloth.

#### Additional cleaning instructions

- Allow the belts to air dry. Ensure that the belts and pads are completely dry before installation.
- Do not expose the belts to direct heat (e.g. sunshine, stove or radiator).
- Do not iron or bleach the belts.

### 8.7.2 Disinfection

- 1) Thoroughly clean the pads before disinfecting.
- 2) Wipe all parts of the product with a disinfectant.

#### Important information about disinfecting

- If the product is used by more than one person, using a conventional disinfectant is required.
- Only use colourless water-based disinfectants. Observe the instructions for use provided by the manufacturer.

## 9 Maintenance and repair

### 9.1 Maintenance

##### WARNING

#### Insufficient maintenance

Severe user injuries, damage to the product due to failure to observe maintenance intervals

- ▶ Only carry out the maintenance tasks described in this section. All other maintenance and service tasks may only be carried out by qualified personnel.
- ▶ The functionality and operating safety of the product must be verified and a service performed at least **once per year**.
- ▶ For users with a changing anatomy (for example body dimensions, weight) or users with a changing clinical picture, have the product inspected, adjusted and serviced at least **once every six months**.

- The function of the product should be checked **before each use**.
- The product may not be used if defects are noted. This applies in particular in case of instability of the product or altered driving characteristics as well as problems with the user's seating position or the stability of the seat. Inform the qualified personnel promptly for the rectification of defects.
- This also applies if loose, worn, bent or damaged components, cracks in the frame or broken frame components are identified.

- Some maintenance tasks can be carried out to a certain extent by the user at home (see the sections "Maintenance intervals" and "Maintenance tasks").
- Failure to maintain the product can lead to injuries for the user of the product.

### 9.1.1 Maintenance intervals

The functions described below must be checked by the user or an attendant at the specified intervals:

Inspection task	Before each use	Monthly	Quarterly
Functional test of the wheel locks	X		
Check of all screw connections		X	
Wear and tear of the folding mechanism		X	
Locking mechanism of the footplates		X	
Check belts for wear		X	
Check that padding and calf band are in proper condition			X
Visual inspection of wear and tear parts (e.g., tyres, bearings)			X
Soiling of bearings			X
Check the legibility of all labels and markings on the product			X

### 9.1.2 Maintenance tasks

To ensure smooth operation at all times, users or attendants with some technical skills can maintain some parts of the product:

- Particles of dirt and hair accumulate on the caster wheel axle and threaded axle of the caster fork over time. This makes the steering stiffer. Regularly remove dirt and oil the axles. See the section "Approach in case of stiffness".
- Dirt particles and hair accumulate on the axle over time. This causes the wheels to turn less easily. Regularly remove dirt from the axle to ensure wheels remain functional. Then lubricate the axle with a few drops of thin, resin-free oil (e.g. sewing machine oil).  
**CAUTION: Do not use pointed or sharp objects that could damage the product when removing dirt. If the wheels still do not rotate smoothly, please contact the qualified personnel.**
- If the rehab buggy gets wet, it should be rubbed dry again.
- Screw connections must be periodically checked for tightness, especially during the initial period of use or after making adjustments to the rehab buggy. If a screw connection loosens repeatedly, contact the qualified personnel promptly.

## 9.2 Repair

### WARNING

#### Prohibited repairs

Severe user injuries, damage to the product due to adjustment and installation errors

- ▶ All repairs may only be carried out by qualified personnel.

## 10 Disposal

### 10.1 Disposal information

Return the product to the qualified personnel for disposal.

All components of the product must be disposed of properly in accordance with the respective national environmental regulations.

### 10.2 Information on re-use

### CAUTION

#### Used seat padding

Functional and/or hygienic risks due to re-use

- ▶ Replace the seat padding if the wheelchair is to be re-used.

The product is suitable for re-use.

Similar to second-hand machines or vehicles, products that are being re-used are subject to increased strain. Features and performance must not change in a way that could impair the safety of users or third parties during the period of use.

The relevant product must be thoroughly cleaned and disinfected before re-use. Then have the product inspected by qualified personnel with respect to its condition, wear and tear, and damage. Worn and damaged parts as well as components that do not fit or are unsuitable for the user must be replaced.

Detailed information on replacing components as well as information on the required tools can be found in the service manual.

## 11 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 11.2 Warranty

Further information on the warranty terms and conditions can be obtained from the qualified personnel that has fitted this product or the manufacturer's service (see inside back cover for addresses).

### 11.3 Service life

Expected lifetime: **4 years**

The design, manufacturing and requirements for the intended use of the product are based on the expected lifetime. These also include the requirements for maintenance, ensuring effectiveness and the safety of the product.

## 12 Technical data

### INFORMATION

- ▶ Much of the technical data below is given in mm. Please note that product settings – unless otherwise specified – cannot be adjusted in the mm range but only in increments of approx. **0.5 cm** or **1 cm**.
- ▶ Note that the values achieved during adjustment may deviate from the values specified below. The deviation can be **±10 mm** and **±2°**.

Dimensions and weights	Size 1	Size 2
Model version	HR32130000	HR32140000
Overall width [mm]	605	655
Min./max. overall length [mm]	1020/1220	1,020/1,220
Push handle height (height of the buggy) [mm]	1,030	1,030
Front wheel diameter [mm]	180	180
Rear wheel diameter [mm]	250	250
Folding size (LxWxH) [mm]	1,030 x 450 x 430	1,030 x 450 x 430
Turning radius [mm]	1,400	1,450
Seat depth [mm]	300/350	350/400
Seat width (with/without padding) [mm]	360/320	410/370
Lower leg length [mm]	270–400	270–400
Seat height, front/rear [mm]	500/380	520/380
Back support height [mm]	650	700
Angle between back support and seat surface (min./max.) [°]	90–120	90–120
Seat angle [°]	22	22
Max. user weight [kg]	50	60
Weight [kg]	13.6	14

**Ambient conditions**

<b>Temperatures and relative humidity</b>	
Temperature during use [°C (°F)]	-10 to +40 (14 to 104)
Transport and storage temperature [°C (°F)]	-10 to +40 (14 to 104)
Relative humidity [%]	10 to 85; non-condensing

**13 Appendices**

**13.1 Threshold values for rehab buggies transportable by train**

**INFORMATION**

► The products in this series fully satisfy the minimum technical requirements of regulation (EU) No. 1300/2014 regarding train accessibility for people with disabilities.

<b>Feature</b>	<b>Threshold value (according to Regulation (EU) No. 1300/2014)</b>
Length [mm]	1200 (plus 50 mm for the feet)
Width [mm]	700 (plus 50 mm on each side for the hands when moving)
Smallest wheels ["]	approx. 3 or greater (according to the regulation, the smallest wheel must be able to overcome a gap measuring 75 mm horizontally and 50 mm vertically)
Height [mm]	max. 1375; including a 1.84 m large male user (95th percentile)
Turning radius [mm]	1500
Maximum weight [kg]	200 (product with user, including luggage)
Maximum obstacle height that can be overcome [mm]	50
Ground clearance [mm]	60 (at an upward slope angle of 10°, ground clearance must measure at least 60 mm under the footrest for going forward at the end of the slope)
Maximum inclination angle on which the product will remain stable [°]	6 (dynamic stability in all directions) 9 (static stability in all directions, also when wheel lock engaged)

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Otto Bock Healthcare Products GmbH  
Brehmstraße 16 · 1110 Wien · Austria  
F +43 1 5267985  
service-admin.vienna@ottobock.com · www.ottobock.at

Otto Bock Adria d.o.o. Sarajevo  
Ramiza Salčina 85  
71000 Sarajevo · Bosnia-Herzegovina  
T +387 33 255-405 · F +387 33 255-401  
obadria@bih.net.ba · www.ottobockadria.com.ba

Otto Bock Bulgaria Ltd.  
41 Tzar Boris III<sup>rd</sup> Blvd. · 1612 Sofia · Bulgaria  
T +359 2 80 57 980 · F +359 2 80 57 982  
info@ottobock.bg · www.ottobock.bg

Otto Bock Suisse AG  
Luzerner Kantonsspital 10 · 6000 Luzern 16 · Suisse  
T +41 41 455 61 71 · F +41 41 455 61 70  
suisse@ottobock.com · www.ottobock.ch

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Otto Bock South Africa (Pty) Ltd  
Building 3 Thornhill Office Park · 94 Bekker Road  
Midrand · Johannesburg · South Africa  
T +27 11 564 9360  
info-southafrica@ottobock.co.za  
www.ottobock.co.za

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Otto Bock HealthCare Canada  
5470 Harvester Road  
Burlington, Ontario, L7L 5N5, Canada  
T +1 800 665 3327 · F +1 800 463 3659  
CACustomerService@ottobock.com  
www.ottobock.ca

Oficina Ottobock Habana  
Calle 3ra entre 78 y 80.  
Edificio Jerusalem · Oficina 112 · Calle 3ra.  
Playa, La Habana. Cuba  
T +53 720 430 69 · +53 720 430 81  
hector.corcho@ottobock.com.br  
www.ottobock.com.br

Otto Bock HealthCare Andina Ltda.  
Calle 138 No 53-38 · Bogotá · Colombia  
T +57 1 8619988 · F +57 1 8619977  
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.  
Prolongación Calle 18 No. 178-A  
Col. San Pedro de los Pinos  
C.P. 01180 México, D.F. · Mexico  
T +52 55 5575 0290 · F +52 55 5575 0234  
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare LP  
11501 Alterra Parkway Suite 600  
Austin, TX 78758 · USA  
T +1 800 328 4058 · F +1 800 962 2549  
USCustomerService@ottobock.com  
www.ottobockus.com

## Asia/Pacific

Otto Bock Australia Pty. Ltd.  
Suite 1.01, Century Corporate Centre  
62 Norwest Boulevard  
Baulkham Hills NSW 2153 · Australia  
T +61 2 8818 2800 · F +61 2 8814 4500  
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.  
B12E, Universal Business Park  
10 Jiuxianqiao Road, Chao Yang District  
Beijing, 100015, P.R. China  
T +8610 8598 6880 · F +8610 8598 0040  
news-service@ottobock.com.cn  
www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.  
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza  
1 Science Museum Road, Tsim Sha Tsui  
Kowloon, Hong Kong · China  
T +852 2598 9772 · F +852 2598 7886  
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India Pvt. Ltd.  
20th Floor, Express Towers  
Nariman Point, Mumbai 400 021 · India  
T +91 22 2274 5500 / 5501 / 5502  
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.  
Yokogawa Building 8F, 4-4-44 Shibaura  
Minato-ku, Tokyo, 108-0023 · Japan  
T +81 3 3798-2111 · F +81 3 3798-2112  
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.  
4F Agaworld Building · 1357-74, Seocho-dong  
Seocho-ku, 137-070 Seoul · Korea  
T +82 2 577-3831 · F +82 2 577-3828  
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.  
1741 Phaholyothin Road  
Kwaeng Chatuchark · Khet Chatuchark  
Bangkok 10900 · Thailand  
T +66 2 930 3030 · F +66 2 930 3311  
obsea@ottobock.co.th · www.ottobock.co.th

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Otto Bock Mobility Solutions GmbH