ottobock.

451F10, 451F11, 451F12, 451F13

DE Gebrauchsanweisung	3
EN Instructions for use	5
FR Instructions d'utilisation	7
IT Istruzioni per l'uso	9
ES Instrucciones de uso	11
PT Manual de utilização	13
SV Bruksanvisning	15
PL Instrukcja użytkowania	17
HU Használati utasítás	19
CS Návod k použití	21
EL Οδηγίες χρήσης	23

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English

Last update: 2019-06-17

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

The 451F10, 451F11, 451F12 and 451F13 Compression Shrinker Socks are to be used solely for compressive treatment following lower limb amputations.

2 Indications

- Reduces post-operative edema
- Prevents edema after removal of the prosthesis

3 Contraindications

- Weeping skin areas
- Open wounds
- Unreduced edema, congestive dermatitis, arterial circulatory disturbance
- Thrombophlebitis acuta, congestive heart failure

4 Indications for use

- Physician should determine the amount of compression used and when product should be fit
- Therapeutic success requires regular check-ups with a physician

5 Recommended application

Compression Shrinker Socks should be worn only when prescribed by a physician and should be given to the patient by trained personnel. Before applying or removing the sock, we recommend removing jewelry and wearing thin protective gloves since rough fingernails might damage the fine material of the sock.

1. Turn the sock inside out (fig. 1).
2. Stretch the double-layer upper part of the sock, so that you can slip your residual limb into the sock (fig. 2).
3. Slowly and evenly slip the sock upwards over the leg without twisting it. For optimal effect, "massage" the sock onto the leg for even distribution and to correct any "overstretched" material (fig. 3).

6 Care & Washing instructions

- Hand or machine wash (gentle cycle) your compression shrinker sock-daily using a standard mild detergent in warm water (max. 40 °C/104 °F).
- If your sock is provided with a hook & loop closure, close it before washing.
- Thoroughly rinse with clear water. Do not use any softeners, spin-dry or wring out.
- After washing, roll the sock in a towel, firmly press out excess water and let air-dry.
- Do not expose the sock to direct heat (radiator, sun, etc.). Do not tumble-dry, iron, chemically clean or use benzine for cleaning.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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