ottobock.

1A30

DE	Gebrauchsanweisung	4
EN	Instructions for use	8
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	17
ES	Instrucciones de uso	21
PT	Manual de utilização	25
NL	Gebruiksaanwijzing	29
SV	Bruksanvisning	33
DA	Brugsanvisning	37
NO	Bruksanvisning	41
PL	Instrukcja użytkowania	45
HU	Használati utasítás	50
CS	Návod k použití	54
RO	Instrucțiuni de utilizare	58
EL	Οδηγίες χρήσης	63
RU	Руководство по применению	67
ZH	使用说明书	72

●	▲	■
(D) Einzelteile-Pack	Mindestmenge	Einzelteile
(GB) Single-Component Pack	Minimum order quantity	Single components
(F) Kit de pièces de rechange	Quantité minimum	Pièces à l'unité
(I) Imballo dei singoli componenti	Quantità minima	Singoli componenti
(E) Kit componentes	Cantidad mínima	Componentes
(P) Conjunto de peças de reposto	Quantidade mínima	Componentes
(NL) Onderdelenpakket	Minimikvantitet	Onderdelen
(S) Service-Set	Minimum aantal	Separata delar
(DK) Komponentpakke	Mindste mængde	Komponenter
(PL) Zestaw naprawczy	Ilość minimalna	Podzespoły
(H) Alkatrész-csomag	Minimum mennyiség	Alkatrészek
(CZ) Sada jednotlivých dílů	Minimální množství pro objednání	Jednotlivé díly
(RO) Set componente individuale	Cantitate minimă	Componente individuale
(HR) Pakiranje pojedinačnih dijelova	Minimalna količina	Pojedinačni dijelovi
(TR) Münferit parça ambalajı	Minimum miktar	Münferit parçalar
(GR) Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα	Μεμονωμένα εξαρτήματα
(RU) Комплект деталей	Минимальное количество	Отдельные детали
(J) 交換パーツセット	最低注文個数あり	個別販売
(中) 维修组件	最低起订量	单个部件

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English

Last update: 2018-08-22

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document in a safe place.

1 Scope of delivery (Fig. 1)

1.1 Single-Component pack (Fig. 2, ●)

2D3 (for size 24–25), *consisting of:*

- 1 of each elastomer soft, medium, hard (1),
- 1 flexible joint retainer (2),

- 1 two hole washer (3),
- 2 attachment bolts (4),
- 1 washer (5),
- 1 lock nut (6).

2D4 (for size 26–29), *consisting of:*

- 1 of each elastomer soft, medium, hard (1),
- 1 flexible joint retainer (2),
- 1 two hole washer (3),
- 2 attachment bolts (4),
- 1 washer (5),
- 1 lock nut (6).

1.2 Minimum order quantity ▲

None

1.3 Single components (Fig. 2, ■)

- **1A31=(side)(size) Greissinger plus** (7)
Shaped foot component without adapter
- **2R86=(size) Foam Connection Cap** (8)

2 Description

2.1 Indications for use

The 1A30 Greissinger plus Foot is to be **exclusively** used for the prosthetic fitting of amputations of the lower limb.

2.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System**:

1A30=24-25

Recommended for amputees with **Mobility Grades 2 and 3**
(restricted outdoor walkers and non-restricted outdoor walkers).

Approved for a patient weight of up to 75 kg/165 lbs.

1A30=26-29

Recommended for amputees with **Mobility Grades 2 and 3**
(restricted outdoor walkers and non-restricted outdoor walkers).

Approved for a patient weight of up to 100 kg/220 lbs.

3 Safety instructions

Attention!

Improper use can weaken components of the prosthetic foot. To avoid the risk of accident and possible injury to the patient, the foot must no longer be used after a sudden loss of function has occurred. This noticeable loss of function can manifest itself as reduced support during roll-over or decreased forefoot resistance.

Measure

Have the prosthetic foot examined and, if necessary, replaced at a specialist workshop.

Please avoid exposing prosthetic components to

- corrosive elements such as fresh water, salt water, and acids, or
- abrasive substances (e.g. sand).

Using this medical product under the mentioned environmental conditions will render all claims against Otto Bock HealthCare null and void.

The material is not resistant to UV rays or hydrolysis.

The surface of this prosthetic foot must not be damaged by grinding, as this would lead to premature wear.

Please inform your patients.

4 Description and functional properties

This prosthetic foot combines the proven functional capabilities of the Greissinger ankle with a new flexible keel foot (similar to 1D10).

The integral titanium adapter, combined with the internal elastomeric connections, permits multi-axial ankle movement. Plantar flexion following heel contact is controlled by compression of the ring-shaped rocking rubber, which also provides an elastic dorsiflexion limit in late stance phase. Controlled inversion, eversion, internal rotation, and external rotation are also provided by this mechanism.

The Greissinger Plus is a complete ankle-foot assembly and includes the 2R86 Foam Connection Cap for finishing.

5 Alignment recommendations

5.1 Bench alignment recommendations

Foot size	Mid-foot anterior to the alignment reference line	Heel height
24-29 cm	30 mm	10 mm ± 5 mm

5.2 Static alignment recommendations

Please use the Ottobock alignment recommendations for the alignment of transtibial and transfemoral amputees.

If available, use the L.A.S.A.R. Posture for static alignment of the prosthesis.

5.3 Dynamic alignment recommendations

Adapt the prosthetic alignment in the frontal plane (ML) and in the sagittal plane (AP) by angle adjustments or sliding adjustments, in order to ensure correct heel contact, easy roll-over, and optimal shifting of weight onto the contralateral side. In case of transtibial fittings make sure that a physiological knee flexion is achieved during the stance phase.

5.4 Maintenance instructions

Note:

As a basic principle, all Ottobock modular foot components are subjected to tests involving two million load cycles. Depending on the amputee's activity this corresponds to a service life of two to three years.

We recommend carrying out regular safety checks once a year.

6 Replacing the rocking rubber

The 1A30 Greissinger plus includes a medium Rocking Rubber (2S2). Both a lower durometer (2S15) and a higher durometer (2S19) Rocking Rubber are also available. To replace the rocking rubber, unscrew the hexagon nut on the plantar surface of the foot and remove the ankle assembly. After installing the desired rocking rubber, the **hexagon nut must be torqued to 30 Nm/22.1 ft.lbs./265.4 in.lbs.** using the **710D4 Torque Wrench (with 701Y5 and 710Y6=17) or equivalent.**

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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