

**SiOCX TF 7T450=1, 7T451=1, 7T431=4, 7T431=5,
7T452=1, 7T452=2**



DE Gebrauchsanweisung	4
EN Instructions for use	9
FR Instructions d'utilisation	13
ES Instrucciones de uso	19
PL Instrukcja użytkowania	24

1

2

INFORMATION

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INFORMATION

Date of last update: 2020-07-28

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

SiOCX inner sockets are made of HTV silicone and, depending on the version, also a carbon prepreg outer socket with adapter. Depending on the version, additional silicone gel pads and flexible socket sections can be integrated as options to provide pressure relief in sensitive or bony areas. The inner sockets are individually fabricated according to the patient data. For the SiOCX-TF with closed outer socket, a fitting step is integrated into the process for better adaptation to the patient. For this purpose, the definitive inner socket is supplied with a ThermoLyn check socket without adapter. With the check socket, the shape and socket trim line can be changed and the adapter position can be set.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

The **SiOCX TF Prox 7T452=1** is suitable for various suspension systems such as shuttle locks and KISS pull-in systems. The **SiOCX TF Seal 7T452=2** is especially well suited for elevated vacuum sockets in combination with the ProSeal SIL liner.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The product must be used in TF prostheses only.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

⚠ CAUTION

Use of the product on damaged skin

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product on healthy skin.

⚠ CAUTION

Volume fluctuations or changes of the affected limb

Injuries, pressure sores, skin irritation, circulatory disorders and paralytic symptoms caused by problems with the fit

- ▶ Do not continue using the product if problems with the fit occur.
- ▶ Please consult your prosthetist in case of problems with the fit or signs of skin damage, even if they are only minor.
- ▶ Check the skin of the affected limb daily for damage and sensitivity disturbances (visually and by palpating).

⚠ CAUTION

Reuse on other patients

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Use the product on only one patient.

NOTICE

Contact with heat, embers or fire

Damage to the product due to high temperatures

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

Fig.	Item	Quantity	Designation	Reference number
-	-	1	Instructions for use	-
see fig. 1, see fig. 2	1	1	Check socket (only with 7T450=1, 7T431=4)	7T455=1

Fig.	Item	Quantity	Designation	Reference number
see fig. 3	1	1	SiOCX TF Pro prosthetic socket	7T451=1
see fig. 1	2		SiOCX TF prosthetic socket	7T450=1
see fig. 2	2		SiOCX TF inner socket	7T431=4
see fig. 4	1		SiOCX TF Basic inner socket	7T431=5
see fig. 5	1		SiOCX TF Seal (illustrated with mandatory outer socket)	7T452=2
see fig. 6	1		SiOCX TF Prox	7T452=1

Spare parts/accessories (not included in the scope of delivery)

Designation	Reference number
Check socket (as accessory only for 7T431=5, 7T452=1, 7T452=2)	7T455=1

5 Preparing the product for use

⚠ CAUTION

Incorrect alignment or assembly

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

NOTICE

Improper thermal treatment of the silicone socket

Damage due to thermal overloading

- ▶ Do not carry out any heat treatment at a temperature above 120 °C (248 °F).

INFORMATION

- ▶ Prevent the silicone socket from being damaged by sharp-edged objects to avoid functional impairment.

5.1 Alignment with check socket

⚠ CAUTION

Lack of stability due to cutting out the seating area

Risk of injury due to breakage of load-bearing components

- ▶ Do not cut out the seating area of the thermoplastic prosthetic socket.

⚠ CAUTION

Heat treatment of the thermoplastic prosthetic socket

Risk of injury due to contact with hot material

- ▶ Do not carry out any heat treatments on the prosthetic socket when the patient is wearing it.
- ▶ Allow the prosthetic socket to cool after heat treatment.

INFORMATION

- ▶ If the check socket needs to be reworked, note the instructions for processing ThermoLyn in the 646D300* technical information.

- 1) Check the fit of the check socket on the patient.
- 2) Perform the bench alignment of the prosthesis with the check socket.
- 3) Position the adapter on the check socket as described in the instructions for the adapter.

- 4) Secure the adapter on the check socket such that the prosthesis can withstand the loads occurring during the dynamic trial fitting and test phase.
- 5) Carry out the static alignment and dynamic trial fitting of the prosthesis. If required, mark a flexible seating area for the definitive socket on the check socket.
- 6) Allow the patient to test the prosthesis for a maximum of 4 weeks.
- 7) After a successful test phase, send the complete check socket with inner socket to the manufacturer of the outer socket. Make sure that the adapter cannot be detached from the check socket in transit.

5.2 Prosthetic alignment with definitive socket

INFORMATION

- ▶ Definitive sockets must not be reworked.

- 1) Perform the bench alignment of the prosthesis with the definitive socket.
- 2) Carry out the static alignment and dynamic trial fitting of the prosthesis.
- 3) Finish the prosthesis.

6 Use

6.1 Caring for the skin

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

INFORMATION

Derma Prevent damages inner sockets made of silicone

- ▶ Do not use Derma Prevent for skin care.

- ▶ For skin care use a pH neutral, water-based and unscented skin lotion (e.g. ProComfort Gel).
- ▶ In consultation with your doctor, you can use e.g. Allpresan® diabetic foam cream intensive to alleviate sore areas.

6.2 Donning the prosthesis

- ▶ Slide your residual limb into the prosthetic socket. Use a donning aid if necessary.

Information for inner sockets that are worn directly on the skin

⚠ CAUTION

Prosthetic socket is not in full contact with residual limb

Skin irritation due to accumulation of liquid in prosthetic socket

- ▶ Wear the prosthetic socket directly on the skin.
- ▶ If you detect any air pockets take the prosthetic socket off and don it again.

7 Cleaning

⚠ CAUTION

Improper cleaning of the inner socket

Skin irritation, formation of eczema or infections due to contamination with germs and damage to the inner socket

- ▶ Clean the inner socket daily according to the instructions given in this section.

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

- > **Required cleaning agents:** pH neutral soap (e.g. 453H10 Derma Clean), scratch and lint-free cloth
- 1) Clean the product with a cloth using fresh water and a pH neutral soap.
Make sure no liquid gets onto the other prosthesis components and between the inner and outer sockets.
- 2) Remove soap residue with a cloth using fresh water.
- 3) Dry the product with a lint-free cloth and allow it to fully air dry.

8 Maintenance

- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Conduct half-yearly safety inspections.
- ▶ During the standard inspection of the prosthesis components, also inspect the product.

Maintenance steps

- > **Required materials:** Allen wrench 2 mm, 646F297* isopropyl alcohol
- 1) Check the prosthetic socket for damage.
- 2) **If present:** Loosen the valve closure.
- 3) Mark the matching screws and insert nuts.
- 4) Remove the screws with the Allen wrench.
- 5) Take out the silicone socket.
- 6) **If present:** Remove the electrodes from the silicone socket.
- 7) **If present:** Clean the myo-contact surfaces with isopropyl alcohol.
- 8) Clean the prosthetic socket (see Page 12).
- 9) **Optional:** Boil the silicone socket.
- 10) Reassemble the prosthetic socket.

9 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

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