1B1-2=* Meridium

DE Gebrauchsanweisung (Benutzer)	3
EN Instructions for use (user)	45
FR Instructions d'utilisation (Utilisateur)	85
IT Istruzioni per l'uso (Utilizzatore)	131
ES Instrucciones de uso (Usuario)	173
PT Manual de utilização (Usuário)	215
NL Gebruiksaanwijzing (Gebruiker)	257
HR Upute za uporabu (Korisnik)	299
TR Kullanma talimatı (Kullanıcı)	339
JA 取扱説明書(ユーザー用)	377

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1	Foreword	48
2	Product description	48
2.1	Design	48
2.2	Function	48
3	Intended use	49
3.1	Indications for use	49
3.2	Conditions of use	49
3.3	Qualification	49
4	Safety	49
4.1	Explanation of warning symbols	49
4.2	Structure of the safety instructions	50
4.3	General safety instructions	50
4.4	Information on the Power Supply/Battery Charging	53
4.5	Battery charger information	53
4.6	Information on Proximity to Certain Areas	54
4.7	Information on Use	55
4.8	Notes on the safety modes	56
4.9	Information on the use of a mobile device with the cockpit app	56
5	Scope of Delivery and Accessories	57
5.1	Scope of delivery	57
5.2	Accessories	57
6	Charging the battery	57
6.1	Connecting the power supply and battery charger	58
6.2	Charging the prosthesis battery	58
6.3	Display of the current charge level	59
7	Cockpit app	59
7.1	Initial connection between cockpit app and prosthesis	60
7.1.1	Starting the cockpit app for the first time	60
7.2	Control elements for cockpit app	61
7.2.1	Cockpit app navigation menu	62
7.3	Managing components	62
7.3.1	Adding component	62
7.3.2	Deleting a component	63
7.3.3	Connecting component with multiple mobile devices	63
8	Use	63
8.1	Setting the heel height	63
8.1.1	Setting heel height using motion pattern	64
8.1.2	Setting the heel height with the cockpit app	64
8.2	Movement patterns in basic mode (mode 1)	64
8.2.1	Standing	64
8.2.2	Walking	65

8.2.3	Sitting down/sitting	65
8.2.4	Standing up	65
8.2.5	Walking up stairs.....	65
8.2.6	Walking down stairs	66
8.2.6.1	Stair function	66
8.2.7	Walking down a ramp	66
8.2.8	Walking backwards	66
8.2.9	Walking up a ramp	67
8.2.10	Standing on inclined surface.....	67
8.2.11	Kneeling.....	67
8.2.12	Relief function	67
8.3	MyModes	67
8.3.1	Switching MyModes with the cockpit app	68
8.3.2	Switching MyModes using motion patterns	68
8.3.3	Switching on ankle lock.....	69
8.3.4	Switching from a MyMode back to basic mode.....	70
8.4	Changing prosthesis settings.....	70
8.4.1	Changing the prosthesis setting using the cockpit app.....	71
8.4.2	Overview of adjustment parameters in basic mode.....	71
8.4.3	Overview of adjustment parameters in MyModes	72
8.5	Turning Bluetooth on the prosthesis on/off.....	72
8.5.1	Switching Bluetooth off/on using the cockpit app	72
8.6	Querying the prosthesis status	73
8.6.1	Query status through cockpit app	73
8.6.2	Status display in the cockpit app	73
9	Additional operating states (modes)	73
9.1	Empty battery mode	73
9.2	Mode for charging the prosthesis	73
9.3	Safety mode.....	73
9.4	Overheating mode	73
10	Cleaning	74
11	Maintenance	74
12	Legal information	74
12.1	Liability	74
12.2	Trademarks	74
12.3	CE conformity	74
12.4	Local Legal Information	75
13	Technical data	76
14	Appendices	78
14.1	Symbols Used.....	78
14.2	Operating states/error signals	79
14.2.1	Signals for operating states	79
14.2.2	Warnings/error signals.....	80
14.2.3	Error messages while establishing a connection with the cockpit app	82

14.2.4	Status signals	82
14.3	Directives and manufacturer's declaration	83
14.3.1	Electromagnetic environment	83

1 Foreword

INFORMATION

Date of last update: 2020-07-28

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The product "1B1-2* Meridium" is referred to as the product/component/prosthesis/prosthetic foot below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Design

The product consists of the following components:

1. Cover/connection plate
2. Cover cap with charging receptacle
3. Battery
4. Ankle spring
5. Heel bracket
6. Ankle axis
7. Heel spring
8. Hydraulic unit
9. Toe axis
10. Toe plate
11. Carbon fibre frame
12. Main electronics
13. Spherical cap with pyramid

2.2 Function

This product features microprocessor-controlled damping of plantar flexion (motion of the foot in the ankle joint in the direction of the sole of the foot) and dorsiflexion (motion of the foot in the ankle joint in the direction of the instep).

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the damping behaviour of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

Thanks to the microprocessor-controlled damping of the plantar flexion and dorsiflexion, the product can be individually adapted to requirements.

For this purpose, the product is adjusted by qualified personnel using adjustment software.

The product features MyModes for special motion types (e.g. golf ...). These have default values configured by the O&P professional using the adjustment software, and can be called up via the Cockpit app or using a special motion pattern (see page 67).

If configured by the prosthetist, an additional mode "Ankle lock" can also be selected, locking the ankle joint of the prosthetic foot in its current position.

In case of a product malfunction, safety mode makes restricted operation possible. Resistance parameters that are predefined by the product are configured for this purpose (see page 73).

Empty battery mode permits safe walking when the battery is drained. Resistance parameters that are predefined by the product are configured for this purpose (see page 73).

The microprocessor-controlled hydraulic unit offers the following advantages

- Approximation of the physiological gait pattern
- Stable standing on even and inclined surfaces
- Adaptation of product characteristics to various surfaces, inclines, gait situations, walking speeds and heel heights

3 Intended use

3.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.

3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive impact load (tennis, basketball, running, ...) or extreme sports (free climbing, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see page 76).

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker).

Foot size [cm]	24 to 25	26 to 29
Max. body weight [kg]	100	125

3.3 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

4 Safety

4.1 Explanation of warning symbols

WARNING	Warning regarding possible serious risks of accident or injury.
CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

WARNING

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

WARNING

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

WARNING

Using the prosthesis while operating a vehicle

Accident due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ All users are required to observe their country's national and state driving laws when operating vehicles with a prosthesis. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.
- ▶ Observe national legal regulations for retrofitting your vehicle in accordance with the type of fitting.
- ▶ The leg on which the prosthesis is worn may not be used to control the vehicle or its peripheral components (e.g. clutch pedal, brake pedal, gas pedal, etc.).

WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

CAUTION

Failure to observe warning/error signals

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 80) and corresponding change in damping settings must be observed.

CAUTION

Independent manipulation of the product and the components

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.

- ▶ The battery may only be handled by authorised, qualified Ottobock personnel (no replacement by the user).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Mechanical stress on the product

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Use of the product when battery charge level is too low

- Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.
- ▶ Check the current charge level before use and charge the prosthesis if required.
 - ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

⚠ CAUTION

Penetration of liquids into the product

Falling due to unexpected product behaviour as the result of a malfunction.

- ▶ When the footshell is intact, the product is protected against splashed water from all directions. However, it is not protected against submersion, jets of water and steam.
- ▶ If water has penetrated the product, have the footshell removed by the O&P professional and allow both components to dry. The product must be inspected by an authorised Ottobock Service Center. The O&P professional is your contact.
- ▶ If salt water has penetrated the product, the footshell must be removed immediately by the O&P professional. The product must be inspected by an authorised Ottobock Service Center.
- ▶ Do not use the product for bathing prostheses.

⚠ CAUTION

Overloading due to unusual activities

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive impact load (tennis, basketball, running, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by a prosthetist. If necessary, the prosthetist will forward the product to an authorised Ottobock Service Centre.

CAUTION

Signs of wear and tear on the product components

Falling due to damage or malfunction of the product.

- ▶ Regular service inspections (maintenance) are mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty.

CAUTION

Use of unapproved components

> Injury due to product malfunction as a result of reduced interference resistance.

> Interference of other electronic devices due to increased emissions.

- ▶ Use the product only in combination with the components, signal converters and cables listed in the section "Scope of delivery" (see page 57).

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage.

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

CAUTION

Use of unapproved accessories

> Falling due to product malfunction as a result of reduced interference resistance.

> Interference of other electronic devices due to increased emissions.

- ▶ Use the product only in combination with the accessories, signal converters and cables listed in the sections "Scope of delivery" (see page 57) and "Accessories" (see page 57).

INFORMATION

When exoprosthetic components are used, the use of control functions which are performed hydraulically or movement of the component in the footshell can cause noise. This kind of noise is normal and unavoidable. It generally does not cause any problems. If movement noise increases noticeably during the lifecycle of the component, the component should be immediately inspected by authorised qualified Ottobock personnel.

Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

4.4 Information on the Power Supply/Battery Charging

⚠ CAUTION

Charging the prosthesis without taking it off

- > Falling due to walking and getting caught on a connected battery charger.
- > Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.
- ▶ For safety reasons, remove the prosthesis prior to charging it.

⚠ CAUTION

Charging the product with damaged power supply unit/charger/charger cable

- Falling due to unexpected behaviour of the product caused by insufficient charging.
- ▶ Check the power supply unit, charger and charger cable for damage before use.
 - ▶ Replace any damaged power supply unit, charger or charger cable.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

4.5 Battery charger information

NOTICE

Penetration of dirt and humidity into the product

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

NOTICE

Mechanical stress on the power supply/battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

NOTICE

Operating the power supply unit/charger outside of the permissible temperature range

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see page 76).

NOTICE

Independent changes or modifications carried out to the battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

4.6 Information on Proximity to Certain Areas

CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

CAUTION

Operating the product in very close proximity to other electronic devices

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for people (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations...).
- ▶ If this cannot be avoided, make sure to at least have a safeguard when walking or standing (e.g. a handrail or the support of another person).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected changes in the damping behaviour of the product.
- ▶ In general, monitor the product for unexpected changes in the damping behaviour when electronic or magnetic devices are in the immediate vicinity.

CAUTION

Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of strong magnetic fields.
- ▶ Take off the product before entering a room or area with strong magnetic fields and store the product outside this room or area.
- ▶ Damage to the product caused by exposure to strong magnetic fields cannot be repaired.

CAUTION

Remaining in areas outside the allowable temperature range

Falling due to malfunction or the breakage of load-bearing product components.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see page 76).

4.7 Information on Use

CAUTION

Walking up stairs

Falling due to foot placed incorrectly on stair.

- ▶ Always use the handrail when walking up stairs and place most of the sole of your foot on the stair surface. The toe plate of the prosthesis may fold away if only the front part of the foot is placed on the edge of the stair.
- ▶ Particular caution is required when carrying children up stairs.

CAUTION

Walking down stairs

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and place most of the sole of your foot on the stair surface. Rolling the foot over the edge of the stair is not necessary.
- ▶ Particular caution is required when carrying children down the stairs.

CAUTION

High temperature of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)

Falling due to unexpected behaviour of the product caused by use in overheating mode.

- ▶ Be sure to pay attention when pulsating vibration signals start. They indicate the risk of overheating.
- ▶ As soon as these pulsating vibration signals begin, you must reduce your level of activity so the hydraulic unit can cool down.
- ▶ Please note that the range of motion in the ankle joint is reduced as the temperature increases, and this may ultimately result in the ankle joint being locked completely. Special caution is therefore required, especially when walking down stairs.
- ▶ Full activity may be resumed after the pulsating vibration signals stop.

CAUTION

Improper mode switching

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.
- ▶ Take the weight off the product and correct the switching, if required.

CAUTION

Using prosthetic foot without footshell

Falling due to slipping when walking on smooth floor surfaces (tiles).

- ▶ Do not use the prosthetic foot without the accompanying footshell.

CAUTION

Using prosthetic foot with a damaged foot shell

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.

- ▶ Do not use the prosthetic foot when the foot shell is damaged. Replace damaged foot shells immediately, before using the product again.

4.8 Notes on the safety modes

CAUTION

Using the product in safety mode

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 80) have to be observed.

CAUTION

Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Do not continue using the defective product.
- ▶ Consult the O&P professional promptly.

CAUTION

Safety mode cannot be deactivated

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ If safety mode cannot be deactivated by recharging the battery, a permanent error has occurred.
- ▶ Do not continue using the defective product.
- ▶ The product must be inspected by an authorised Ottobock Service Center. The O&P professional is your contact.

CAUTION

Safety signal occurs (ongoing vibration)

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Observe the warnings/error signals (see page 80).
- ▶ Do not continue using the product after the safety signal has been emitted.
- ▶ The product must be inspected by an authorised Ottobock Service Center. The O&P professional is your contact.

4.9 Information on the use of a mobile device with the cockpit app

CAUTION

Improper use of the mobile device

Falling due to changed damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Make sure you have been instructed in the proper use of the mobile device with the Cockpit app.

CAUTION

Independently applied changes or modifications made to the mobile device

Falling due to altered damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Do not make any independent changes to the hardware of the mobile device on which the app is installed.
- ▶ Do not make any independent changes to the software/firmware of the mobile device that are not included in the update function of the software/firmware.

⚠ CAUTION

Improper mode switching with the device

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter (beeper) and the display on the device.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.

NOTICE

Failure to observe the system requirements for the installation of the Cockpit app

Mobile device malfunction.

- ▶ The Cockpit App should only be installed on mobile devices and versions which comply with the specifications in the respective online stores (e.g. Apple App Store, Google Play Store, ...)

INFORMATION

The illustrations in these instructions for use are only examples and may deviate from the respective mobile device being used and the version.

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 1 pc. 1B1-2 Meridium
- 1 pc. 757L16-4 power supply
- 1 pc. 4E50-2 C-Leg battery charger
- 1 pc. 646C107 Bluetooth PIN card
- 1 pc. 647F542 prosthesis passport
- 1 pc. Instructions for use (user)

Cockpit app for download from the website: <https://www.ottobock.com/cockpitapp>

- "Cockpit 4X441-IOS=V*" app for iOS
- "Cockpit 4X441-ANDR=V*" app for Android

5.2 Accessories

The following components are not included in the scope of delivery and may be ordered separately:

- 1 pc. 757P48 Y adapter cable
This is for charging more than one product at the same time (e.g. 1B1-2; 1B1; 3B1/3B1=ST; 3B1-2/3B1-2=ST; 3B5-X3/3B5-X3=ST; 3C98-2/3C88-2; 3C98-3/3C88-3; 3C96-1/3C86-1) with the 757L16-4 power supply.

6 Charging the battery

The following points must be observed when charging the battery:

- Use the 757L16-4 power supply and 4E50-2 battery charger to charge the battery.
- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- For the maximum operating time with one battery charge, disconnecting the battery charger from the product only immediately before using the product is recommended.

- The battery should be charged until the yellow LED on the battery charger turns off prior to initial use, and for at least 4 hours. This calibrates the charge level indicator via the Cockpit app and by turning over the prosthesis.
If the battery charger is disconnected from the prosthesis too soon, the charge level indicator via the Cockpit app and by turning over the prosthesis may not correspond to the actual charge level.
- The ankle joint of the prosthetic foot is locked during the charging process.
- The battery may discharge while the product is not being used.

6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
- 2) Connect the round, **four-pin** plug of the charging cable to the **OUT** receptacle on the battery charger so that the plug locks into place (see fig. 2).
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
- 3) Connect the round, **three-pin** plug of the power supply to the **12 V** receptacle on the battery charger so that the plug locks into place (see fig. 2).
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
- 4) Plug the power supply unit into the wall socket.
→ The green LED on the back of the power supply and the green LED on the battery charger light up (see fig. 3).
→ If the green LED on the power supply and the green LED on the battery charger do not light up, there is an error (see page 80).

6.2 Charging the prosthesis battery

- 1) Open the charging receptacle cover.
- 2) Connect the charging plug to the charging receptacle of the product.
INFORMATION: Take note of the plug direction!
→ The correct connection of the battery charger to the product is indicated by feedback (see page 82).
- 3) The charging process starts.
→ Once the product battery is fully charged, the yellow LED on the battery charger turns off.
- 4) Disconnect the product after the charging process is complete.
→ This is followed by an electronics self-test which is confirmed by feedback (see page 82).
- 5) Close the charging receptacle cover.

6.3 Display of the current charge level

INFORMATION

The charge level cannot be displayed during the charging process.

- 1) Turn the prosthesis by 180° (the sole of the foot must face up).
- 2) Hold the prosthesis still and wait for beep signals.

Prosthetic foot with knee joint:

The beep signal for the knee joint sounds after approx. 2 seconds.

The beep signal for the prosthetic foot sounds after approx. 4 seconds.

Prosthetic foot without knee joint:

The beep signal for the prosthetic foot sounds after approx. 2 seconds.

Beep signal	Vibration signal	Battery charge level
5x short	–	More than 80%
4x short	–	66% to 80%
3x short	–	51% to 65%
2x short	–	36% to 50%
1x short	3x long	20% to 35%
1x short	5x long	Less than 20%

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 70).

Display of the current charge level using the cockpit app:

Once the Cockpit app has been started, the current charge level is displayed in the bottom line of the screen:

1.  38% – Charge level of battery for currently connected component

7 Cockpit app

The Cockpit app enables switching from basic mode to the pre-configured MyModes. In addition, information about the product (step counter, battery charge level, etc.) can be retrieved.

The everyday behaviour of the product can be changed to a certain extent using the app (e.g. while becoming accustomed to the product). The O&P professional can use the adjustment software to track these changes at the next appointment.

Information on the Cockpit app

- The Cockpit app can be downloaded free of charge from the respective online store. For more information, please visit the following website: <https://www.ottobock.com/cockpitapp>. To download the Cockpit app, the QR code on the supplied Bluetooth PIN card can also be read with the mobile device (requirement: QR code reader and camera).
- The language of the Cockpit app user interface can only be changed by the O&P professional using the adjustment software.
- The serial number of the component to be connected has to be registered with Ottobock the first time it is connected. If the registration is not accepted, use of the Cockpit app for this component will be limited.
- Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (sole of the foot must point up) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 72).
- Keep the mobile app up to date at all times.
- Please contact the manufacturer if you suspect cybersecurity problems.

7.1 Initial connection between cockpit app and prosthesis

The following points need to be observed before establishing the connection:

- Bluetooth of the component must be switched on (see page 72).
- Bluetooth on the mobile device must be switched on.
- The mobile device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- **The mobile device must be connected to the Internet.**
- The serial number and Bluetooth PIN of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The serial number begins with the letters "SN".

INFORMATION

If the Bluetooth PIN card with the Bluetooth PIN and serial number of the component is lost, contact your O&P professional.

7.1.1 Starting the cockpit app for the first time

- 1) Tap the symbol of the Cockpit app ().
→ The end user license agreement (EULA) is displayed.
- 2) Accept the end user license agreement (EULA) by tapping the **Accept** button. If the end user license agreement (EULA) is not accepted, the Cockpit app cannot be used.
→ The welcome screen appears.
- 3) Hold the prosthesis with the sole of the foot facing up, or connect and then disconnect the battery charger, in order to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
- 4) Tap the **Add component** button.
→ The Connection Wizard opens and guides you through the process of establishing a connection.
- 5) Follow the subsequent instructions on the screen.
- 6) After the Bluetooth PIN is entered, a connection to the component is established.
→ While the connection is being established, 3 beep signals sound and the  symbol appears.
The  symbol is displayed when the connection has been established.

→ Once the connection has been established, the data are read from the component. This process may take up to a minute.

Then the main menu appears with the name of the connected component.


INFORMATION

After the initial connection to the component has been established successfully, the app will connect automatically each time it is started. No further steps are required.

INFORMATION


After activating the "visibility" of the component (holding the component with the sole of the foot facing up, or connecting and then disconnecting the battery charger), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case, hold the component with the sole of the foot facing up again, or connect and then disconnect the battery charger.

7.2 Control elements for cockpit app

- ☰ Access the navigation menu (see page 62)
- Product
The component name can only be changed with the adjustment software.
- If connections to more than one component have been saved, you can switch between the saved components by tapping the **change** item (see page 62).
- MyModes configured with the adjustment software.
Switch the mode by tapping the corresponding symbol and confirming by tapping "**OK**".
- Currently selected mode
- Charge level of the component.
 - 🔋 Component battery fully charged
 - 🔌 Component battery empty
 - 🔌 Component battery chargingThe current charge level is also displayed in %.
- Display and name of the currently selected mode (e.g. **1. Basic Mode**)
- Connection to component has been established
(🔌) Connection to component has been interrupted. The app is attempting to re-establish the connection automatically.
(🔌) No existing connection to the component.

7.2.1 Cockpit app navigation menu

Tap the ☰ symbol in the menus to display the navigation menu. Additional settings for the connected component can be configured in this menu.

Product

Name of the connected component

MyModes

Return to the main menu to switch MyModes

Heel Height

Setting the heel height (see page 64)

Functions

Access additional functions of the component (e.g. turn off Bluetooth) (see page 72)

Settings

Change settings of the selected mode (see page 70)

Status

Query status of the connected component (see page 73)

Manage components

Add or delete components (see page 62)

Imprint/Info

Display information/legal notices for the Cockpit app

7.3 Managing components

Connections with up to four different components can be stored in the app. However, a component can only be connected to one mobile device at a time.

INFORMATION

Before establishing the connection, observe the points in the section "Initial connection between Cockpit app and component" (see page 60).

7.3.1 Adding component

- 1) Tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the "**Manage components**" item.
- 3) Hold the prosthesis with the sole of the foot facing up, or connect and then disconnect the battery charger, in order to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
- 4) Tap the "+" button.
→ The Connection Wizard opens and guides you through the process of establishing a connection.
- 5) Follow the subsequent instructions on the screen.
- 6) After the Bluetooth PIN is entered, a connection to the component is established.
→ While the connection is being established, three beep signals sound and the (📶) symbol appears.
The (📶) symbol is displayed when the connection has been established.

- Once the connection has been established, the data are read from the component. This process may take up to a minute.

The main menu will then appear with the name of the connected component.

INFORMATION

If establishing a connection to a component is not possible, perform the following steps:

- ▶ Delete the component from the Cockpit app if applicable (see the section "Deleting a component")
- ▶ Add the component again in the Cockpit app (see the section "Adding a component")

INFORMATION

After activating the "visibility" of the component (holding the component with the sole of the foot facing up, or connecting and then disconnecting the battery charger), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case, hold the component with the sole of the foot facing up again, or connect and then disconnect the battery charger.


7.3.2 Deleting a component

- 1) Tap the ☰ symbol in the main menu.
 - The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Tap the "**Edit**" button.
- 4) Tap the 🗑 symbol under the component you want to delete.
 - The component is deleted.

7.3.3 Connecting component with multiple mobile devices

The connection for a component can be stored on more than one mobile device. However, only one mobile device can be connected to the component at one time.

If there is an existing connection between the component and a different mobile device, the following information appears while the connection is being established with the current mobile device:

- ▶ Tap the "**OK**" button.
 - The connection to the last connected mobile device is broken off and established with the current mobile device.

8 Use

8.1 Setting the heel height

The heel height of the prosthesis must be set on a level surface. A sloped ground surface distorts the heel height measurement and leads to improper control of damping behaviour.


A heel height that is too high may cause control of the prosthetic foot to function incorrectly because of too little motion in the ankle joint. This occurs especially in case of small feet, heels that are moved forward, when walking down stairs and ramps, and standing on surfaces that are sloping downwards. The maximum heel height in the section "Technical data" must therefore be taken into account (see page 76).

8.1.1 Setting heel height using motion pattern

- 1) Put on shoes with the new heel height.
 - 2) Stretch out foot with the prosthetic foot sideways.
 - 3) Swing the foot to the side three times.
→ A beep signal will sound to confirm that the motion pattern has been recognised.
 - 4) Put feet at the same level and make sure the heel and toes are resting on the ground surface.
 - 5) Distribute weight on feet evenly.
- A confirmation signal will sound to indicate that the new heel height was successfully stored.

INFORMATION: If there is no feedback (e.g. beep signal), this indicates that the new heel height could not be stored. Repeat the heel height measurement.


8.1.2 Setting the heel height with the cockpit app

- 1) Once the component is connected and in the desired mode, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) Tap the **"Heel Height"** menu option.
- 3) Follow the on-screen instructions.
- 4) Tap the **"Set the heel height"** option.
- 5) Follow the subsequent instructions on the screen.

8.2 Movement patterns in basic mode (mode 1)

8.2.1 Standing


The intuitive stance function recognises any situation where the prosthesis is kept still while standing. The prosthesis prevents tipping forward, stabilising the user.

During rollover to the front or when the prosthesis is lifted off the floor, the prosthesis leaves the intuitive stance function automatically and switches back to the walking function.

When coming to a stop while walking, the leg must be set under the body and extended or weight has to be placed on the heel.

The relief function can be used when standing (see page 67).

8.2.2 Walking

Initial attempts at walking with the prosthesis always require the instruction of trained, qualified personnel. The prosthesis stabilises the user in the stance phase. The rollover behaviour automatically adjusts itself to the walking speed. Dropping the tip of the foot is prevented in the swing phase to maintain ground clearance. Even before ground contact, damping of the prosthesis is adjusted to make stepping down comfortable and quickly achieve full ground contact.

8.2.3 Sitting down/sitting

Sitting down

- 1) Place both feet side by side at the same level.
- 2) While sitting down, distribute weight evenly between both legs and use armrests, if available.
- 3) Move the buttocks in the direction of the backrest and lean the upper body forward.

Sitting

The relief function can be used when sitting; the tip of the foot lowers to reach a more natural foot position. (see page 67).

8.2.4 Standing up


- 1) Place the feet at the same level. Make sure the foot is positioned perpendicularly under the knee or shifted further forward, and that weight is distributed evenly over the feet.

INFORMATION: The ankle joint may lock if the prosthetic foot is set further back than at a perpendicular angle under the knee.

- 2) Lean the upper body forward.
- 3) Put the hands on armrests, if available.
- 4) Stand up with support from the hands, while keeping weight evenly distributed on the feet.


8.2.5 Walking up stairs

When the lower leg reaches a vertical position, the prosthesis prevents tipping forward in order to stabilise. Walking up stairs step-over-step is only possible if certain physical conditions are met. This function must be practised and executed consciously.

- 1) Hold the handrail with one hand.
- 2) Place the foot of the sound leg on the first step.
- 3) Pull up the leg with the prosthetic foot and position it on the step with the full surface of the sole.

8.2.6 Walking down stairs

This function must be practised and executed consciously. Only by properly stepping down with the sole can the system switch correctly and permit controlled rollover. The motion must be carried out in a continuous pattern in order to allow the motion to proceed in a fluid manner.

A stair function can be activated by using the adjustment software. See the following section for further information regarding the stair function.

- 1) Hold on to the handrail with one hand.
- 2) Position the leg with the prosthetic foot on the step so that as much of the foot as possible is on the step.


INFORMATION: Rolling the foot over the edge of the stair is not necessary.

- 3) Place the contralateral side onto the next step.
While doing so, check whether the knee joint and prosthetic foot will allow this motion.
- 4) Place the leg with the prosthetic foot on one step beyond the next step.
- 5) At the end of the stairs, take a bigger step when transitioning to level ground, so that the prosthetic foot switches over correctly from walking down the stairs to normal walking phase.

8.2.6.1 Stair function

The stair function increases the rollover angle when walking down stairs. This function should be switched on when walking down stairs step-over-step. If walking down stairs step-over-step is not required, this function can be switched off. Further information on switching on/off see page 71.


8.2.7 Walking down a ramp

The foot is already adjusted to the inclination of the ramp with the first step. This makes it possible for the heel to strike as you lower the tip of the foot so the entire sole of the foot is in contact with the surface during the rollover. Lowering the tip of the foot is restricted for walking with a prosthetic knee joint.


After the prosthetic foot is put down on the ramp, you should not use the knee to work against it, but should allow movement in the knee joint during heel strike (yielding). The prosthesis will then recognise the movement as walking.

8.2.8 Walking backwards

The foot enables plantar flexion from the stance phase when walking backwards. When the toe is subsequently put down, the ankle joint drops to the neutral position in the dorsiflexion direction.


8.2.9 Walking up a ramp

With the first step, the foot is already adjusted to the inclination of the ramp, making a rollover possible if the step is on the heel or middle of the foot. In order to do so, the lower leg should be positioned almost perpendicular to the surface of the ramp and the full area of the foot should be placed on the surface.

When stepping on the forefoot with the lower leg at a steep angle (e.g. on very steep ramps), the foot secures dorsiflexion, which enables the body to be lifted in a stable manner.

8.2.10 Standing on inclined surface

Standing on an inclined ground surface is no different to standing on a level surface. The foot stabilises in dorsiflexion when the lower leg is perpendicular. Put weight on the heel in order to lower the forefoot (e.g. when standing in a downward direction).

Perform one of the following motions in order to continue to walk on a downward-sloping ground surface from the standing position:

- Start the first step with the prosthesis side.
- Deliberately produce a rollover motion with the prosthesis side. The prosthetic foot then yields in dorsiflexion, in order to allow the body's centre of gravity to drop before the heel of the other leg strikes.

The relief function can be used when standing on an inclined surface (see page 67).


Wearing shoes with heels limits the range of incline, so it may not be possible to achieve a perpendicular lower leg.

8.2.11 Kneeling

If the leg with the joint is tilted backwards, plantar flexion damping is reduced, making it possible for the foot to bend so the lower leg can lie flatter to the ground surface.





8.2.12 Relief function

When weight is evenly placed on the heel without movement for more than 2 seconds, the tip of the foot lowers to reach a more natural foot position. Possible applications include sitting with the heel in front of the knee axis, standing while leaning against something, and standing on a downward sloping surface.

8.3 MyModes

With the help of adjustment software, the O&P professional can activate and configure MyModes in addition to the basic mode. These can be selected by using the Cockpit app or movement patterns. Switching by using movement patterns has to be activated in the adjustment software by the O&P professional.

The MyModes are intended for specific motion patterns or postures (e.g. golfing, ...). Settings can be adjusted using the Cockpit app (see page 72).

8.3.1 Switching MyModes with the cockpit app


INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (function only available in basic mode) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 72).

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 70).

Once a connection to a prosthesis has been established, the cockpit app can be used to switch between the MyModes.

- 1) Tap the symbol of the MyMode (1) you want in the main menu of the app.
→ A security question for changing the MyMode appears.
- 2) If you want to change the mode, tap the "OK" button.
→ A beep signal sounds to confirm the switch.
- 3) After switching, a symbol (2) is displayed to identify the active mode.
→ The current mode is also indicated by the name on the lower edge of the screen (3).

8.3.2 Switching MyModes using motion patterns

Information on switching

- Switching and the number of motion patterns must be activated in the adjustment software.
- Before further activities, always check whether the selected mode corresponds to the required motion type.
- If the **Volume** parameter is set to "0" in the Cockpit app, there are no beep signals (see page 70).

Switching process

- 1) Hold the prosthetic leg under the body.
- 2) Using the heel of the prosthetic foot, knock on a solid obstacle behind you (e.g. a wall) according to the number of the configured MyMode (MyMode 1 = three times, MyMode 2 = four times, MyMode 3 = five times). You can also knock against the tip of the shoe on the contralateral leg.
→ A beep and vibration signal will occur to confirm that the movement pattern has been recognised.
INFORMATION: If this beep and vibration signal is not emitted, the knocking was not recognised.
- 3) Lean the prosthetic foot slightly to the back and place weight on the forefoot.
INFORMATION: If the prosthetic foot is strongly dorsiflexed, weight can be placed on the heel.
→ A confirmation signal will sound to indicate that the foot successfully switched to the corresponding mode (two times = MyMode 1, three times = MyMode 2, four times = MyMode 3).
INFORMATION: If this confirmation signal does not occur, weight was not placed on the prosthetic foot correctly or not for enough time. Repeat the process for correct switching.
- 4) Remove the load from the prosthetic leg.
→ The mode has been changed.

8.3.3 Switching on ankle lock

Information on switching

- The ankle lock must be selected in the adjustment software as MyMode. In addition, the number of motion patterns used to switch it on must be activated in the adjustment software.
- Before further activities, always check whether the selected mode corresponds to the required motion type.
- If the **Volume** parameter is set to "0" in the Cockpit app, there are no beep signals (see page 70).

Switching process

- 1) Hold the prosthetic leg under the body.
- 2) Using the heel of the prosthetic foot, knock on a solid obstacle behind you (e.g. a wall) according to the number of the configured MyMode (MyMode 1 = 3 times, MyMode 2 = 4 times, MyMode 3 = 5 times). It is also possible to knock against the tip of the shoe of the contralateral leg.
→ A beep and vibration signal will occur to confirm that the movement pattern has been recognised.
- 3) Lean the prosthetic foot slightly to the back and place weight on the forefoot.
INFORMATION: If the prosthetic foot is strongly dorsiflexed, weight can be placed on the heel.
→ A confirmation signal will sound to indicate that the prosthesis has successfully switched to the corresponding mode (2 times = MyMode 1, 3 times = MyMode 2, 4 times = MyMode 3).
INFORMATION: If this confirmation signal does not occur, weight was not placed on the prosthetic foot correctly or not for enough time. Repeat the process to correctly switch to the required mode.
- 4) Remove the load from the prosthetic leg.
→ The mode has been changed.
- 5) Within 2 seconds, lower the prosthetic leg and assume the required position of the ankle angle.
→ At the end of this time, a signal will occur to indicate the locking of the ankle joint.

8.3.4 Switching from a MyMode back to basic mode

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 70).

Information on switching

- Regardless of the configuration of additional MyModes in the adjustment software, it is always possible to switch back to basic mode (mode 1) with a motion pattern.
- It is always possible to switch back to basic mode (mode 1) by connecting/disconnecting the battery charger.
- Before further activities, always check whether the selected mode corresponds to the required motion type.
- If the **Volume** parameter is set to "0" in the Cockpit app, there are no beep signals (see page 70).

Switching process

- 1) Hold the prosthetic leg under the body.
- 2) Using the heel of the prosthetic foot, knock on a solid obstacle behind you at least 3 times but no more than 5 times.
→ A beep and vibration signal will occur to confirm that the motion pattern has been recognised.
- 3) Lean the prosthetic foot slightly to the back and place weight on the forefoot.

INFORMATION: If the prosthetic foot is strongly dorsiflexed, weight can be placed on the heel.

→ A confirmation signal will sound to indicate that the prosthesis has successfully switched over to basic mode.

INFORMATION: If this confirmation signal does not occur, weight was not placed on the prosthetic foot correctly or not for enough time. Repeat the process to correctly switch to the required mode.

- 4) Remove the load from the prosthetic leg.
→ The mode has been changed.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

8.4 Changing prosthesis settings


Once an active connection to a component has been established, the settings **of the respective active mode** can be changed using the Cockpit app.

INFORMATION

Bluetooth on the prosthesis must be switched on to change the prosthesis settings.


If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. The connection must be established during this period.

Information for changing the prosthesis settings

- Before changing settings, always check the main menu of the Cockpit app to make sure the correct component has been selected. Otherwise parameters could be changed for the wrong component.
- It is not possible to change prosthesis settings nor to switch to a different mode while the prosthesis battery is being charged. Only the status of the prosthesis can be called up. Instead of the  symbol, the  symbol appears in the bottom row of the screen in the cockpit app.

- The O&P professional's setting is in the middle of the scale. After making adjustments, this setting can be restored by tapping the **"Standard"** button in the Cockpit app.
- Prosthesis settings should be optimised using the adjustment software. The Cockpit app is not intended for use by the O&P professional to set up the prosthesis. The everyday behaviour of the prosthesis can be changed to a certain extent using the app (e.g. while becoming accustomed to the prosthesis). The O&P professional can use the adjustment software to track these changes at the next appointment.
- If the settings of a MyMode are to be modified, one must first switch to this MyMode.

8.4.1 Changing the prosthesis setting using the cockpit app

- 1) Once the component is connected and in the desired mode, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) Tap the **"Settings"** menu option.
→ A list appears with the parameters for the currently selected mode.
- 3) Change the setting of the desired parameter by tapping the "<", ">" symbols.

INFORMATION: The O&P professional's setting is marked and, after the setting has been changed, can be restored by tapping the "Standard" button.

8.4.2 Overview of adjustment parameters in basic mode

The parameters in basic mode describe the dynamic behaviour of the prosthesis in a normal gait cycle. These parameters act as basic settings for automatically adjusting the damping behaviour to the current motion situation (e.g. ramps, slow walking speed, etc.).

The following parameters can be modified:

Parameter	Adjustment software range	Setting range, app	Meaning
Pitch	1000 Hz–4000 Hz	1000 Hz–4000 Hz	Pitch (frequency) of beep signal for confirmation tones
Volume	0–4	0–4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals are still generated when errors occur (see page 80).

Parameter	Adjustment software range	Setting range, app	Meaning
Heel Resist.	10–60	± 20	Plantar flexion damping. How fast the forefoot is lowered when weight is placed on the heel.
Rollover Resist.	110–170	± 10	This parameter defines how easy the rollover is.
Stair Function	ON–OFF	ON–OFF	Switching on this function increases the rollover angle when walking down stairs. In order to do so, this function needs to be enabled in the adjustment software.

8.4.3 Overview of adjustment parameters in MyModes

The parameters in the MyModes describe the static behaviour of the prosthesis for a specific motion pattern such as golf. Damping behaviour is not automatically controlled and adjusted in MyModes.

The following parameters can be modified in MyModes:

Parameter	Adjustment software range	Setting range, app	Meaning
Heel Resist.	0–195	± 20	Plantar flexion damping. How fast the forefoot is lowered when weight is placed on the heel.
Rollover Resist.	0–195	± 10	Dorsiflexion damping. How easily the value of the "Stop angle" parameter can be reached, in other words how strong the resistance is to reaching the value of the "Stop angle" parameter.
Stop angle	-200–200	± 10 displayed in 0.1°	Ankle angle from which motion in the rollover direction (in direction of dorsiflexion) is locked.

8.5 Turning Bluetooth on the prosthesis on/off

INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the Cockpit app.

If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (function only available in basic mode) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 72).

8.5.1 Switching Bluetooth off/on using the cockpit app

Switching off Bluetooth

- 1) When the component is connected, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".

- 3) Tap the entry "**Deactivate Bluetooth**".
- 4) Follow the on-screen instructions.

Switching on Bluetooth

- 1) Rotate the component or connect/disconnect the battery charger.
 - Bluetooth is switched on for approx. 2 minutes. During this time, the app must be started in order to establish a connection to the component.
- 2) Follow the on-screen instructions.
 - If Bluetooth is switched on, the (📶) symbol appears on the screen.

8.6 Querying the prosthesis status

8.6.1 Query status through cockpit app

- 1) When the component is connected, tap the ☰ symbol in the main menu.
- 2) In the navigation menu, tap the entry "**Status**".

8.6.2 Status display in the cockpit app

Menu option	Description	Possible actions
Trip: 1747	Daily step counter	Reset the counter by tapping the " Reset " button.
Step: 1747	Total step counter	For informational purposes only
Batt.: 68	Current prosthesis charge level, as a percentage	For informational purposes only

9 Additional operating states (modes)

9.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level drops to 0% (see page 80). During this time, damping settings are set to their safety mode values. The prosthesis is then switched off. You can switch back to basic mode (mode 1) from empty battery mode by charging the product.

9.2 Mode for charging the prosthesis

The ankle joint of the prosthetic foot is locked during the charging process.

9.3 Safety mode

The product automatically switches to safety mode if a critical error occurs (e.g. failure of a sensor signal) or if the battery is drained. Safety mode remains in effect until the error has been rectified.

Default damping values are activated in safety mode. This makes limited walking possible for the user even though the product is not active.

The switch to safety mode is indicated by beeps and vibration signals immediately prior to switching (see page 80).

Safety mode can be disabled by connecting then disconnecting the battery charger. If the product switches into safety mode again, this means a permanent error exists. The product must be inspected by an authorised Ottobock Service Centre.

9.4 Overheating mode

To prevent overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill), the range of motion in the ankle joint is restricted as the temperature rises. Depending on the temperature, this restriction may also result in the ankle joint being locked completely. When the hydraulic unit cools down, the product switches back to the settings that existed before the overheating mode.

Overheating mode is indicated by a brief vibration every 5 seconds.

10 Cleaning

- 1) Switch off the product before cleaning.
- 2) Clean the product with a damp cloth and mild soap when needed.
Make sure that no liquids get into the product and product components.
- 3) Dry the product with a lint-free cloth and allow it to air dry fully.

11 Maintenance

INFORMATION

The foot shell for the prosthetic foot has an expected working life of approx. one year when professionally mounted and properly used. Damaged foot shells must be replaced immediately, before the prosthetic foot is used again.

Regular maintenance (service inspections) at 24-month intervals is mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty, maintain basic safety and the essential performance characteristics, and ensure safety in regards to EMC. When maintenance is due, this is indicated by feedback after disconnecting the battery charger (see the section "Operating states/error signals", see page 79). The manufacturer grants a grace period of no more than one month before, or two months after, the due date.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components always have to be submitted to the O&P professional for maintenance and repairs:

The prosthesis, battery charger and power supply unit.

12 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

12.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

12.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

12.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address:
<http://www.ottobock.com/conformity>

12.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
3820 West Great Lakes Drive
Salt Lake City, Utah 84120-7205 USA
Phone + 1-801-956-2400
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

http://www.hc-sc.gc.ca/rpb.
 Responsible party:
 Otto Bock Healthcare Canada Ltd.
 5470 Harvester Road
 L7L 5N5 Burlington, Ontario
 Canada
 Phone + 1-800-665-3327

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

13 Technical data

Environmental conditions	
Storage and transport in original packaging (≤3 months)	-20 °C/-4 °F to +40 °C/+104 °F
Storage and transport without packaging (<48 hours)	-25 °C/-13 °F to +70 °C/+122 °F Max. 93% relative humidity, non-condensing
Long-term storage (>3 months)	-20°C/-4°F to +20°C/+68°F Max. 93% relative humidity, non-condensing
Operation	-10°C/+14°F to +40°C/+104°F Max. 93% relative humidity, non-condensing
Charging the battery	+10 °C/+50 °F to +45 °C/+113 °F

Product	
Reference number	1B1-2
Maximum adjustable heel height	50 mm/2 inches
Dorsiflexion with heel height of 1 cm/0.39 inches	14.5°
Plantar flexion with heel height of 1 cm/0.39 inches	22°
Mobility grade (MOBIS)	2–3
Footshell colours	Translucent, beige, brown
Max. system height with heel height of 2 cm/0.79 inches	18.5 cm/7.28 inches
Protection rating	IP54
Water resistance	Weatherproof but not corrosion-resistant Not designed for underwater use or submer- sion
Range of Bluetooth connection	Max. 10 m
Information on the product's ruleset and firm- ware version	Accessible via the Cockpit app navigation menu and the menu item "Imprint/Info"
Expected lifetime given compliance with the prescribed maintenance intervals	6 years
Test procedure (foot sizes 24 and 25)	ISO 22675-P5-100 kg/2 million load cycles
Test procedure (foot sizes 26 through 29)	ISO 22675-P6-125 kg/2 million load cycles

Foot size [cm]	24	25	26	27	28	29
Max. body weight	100 kg/220 lbs	125 kg/275 lbs	125 kg/275 lbs			

Foot size [cm]	24	25	26	27	28	29
Max. weight incl. footshell	approx. 1275 g/45 oz		approx. 1485 g/52 oz		approx. 1555 g/55 oz	

Data transfer	
Wireless technology	Bluetooth Smart Ready
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK, $\pi/4$ DQPSK, 8DPSK
Data rate (over the air)	2178 kbps (asymmetrical)
Maximum output power (EIRP):	+8.5 dBm

Prosthesis battery	
Battery type	Li-Ion
Charging cycles (charging and discharging cycles) after which at least 80% of the original battery capacity remains available	500
Charging time until battery is fully charged	8 hours
Behaviour of prosthetic foot while being charged	The ankle joint of the prosthetic foot is locked
Operating time of prosthesis with fully charged battery	1 day with average use

Power supply unit	
Reference number	757L16-4
Type	FW8001M/12
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Operation	0 °C/+32 °F to +50 °C/+122 °F Max. 95% relative humidity Air pressure: 70–106 kPa (up to 3,000 m without pressure equalisation)
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ===

Battery charger	
Reference number	4E50-2
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F Max. 93% relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F Max. 93% relative humidity, non-condensing
Input voltage	12 V ===
Lifetime	8 years

Cockpit app	
Reference number	4X441-IOS=*/4X441-Andr=V* Cockpit
Supported operating system	See the information in the respective online store (e.g. Apple App Store, Google Play Store, etc.) regarding compatibility with mobile devices and versions.
Website for download	https://www.ottobock.com/cockpitapp

14 Appendices

14.1 Symbols Used

Manufacturer

Compliance with the requirements according to "FCC Part 15" (USA)

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Non-ionising radiation

IP54

Protected against dust, protected against splashing water

DUAL

The product's Bluetooth wireless module can establish a connection to mobile devices with the following operating systems: iOS (iPhone, iPad, iPod...) and Android

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)
 YYYY – year of manufacture
 WW – week of manufacture
 NNN – sequential number

LOT Lot number (PPPP YYYY WW)
 PPPP – plant
 YYYY – year of manufacture
 WW – week of manufacture

MD Medical device

REF Article number

Please note the instructions for use

14.2 Operating states/error signals

The prosthesis indicates operating states and error messages through beeps and vibration signals.

14.2.1 Signals for operating states

Battery charger connected/disconnected

Beep signal	Vibration signal	Event
1x short		Battery charger is connected or battery charger already disconnected prior to start of charging mode
	3x short	Charging mode started (3 sec. after connection of battery charger)
1x short	1x before beep signal	Battery charger disconnected after start of charging mode

Mode switching

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 70).

Beep signal	Vibration signal	Additional action performed	Event
1 x short	1 x short	Mode switching using the Cockpit app	Mode switching is performed using the Cockpit app
1 x short	1 x short	Tapping with heel to switch modes or swing sideways three times to set heel height	Motion pattern recognised
1 x short	1 x short	Weight placed on prosthetic leg and leg kept still for 1 second to switch mode or feet placed on same level and weight distributed evenly on both feet to set heel height	Switching to basic mode (mode 1) carried out.

Beep signal	Vibration signal	Additional action performed	Event
2 x short	2 x short	Weight placed on prosthetic leg and leg kept still for 1 second	Switching to MyMode 1 (mode 2) carried out.
3 x short	3 x short	Weight placed on prosthetic leg and leg kept still for 1 second	Switching to MyMode 2 (mode 3) carried out.

14.2.2 Warnings/error signals

Error during use

Beep signal	Vibration signal	Event	Required action
-	1 x long at interval of approx. 5 seconds	Overheated hydraulic unit	Reduce activity.
-	3 x long	Battery charge level under 25%	Charge battery soon.
-	5 x long	Battery charge level under 15%	Charge battery immediately; the product will be switched off after the next warning signal.
10 x short	10 x long	Charge level 0% After the beep and vibration signals, the product switches to empty battery mode and then switches off.	Charge the battery.
30 x long	1 x long, 1 x short repeated every 3 seconds	Severe error/indication of safety mode activation e.g. one or more sensors are not operational.	Walking possible with restrictions. Please note the possible change in flexion/extension resistance. Attempt to reset this error by connecting/disconnecting the battery charger. The battery charger must be connected for at least 5 seconds before it is disconnected. If the error persists, use of the product is prohibited. The product must be inspected by an O&P professional.
-	Continuous	Total failure Electronic control no longer possible. Safety mode active or undetermined valve	Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohib-

Beep signal	Vibration signal	Event	Required action
		state. Undetermined product behaviour.	ited. The product must be inspected by an O&P professional.

Error while charging the product

LED on power supply	LED on battery charger	Error	Resolution
○	 ○ ○ 	Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.
		Non-functional socket	Check socket with another electric device.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●	 ○ ○ 	No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●	 ○ ● 	Battery is fully charged (or connection with product is interrupted).	<p>Take note of the confirmation signal for differentiation.</p> <p>When the battery charger is connected or disconnected, a self-test is conducted and confirmed by a beep and vibration signal. The battery is fully charged if this signal is heard.</p> <p>If no signal is emitted, the connection to the product is interrupted.</p> <p>If the connection to the product is interrupted, an authorised Ottobock Service Centre must inspect the product, battery charger and power supply.</p>

Beep signal	Error	Resolution
4 x short at intervals of approx. 20 sec. (continuously)	Charging the battery outside the allowable temperature range	Check whether the specified ambient conditions for charging the battery are met (see page 76).

14.2.3 Error messages while establishing a connection with the cockpit app

Error message	Cause	Correction
Component was connected to another device. Establish connection?	The component was connected to another mobile device	To disconnect the original connection, tap the "OK" button. If the original connection is not to be disconnected, tap the "Cancel" button.
Mode change failed	An attempt was made to switch to a different MyMode while the component was in motion (e.g. while walking)	For safety reasons, switching MyModes is only permitted when components are at rest, e.g. while standing or sitting.
	A current connection to the prosthesis was interrupted	Check the following points: <ul style="list-style-type: none"> • Distance from the prosthesis to the mobile device • Charge level of the prosthesis • Bluetooth on the prosthesis switched on? (see page 72) • Hold the component with the sole of the foot facing up to make the component "visible" for 2 minutes. • Prosthesis switched on? (Switching off the product) • If multiple prostheses were stored, was the correct prosthesis selected?

14.2.4 Status signals

Battery charger is connected

LED on power supply	LED on battery charger	Event
		Power supply and battery charger operational

Battery charger disconnected

Beep signal	Vibration signal	Event
1x short	1x short	Self-test completed successfully. Product is operational.

Battery charge level

Battery charger	
	Battery is being charged, battery charge level is less than 50%
	Battery is being charged, battery charge level is over 50%
	Battery is fully charged (or connection with product is interrupted). Take note of the confirmation signal for differentiation. When the battery charger is connected or disconnected, a self-test is conducted and confirmed by a beep and vibration signal.

Battery charger	
	The battery is fully charged if this signal is heard. If no signal is emitted, the connection to the product is interrupted.

14.3 Directives and manufacturer's declaration

14.3.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (see page 54).

Electromagnetic emissions

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
Harmonics according to IEC 61000-3-2	Not applicable – power below 75 W	–
Voltage fluctuations/flicker according to IEC 61000-3-3	Product meets the requirements of the standard.	–

Electromagnetic interference immunity

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air,
High-frequency electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate
Surges Line against line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Voltage drops	IEC 61000-4-11	0% U _T ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0% U _T ; 1 period and 70% U _T ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% U _T ; 250/300 periods

Interference resistance against wireless communication devices

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/90- 0, TETRA 800, iDEN 820, CDMA 850, GSM 800/90- 0, LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1,720	1,700 to 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400 to 2,570	Bluetooth WLAN 802.1- 1 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5,100 to 5,800	WLAN 802.1- 1 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

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