10V51

DE	Gebrauchsanweisung	5
EN	Instructions for use	11
FR	Instructions d'utilisation	16
IT	Istruzioni per l'uso	22
ES	Instrucciones de uso	28
PT	Manual de utilização.....	34
NL	Gebruiksaanwijzing.....	40
SV	Bruksanvisning.....	46
DA	Brugsanvisning.....	51
NO	Bruksanvisning.....	57
TR	Kullanma talimatı.....	62
EL	Οδηγίες χρήσης	68
JA	取扱説明書.....	74

INFORMATION

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INFORMATION

Last update: 2020-09-15

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

Before using the product, you should become familiar with the handling, function and use of the product and **then practice its use with child.**

Explanation of symbols

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting/use.

The scope of delivery is shown on the front page.

1 Single components

1.1 Single components

For available single components refer catalogue.

2 Description

2.1 Indications for use

The 10V51 MyolinoWrist 2000 is intended exclusively for the exoprosthetic fitting of the upper limbs in combination with the 8E51 Ottobock electric hand 2000 and the 9E420=L/9E420=R 7in1 Controller or 9E369/9E370 4in1 Controller LS.

2.2 Field of application

The MyolinoWrist 2000 is used as flexible wrist joint with an Ottobock Electrohand 2000. It is suitable mainly for children aged from 3 to 13 years. It supports the physiologically correct posture when handling the prosthesis and particularly reduces the need for unnatural compensating movements when performing longer monotonous activities.

2.3 Safety instructions

Please forward the following safety instructions to your patients:

CAUTION

Risk of injuries caused by overuse. The Ottobock prosthetic components described in this manual have been developed for everyday life activities and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service Centre.

NOTICE

Risk of corrosion. Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids.

Using this medical product in such environmental conditions will render all claims against Ottobock HealthCare null and void.

NOTICE

Risk of damage due to improper environmental conditions. The prosthesis and its components should not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthesis and its components. Failure to follow these instructions can lead to malfunctioning and damage of the prosthesis.

NOTICE

Risk of damage due to use of improper cleansing agents. Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

2.4 Function

The MyolinoWrist 2000 is available in the following designs:

10V51=1 for lamination into a forearm socket.

10V51=2 for bonding into a forearm component with a \varnothing 40 mm wrist connection.

The MyolinoWrist 2000 is a ball joint which allows a continuous flexion of the wrist joint up to 30° in all directions. A passive pronation and supination is provided by the 8E51 Ottobock Electrohand 2000.

3 Technical data

Article number	10V51=1	10V51=2
Flexion	continuous up to 30°	
Weight [g]	47	50
Total height [mm]	32	40
Wrist joint- \varnothing [mm]	40	
Maximum load capacity	32.2 Nm $\hat{=}$ ~27 kg	
System length [mm]	22	24

4 Handling

4.1 Assembly preparations

Disassemble the upper part of the 10V51 MyolinoWrist 2000 with the 709S42 Allen Wrench included (fig. 1). Remove the O-ring from the lamination ring guide (fig. 2).

4.2 Assembly

4.2.1 Assembly 10V51=1

1) Close threaded holes on distal end of the 11D3 Lamination Ring with 636K8 Plasta Tape. Apply 633F11 Silicone Grease (fig. 3) and snap

- the lamination dummy in place audibly. Remove excess silicone grease.
- 2) Close threaded holes on proximal end of the lamination ring with 636K8 Plasta Tape (fig. 4). Fill inner contours with 636K8 Plasta Tape (fig. 5).
 - 3) Place lamination ring on the prepared contoured foam part isolated by 99B81 PVA Bag (fig. 6). Use the circumferential groove to tie off the 623T3 Perlon Stockinette (fig. 7). Additionally apply the carbon fibre cloth for reinforcement of the outer socket.
 - 4) Laminate the outer socket while observing the manufacturer's processing instructions.
 - 5) After the resin has cured, take off the PVA bag and cut the lamination ring free (fig. 8). Remove the resin from lamination dummy and remove the lamination dummy from the socket.
 - 6) Cut off the excess laminating resin. For proper function of the Myolino Wrist 2000, the laminating resin on the outer brim of the lamination ring (fig. 9, arrow) should be removed.

4.2.2 Assembly 10V51=2

CAUTION

Risk of injury due to safety-relevant ingredients. Before processing the 636K18 Orthocryl® Sealing Resin Compact Glue and the 617H14 Hardener Paste, read the relevant safety data sheet.

INFORMATION

Before bonding the lamination ring, shorten the forearm component to the required measurement.

- 1) Mix 636K18=1 Orthocryl® Sealing Resin Compact Glue and 617P14 Hardener Paste in accordance with the manufacturer's specifications.
- 2) Roughen the socket inner side and clean with 634A3 Acetone. Allow the acetone evaporate and apply Sealing Resin Compact Glue.
- 3) Clean the knurled edge of the lamination ring with acetone. Allow the acetone evaporate and apply Sealing Resin Compact Glue (fig. 10).
- 4) Insert the lamination ring into the socket. Remove excess Sealing Resin Compact Glue. Let the adhesive joint harden.

4.2.3 Assembling the contact plate with cable for the 7in1 and 4in1 Controllers

- 1) Loosen the three screws connecting the adapter ring with the attachment using the Phillips screwdriver (fig. 11).
- 2) Place the 711M50-1 / 711M50 assembly adapter into the 7in1 Controller / 4in1 Controller LS.
- 3) The contact plate should audibly engage in the adapter ring (fig. 12).
- 4) Guide the cable of the 7in1 Controller / 4in1 Controller LS through the hole in the attachment (Fig. 13).
- 5) Screw the adapter ring to the attachment again (fig. 11).

4.2.4 Wrist unit assembly

- 1) Re-insert the O-ring into the guide (fig. 2).
- 2) Insert the cable centrally through the lamination ring and put on the upper part of the MyolinoWrist 2000. Mount the wrist unit with the oval countersunk head screws (fig. 14).

The friction force of the MyolinoWrist 2000 can be adjusted by evenly tightening or loosening the screws (fig. 14).

Please provide your patient or his guardian / attendant with the Allen wrench included so they can readjust the friction force according to the situation.

NOTICE

Risk of damage due to loosened screws. The countersunk head screws are designed so they should not loosen during use. Nevertheless, when this happens, the screws must be retightened. Therefore, the friction force of the wrist unit must be checked on regular basis.

5 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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