
## 21A47=\*<sup>\*</sup>

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b> Instructions for use (qualified personnel) .....	12
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	21
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	31
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	40
<b>PT</b> Manual de utilização (Pessoal técnico) .....	50
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	60
<b>SV</b> Bruksanvisning (Fackpersonal) .....	69
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	78
<b>NO</b> Bruksanvisning (Fagpersonell) .....	88
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	97
<b>CS</b> Návod k použití (Odborný personál) .....	106
<b>HR</b> Upute za uporabu (Stručno osoblje) .....	115
<b>JA</b> 取扱説明書 (有資格担当者) .....	125

	Größe / Size	Armlochumfang / armhole size	
		cm	inch
	S	38 - 45,3	15 - 17.8
	M	45,3 - 52,6	17.8 - 20.7
	L	52,6 - 60	20.7 - 23.6

	Größe / Size	Rückenbreite / back width	
		cm	inch
	S	30 - 36,6	11.8 - 14.4
	M	36,6 - 43,2	14.4 - 17
	L	43,2 - 50	17 - 19.7

## INFORMATION

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### INFORMATION

Date of last update: 2020-08-20

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The product "transhumeral soft harness" is referred to simply as the product/support in the following.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.


## 2 Product description

### 2.1 Function

The product is for additional fixation of the prosthetic socket.

### 2.2 Design

The product consists of the following components:

1. Strap
2. Shoulder flap
3. Hook-and-loop shoulder closure
4. Sleeve
5. Hook-and-loop sleeve closure
6. Shoulder loop

### 3 Intended use

#### 3.1 Indications for use

The product is intended for upper limb exoprosthesis fittings on youths and adults.

#### 3.2 Conditions of use

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

The product has been developed for use in everyday activities. It must not be used for activities where it would be exposed to unallowable environmental conditions or strains (e.g. extreme sports).

The product can be used with myoelectric or passive prostheses.

#### 3.3 Indications

- For patients with unilateral transhumeral amputation.

#### 3.4 Contraindications

- Transradial amputation
- Amputation in the shoulder area
- Bilateral amputation
- Fitting with a body-powered prosthesis

#### 3.5 Further Usage Restrictions

- Patients with measurements outside the planned sizes S, M and L

#### 3.6 Qualification

The fitting of a patient with the product may only be carried out by an O&P professional.

## 4 Safety

### 4.1 Explanation of warning symbols

Caution indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

Warning regarding possible technical damage.

## 4.2 Structure of the safety instructions

### CAUTION

#### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

## 4.3 General safety instructions

### CAUTION

#### **Non-observance of safety instructions**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety instructions and the stated precautions in this accompanying document.

### CAUTION

#### **Use of the product on damaged skin**

Possible skin irritation (reddening) and pressure points due to friction of the product.

- ▶ Only use the product on healthy skin.
- ▶ Increased attention is required when using the product in areas with sensitive skin.

### CAUTION

#### **Mechanical damage to the product**

Injury due to faulty control or malfunction of the prosthesis system resulting from insufficient contact between the electrodes in the prosthetic socket and the user's skin.

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use (e.g., hook-and-loop closures, etc.)
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### NOTICE

#### **Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product only as described in the section "Cleaning and care" (Cleaning and care).
- ▶ The hook-and-loop closure must always be fastened and the product placed in a laundry bag before washing.
- ▶ Be sure to clean the product in order to avoid possible skin irritation.

### NOTICE

#### **Improper assembly/disassembly of the product**

Damage to the product.

- ▶ Observe the process steps, instructions and required tools in the section "Preparation for use".
- ▶ Instruct the patient in the proper disassembly/ assembly of the product.

## NOTICE

### Independent user changes or modifications made to the product

Damage to the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.

## 5 Scope of Delivery and Accessories

### Scope of Delivery

- 1 pc. Transhumeral soft harness L/R
- 1 pc. Sewing silk (100 m)
- 1 pc. Instructions for use (qualified personnel)
- 1 pc. Instructions for use (user)

### Accessories

- 1 pc. 21A48\* sling

## 6 Preparing the product for use

### 6.1 Preparation


- 1) Lay the product flat on the table.

- 2) Thread the loose end of the strap through the eyelet of the shoulder loop from above.

**INFORMATION: Make sure the strap is not twisted.**

- 3) Lay the loose end of the strap onto the back of the sleeve.

- 4) Secure the loose end to the sleeve with pins.

**INFORMATION: The strap should overlap the sleeve by about 7 cm.**

### 6.2 Application


Complete the following steps to put on the product:

- 1) Open the hook-and-loop closure on the sleeve.
- 2) Hold the prosthesis between your legs (the fingers of the prosthesis face up).
- 3) Pull the sleeve of the product over the prosthesis.

- 4) Turn the prosthesis over (the fingers of the prosthesis face down).
- 5) Pull the sleeve of the product over the socket.

- 6) Note correct positioning of the product.

- 7) Fold up the sleeve of the product and keep the donning hole clear.
- 8) Fasten the hook-and-loop closure on the sleeve.

- 9) Put on the prosthesis.  
**INFORMATION: Make sure that no part of the product is drawn in during donning.**

- 10) Lay the shoulder flaps onto the shoulder.  
11) Open the hook-and-loop closure on the shoulder.

- 12) With the other hand, put on the shoulder loop.  
13) Fasten the hook-and-loop closure on the shoulder.

### 6.3 Fitting

The correct size is determined using the table (see page 2).

If a size cannot be clearly established (arm hole circumference, width of the back), the arm hole circumference is the decisive measurement. The back width can be adjusted via the strap.

Complete the following steps to adjust the product:

- 1) Have the patient stand.  
→ The patient has to stand.
- 2) Put the product on the patient.  
**INFORMATION: The hook-and-loop closure on the shoulder must be fastened so it is centred in the area of the sleeve.**
- 3) Check the strap position again. The strap should overlap the shoulder flap by at least 7 cm.
- 4) Tighten the strap and secure it in coordination with the patient. Shorten the hook-and-loop sleeve closure if necessary.

**INFORMATION: Make sure the strap is not twisted.**

- 5) Trim any excess strap.
- 6) Remove the product.

- 7) Melt the cut edge of the strap.  
**INFORMATION: All fibres must be melted to prevent fraying of the strap.**

- 8) Sew on the strap in a rectangle. Thread the sewing machine with the supplied silk thread.
- 9) Repeat this process twice. The start and end of the seam are secured by double sewing.

**INFORMATION: On a sewing machine, select the "straight" stitch with a stitch length of "2.5".**

- 10) Check the seam on the back of the shoulder flap to make sure there is an even stitch pattern.
- 11) Trim the threads as short as possible without damaging the seam.

- 12) Check the strength of the seam by pulling on the strap.
- 13) Put the product on the patient. Check that the product fits firmly and comfortably.

## 6.4 Removal

Complete the following steps to take off the product:

- 1) Open the hook-and-loop closure on the shoulder.
- 2) Pull the arm out of the shoulder loop.
- 3) Fasten the hook-and-loop closure on the shoulder.
- 4) Pull the residual limb out of the socket.

### **INFORMATION**

For putting it on the next time, the sleeve of the product does not have to be pulled off the prosthesis.

## 7 Cleaning and Care

- **INFORMATION: Fasten all hook and loop closures before washing.**

- 1) Wash the product at 30 °C/86 °F in a washing machine (setting: easy care or delicates).
- 2) Wash the product as recommended with a mild detergent.
- 3) Allow the product to fully air dry.

**INFORMATION**

Note the washing symbols on the label.

**INFORMATION**

Drying may take between 12 and 24 hours. A second product for changing is recommended.

## 8 Legal information

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

### 8.3 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 9 Technical data

<b>Ambient conditions</b>	
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F Max. 93% relative humidity, non-condensing
Operation	-10 °C/+14 °F to +70 °C/+158 °F Max. 93% relative humidity, non-condensing

<b>General information</b>	
Reference number	21A47=*
Weight	190 g
Lifetime	6 months (wear and tear part)

<b>Load limits</b>	
Maximum vertical load (weight of the prosthesis + permitted loading)	6.5 kg (14.3 lbs) Corresponds to 65 N

## 10 Appendix

### 10.1 Symbols on the product

Declaration of conformity according to the applicable European directives

Lot number (PPPP YYYY WW)  
PPPP – plant  
YYYY – year of manufacture  
WW – week of manufacture

Please note the instructions for use

Do not dry clean

Do not tumble dry

Do not bleach

Gentle cycle at 30°

Iron at medium temperature, max. 150°

Not made with natural rubber latex.

Natural rubber latex was not used as a material

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