## Genium X3 - Bionic Prosthetic System

3B5-X3/3B5-X3=ST,

Genium X3 3B5-2/3B5-2=ST

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	57
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	111
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	169
<b>ES</b> Instrucciones de uso (Usuario) .....	225
<b>PT</b> Manual de utilização (Usuário) .....	281
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	337
<b>HR</b> Upute za uporabu (Korisnik) .....	393
<b>TR</b> Kullanma talimatı (Kullanıcı) .....	445
<b>JA</b> 取扱説明書(ユーザー用) .....	497

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
	<b>Important information for the products 3B5-X3/3B5-X3=ST and 3B5-2/3B5-2=ST</b> .....	<b>60</b>
<b>1</b>	<b>Foreword</b> .....	<b>61</b>
<b>2</b>	<b>Product description</b> .....	<b>61</b>
2.1	Design .....	61
2.2	Function .....	61
<b>3</b>	<b>Intended use</b> .....	<b>62</b>
3.1	Indications for use .....	62
3.2	Conditions of use .....	62
3.3	Indications .....	62
3.4	Contraindications .....	63
3.5	Qualification .....	63
<b>4</b>	<b>Safety</b> .....	<b>63</b>
4.1	Explanation of warning symbols .....	63
4.2	Structure of the safety instructions .....	63
4.3	General safety instructions .....	63
4.4	Information on the Power Supply/Battery Charging .....	66
4.5	Battery charger information .....	66
4.6	Information on Proximity to Certain Areas .....	67
4.7	Information on Use .....	68
4.8	Notes on the safety modes .....	70
4.9	Instructions for use with an osseointegrated implant system .....	71
4.10	Information on the use of a mobile device with the cockpit app .....	71
<b>5</b>	<b>Scope of Delivery and Accessories</b> .....	<b>72</b>
5.1	Scope of delivery .....	72
5.2	Accessories .....	72
<b>6</b>	<b>Charging the battery</b> .....	<b>72</b>
6.1	Connecting the power supply and battery charger .....	73
6.2	Charging the prosthesis battery .....	73
6.3	Display of the current charge level .....	74
6.3.1	Display of battery charge level without additional devices .....	74
6.3.2	Display of the current charge level using the Cockpit app .....	74
6.3.3	Display of the current charge level during the charging process .....	75
6.3.4	Display of current charge level using remote control (optional accessory) .....	75
<b>7</b>	<b>Cockpit app</b> .....	<b>76</b>
7.1	System Requirements .....	76
7.2	Initial connection between cockpit app and prosthesis .....	76
7.2.1	Starting the cockpit app for the first time .....	77
7.3	Control elements for cockpit app .....	78
7.3.1	Cockpit app navigation menu .....	78
7.4	Managing components .....	79
7.4.1	Adding component .....	79

7.4.2	Deleting a component .....	79
7.4.3	Connecting component with multiple mobile devices .....	80
<b>8</b>	<b>Use .....</b>	<b>80</b>
8.1	Movement patterns in basic mode (mode 1) .....	80
8.1.1	Standing .....	80
8.1.1.1	Stance function .....	80
8.1.2	Walking.....	81
8.1.3	Running short distances ("walk-to-run" function) .....	81
8.1.4	Sitting down .....	81
8.1.5	Sitting .....	81
8.1.5.1	Sitting function.....	82
8.1.6	Standing up.....	82
8.1.7	Walking up stairs step-over-step.....	82
8.1.8	Overcoming obstacles .....	83
8.1.9	Walking down stairs .....	84
8.1.10	Walking down a ramp .....	84
8.2	Changing prosthesis settings.....	84
8.2.1	Changing the prosthesis setting using the cockpit app.....	85
8.2.2	Overview of adjustment parameters in basic mode.....	85
8.2.3	Overview of adjustment parameters in MyModes .....	86
8.3	Turning Bluetooth on the prosthesis on/off.....	87
8.3.1	Switching Bluetooth off/on using the cockpit app .....	87
8.4	Querying the prosthesis status .....	88
8.4.1	Query status through cockpit app.....	88
8.4.2	Status display in the cockpit app .....	88
8.5	Mute mode (silent mode) .....	88
8.5.1	Turning mute mode on/off using the Cockpit app .....	88
8.6	Deep sleep mode .....	88
8.6.1	Turning deep sleep mode on/off using the Cockpit app.....	89
8.7	OPG function (Optimised Physiological Gait).....	89
<b>9</b>	<b>MyModes .....</b>	<b>90</b>
9.1	Running function as configured MyMode.....	90
9.2	Switching MyModes with the cockpit app .....	90
9.3	Switching MyModes using motion patterns.....	91
9.4	Switching from a MyMode back to basic mode.....	92
<b>10</b>	<b>Additional operating states (modes) .....</b>	<b>93</b>
10.1	Empty battery mode .....	93
10.2	Mode for charging the prosthesis .....	93
10.3	Safety mode.....	93
10.4	Overheating mode .....	94
<b>11</b>	<b>Storage and bleeding .....</b>	<b>94</b>
<b>12</b>	<b>Cleaning .....</b>	<b>94</b>
<b>13</b>	<b>Maintenance .....</b>	<b>94</b>

<b>14</b>	<b>Legal information .....</b>	<b>95</b>
14.1	Liability .....	95
14.2	Trademarks .....	95
14.3	CE conformity .....	95
14.4	Local Legal Information .....	95
<b>15</b>	<b>Technical data .....</b>	<b>97</b>
<b>16</b>	<b>Appendices .....</b>	<b>99</b>
16.1	Symbols Used.....	99
16.2	Operating states/error signals .....	100
16.2.1	Signals for operating states.....	100
16.2.2	Warnings/error signals .....	101
16.2.3	Error messages while establishing a connection with the cockpit app.....	103
16.2.4	Status signals.....	104
16.3	Directives and manufacturer's declaration.....	105
16.3.1	Electromagnetic environment .....	105

## Important information for the products 3B5-X3/3B5-X3=ST and 3B5-2/3B5-2=ST

Using the 4X350-X3 remote control is no longer possible on products with the marking "DUAL"

# 1 Foreword

## INFORMATION

Date of last update: 2020-09-30

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The product "Genium X3 - Bionic Prosthetic System" or "Genium X3" is called the product/prosthesis/knee joint/component in the following.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

The product consists of the following components:

1. Proximal pyramid adapter
2. Optional flexion stops
3. Rechargeable battery
4. Hydraulic unit
5. LED (blue) as indicator for the Bluetooth connection
6. Receiver of the inductive charging unit

### 2.2 Function

This product features microprocessor control of the stance and swing phase.

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the damping behaviour of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

The product can be individually adapted to your needs with adjustment software.

The product features MyModes for special motion types (e.g. cycling etc.). These are pre-configured by the orthotics and prosthetics professional (O&P professional) using the adjustment software and can be activated with special movement patterns and the Cockpit app (see page 90).

In case of a product malfunction, safety mode makes restricted operation possible. Resistance parameters that are predefined by the product are configured for this purpose (see page 93).

Empty battery mode permits safe walking when the battery is drained. Resistance parameters that are predefined by the product are configured for this purpose (see page 93).

**The microprocessor-controlled hydraulic unit offers the following advantages**

- Approximation of the physiological gait pattern
- Stability while standing and walking
- Adaptation of product characteristics to various surfaces, inclines, gait situations and walking speeds

**Essential performance of the product**

- Stability in the stance phase
- Adjustable swing phase extension resistance

**3 Intended use**

**3.1 Indications for use**

The product is to be used **solely** for lower limb exoprosthetic fittings.

**3.2 Conditions of use**

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see page 97).

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of **up to 125 kg (275 lbs)**.

**3.3 Indications**

- For patients with knee disarticulation, transfemoral amputation and hip disarticulation (patients with hip disarticulation or hemipelvectomy must be fitted with the 7E10=\* Helix<sup>3D</sup> hip joint).
- For unilateral or bilateral amputation
- Dysmelia patients with residual limb characteristics corresponding to knee disarticulation, transfemoral amputation or hip disarticulation
- The patient must fulfil the physical and mental requirements for perceiving visual/acoustic signals and/or mechanical vibrations.

### 3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.5 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

#### Only valid for the products 3B5-2/3B5-2=ST

If the product is to be connected to an osseointegrated implant system, the qualified personnel must also be authorised for the connection to the osseointegrated implant system.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

 <b>CAUTION</b>
<p><b>The heading describes the source and/or the type of hazard</b></p> <p>The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:</p> <ul style="list-style-type: none"><li>&gt; E.g.: Consequence 1 in case of failure to observe the hazard</li><li>&gt; E.g.: Consequence 2 in case of failure to observe the hazard</li></ul> <p>▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.</p>

### 4.3 General safety instructions

 <b>WARNING</b>
<p><b>Non-observance of safety notices</b></p> <p>Personal injury/damage to the product due to using the product in certain situations.</p> <p>▶ Observe the safety notices and the stated precautions in this accompanying document.</p>

 <b>WARNING</b>
<p><b>Using the prosthesis while operating a vehicle</b></p> <p>Accident due to unexpected behaviour of the prosthesis because of changed damping behaviour.</p> <ul style="list-style-type: none"><li>▶ All users are required to observe their country's national and state driving laws when operating vehicles with a prosthesis. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.</li><li>▶ Observe national legal regulations for retrofitting your vehicle in accordance with the type of fitting.</li><li>▶ The leg on which the prosthesis is worn may not be used to control the vehicle or its peripheral components (e.g. clutch pedal, brake pedal, gas pedal, etc.).</li></ul>

**⚠ WARNING**

**Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

**⚠ CAUTION**

**Failure to observe warning/error signals**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 101) and corresponding change in damping settings must be observed.

**⚠ CAUTION**

**Failure to observe activated mute mode (silent mode)**

Falling due to unexpected product behaviour because of changed damping behaviour.

The following feedback signals are deactivated when mute mode is activated:

- > Long vibration signal if the hydraulic unit overheats.
- > Beep and vibration signal to confirm that the movement pattern has been recognised (switching to a MyMode/basic mode with movement pattern).
- > Beep and vibration signal to indicate successful switching to a MyMode/basic mode.
- > Beep and vibration signal upon successfully switching to deep sleep mode.
- ▶ Before activating mute mode, note that these feedback signals will be deactivated. For more information about mute mode, see the section "Mute mode" (see page 88).
- ▶ After switching to a MyMode/basic mode, verify the changed damping settings.
- ▶ Ensure that you stand securely during all switching processes.
- ▶ To deactivate mute mode, connect and then disconnect the battery charger.

**⚠ CAUTION**

**Independent manipulation of the product and the components**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by authorised, qualified Ottobock personnel (no replacement by the user).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

**⚠ CAUTION**

**Mechanical stress on the product**

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

**⚠ CAUTION**

**Use of the product when battery charge level is too low**

Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the prosthesis if required.
- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

**⚠ CAUTION**

**Risk of pinching in the joint flexion area**

Injuries due to pinching of body parts.

- ▶ Ensure that fingers/body parts or soft tissue of the residual limb are not in this area when bending the joint.

**⚠ CAUTION**

**Penetration of dirt and humidity into the product**

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that no solid particles or foreign objects can penetrate into the product.
- ▶ The knee joint and the AXON tube adapter are waterproof, corrosion-resistant and protected against the penetration of jets of water. The knee joint and AXON tube adapter can be operated in fresh and salt water. Do not use the knee joint under extreme conditions like diving or jumping into water. The knee joint and AXON tube adapter are designed for underwater use (for the maximum duration and water depth, see the section "Technical data" (see page 97).
- ▶ After contact with water, hold the prosthesis with the sole of the foot facing up until the water has drained from the knee joint and AXON tube adapter.
- ▶ After using the knee joint in a salt water environment, remove the Protective Cover and rinse knee joint, AXON tube adapter and Protective Cover with fresh water. Dry the knee joint and components with a lint-free cloth and allow the components to fully air dry.
- ▶ Should the knee joint or AXON tube adapter come into contact with **solutions other than fresh or salt water**, **promptly** remove the Protective Cover and clean the **knee joint**. In order to do so, rinse knee joint, AXON tube adapter and Protective Cover with fresh water and let them dry.
- ▶ In case of a malfunction after drying, the knee joint and AXON tube adapter must be inspected by an authorised Ottobock Service Center.
- ▶ The knee joint and the AXON tube adapter are not protected against the penetration of steam.

**⚠ CAUTION**

**Use of the product without protector or with damaged protector**

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ If the protector has been removed, ensure prior to the next use of the product that the protector has been properly installed.
- ▶ Using the product with a damaged protector or without the protector is not permitted.
- ▶ Using the product with a foam cover is not possible, since the protector would have to be removed for this purpose.

### CAUTION

#### **Signs of wear and tear on the product components**

Falling due to damage or malfunction of the product.

- ▶ Regular service inspections (maintenance) are mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty.

### **NOTICE**

#### **Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product with a damp cloth only (fresh water).
- ▶ Only use fresh water with a temperature of no more than 65 °C for cleaning.
- ▶ If dirt cannot be removed, the product must be sent to an authorised Ottobock Service Center. Your contact is the O&P professional.

### **INFORMATION**

#### **Knee joint movement noise**

When using exoprosthetic knee joints, servomotor, hydraulic, pneumatic or brake load dependent control functions can cause movement noises. This kind of noise is normal and unavoidable. It generally does not indicate any problems. If movement noise increases noticeably during the lifecycle of the knee joint, the knee joint should be inspected by the O&P professional immediately.

## **4.4 Information on the Power Supply/Battery Charging**

### CAUTION

#### **Charging the product without taking it off**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ For safety reasons, take the product off before charging the battery.

### CAUTION

#### **Charging the product with damaged power supply unit/charger/charger cable**

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

### **NOTICE**

#### **Use of incorrect power supply unit/battery charger**

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

## **4.5 Battery charger information**

### WARNING

#### **Storing/transporting the product near active implanted systems**

Interference with active implantable systems (e.g. pacemaker, defibrillator, etc.) due to the product's magnetic field.

- ▶ When storing/transporting the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety notices stipulated by the manufacturer of the implant.

#### NOTICE

##### **Penetration of dirt and humidity into the product**

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

#### NOTICE

##### **Mechanical stress on the power supply/battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

#### NOTICE

##### **Operating the power supply unit/charger outside of the permissible temperature range**

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see page 97).

#### NOTICE

##### **Independent changes or modifications carried out to the battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

#### NOTICE

##### **Contact of the battery charger with magnetic data storage devices**

Wiping of the data storage device.

- ▶ Do not place the battery charger on credit cards, diskettes, audio or video cassettes.

## 4.6 Information on Proximity to Certain Areas

#### **⚠ CAUTION**

##### **Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.50 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.35 m
  - DECT cordless phones incl. base station: 0.18 m
  - WiFi (routers, access points,...): 0.11 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.11 m

**⚠ CAUTION**

**Operating the product in very close proximity to other electronic devices**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

**⚠ CAUTION**

**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for people (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations...).
- ▶ If this cannot be avoided, make sure to at least have a safeguard when walking or standing (e.g. a handrail or the support of another person).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected changes in the damping behaviour of the product.
- ▶ In general, monitor the product for unexpected changes in the damping behaviour when electronic or magnetic devices are in the immediate vicinity.

**⚠ CAUTION**

**Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)**

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of strong magnetic fields.
- ▶ Take off the product before entering a room or area with strong magnetic fields and store the product outside this room or area.
- ▶ Damage to the product caused by exposure to strong magnetic fields cannot be repaired.

**⚠ CAUTION**

**Remaining in areas outside the allowable temperature range**

Falling due to malfunction or the breakage of load-bearing product components.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see page 97).

## 4.7 Information on Use

**⚠ CAUTION**

**Walking up stairs**

Falling due to foot placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking up stairs and place most of the area of the sole of your foot on the stair surface.
- ▶ Particular caution is required when carrying children up stairs.

**⚠ CAUTION**

**Walking down stairs**

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and roll over the edge of the step with the middle of the shoe.
- ▶ Observe the warning/error signals (see page 101).
- ▶ Be aware that resistance in the flexion and extension direction can change in case of warning and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

**⚠ CAUTION**

**Overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)**

- > Falling due to unexpected behaviour of the product because of switching into overheating mode.
- > Burns due to touching overheated components.
- ▶ Be sure to pay attention when pulsating vibration signals start. They indicate the risk of overheating.
- ▶ As soon as these pulsating vibration signals begin, you must reduce your level of activity so the hydraulic unit can cool down.
- ▶ Full activity may be resumed after the pulsating vibration signals stop.
- ▶ If the activity level is not reduced in spite of the pulsating vibration signals, this could lead to the hydraulic element overheating and, in extreme cases, cause damage to the product. In this case, the product should be inspected for damage by an O&P professional. If necessary, they will forward the product to an authorised Ottobock Service Center.

**⚠ CAUTION**

**Overloading due to activities with unusual loads**

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ The product was developed for everyday use and must not be used for activities with unusual loads. These unusual activities include, for example, extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional. If necessary, he or she will forward the product to an authorised Ottobock Service Center.

**⚠ CAUTION**

**Improper mode switching**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.
- ▶ Take the weight off the product and correct the switching, if required.

**⚠ CAUTION**

**Improper use of the stance function**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Make sure that you are standing safely when using the stance function and check the lock of the knee joint before placing your full weight on the prosthesis.
- ▶ Make sure you have been instructed in the correct use of the stance function by the O&P professional and/or therapist. Information on the stance function see page 80.

**⚠ CAUTION**

**Quickly pushing the hip forward with the prosthesis extended (e.g. serve while playing tennis)**

- > Falling due to unexpected activation of the swing phase.
- ▶ Note that the knee joint may flex unexpectedly when the hip is pushed forward quickly while the prosthesis is extended.
- ▶ Therefore, familiarise yourself with swing phase activation in such situations under secure conditions (e.g. while holding on to parallel bars) and with the instruction of trained, qualified personnel.
- ▶ For sports where this movement pattern may occur, use a corresponding pre-configured MyMode. For further information about the MyModes, see the section 'MyModes' (see page 90).

#### 4.8 Notes on the safety modes

**⚠ CAUTION**

**Using the product in safety mode**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 101) have to be observed.
- ▶ Particular caution is necessary when using a bicycle without a freewheel (with a fixed gear).

**⚠ CAUTION**

**Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Do not continue using the defective product.
- ▶ Consult the O&P professional promptly.

**⚠ CAUTION**

**Safety mode cannot be deactivated**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ If safety mode cannot be deactivated by recharging the battery, a permanent error has occurred.
- ▶ Do not continue using the defective product.
- ▶ The product must be inspected by an authorised Ottobock Service Center. The O&P professional is your contact.

**⚠ CAUTION**

**Safety signal occurs (ongoing vibration)**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Observe the warnings/error signals (see page 101).

- ▶ Do not continue using the product after the safety signal has been emitted.
- ▶ The product must be inspected by an authorised Ottobock Service Center. The O&P professional is your contact.

## 4.9 Instructions for use with an osseointegrated implant system

### Only valid for the products 3B5-2/3B5-2=ST

#### **WARNING**

#### **High mechanical loads due to normal or unusual situations, such as falling**

- > Overloading of the bone, which can lead to pain, loosening of the implant, death of bone tissue or bone fracture, among other things.
- > Damage or breakage of the implant system or its components (safety components...).
- ▶ Verify compliance with the fields of application, conditions of use and indications according to the information of the manufacturers, both for the knee joint and for the implant system.
- ▶ Note the instructions of the clinical personnel that indicated the use of the osseointegrated implant system.
- ▶ Note changes in your state of health that result in restrictions or doubt regarding the use of the osseointegrated connection.

## 4.10 Information on the use of a mobile device with the cockpit app

#### **CAUTION**

#### **Improper use of the mobile device**

Falling due to changed damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Make sure you have been instructed in the proper use of the mobile device with the Cockpit app.

#### **CAUTION**

#### **Improper use of the setting parameters in the MyModes**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ask the O&P professional and/or therapist to instruct you regarding the functionality and adjustment options for **all parameters** of the MyModes.

#### **CAUTION**

#### **Independently applied changes or modifications made to the mobile device**

Falling due to altered damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Do not make any independent changes to the hardware of the mobile device on which the app is installed.
- ▶ Do not make any independent changes to the software/firmware of the mobile device that are not included in the update function of the software/firmware.

#### **CAUTION**

#### **Improper mode switching with the device**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter (beeper) and the display on the device.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.

## NOTICE

### **Destruction of the mobile device due to falling or penetration of water**

Malfunction of the mobile device.

- ▶ Note the instructions for use for the mobile device.
- ▶ If it should no longer be possible to switch back from a MyMode to basic mode, the component can only be switched back to basic mode by using a movement pattern (see page 92) or by connecting/disconnecting the inductive battery charger.

## NOTICE

### **Failure to observe the system requirements for the installation of the Cockpit app**

Mobile device malfunction.

- ▶ The Cockpit App should only be installed on mobile devices and versions which comply with the specifications in the respective online stores (e.g. Apple App Store, Google Play Store, ...)

## INFORMATION

The illustrations in these instructions for use are only examples and may deviate from the respective mobile device being used and the version.

## 5 Scope of Delivery and Accessories

### 5.1 Scope of delivery

- 1 pc. 3B5-X3=ST Genium X3 - Bionic Prosthetic System (with threaded connector) or
- 1 pc. 3B5-X3 Genium X3 - Bionic Prosthetic System (with pyramid connector) or
- 1 pc. 3B5-2=ST Genium X3 (with threaded connector) or
- 1 pc. 3B5-2 Genium X3 (with pyramid connector)
- 1 pc. 2R19 AXON tube adapter
- 1 pc. 757L16\* power supply
- 1 pc. 4E60\* inductive charger
- 1 pc. cosmetic case for battery charger and power supply
- 1 pc. 646C107 Bluetooth PIN card
- 1 pc. 647F542 prosthesis passport
- 1 pc. Instructions for use (user)

Cockpit app for download from the website: <https://www.ottobock.com/cockpitapp>

- "Cockpit 4X441-ANDR=V\*" app for Android

**Only when "DUAL" marking is present** (see page 60)

- "Cockpit 4X441-IOS=V\*" app for iOS

### 5.2 Accessories

The following components are not included in the scope of delivery and may be ordered separately:

- Genium X3 Protective Cover: 4X900
- Genium X3 Protective Cover: 4X193-1

**Only when "DUAL" marking is not present** (see page 60)

- 4X350-X3 remote control  
The user interface is available in the following languages: German, English, French, Italian, Spanish, Dutch, Swedish
- 4X194 wrist band for 4X350-X3 remote control

## 6 Charging the battery

The following points must be observed when charging the battery:

- Use the 757L16\* power supply and 4E60\* battery charger to charge the battery.
- With average use, the capacity of the fully charged battery is sufficient for about 5 days.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.
- Note the permissible temperature range for charging the battery (see page 97).
- The distance between the battery charger and the receiver on the product must not exceed 2 mm.

### 6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
- 2) Connect the round, **three-pin** plug of the power supply to the receptacle on the battery charger (see fig. 2) so that the plug locks into place.

**INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.**

- 3) Plug the power supply unit into the outlet (see fig. 3).
    - The green LED on the back of the power supply lights up.
    - The LED ring (status indicator) on the rear of the charger lights up green to indicate the correct connection to the power supply.
- If the green LED on the power supply and the LED ring on the battery charger do not light up, there is an error (see page 101).

### 6.2 Charging the prosthesis battery

#### INFORMATION

When the Protective Cover is installed, the battery charger cable has to point to the upper closure. A correct knee joint charging process is only ensured with this alignment.

- 1) Connect the inductive charger to the receiver of the charging unit on the rear of the product. The charger is held in place by a magnet.
  - The LED ring on the rear of the charger pulsates purple (4-second cycle).
  - If the LED ring lights up in a different colour, this indicates an error (see page 101).
- 2) The charging process starts.
  - Once the product battery is fully charged, all LEDs on the side of the battery charger light up.
- 3) After the charging process is complete, hold the prosthesis still and remove the inductive charger from the receiver.
  - A self-test is performed. The joint is operational only after corresponding feedback (see page 104).

## 6.3 Display of the current charge level

### 6.3.1 Display of battery charge level without additional devices

#### INFORMATION

The charge level cannot be displayed during the charging process, e.g. by turning the prosthesis over. The product is in charging mode.

- 1) Turn the prosthesis by 180° (the sole of the foot must face up).
- 2) Hold still for 2 seconds and wait for beep signals.

Beep signal	Battery charge level
5x short	more than 80%
4x short	60% to 80%
3x short	40% to 60%
2x short	20% to 40%
1x short	less than 20%

#### INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app (see page 84) or if mute mode (silent mode) is activated, there are no beep signals.

### 6.3.2 Display of the current charge level using the Cockpit app

Once the Cockpit app has been started, the current charge level is displayed in the bottom line of the screen:

1.  38% – Charge level of battery for currently connected component

### 6.3.3 Display of the current charge level during the charging process


During the charging process, the current battery charge level is indicated by the number of LEDs lit on the side of the charger.

	Quantity	Battery charge level
	0	0%-10%
	1	10%-30%
	2	30%-50%
	3	50%-70%
	4	70%-90%
	5	> 90%

### 6.3.4 Display of current charge level using remote control (optional accessory)

**Only when "DUAL" marking is not present** (see page 60)

The current charge level is displayed in the status bar of the remote control when it is connected and switched on:

1.  – Charge level of battery for currently connected prosthesis
2.  – Prosthesis is being charged
3. 68% – Charge level of battery for currently connected prosthesis, as a percentage

## 7 Cockpit app

The Cockpit app enables switching from basic mode to the pre-configured MyModes. In addition, information about the product (step counter, battery charge level, etc.) can be retrieved.

The everyday behaviour of the product can be changed to a certain extent using the app (e.g. while becoming accustomed to the product). The O&P professional can use the adjustment software to track these changes at the next appointment.

### Information on the Cockpit app

- The Cockpit app can be downloaded free of charge from the respective online store. For more information, please visit the following website: <https://www.ottobock.com/cockpitapp>. To download the Cockpit app, the QR code on the supplied Bluetooth PIN card can also be read with the mobile device (requirement: QR code reader and camera).
- The language of the Cockpit app user interface can only be changed using the adjustment software.
- The serial number of the component to be connected has to be registered with Ottobock the first time it is connected. If the registration is not accepted, use of the Cockpit app for this component will be limited.
- Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (sole of the foot must point up) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 87).
- Keep the mobile app up to date at all times.
- Please contact the manufacturer if you suspect cybersecurity problems.

### 7.1 System Requirements

The functioning of the Cockpit app is assured on mobile devices that support the following operating systems:

**Only when "DUAL" marking is not present** (see page 60)

- **Android:** version 4.0.3 or higher

**Only when "DUAL" marking is present** (see page 60)

- The 4X441=\* Cockpit app must be installed in version 2.0 or higher
- **iOS (for iPhone, iPad, iPod):** version 10.0 or higher
- **Android:** version 5.0 or higher

### 7.2 Initial connection between cockpit app and prosthesis

**The following points should be observed before the initial connection:**

- Bluetooth of the component must be switched on (see page 87).
- Bluetooth on the mobile device must be switched on.
- The mobile device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- **The mobile device must be connected to the Internet.**
- The serial number and Bluetooth PIN of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The serial number begins with the letters "SN".

## INFORMATION

If the Bluetooth PIN card with the Bluetooth PIN and serial number of the component is lost, contact your O&P professional.

### 7.2.1 Starting the cockpit app for the first time

- 1) Tap the symbol of the Cockpit app (  ).
  - The end user license agreement (EULA) is displayed.
- 2) Accept the end user license agreement (EULA) by tapping the **Accept** button. If the end user license agreement (EULA) is not accepted, the Cockpit app cannot be used.
  - The welcome screen appears.
- 3) Hold the prosthesis with the sole of the foot facing up, or connect and then disconnect the battery charger, in order to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
- 4) Tap the **Add component** button.
  - The Connection Wizard opens and guides you through the process of establishing a connection.
- 5) Follow the subsequent instructions on the screen.
- 6) After the Bluetooth PIN is entered, a connection to the component is established.
  - While the connection is being established, 3 beep signals sound and the  symbol appears.
    - The  symbol is displayed when the connection has been established.
  - Once the connection has been established, the data are read from the component. This process may take up to a minute.
    - Then the main menu appears with the name of the connected component.

## INFORMATION


After the initial connection to the component has been established successfully, the app will connect automatically each time it is started. No further steps are required.

**Only when "DUAL" marking is present** (see page 60)

## INFORMATION


After activating the "visibility" of the component (holding the component with the sole of the foot facing up, or connecting and then disconnecting the battery charger), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case, hold the component with the sole of the foot facing up again, or connect and then disconnect the battery charger.

## 7.3 Control elements for cockpit app

1. ☰ Access the navigation menu (see page 78)
2. Product  
The component name can only be changed with the adjustment software.
3. If connections to more than one component have been saved, you can switch between the saved components by tapping the entry **change**.
4. MyModes configured with the adjustment software.  
Switching the mode by tapping the corresponding symbol and confirming by tapping "OK".
5. Currently selected mode
6. Charge level of the component.  
  - 🔋 Component battery fully charged
  - 🔌 Component battery empty
  - 🔌🔋 Component battery is being charged
 The current charge level is also displayed in %.
7. Display and name of the currently selected mode (e.g. **1. Basic Mode**)
8. 🔇 Mute mode is activated
9. (i+)  
  - (i+) Connection to component has been established
  - (i) Connection to component has been interrupted. The app is attempting to re-establish the connection automatically.
  - (i-) No existing connection to the component.

### 7.3.1 Cockpit app navigation menu

Tap the ☰ symbol in the menus to display the navigation menu. Additional settings for the connected component can be configured in this menu.

#### Product

Name of the connected component

#### MyModes

Return to the main menu to switch MyModes

#### Functions

Call up additional functions of the component (e.g. turn off Bluetooth) (see page 87)

#### Settings

Change settings of the currently selected mode (see page 84)

#### Status

Query status of the connected component (see page 88)

#### Manage components

Add or delete components (see page 79)

#### Imprint/Info

Display information/legal notices for the cockpit app

## 7.4 Managing components

Connections with up to four different components can be stored in the app. However, a component can only be connected to one mobile device at a time.

### INFORMATION

Before establishing the connection, observe the points in the section "Initial connection between Cockpit app and component" (see page 76).

### 7.4.1 Adding component

- 1) Tap the ☰ symbol in the main menu.  
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Hold the prosthesis with the sole of the foot facing up, or connect and then disconnect the battery charger, in order to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
- 4) Tap the "+" button.  
→ The Connection Wizard opens and guides you through the process of establishing a connection.
- 5) Follow the subsequent instructions on the screen.
- 6) After the Bluetooth PIN is entered, a connection to the component is established.  
→ While the connection is being established, 3 beep signals sound and the (📶) symbol appears.  
The (📶) symbol is displayed when the connection has been established.  
→ Once the connection has been established, the data are read from the component. This process may take up to a minute.  
Then the main menu appears with the name of the connected component.

### INFORMATION

If establishing a connection to a component is not possible, perform the following steps:

- ▶ Delete the component from the Cockpit app if applicable (see the section "Deleting a component")
- ▶ Add the component again in the Cockpit app (see the section "Adding a component")

**Only when "DUAL" marking is present** (see page 60)

### INFORMATION

After activating the "visibility" of the component (holding the component with the sole of the foot facing up, or connecting and then disconnecting the battery charger), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case, hold the component with the sole of the foot facing up again, or connect and then disconnect the battery charger.


### 7.4.2 Deleting a component

- 1) Tap the ☰ symbol in the main menu.  
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Tap the "**Edit**" button.
- 4) Tap the 🗑️ symbol under the component you want to delete.  
→ The component is deleted.

### 7.4.3 Connecting component with multiple mobile devices

The connection for a component can be stored on more than one mobile device. However, only one mobile device can be connected to the component at one time.

If there is an existing connection between the component and a different mobile device, the following information appears while the connection is being established with the current mobile device:


▶ Tap the "OK" button.

→ The connection to the last connected mobile device is broken off and established with the current mobile device.

## 8 Use

### 8.1 Movement patterns in basic mode (mode 1)

#### 8.1.1 Standing

Knee control through high hydraulic resistance and static alignment.

A stance function can be enabled by the O&P professional. Please see the following section for further information on the stance function.

##### 8.1.1.1 Stance function

###### INFORMATION

To use this function, it needs to be enabled in the adjustment software. It also has to be activated using the Cockpit app (see page 85).

The stance function (standing mode) is a functional supplement to the basic mode (mode 1). This function makes it easier, for example, to stand on an inclined surface for a longer time. The joint is fixed in the flexion direction.

The stance function must be enabled by the O&P professional. In addition, the type of joint locking (conscious/intuitive) has to be established by the O&P professional. Changing the locking type using the Cockpit app is not possible.

###### Intuitive locking of the joint

The intuitive stance function recognises any situation that puts strain on the prosthesis in the flexion direction but where flexion is not permitted. Examples of this include standing on uneven or sloping surfaces. The knee joint is always locked in the flexion direction when the prosthetic leg is not fully extended, under some amount of load and at rest. Upon forward or backward rollover or extension, the level of resistance is immediately reduced to stance phase resistance again.

The knee joint is not locked when the above conditions are met and a sitting position is assumed (for example while driving).

###### Deliberate locking of the joint

- 1) Assume the desired knee angle.
- 2) Do not remove the entire load from the prosthesis.


3) Do not change the knee angle for a brief period (1/8 second). This time period prevents unintentional activation of the stance function while walking.

→ The blocked joint can now be loaded in the flexion direction.

### Deliberate unlocking of the joint

► By deliberately extending or unloading the knee joint, it is unlocked again.

### 8.1.2 Walking


Initial attempts at walking with the prosthesis always require the instruction of trained, qualified personnel.

The hydraulics stabilise the knee joint in the stance phase and release the knee joint in the swing phase so that the leg can swing forward freely.

Switching to the swing phase requires that the prosthesis roll over to the front out of the stride position.

### 8.1.3 Running short distances ("walk-to-run" function)

For covering short distances quickly, the knee joint detects a transition from walking to running in basic mode and automatically changes the following settings:

- The swing phase angle is increased
- Preflexion of 4° at heel strike (PreFlex) is reduced to 0°

The requirements to automatically switch to the running motion are fast forward movement of the prosthetic leg and high dynamic load on the knee joint. When stopping from the running motion, the changed settings are set back to the standard values.

#### INFORMATION

For covering longer distances, the prosthetist can configure a "Running" MyMode (see page 90).

### 8.1.4 Sitting down

The resistance in the prosthetic knee joint while sitting down ensures even bending into the sitting position.

The O&P professional can use the adjustment software to configure whether the sitting process is to be supported or not.


- 1) Place both feet side by side at the same level.
- 2) While sitting down, weight should be distributed evenly between both legs and the arm supports used where applicable.
- 3) Move the buttocks in the direction of the back support and lean the upper body forward.

**INFORMATION: Resistance while sitting down can be changed with the Cockpit app via the parameter "Resistance" (see page 85).**

### 8.1.5 Sitting

#### INFORMATION

While sitting, the knee joint also switches to energy saving mode. This energy saving mode is activated regardless of whether the sitting function is activated or not.

If the patient is in a sitting position for more than two seconds (i.e. the thigh is close to horizontal and there is no load on the leg), the knee joint switches the resistance to a minimum in the extension direction. A sitting function can be enabled by the O&P professional. For more information about the sitting function, see the following section.

### 8.1.5.1 Sitting function

#### INFORMATION

To use this function, it needs to be enabled in the adjustment software. It also has to be activated using the Cockpit app (see page 85).

In the sitting position, the resistance in the flexion direction is reduced in addition to the reduction of resistance in the extension direction. This makes it possible to swing the prosthetic leg freely.

### 8.1.6 Standing up

Flexion resistance is increased steadily while standing up.


- 1) Place the feet at the same level.
- 2) Lean the upper body forward.
- 3) Put the hands on armrests, if available.
- 4) Stand up with support from the hands while keeping weight evenly distributed on feet.

### 8.1.7 Walking up stairs step-over-step

#### INFORMATION

To use this function, it needs to be enabled in the adjustment software. It also has to be activated using the Cockpit app (see page 85).

Although the knee joint is passive, which means it cannot actively initiate movements, negotiating stairs step-over-step is possible. This function must be practised and executed consciously.

- 1) Lift the extended prosthesis off the floor.
- 2) Immediately after lifting the extended leg off the floor, extend the hip briefly and then abruptly flex it. This requires a sufficiently secure hold in the socket and a certain level of residual limb strength.  
→ This whip motion flexes the knee, because the knee joint automatically recognises the movement and sets the flexion resistance to minimum.


**INFORMATION: Take note of people behind you before executing the whip motion.**

- 3) When sufficient knee flexion has been achieved, the knee joint increases extension resistance so that there is enough time to position the foot on the next step before the knee joint is extended again.
- 4) Set the foot onto the next step.  
The support area for the foot on the step must be large enough that the heel does not extend back too far over the edge. With too little support area, the lower leg would extend too early and position the leg too far backwards. In this phase, the knee joint has already set the flexion resistance to maximum (blocked). The knee joint can no longer be flexed but only extended. This ensures that the leg does not buckle if the hip strength is not sufficient for the extending motion.
- 5) Support yourself with your hand on the contralateral side. A smooth wall will also work. This lateral support is intended to prevent the residual limb from twisting in the socket. Twisting can lead to unpleasant surface tension between the skin and the socket. Lateral support also improves balance.
- 6) Bring the knee into extended position. When the knee joint is fully extended, the initial position has been reached.
- 7) You can climb the next step or continue walking normally.

### 8.1.8 Overcoming obstacles

**INFORMATION**


To use this function, it needs to be enabled in the adjustment software. It also has to be activated using the Cockpit app (see page 85).

The stair function can also be used to cross obstacles:

- 1) Lift the extended prosthesis off the floor.
- 2) Briefly extend the hip.
- 3) Quickly flex the hip. This causes the knee to flex.
- 4) With the knee flexed, step over the obstacle.  
With sufficient knee flexion, the extension resistance is increased to allow enough time for crossing the obstacle.

### 8.1.9 Walking down stairs


This function must be practised and executed consciously. Only when the sole is properly positioned can the knee joint react correctly and permit controlled flexion.

- 1) Hold the handrail with one hand.
- 2) Position the leg with the prosthesis on the step so that the foot projects halfway over the edge of the step.  
→ This is the only way to ensure a secure rollover.
- 3) Roll the foot over the edge of the step.  
→ This flexes the prosthesis slowly and evenly in the knee joint.
- 4) Place the foot of the other leg onto the next step.
- 5) Place the foot of the prosthetic leg on the next step after that.

**INFORMATION: The flexion speed of the knee joint can be changed using the Cockpit app via the parameter "Resistance" (see page 85).**

### 8.1.10 Walking down a ramp

Under increased flexion resistance, permit controlled flexion of the knee joint which lowers the body's centre of gravity.

**INFORMATION: The flexion resistance at which the knee bends can be changed using the Cockpit app via the parameter "Resistance" (see page 85).**

## 8.2 Changing prosthesis settings


Once an active connection to a component has been established, the settings **of the respective active mode** can be changed using the Cockpit app.

### INFORMATION


Bluetooth on the prosthesis must be switched on to change the prosthesis settings.

If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. The connection must be established during this period.

### Information for changing the prosthesis settings

- Before changing settings, always check the main menu of the Cockpit app to make sure the correct component has been selected. Otherwise parameters could be changed for the wrong component.
- It is not possible to change prosthesis settings nor to switch to a different mode while the prosthesis battery is being charged. Only the status of the prosthesis can be called up. Instead of the  symbol, the  symbol appears in the bottom row of the screen in the cockpit app.
- The O&P professional's setting is in the middle of the scale. After making adjustments, this setting can be restored by tapping the "**Standard**" button in the Cockpit app.
- Prosthesis settings should be optimised using the adjustment software. The Cockpit app is not intended for use by the O&P professional to set up the prosthesis. The everyday behaviour of the prosthesis can be changed to a certain extent using the app (e.g. while becoming accustomed to the prosthesis). The O&P professional can use the adjustment software to track these changes at the next appointment.
- If the settings of a MyMode are to be modified, one must first switch to this MyMode.

### 8.2.1 Changing the prosthesis setting using the cockpit app

- 1) Once the component is connected and in the desired mode, tap the ☰ icon in the main menu.  
→ The navigation menu opens.
- 2) Tap the **"Settings"** menu option.  
→ A list appears with the parameters for the currently selected mode.
- 3) Change the setting of the desired parameter by tapping the "<" , ">" icons.

**INFORMATION: The O&P professional's setting is marked and, after the setting has been changed, can be restored by tapping the "Standard" button.**

### 8.2.2 Overview of adjustment parameters in basic mode

#### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

The parameters in basic mode describe the dynamic behaviour of the prosthesis in a normal gait cycle. These parameters act as basic settings for automatically adjusting the damping behaviour to the current motion situation (e.g. ramps, slow walking speed, etc.).

The stance function, sitting function and/or stairs and obstacles function can also be activated/deactivated. Further information on the stance function (see page 80), sitting function (see page 81), stairs and obstacles function (see page 82).

**The following parameters can be modified:**

Parameter	Adjustment software range	Setting range, app	Meaning
Resistance	120 – 180	+/- 10	Resistance against flexion motion e.g. when walking down stairs or when sitting down
Angle	55° – 70°	+/- 3°	Maximum flexion angle in the swing phase
Stance function	Deactivated Activated	0 - deactivated 1 - activated	Activation/deactivation of the stance function. For switching with the Cockpit app, this function has to be activated in the adjustment software. Further information (see page 80).

Parameter	Adjustment software range	Setting range, app	Meaning
Sitting function	Deactivated Activated	0 - deactivated 1 - activated	Activation/deactivation of the sitting function. For switching with the Cockpit app, this function has to be activated in the adjustment software. Further information (see page 82).
Stair Function	Deactivated Activated	0 - deactivated 1 - activated	Activating/deactivating the stair and obstacle function. For switching with the Cockpit app, this function has to be activated in the adjustment software. Further information (see page 82).
Pitch	1000 Hz – 4000 Hz	1000 Hz – 4000 Hz	Pitch of beep signal for confirmation tones
Volume	0 – 4	0 – 4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals in case of errors are still generated.

### 8.2.3 Overview of adjustment parameters in MyModes

#### CAUTION

##### Improper use of the setting parameters in the MyModes

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- Ask the O&P professional and/or therapist to instruct you regarding the functionality and adjustment options for **all parameters** of the MyModes.

#### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

The parameters in the MyModes describe the static behaviour of the prosthesis for a specific motion pattern such as cross-country skiing. Damping behaviour is not automatically controlled and adjusted in MyModes.

Parameter	Adjustment software range	Setting range, app	Meaning
Basic flex.	0–200	+/- 20	Level of flexion resistance at the start of flexing the knee joint
Gain	0–100	+/- 10	Increase in flexion resistance (starting with the parameter " <b>Basic flex.</b> ") when flexing the knee joint. The knee joint locks at a certain flexion angle, which depends on the setting of the parameters " <b>Basic flex.</b> " and " <b>Gain</b> ".
Basic ext.	0–60	+/- 20	Level of extension resistance

Parameter	Adjustment software range	Setting range, app	Meaning
Locking angle	0–90	+/- 10	Angle up to which the knee joint can be extended. <b>Information:</b> If this parameter is >0, the knee joint is locked in a flexed position in the extension direction. To unlock it, unload the prosthesis and tilt it back for at least 1.5 seconds. This makes extension of the joint possible independently of the settings for the parameters " <b>Basic ext.</b> " and " <b>Locking angle</b> ". This may be necessary to switch to basic mode using a movement pattern.
Pitch	1000 Hz–4000 Hz	1000 Hz–4000 Hz	Pitch of beep signal for confirmation tones
Volume	0–4	0–4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals in case of errors are still generated.

### 8.3 Turning Bluetooth on the prosthesis on/off

#### INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (function only available in basic mode) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 87).

#### INFORMATION

Basic mode (mode 1) has to be active to turn off Bluetooth. If a MyMode is activated, one has to switch to basic mode to turn off Bluetooth.

#### 8.3.1 Switching Bluetooth off/on using the cockpit app

##### Switching off Bluetooth

- 1) When the component is connected, tap the ☰ symbol in the main menu.  
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".
- 3) Tap the entry "**Deactivate Bluetooth**".
- 4) Follow the on-screen instructions.

##### Switching on Bluetooth

- 1) Turn the component over or connect/disconnect the battery charger.  
→ Bluetooth is switched on for approx. 2 minutes. During this time, the app must be started in order to establish a connection to the component.

- 2) Follow the on-screen instructions.  
→ If Bluetooth is switched on, the  symbol appears on the screen.

## 8.4 Querying the prosthesis status

### 8.4.1 Query status through cockpit app

- 1) When the component is connected, tap the  symbol in the main menu.
- 2) In the navigation menu, tap the entry "**Status**".

### 8.4.2 Status display in the cockpit app

Menu option	Description	Possible actions
Trip	Daily step counter (steps with the prosthesis side)	Reset the counter by tapping the " <b>Reset</b> " button.
Step	Total step counter (steps with the prosthesis side)	For informational purposes only
Service	Display of the next maintenance date	For informational purposes only
Batt.	Current prosthesis charge level, as a percentage	For informational purposes only
Stb/Act: 58/29	Estimated remaining operating time of prosthesis in hours. Rest mode (Stb) e.g. 58 hours, active use (Act) e.g. 29 hours	For informational purposes only

## 8.5 Mute mode (silent mode)

Activating mute mode (silent mode) turns off the audible feedback signals and the vibration signals. However, warnings in case of component errors are still generated (see page 101). Mute mode can be activated/deactivated using the Cockpit app.

### INFORMATION

Connecting the battery charger automatically deactivates mute mode again.

### 8.5.1 Turning mute mode on/off using the Cockpit app

- 1) When the component is connected, tap the  symbol in the main menu.  
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".
- 3) Tap the entry "**Mute mode**".
- 4) Follow the on-screen instructions.

## 8.6 Deep sleep mode

### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

### INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 84).

The Cockpit app can be used to place the knee joint into a deep sleep mode, in which power consumption is minimised. The knee joint offers no functionality in this mode. The safety mode damping values are activated.

It can be awakened from deep sleep mode with the Cockpit app or by connecting the battery charger. Waking from deep sleep mode using the Cockpit app can take up to 30 seconds. After ending deep sleep mode, the knee joint is in basic mode again.

### 8.6.1 Turning deep sleep mode on/off using the Cockpit app

#### Activating deep sleep mode

- 1) When the component is connected, tap the ☰ symbol in the main menu.  
→ The navigation menu opens.
  - 2) In the navigation menu, tap the entry "**Functions**".
  - 3) Tap the entry "**Activate deep sleep mode**".
  - 4) Follow the on-screen instructions.
- The activation of deep sleep mode is indicated by a short beep signal and a short vibration signal, assuming that mute mode (silent mode) is not active.

#### Deactivating deep sleep mode

- 1) When deep sleep mode is active for the currently connected prosthesis, the **Exit deep sleep mode** button automatically appears when the Cockpit app is started.
- 2) Tapping this button establishes a connection to the prosthesis and deactivates deep sleep mode.

**INFORMATION: Establishing a connection in deep sleep mode can take up to 30 seconds.**

If a prosthesis is in deep sleep mode but not connected to the Cockpit app, a connection to the prosthesis has to be established (see page 79).

## 8.7 OPG function (Optimised Physiological Gait)

### INFORMATION

The O&P professional can turn the "PreFlex" function on or off using the adjustment software. All other parameters of the OPG function are always active and cannot be influenced.

The OPG function minimises the prosthesis wearer's prosthetic deviations from a harmonious gait pattern and promotes more biomechanically correct walking. This function enables the following features:

#### PreFlex

PreFlex ensures the knee is in 4° of flexion at the end of swing phase in preparation for loading response. This makes initiating stance phase flexion easier and forward movement is less restricted.

#### Adaptive yielding control

The knee joint has auto-adaptive stance and swing extension resistance. The stance flexion resistance experienced by the user is dependent on the slope or incline when walking downhill. When walking on a ramp, adaptive yielding control manages flexion depending on the angle of the ramp. The knee joint flexes slowly if the ramp is flat, and flexes quickly if the ramp is steep.

#### Dynamic stability control (DSC)

DSC ensures the knee will not release stance resistance during biomechanically unstable static and dynamic conditions. Constantly checking multiple parameters, DSC ensures the optimally timed decision for the knee to safely switch from stance to swing. Because DSC is always monitoring knee function, multi-directional movement and walking backward are also possible without risk of stance resistance releasing.


#### Adaptive swing phase control

Instantaneous adaptation to varied walking cadences and to changes of the pendular mass (e.g. varying footwear) ensures the knee always achieves the swing flexion target angle within (+/-) one degree. The swing phase extension and flexion resistance experienced by the user are auto-adaptive.

The flexed and partially loaded knee will also disable the stance phase on slopes and ramps to allow for greater knee flexion and more ground clearance in the swing phase.

## 9 MyModes

The O&P professional can activate and configure up to five MyModes in addition to basic mode using adjustment software. These can be selected by using the Cockpit app. Only the first three MyModes can be selected using movement patterns. Switching by using movement patterns has to be activated in the adjustment software by the O&P professional.


These modes are intended for specific motion patterns or postures (e.g. inline skating, running (jogging) ...). Settings can be adjusted using the Cockpit app (see page 86).

**Only when "DUAL" marking is not present** (see page 60)

For information on switching MyModes with the remote control (optional accessory), see the remote control instructions for use.

### 9.1 Running function as configured MyMode

For running over longer periods of time, the O&P professional can configure a "Running" MyMode which can be activated using the cockpit app or a movement pattern.

In this mode, every step will be performed as a running step with larger swing phase angle and no preflexion at heel strike (PreFlex) (see page 89).

#### INFORMATION

The running function will work with specialised running feet such as the 1E95 Challenger as well as with prosthetic feet with axial compression such as the 1C61 Triton Vertical Shock. For details on assembly and alignment please refer to the instructions for use of the foot. Feet without axial compression are generally not suited for running.

### 9.2 Switching MyModes with the cockpit app


#### INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (function only available in basic mode) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 87).

## INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app (see page 84) or if mute mode (silent mode) is activated, there are no beep signals.

Once a connection to a prosthesis has been established, the cockpit app can be used to switch between the MyModes.

- 1) Tap the symbol of the MyMode (1) you want in the main menu of the app.  
→ A security question for changing the MyMode appears.
- 2) If you want to change the mode, tap the "OK" button.  
→ A beep signal sounds to confirm the switch.
- 3) After switching, a symbol (2) is displayed to identify the active mode.  
→ The current mode is also indicated by the name on the lower edge of the screen (3).

## 9.3 Switching MyModes using motion patterns

### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

### INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app (see page 84) or if mute mode (silent mode) is activated, there are no beep signals.

#### Information on switching


- Switching and the number of movement patterns must be activated by the O&P professional in the adjustment software.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

#### Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Switching using movement patterns has to be enabled by the O&P professional.
- Position the prosthetic leg back slightly and bounce on the forefoot with the leg extended while maintaining constant contact with the floor.
- Weight must be placed on the forefoot while bouncing.
- The load may not be taken off the foot fully when relieving the load.

## Switching process

- 1) Position the prosthetic leg back slightly.
- 2) While maintaining constant contact with the floor, bounce on the forefoot with the leg extended a number of times in one second depending on the desired MyMode (MyMode 1 = 3 times, MyMode 2 = 4 times).
- 3) Keep the prosthetic leg still in this position with no load.

→ A beep and vibration signal will sound to confirm that the movement pattern has been recognised.

**INFORMATION: If this beep and vibration signal does not sound, the requirements were not observed when bouncing the foot or mute mode (silent mode) is activated. For more information about mute mode, see the section "Mute mode" (see page 88).**

- 4) After the beep and vibration signal sounds, keep the prosthetic leg extended and still for 1 second.

→ A confirmation signal will sound to indicate that the prosthesis has successfully switched to the corresponding MyMode (2 times = MyMode 1, 3 times = MyMode 2).

**INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly kept still or mute mode (silent mode) is activated. Repeat the process to correctly switch to the required mode. For more information about mute mode, see the section "Mute mode" (see page 88).**

## 9.4 Switching from a MyMode back to basic mode

### Information on switching


- Regardless of the configuration of additional MyModes in the adjustment software, it is always possible to switch back to basic mode (mode 1) with a motion pattern.
- It is always possible to switch back to basic mode (mode 1) by connecting/disconnecting the battery charger.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

### Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Position the prosthetic leg back slightly and bounce on the forefoot with the leg extended while maintaining constant contact with the floor.
- Weight must be placed on the forefoot while bouncing.
- The load may not be taken off the foot fully when relieving the load.

## Switching process

- 1) Position the prosthetic leg back slightly.

- 2) While maintaining constant contact with the floor and with the leg extended, bounce on the forefoot at 3 times or more.
- 3) Keep the prosthetic leg still in this position with no load.  
→ A beep and vibration signal will sound to confirm that the movement pattern has been recognised.

**INFORMATION: If this beep and vibration signal does not sound, the requirements were not observed when bouncing the foot or mute mode (silent mode) is activated. For more information about mute mode, see the section "Mute mode" (see page 88).**

- 4) After the beep and vibration signal sounds, keep the prosthetic leg extended and still for 1 second.  
→ A confirmation signal will sound to indicate that the prosthesis has successfully switched over to basic mode.

**INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly kept still or mute mode (silent mode) is activated. Repeat the process to correctly switch to the required mode. For more information about mute mode, see the section "Mute mode" (see page 88).**

## 10 Additional operating states (modes)

### 10.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level is 5% (see page 101). During this time, damping settings are set to their safety mode values. This may be low or high depending on the setting in the adjustment software. The prosthesis is then switched off. You can switch back to basic mode (mode 1) from empty battery mode by charging the product.

### 10.2 Mode for charging the prosthesis

The product is non-functional during charging.

The product is set to the flexion resistance of safety mode. This may be low or high depending on the setting configured by the O&P professional.

### 10.3 Safety mode

The product automatically switches to safety mode if a critical fault occurs (e.g. failure of a sensor signal). Safety mode remains in effect until the error has been rectified.

The switch to safety mode is indicated by beeps and vibration signals immediately prior to switching (see page 101).

Safety mode can be disabled by connecting and disconnecting the battery charger. If the product switches into safety mode again, this means a permanent error exists. The product must be inspected by an authorised Ottobock Service Centre.

Depending on the type and severity of the error, different remaining functionality is offered in safety mode. This makes limited walking possible for the user depending on the type of error.

#### **The following remaining functionality is available:**

- **Minor error:** Continuous stance phase flexion resistance is set, with the ability to initiate the swing phase.
- **Moderate error:** Continuous stance phase flexion resistance is set, with the ability to initiate the swing phase. The swing phase control and stance phase extension resistance may be available or not depending on the type of error.
- Safety mode flexion resistance is configured. This may be low or high depending on the setting configured by the O&P professional.

#### **The following functions are deactivated in safety mode:**

- OPG function
- Stairs and obstacles function
- Stance function

- Sitting function

## 10.4 Overheating mode

### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

When the hydraulic unit overheats due to uninterrupted, increased activity (e.g. extended walking downhill), damping is increased along with the rising temperature in order to counteract the overheating. When the hydraulic unit cools down, the product switches back to the damping settings that existed before the overheating mode.

Overheating mode is not activated in the MyModes.

Overheating mode is indicated by a long vibration every 5 seconds.

### The following functions are deactivated in overheating mode:

- Sitting function
- Display of the battery charge level without additional equipment
- Switching to a MyMode

## 11 Storage and bleeding

Air may accumulate in the hydraulic unit if the product is stored for longer periods and not in an upright position. This is noticeable through sounds and irregular damping behaviour.

The automatic bleeding mechanism ensures that all functions of the product are again intact after approximately 10 - 20 steps.

### Storage

- Before storing the knee joint, the knee head has to be extended. The knee head must not be flexed!
- Avoid extended disuse of the product (use the product regularly).

## 12 Cleaning

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

### INFORMATION

Please note that the weight of dirt adhering to the prosthesis can affect the gait pattern.

## 13 Maintenance

### INFORMATION

This component was tested for three million load cycles in accordance with ISO 10328.

Depending on the activity level, this corresponds to a service life of three to five years.

The duration of use can be individually extended depending on the intensity of use by performing regular service inspections.

Regular maintenance (service inspections) at 12-monthly intervals is mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty, to maintain basic safety and the essential performance characteristics, and to ensure safety in regards to EMC.

When maintenance is due, this is indicated by feedback after disconnecting the battery charger (see the section "Operating states/error signals", see page 100). The manufacturer grants a grace period of no more than two months before, or three months after, the due date.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components always have to be submitted to the O&P professional for maintenance and repairs:

The prosthesis, battery charger and power supply unit.

**INFORMATION**

If a remote control was included in the scope of delivery as an optional product accessory, it must also be sent along with the product for service inspections.

**14 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

**14.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

**14.2 Trademarks**

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

**14.3 CE conformity**

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

**Only valid for the products 3B5-X3/3B5-X3=ST**

This product meets the requirements of the European Directive 1999/5/EC for radio equipment and telecommunications terminal equipment. The conformity assessment was drawn up by the manufacturer in accordance with Annex IV of the directive.

**Only valid for the products 3B5-2/3B5-2=ST**

This product meets the requirements of the 2014/53/EU directive.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

**14.4 Local Legal Information**

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and

2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Caution: Exposure to Radio Frequency Radiation.**

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP  
3820 West Great Lakes Drive  
Salt Lake City, Utah 84120-7205 USA  
Phone + 1-801-956-2400  
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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**Caution: Exposure to Radio Frequency Radiation.**

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.  
5470 Harvester Road  
L7L 5N5 Burlington, Ontario  
Canada  
Phone + 1-800-665-3327

---

**Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.**

## 15 Technical data

<b>Environmental conditions</b>	
Transportation in original packaging	-25°C/-13°F to +70°C/+158°F
Transportation without packaging	-25°C/-13°F to +70°C/+158°F Max. 93% relative humidity, non-condensing
Storage (≤3 months)	-20°C/-4°F to +40°C/+104°F Max. 93% relative humidity, non-condensing
Long-term storage (>3 months)	-20°C/-4°F to +20°C/+68°F Max. 93% relative humidity, non-condensing
Operation	-10°C/+14°F to +60°C/+140°F Max. 93% relative humidity, non-condensing
Charging the battery	+10°C/+50°F to +45°C/+113°F

<b>Product</b>	
Reference number	3B5-X3/3B5-X3=ST/3B5-2/3B1-2=ST
Mobility grade according to MOBIS	3 and 4
Maximum body weight	125 kg
Protection rating	IP68 Maximum water depth: 3 m Maximum time: 1 hour
Water resistance	Waterproof, corrosion-resistant
Range of Bluetooth connection to mobile device	max. 10 m
Weight of the prosthesis without tube adapter and Protective Cover	approx. 1,710 g

<b>Tube adapter</b>	
Reference number	2R19
Weight	190 g -300 g
Material	Aluminium
Max. body weight	125 kg (275 lbs)
Protection rating	IP68 Maximum water depth: 3 m Maximum time: 1 hour

<b>Prosthesis battery</b>	
Battery type	Li-Ion
Charging cycles (charging and discharging cycles) after which at least 80% of the original battery capacity remains available	500
Charge level after 1 hour charging time	30 %
Charge level after 2 hours charging time	50 %
Charge level after 4 hours charging time	80 %
Charge level after 8 hours charging time	Fully charged
Product behaviour during the charging process	The product is non-functional
Operating time of the prosthesis with new, fully charged battery at room temperature	Approx. 5 days with average use

<b>Power supply</b>	
Reference number	757L16*
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10 % to 93 % relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F max. 90 % relative humidity, non-condensing
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ===

<b>Battery charger</b>	
Reference number	4E60*
Storage and transport in original packaging	-25 °C to 70 °C/-13 °F to 158 °F
Storage and transport without packaging	-25 °C to 70 °C/-13 °F to 158 °F Max. 93% relative humidity, non-condensing
Operation	5 °C to 40 °C/41 °F to 104 °F Max. 93% relative humidity, non-condensing
Protection rating	IP40
Input voltage	12 V ===
Wireless technology	Proprietary protocol
Frequency range	270 kHz to 450 kHz
Modulation	ASK, load modulation
Maximum output power (EIRP)	-12.7 dBμA/m @ 10 m

**Only when "DUAL" marking is not present** (see page 60)

<b>Cockpit app</b>	
Article number	Cockpit 4X441-Andr=V*
Supported operating system	Android 4.0.3 or higher
Website for download	<a href="https://www.ottobock.com/cockpitapp">https://www.ottobock.com/cockpitapp</a>

<b>Data transfer</b>	
Wireless technology	Bluetooth 2.0
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK
Data rate (over the air)	704 kbps
Maximum output power (EIRP)	-1.78 dBm (0.66 mW)

**Only when "DUAL" marking is present** (see page 60)

<b>Cockpit app</b>	
Reference number	4X441-IO=* / 4X441-Andr=V* Cockpit
Supported operating system	iOS 10.0/Android 5.0 or higher
Website for download	<a href="https://www.ottobock.com/cockpitapp">https://www.ottobock.com/cockpitapp</a>

<b>Data transfer</b>	
Wireless technology	Bluetooth Smart Ready
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK, $\pi/4$ DQPSK, 8DPSK
Data rate (over the air)	2178 kbps (asymmetrical)
Maximum output power (EIRP):	+8.5 dBm

## 16 Appendices

### 16.1 Symbols Used

Manufacturer

Type BF applied part

Compliance with the requirements according to "FCC Part 15" (USA)

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Non-ionising radiation

**IP40**

Protection against penetration of solid foreign objects with a diameter greater than 1 mm, no protection against water

**IP68**

Dustproof, protection against permanent submersion.  
Maximum depth: 3 m  
Maximum time: 1 hour

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

**DUAL**

The product's Bluetooth wireless module can establish a connection to mobile devices with the following operating systems: iOS (iPhone, iPad, iPod...) and Android

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)

YYYY – year of manufacture

WW – week of manufacture

NNN – sequential number

Lot number (PPPP YYYY WW)

PPPP – plant

YYYY – year of manufacture

WW – week of manufacture

Caution, hot surface

Protect from moisture

## 16.2 Operating states/error signals

The prosthesis indicates operating states and error messages through beeps and vibration signals.

### 16.2.1 Signals for operating states

#### Battery charger connected/disconnected

Beep signal	Vibration signal	Event
–	3x long	Charging mode started (3 sec. after connecting the battery charger)
1x short	1x short	Self-test completed successfully, product is operational

#### Mode switching

##### INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

##### INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 84).

Beep signal	Vibration signal	Additional action performed	Event
1x short	1x short	Mode switching using the Cockpit app	Mode switching is performed using the Cockpit app.
1x long	1x long	Bouncing on the forefoot followed by unloading the prosthetic leg	Bouncing pattern recognised.

Beep signal	Vibration signal	Additional action performed	Event
1x short	1x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to basic mode (mode 1) carried out.
2x short	2x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to MyMode 1 (mode 2) carried out.
3x short	3x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to MyMode 2 (mode 3) carried out.



## 16.2.2 Warnings/error signals

### Error during use



Beep signal	Vibration signal	Event	Required action
-	1x long at intervals of approx. 5 seconds (if mute mode (silent mode) is activated, this signal is not generated)	Overheated hydraulics	Reduce activity.
-	3x long	Battery charge level under 25%	Charge battery soon. Remaining operating time approx. 24 hours
-	5x long	Battery charge level under 10%	Charge battery soon. Remaining operating time approx. 6 hours
5x long	5x long repeated every 60 seconds	<b>Error of moderate severity (see page 93)</b> e.g. a sensor is not operational	Walking possible with restrictions. Note the change in flexion resistance. The product must be inspected by an O&P professional immediately.
10x long	10x long	Charge level 5% After the beep and vibration signals, the product switches to empty battery mode and then switches off.	Charge the battery.
30x long	1x long, 1x short repeated every 3 seconds	<b>Severe error/indication of safety mode activation (see page 93)</b> e.g. one or more sensors are not operational	Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohibited. The product must be inspected by an O&P professional immediately.

Beep signal	Vibration signal	Event	Required action
-	Continuous	<b>Total failure</b> Electronic control no longer possible. Safety mode active or undetermined valve state. Undetermined product behaviour.	Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohibited. The product must be inspected by an O&P professional immediately.

### Error while charging the product

LED on power supply	Status LED on battery charger	Error	Resolution
○		Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.
		Non-functional socket	Check socket with another electric appliance.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●		No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.

	Status LED	Charging status indicator (5 LEDs)	Error	Resolution
	The LED ring is lit in weak purple	No LED is lit	Distance between battery charger and receiver of the charging unit on the prosthesis too great. If the distance is more than 2 mm, the prosthesis cannot be charged.	Reduce distance between battery charger and receiver of the charging unit.

	Status LED	Charging status indicator (5 LEDs)	Error	Resolution
	The LED ring is lit up yellow	LED 2 and 4 are lit up	Charger excess temperature	Check whether the specified ambient conditions for charging the battery are met (see page 97).
		LED 1, 3 and 5 are lit up	Excessively high or low temperature of the prosthesis	
		LED 3 is lit up	The prosthesis is not being charged Distance between battery charger and receiver of the charging unit too great.	Connection may be improved by reducing the distance between the charger and receiver of the charging unit.
	The LED ring is lit up green		The battery charger is functional but not connected to the receiver, or the distance between the battery charger and receiver of the charging unit is too great.	Connect the battery charger or reduce the distance between the charger and the receiver of the charging unit on the prosthesis.
	The LED ring is flashing red		The prosthesis is not being charged Defective battery charger.	Reset the error by disconnecting and connecting the power supply. If the error persists, the battery charger and power supply must be inspected by an authorised Ottobock Service Center.

### 16.2.3 Error messages while establishing a connection with the cockpit app

Error message	Cause	Correction
<b>Component was connected to another device. Establish connection?</b>	The component was connected to another device	To disconnect the original connection, tap the <b>"OK"</b> button. If the original connection is not to be disconnected, tap the <b>"Cancel"</b> button.
<b>Mode change failed</b>	An attempt was made to switch to a different MyMode while the component was in motion (e. g. while walking)	For safety reasons, switching MyModes is only permitted when components are at rest, e. g. while standing or sitting.

Error message	Cause	Correction
	A current connection to the component was interrupted	Check the following points: <ul style="list-style-type: none"> <li>• Distance from the component to the device</li> <li>• Charge level of the component's battery</li> <li>• Bluetooth of the component switched on? (Switching Bluetooth of the component on/off)</li> <li>• Hold the component with the sole of the foot facing up to make the component "visible" for 2 minutes.</li> <li>• If multiple components were stored, was the correct component selected?</li> </ul>

### 16.2.4 Status signals

#### Battery charger connected

LED on power supply	Status LED on battery charger	Event
		Power supply and battery charger operational. Battery charger not yet connected to receiver.
		Battery charger connected to receiver, good connection. This display turns off automatically after approximately one minute to avoid bothersome light at night. The charging process is not affected.

#### Battery charger disconnected

Beep signal	Vibration signal	Event	Resolution
1x short	1x short	Self-test completed successfully. Product is operational.	

Beep signal	Vibration signal	Event	Resolution
3x short	3x short	Maintenance note: E.g.: Maintenance interval has been exceeded, temporary disruption of a sensor signal	<ul style="list-style-type: none"> <li>Use the Cockpit app to check the next maintenance date for the prosthesis (see page 88). If the date is within the next month, make a maintenance appointment with the O&amp;P professional. For this appointment, the O&amp;P professional has to be provided with the battery charger and power supply in addition to the prosthesis and tube adapter.</li> <li>Conduct the self-test again by connecting/disconnecting the battery charger.</li> <li>If the beep is emitted again and the maintenance date has not been reached or exceeded, you should visit the O&amp;P professional soon. If necessary, the O&amp;P professional will forward the prosthesis to an authorised Ottobock Service Center.</li> <li>The product can be used without restrictions. However, vibration signals may not be generated.</li> </ul>

### Battery charge level

During the charging process, the current battery charge level is indicated by the number of LEDs lit on the side of the charger.

LEDs	0	1	2	3	4	5
Battery charge level	0%-10%	10%-30%	30%-50%	50%-70%	70%-90%	>90%

## 16.3 Directives and manufacturer's declaration

### 16.3.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (see page 67).

**Table 1 – Electromagnetic emissions for all devices and systems**

#### Directives and manufacturer's declaration – electromagnetic interference

The product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' must ensure that the device is operated in such an environment.

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1	The product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' uses HF energy exclusively for its internal functioning. Therefore its HF emissions are very low and interference with neighbouring electronic devices is unlikely.
HF emissions according to CISPR 11	Class B	The product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' is intended for use in other than residential facilities and in facilities connected directly to a public low-voltage network that also supplies residential buildings.
Harmonics according to IEC 61000-3-2	N/A	
Voltage fluctuations/flicker according to IEC 61000-3-3	N/A	

**Table 2 – Electromagnetic interference immunity for all devices and systems**

**Directives and manufacturer's declaration – electromagnetic interference immunity**  
The product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' must ensure that the device is operated in such an environment.

Interference immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 6 kV discharge on contact ± 8 kV discharge through air	± 6 kV discharge on contact ± 8 kV discharge through air	Floors should be wood, concrete or ceramic tile. If the floor covering is made of synthetic material, the relative humidity must be at least 30 %.
Electrical fast transient/bursts according to IEC 61000-4-4	± 2 kV for power lines ±1 kV for input and output lines	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply.
Surges according to IEC 61000-4-5	±1 kV phase-to-phase voltage ± 2 kV even voltage	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply.

Interference immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Voltage drops, short interruptions and fluctuations of the supply voltage according to IEC 61000-4-11	$< 5\% U_T$ ( $> 95\%$ drop of $U_T$ ) for 1/2 period $40\% U_T$ ( $60\%$ drop of $U_T$ ) for 5 periods $70\% U_T$ ( $30\%$ drop of $U_T$ ) for 25 periods $< 5\% U_T$ ( $> 95\%$ drop of $U_T$ ) for 5 s	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply. If the user of the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' demands continuous operation even in case of an interruption in the supply voltage, we recommend supplying the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' with an uninterruptible power supply or battery.
Magnetic field at mains frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should be equivalent to the typical levels for commercial or hospital environments.

Note:  $U_T$  is the rated voltage before applying the test levels.

**Table 4 – Electromagnetic interference immunity for non-life-sustaining devices and systems**

**Directives and manufacturer's declaration – electromagnetic interference immunity**  
 The product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' must ensure that the device is operated in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Conducted HF interference according to IEC 61000-4-6	3 V <sub>eff</sub> 150 kHz to 80 MHz	Not applicable	Portable and mobile wireless equipment should not be used at a lesser distance from the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' including the lines than the recommended safety distance calculated using the equation applicable for the transmission frequency. <b>Recommended safety distance:</b> $d = 0.3 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ for 800 MHz to 2.5 GHz where P is the nominal output of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended safety distance in metres (m). The field strength of stationary radio transmitters at all frequencies according to an on-site investigation <sup>a</sup> should be less than the compliance level. Interference is possible in the vicinity of devices that bear the following symbol. 
Radiated HF interference according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	12 V/m 80 MHz to 1 GHz 3 V/m 1 GHz to 2.5 GHz	
<b>Note 1:</b> The higher frequency range applies at 80 MHz and 800 MHz. <b>Note 2:</b> These directives may not be applicable in all cases. The propagation of electromagnetic factors is influenced by absorption and reflection by buildings, objects and people.			

<sup>a</sup> The field strength of stationary transmitters such as base stations of radio telephones and mobile land radio equipment, amateur radio stations, AM and FM radio and television stations cannot be precisely determined theoretically in advance. A site survey is recommended to establish the electromagnetic conditions as a result of stationary HF transmitters. If the field strength determined at the site of the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' exceeds the compliance level specified above, the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' has to be monitored in regards to normal operation at each place of use. If unusual performance characteristics are noted, taking additional steps may be necessary, for example changing the orientation or moving the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2'

**Table 6 – Recommended safe distances between portable and mobile HF telecommunication equipment and the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' – for non-life-sustaining devices and systems**

**Directives and manufacturer's declaration – Recommended safe distances between portable and mobile HF telecommunication equipment and the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2'**

The product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' is designed for operation in an electromagnetic environment where HF interference is monitored. The customer or user of the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2' can help prevent electromagnetic interference by complying with the minimum distances between wearable and mobile HF communication equipment (transmitters) and the product 'Genium X3 - Bionic Prosthetic System 3B5-X3 / Genium X3 3B5-2', as below according to the maximum output of the communication equipment.

Nominal output of transmitter [W]	Safety distance [m] by transmission frequency		
	150 kHz to 80 MHz	80 MHz to 800 MHz $d = 0.3 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	N/A	0.03	0.23
0.1	N/A	0.09	0.73
1	N/A	0.3	2.3
10	N/A	0.95	7.3
100	N/A	3	23

For transmitters with no maximum nominal output specified in the table above, the recommended safety distance  $d$  in metres (m) can be determined using the equation in the respective column, where  $P$  stands for the maximum nominal output of the transmitter in watts (W) according to the transmitter manufacturer's information.

**Note 1:** The higher frequency range applies at 80 MHz and 800 MHz.

**Note 2:** These directives may not be applicable in all cases. The propagation of electromagnetic factors is influenced by absorption and reflection by buildings, objects and people.

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- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

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### **Caution: Exposure to Radio Frequency Radiation.**

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- (1) this device may not cause interference, and

(2) this device must accept any interference, including interference that may cause undesired operation of this device.

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